Spirosart® Wickelmodul

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Sartorius Corporation

131 Heartland Boulevard,
Edgewood, New York 11717, USA
Phone +1.631.254.4274 | Fax +1.631.254.4253
Toll-Free +1.800.368.7178

Specifications subject to change without notice.
Printed and copyrighted by Sartorius · W102.06 · G
Publication No.: SGC6001-p03084

Spirosart® Spiral-Wound Module

Directions for Use of the Sartorius Spiral-Wound Module

Order number code:

310 15358 50 -- A
 ↓ ↓ ↓ ↓ ↓
 Model code Membrane code Length code Type code

Notes for Your Safety

Be careful when you cut open the aluminum bag to avoid damaging the module. The bag contains alcohol: Do not smoke, do not expose to an open flame and avoid eye contact. For further information, please refer to the safety data sheet on our home page at www.sartorius.com

1. Unpacking the Module

The box contains an aluminum bag in which the module and baffle are packaged. Remove the aluminum bag containing the module from the outer box. Place the aluminum bag containing the module in an upright position and open the bag carefully along the heat-sealed seam. Remove the module from the bag and place it in an upright position on a protective cap.

2. Installing the Module in the System

Important Note


Before installing the modules in the system, check the O-rings on the anti-telescoping devices (ATDs) and on the stopper and replace them, if necessary. Use only current ATD versions. Insert module in the direction of the arrow (see label) into pipe. Do not drop module onto the ATDs!

Remove the red protective caps and insert the modules into the module pipes in the following order:

Wet the O-rings on the stopper and ATDs and insert into the module core.

When installing the modules, be sure to position them as shown by the arrows on the outside.

Direction of the arrow(s) on the module = Feed direction

3. Getting Started

Integrity testing: The modules can be integrity tested individually and "on their own," as well as after they have been installed in the module pipe. Please contact your Sartorius Service Center or representative for assistance. The system must be cleaned every time it is completely or partially re-equipped. (Please refer to the "System" section in the Operating Instructions).

4. Technical Specifications

Materials:	PESU, PP, silicone
Dimensions:	500 × 148 mm
Filter area:	5 m ²
Pressure values:	
Filtrate pressure P_f :	$P_f < P_{outlet}$
Inlet pressure:	max. 4 bar (~ 58 psi) per pipe with 2 modules
Transmembrane pressure:	max. 2 bar (29 psi)
Temperature:	
Operating temperature:	0–40 °C
Cleaning temperature:	50 °C min., 65 °C max.
pH:	2–14

sartorius

5. Filtration and Chemical Cleaning

Please follow these instructions to achieve good and uniform flow rates during filtration:

- Make sure that you always maintain crossflow conditions
- Some fining agents may reduce filtration performance (bentonite, activated carbon, caseins)
- Optimize your cleaning cycles and concentrate drainage
- Maintain the prescribed cleaning temperatures

Please refer to the Operating Instructions for more details.

Recommended values

Filtration:	Differential pressure ($P_{inlet} - P_{outlet}$): 2 bar min. (29 psi)
Concentrate drainage:	At approx. 50–70% of the initial flow rate
Intermediate cleaning:	At least after 4–8 hours' operation; 50 °C min., 65 °C max.
Chemical cleaning:	At least after 40 hours' operation 50 °C min., 65 °C max.
Recommended cleaning agents:	Alkaline: Divos 123 P3-Ultrasil 130 or 132 (country-specific) caustic Acidic: Divos 35 Citric acid

For more information and advice, please contact your local field sales engineer.

6. Storage and Preservation

If the system is not operated for an extended period (2 days or longer), the modules must be preserved. It is best to leave them in the system for this purpose. Adequate preservation is achieved by filling the **cleaned** system with NaOH at a concentration of 4 %. Then visually inspect the modules once a month.

The modules must be stored in a frost-free place.

7. Disposal (status: April 2003)

If the used Spirosart contains non-hazardous residual impurities such as those retained during wine filtration, you can use waste code EAK 150203 (European waste code according to the European Waste Catalogue "Europäischer Abfall-Katalog") for adsorbents, wiping cloths, filter materials and protective clothing. If the spiral-wound modules are contaminated with hazardous waste, you should use EAK 150202* for adsorbents, wiping cloths and filter materials contaminated by dangerous substances. If you are unsure about which code to use, please consult your waste management authorities.

* Hazardous waste that requires particular monitoring in line with the recycling management and waste legislation.

8. Service

If you have any further questions, please contact our Service Center:

Sartorius Service Center
Sartorius AG
Weender Landstrasse 94–108
37075 Goettingen, Germany
Phone +49.551.308.3381
Fax +49.551.308.3289
www.sartorius.com

Important Note

Before sending in used modules to be tested, do not forget to clean, preserve (see item 6.) and seal them in **leak-proof** packaging. Otherwise, it will not be possible for us to process them properly. Please fill out and sign the Sartorius Return Shipment form and include it the package of modules you wish to return.

Sartorius AG

Weender Landstrasse 94–108
37075 Goettingen, Germany
Phone +49.551.308.0 | Fax +49.551.308.3289
www.sartorius.com

Sartorius Corporation

131 Heartland Boulevard,
Edgewood, New York 11717, USA
Phone +1.631.254.4274 | Fax +1.631.254.4253
Toll-Free +1.800.368.7178

Specifications subject to change without notice.
Printed and copyrighted by Sartorius · W102.06 · G
Publication No.: SGC6001-p03084

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Sartorius S.p.A.