Operation Manual

Important: Do not operate the John Bunn DigiOx Finger Pulse Oximeter without first reading and understanding this manual! Save this manual for future use.

Info: The most current version of this manual can be found online at www.grahamfield.com.

JB02008-INS-LAB-RevB14


CONTENTS

INTRODUCTION	3
INTENDED USE	3
CONTRAINDICATIONS	3
IMPORTANT SAFETY GUIDELINES - PLEASE READ BEFORE USE	3
OVERVIEW	7
APPEARANCE	7
FEATURES	7
NAME AND MODEL	8
BATTERY INSTALLATION	8
OPERATION INSTRUCTIONS	8
ACCESSORIES (INCLUDED)	11
REPAIR, MAINTENANCE, AND DISPOSAL	11
MAINTENANCE	11
CLEANING AND DISINFECTION	12
DISPOSAL	13
TROUBLESHOOTING	13
SYMBOL KEY	15
TECHNICAL SPECIFICATIONS	16
LIMITED WARRANTY	18
MANUAL DE OPERACIÓN (ESPAÑOL)	21

2

OVERVIEW

Appearance

The DigiOx Finger Pulse Oximeter is shown above.

Features

The John Bunn DigiOx Finger Pulse Oximeter provides a simple way to spot-check patients by combining the sensor and monitor into one integrated, compact, easy-to-use device with single-button operation. The oximeter can accurately measure oxygen saturation and pulse rate. When a finger is inserted into the sensor's rubber cushion, and the Power/Display Button is pressed to turn the oximeter on, then the SpO₂ value and pulse rate value display on the screen. The oximeter features a bright-color LED screen with adjustable brightness and multi-directional displays that enable the user to read results from any angle. The oximeter features a low voltage

INTRODUCTION

Please read and understand all instructions before using the John Bunn DigiOx Finger Pulse Oximeter, JB02008.

Intended use

The John Bunn DigiOx Finger Pulse Oximeter is used to measure oxygen saturation (SpO₂) and pulse rate in home or clinical settings. It is designed for adult and pediatric patients.

Contraindications

⚠ WARNING: The oximeter is not intended for continuous patient monitoring. DO NOT use the oximeter while patient is undergoing MRI or CT scanning.

Important Safety Guidelines - Please Read Before Use

The safety statements presented in this chapter refer to the basic safety information that the operator of the John Bunn DigiOx Finger Pulse Oximeter shall pay attention to and abide by. There are additional safety statements in other chapters or sections, which may be the same as or similar to the following, or specific to the operations.

⚠ WARNING: Indicates a potential hazard situation or unsafe practice that, if not avoided, could result in death or serious injury.

⚠ CAUTION: Indicates a potential hazard or unsafe practice that, if not avoided, could result in minor personal injury.

⚠ NOTICE: Indicates a potential hazard or unsafe practice that, if not avoided, could result in product or property damage.

Info: Provides application recommendations or other useful information to ensure that you get the most from your product.

⚠ WARNINGS:

Important! Read and understand this manual before using the oximeter. Do not use this device without proper instruction from a Healthcare Professional. If the oximeter is not properly operated, personal injury and/or damage to the oximeter could result. If components are damaged or missing, contact your dealer immediately. DO NOT use substitute parts. Inspect the device before use to ensure that there is no visible damage that may affect patient or user's safety or measurement performance. Thereafter, inspect the device at least once a week.

GF Health Products, Inc. assumes no responsibility for any damage or injury caused by improper installation or use of this product.

DO NOT open the oximeter housing. It contains no serviceable parts. Service must be performed by trained, authorized personnel only. Otherwise, device failure and health hazard may occur.

EXPLOSION HAZARD: Do not use the oximeter in the presence of flammable anesthetics, explosive substances, vapors or liquids.

Use only parts and accessories specified in this manual.

The oximeter is designed for real-time and rapid measurement. It is not suitable for long-time continuous patient monitoring. DO NOT apply the sensor to the same finger for more than two hours. If any readings are found to be abnormal, please change the position of the sensor.

The oximeter is intended only as an adjunct in patient assessment. It must be used in conjunction with clinical signs and symptoms. It is intended as a measuring, not a treatment, device.

Magnetic and electrical fields are capable of interfering with the proper performance of the oximeter. For this reason, ensure that all external devices operated in the vicinity of the oximeter comply with the relevant EMC requirements. Mobile phones, X-ray equipment, and/or MRI devices are a possible source of interference, as they may emit high levels of electromagnetic radiation. This device complies with IEC60601-1-2.

DO NOT clip this device on edema or tender tissue. The infrared light emitted from the device (it is invisible) is harmful to the eyes. DO NOT stare into the light.

Keep this device out of the reach of children.


Notice for California Customers- California Proposition 65 WARNING: This product contains a chemical known to the State of California to cause cancer and reproductive or developmental harm.

5

a. Press the Power/Display Button once more; the screen will rotate 180° and display as shown in Figure 6.


b. Press the Power/Display Button once more; the screen will rotate 90° and display as shown in Figure 7.


c. Press the Power/Display Button once more; the screen will rotate 180° and display as shown in Figure 8.

d. Press the Power/Display Button once more; the screen will rotate -90° and also show pulse as a filled wave-form bar graph along bottom of screen, as shown in Figure 9.


e. Press the Power/Display Button once more; the screen will now show pulse as a hollow wave-form bar graph along bottom of screen, as shown in Figure 10.

f. Press the Power/Display Button once more; the screen will return to the original screen as shown in Figure 5 on previous page.

Ensure that no obstruction (contamination, scar tissue, fingernail polish, acrylic fingernail, etc.) exists at the site where the sensor is placed; such obstruction could affect the signal received by the sensor and result in an incorrect measurement.

Vigorous exercise and electrosurgical device interference may affect the measuring accuracy.

1. Open the clip as shown in Figure 4.
2. Insert finger, nail up, into the clip's rubber cushion (ensure the finger is correctly positioned), and then gently clip the finger.
3. Press the Power/Display Button to turn the oximeter on.
4. The screen will then display the data as shown in Figure 5. After the values become stable, the user can read the information on the display screen.
5. When the screen displays as shown in Figure 5, pressing the Power/Display Button will change the screen's orientation:

9

Troubleshooting continued

Symbol	Description
	Type BF
	Warning: See operation manual
	The pulse oxygen saturation
	Pulse rate (beats per minute)
	Low battery voltage
	Serial number
	Not for continuous monitoring

PURPOSE. IF ANY MODEL OR SAMPLE WAS SHOWN TO THE CUSTOMER, SUCH MODEL OR SAMPLE WAS USED MERELY TO ILLUSTRATE THE GENERAL TYPE AND QUALITY OF THE PRODUCT AND NOT TO REPRESENT THAT THE PRODUCT WOULD NECESSARILY CONFORM TO THE MODEL OR SAMPLE IN ALL RESPECTS.

GF'S TOTAL LIABILITY FOR ANY PRODUCT OR SERVICE PROVIDED IS LIMITED TO THE COST OF THE PRODUCT GIVING RISE TO THE CLAIM. IN NO EVENT WHETHER IN CONTRACT, INDEMNITY, WARRANTY, NEGLIGENCE, STRICT LIABILITY OR OTHERWISE WILL GF BE LIABLE FOR ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY OR CONSEQUENTIAL DAMAGES, INCLUDING, BUT NOT LIMITED TO: DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF USE, DOWNTIME, COVER, OR EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS AND BENEFITS.

This warranty gives you specific legal rights. You may have additional rights which vary from state to state.

Manufactured for:
GF Health Products, Inc.
2935 Northeast Parkway
Atlanta, Georgia 30360
© 2009 GF Health Products, Inc.
Made in China
tel: 770-368-4700
fax: 770-368-2386

www.grahamfield.com

GF Health Products, Inc. is not responsible for typographical errors. All illustrations, specifications, packaging and warranties contained in this literature are based on the latest product information available at the time of printing. The most current product information, including the most current version of these instructions, can be found online at www.grahamfield.com. Graham-Field and John Bunn are registered trademarks of GF Health Products, Inc.

© 2009 GF Health Products, Inc., Atlanta GA 30360, 770-368-4700

6

Change the batteries when the low-voltage indicator displays.

Clean the surface of the device before using: Wipe the device with a clean, soft cloth saturated with a solution of 70% isopropyl alcohol, then let it air dry or wipe it dry. Repeat as necessary.

Take out the batteries if the oximeter will not be used for a long period of time.

Environment (NOTICE: Humidity may shorten the life of the device, or even damage it)	Operation: 41°F - 104°F (5°C - 40°C), 15% - 80% RH Storage: -4°F - 158°F (-20°C - 70°C), < 93% RH
--	--

Cleaning and Disinfection

1. Do not immerse the device in liquid. Clean the device only by using a surface-clean technique: Wipe with a clean, soft cloth saturated with a solution of 70% isopropyl alcohol.

2. Wipe with a clean, soft cloth saturated with clean water and air dry, or dry with a clean, soft cloth.

⚠ WARNING: Do not use the sensor if it is damaged.

⚠ NOTICE: Do not sterilize this device. Do not subject this device to irradiation, steam, ethylene oxide, high temperature, or immersion in liquid.

TECHNICAL SPECIFICATIONS

Display screen	Organic LED Display
Power supply requirement	Two 1.5V AAA alkaline batteries
SpO ₂ Parameter Specifications	Measuring range: 0%-100% Accuracy: 70%-100%, ±3% < 70%, no definition
Pulse Rate Parameter Specifications	Measuring range: 30-235 bpm Accuracy: 30-99bpm, ±2 bpm; 100-235 bpm, ±2%
Dimensions (L x W x H)	2.28" x 1.26" x 1.34" (58 mm x 32 mm x 34 mm)
Environment	Operation: 41°F - 104°F (5°C - 40°C), 15% - 80% RH Storage: -4°F - 158°F (-20°C - 70°C), < 93% RH
Net weight with batteries	2 oz. (50g)

12

TECHNICAL SPECIFICATIONS continued

Classification	The type of protection against electric shock: Internally powered equipment
	The degree of protection against electric shock: Type BF applied part
	The degree of protection against harmful ingress of liquids: Ordinary equipment without protection against ingress of water

Electro-Magnetic Compatibility: Group I, Class B

13

LIMITED WARRANTY

SCOPE OF WARRANTY

GF Health Products, Inc. ("GF") warrants to the original purchaser ("Customer") only, that it will replace or repair components, at GF's sole discretion, that are defective in material or workmanship under normal use for a period of one year after the purchase date unless there is an expiration date on the component in which case the warranty shall expire on the earlier of the warranty period or the expiration date. The warranty does not extend to non-durable parts and does not include labor or costs of shipping. This limited warranty is not transferable. All warranties are conditioned upon the proper use of the product strictly in accordance with good commercial practice and applicable GF instructions and manuals, including proper use and maintenance. The warranty is void if the defect is caused by any other reason related to defects in materials or workmanship.

OBTAINING WARRANTY SERVICE

GF's customer service team must be notified of any warranty claim within the applicable warranty period. Call 678-291-3207, or fax 770-368-2386 or e-mail cs@grahamfield.com. Failure to follow the specific directions provided by the GF customer service team will result in denial of the warranty claim.

ENTIRE WARRANTY, EXCLUSIVE REMEDY AND CONSEQUENTIAL DAMAGES DISCLAIMER

THIS WARRANTY IS GF'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. GF MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR

14

11

DISPOSAL

Do not dispose of the device in regular trash.

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