QUÍMICA EXPERIMENTAL III

PROGRAMA DE QUÍMICA INORGÁNICA

ÍNDICE

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
Bibliografía: Nakamoto, K., *Infrared and Raman Spectra of Inorganic and Coordination Compounds, Part A, Theory and Applications in Inorganic Chemistry*, Wiley, 1997. **Ver páginas siguientes.** (Extraídas del texto de K. Nakamoto, *Infrared and raman spectra of inorganic and coordination compounds*, 4th edition, John Wiley and Sons: New York, 1986).

The Nitrate Ion

This ion has several structural roles,¹³⁴ as follows:

The most common forms of the nitrate ion are symmetrical bidentate followed by η^1 . The free NO_3^- ion has relatively high symmetry (D_{3h}); thus its ir spectrum is fairly simple. The totally symmetric N—O stretching mode is not ir active, but the doubly degenerate N—O stretching mode gives rise to a strong band at $\sim 1390 \text{ cm}^{-1}$. There are also two ir-active deformation modes, one of which is doubly degenerate, at 830 and 720 cm^{-1} . When NO_3^- is coordinated, its effective symmetry is reduced, causing the degeneracies to split and all modes (six) to be ir active. Hence, it is possible to distinguish between ionic and coordinated NO_3^- groups.

Because the two commonest forms of coordinated NO_3^- have the same effective symmetry, hence the same number of ir-active vibrational modes, criteria for distinguishing between them must be based on the positions of the bands rather than their number. In practice, the situation is quite complex and there are no entirely straightforward criteria. This is because the array of frequencies depends on both the geometry and strength of coordination.

¹³⁴C. Orvig *et al.*, *Inorg. Chem.* **1995**, *34*, 4921; P. N. V. P. Kumar and D. S. Marynick, *Inorg. Chem.* **1993**, *32*, 1857; R. Han and G. Parkin, *J. Am. Chem. Soc.* **1991**, *113*, 9707; O. Yamauchi *et al.*, *Inorg. Chem.* **1996**, *35*, 7148.

III-6. THIOCYANATO AND OTHER PSEUDOHALOGENO COMPLEXES

The CN^- , OCN^- , SCN^- , SeCN^- , CNO^- , and N_3^- ions are called "pseudo-halide ions," since they resemble halide ions in their chemical properties. These ions may coordinate to a metal through either one of the end atoms. As a result, the following linkage isomers are possible:

M—CN, cyano complex	M—NC, isocyano complex
M—OCN, cyanato complex	M—NCO, isocyanato complex
M—SCN, thiocyanato complex	M—NCS, isothiocyanato complex
M—SeCN, selenocyanato complex	M—NCSe, isoselenocyanato complex
M—CNO, fulminato complex	M—ONC, isofulminato complex



Two compounds are called true linkage isomers if they have exactly the same composition and different linkages mentioned above. A well-known example is nitro (and nitrito) pentammine Co(III) chloride discussed in Sec. III-2. A pair of true linkage isomers is difficult to obtain since, in general, one form is much more stable than the other. As will be shown later, several pairs of new linkage isomers have been isolated, and infrared spectroscopy has proved to be very useful in distinguishing them.

(1) Thiocyanato (SCN) Complexes

The SCN group may coordinate to a metal through the nitrogen or the sulfur or both (M—NCS—M'). In general, Class A metals (first transition series such as Cr, Mn, Fe, Co, Ni, Cu, and Zn) form the M—N bonds, whereas Class B metals (second half of the second and third transition series such as Rh, Pd, Ag, Cd, Ir, Pt, Au, and Hg) form the M—S Bonds.¹⁹¹ However, other factors such as the oxidation state of the metal, the nature of other ligands in a complex, and steric consideration also influence the mode of coordination. Mitchell and Williams¹⁹² have shown that the CN stretching frequencies are generally lower in the M—NCS complexes than in the M—SCN complexes. The C—S stretching frequency is more useful in distinguishing these two isomers: 780–860 cm^{-1} for the M—NCS and 690–720 cm^{-1} for the M—SCN group.^{193–195} The NCS bending frequency is also different between two isomers; 450–490 cm^{-1} for the M—NCS and 400–440 cm^{-1} (often accompanied by weaker bands at higher frequency side) for the M—SCN group.^{194,195} Table III-24 lists some results obtained by Sabatini and Bertini.¹⁹⁵ The M—NCS group is linear or bent, while the M—SCN group is always bent. This seems to suggest that the following resonance structures are predominant in each case:

M—NCS Group

M—SCN Group

Pecile¹⁹⁶ has noted that the integrated intensity of the CN stretching band can also be used to distinguish two linkage isomers.

Clark and Williams¹⁹⁷ have obtained the infrared spectra of tetrahedral $[\text{M}(\text{NCS})_2\text{L}_2]$ and octahedral $[\text{M}(\text{NCS})_2\text{L}_4]$ -type complexes (M = Fe, Co, Ni, etc.; L = pyridine, etc.) and discussed stereochemistry from infrared spectra. A similar study has been made by Nelson and Shepherd.¹⁹⁸ Bennett et al.¹⁹⁹ have obtained the infrared spectra of $[\text{Cr}(\text{NCS})_4\text{L}_2]^{-}$, where L is an amine or phosphine, and have assigned the Cr—N and Cr—L stretching bands in addition to the NCS group vibrations. In $[\text{M}(\text{NCS})_2\text{L}_2]$ -type complexes in which M is Pt(II) or Pd(II), the NCS group is S-bonded if L is an amine and N-bonded if L is a phosphine.¹⁹³ This is because a strong π -acceptor such as a tertiary phosphine makes the d -orbitals of the metal less available for bonding with the π -orbitals of sulfur. Jennings and Wojcicki²⁰⁰ obtained the infrared spectra of square-planar $[\text{Rh}(\text{CO})(\text{NCS})\text{L}_2]$, $[\text{Rh}(\text{NCS})\text{L}_3]$, $[\text{Rh}(\text{CO})_2(\text{NCS})_2]^{-}$, and $[\text{Rh}_2(\text{NCS})_2\text{L}_4]$ -type complexes in

TABLE III-24. VIBRATIONAL FREQUENCIES OF ISOTHIOCYANATO AND THIOCYANATO COMPLEXES (cm^{-1})¹⁹⁵


Compound ^a	$\nu(\text{CN})$	$\nu(\text{CS})^b$	$\delta(\text{NCS})$
K[NCS]	2053	748	486, 471
(Et ₄ N) ₃ [Cr(—NCS) ₆]	2078	—	483
(Et ₄ N) ₂ [Co(—NCS) ₄]	2065	844 } 838 }	481
(Et ₄ N) ₄ [Ni(—NCS) ₆]	2112 } 2103 }	828	470
(Et ₄ N) ₂ [Zn(—NCS) ₄]	2072	837	482
(Et ₄ N) ₂ [Pd(—SCN) ₄]	2112 } 2109 }	698 } 694 }	465, 433 } 429, 418 }
(Et ₄ N)[Au(—SCN) ₄]	2127	695	454, 415
(Et ₄ N) ₂ [Pt(—SCN) ₆]	2120	692	461, 457 } 418, 415 }

^aEt, C₂H₅.


^bThe S-bonded complex may exhibit a band at 880–810 cm^{-1} which is the first overtone of $\delta(\text{NCS})$. Care must be taken, therefore, in distinguishing this band from $\nu(\text{CS})$.²⁰⁵

which L is a phosphine, an arsine, a stibine, or a phosphite. Infrared spectra show that all these Rh(I) complexes are N-bonded. The *cis* and *trans* isomers of Co(III)-NCS complexes can be distinguished from infrared spectra.^{201,202}

As stated before, the distinction of true linkage isomers by infrared spectroscopy has been made for $[\text{Co}(\text{NH}_3)_5(-\text{NO}_2 \text{ or } -\text{ONO})]\text{Cl}_2$ ¹⁰ and $[\text{KFeCr}(\text{CN})_6]$ ¹⁸⁴. Basolo et al.²⁰³ succeeded in isolating true linkage isomers of


The bright-yellow, N-bonded complex absorbs at 2089 (CN stretching) and 854 cm^{-1} (CS stretching); the yellow-orange, S-bonded complex absorbs at 2119 cm^{-1} (CN stretching), although the phenyl absorption near 700 cm^{-1} obscures its C—S stretching band. Infrared spectra of linkage isomers have also been reported for $[\text{Pd}(\text{bipy})(\text{NCS}/\text{SCN})_2]$ ^{204,205} and $[\text{Pd}(\text{As}(n\text{-C}_4\text{H}_9)_3)_2(\text{NCS}/\text{SCN})_2]$.²⁰⁵ A complex of the composition $[\text{Cu}(\text{tren})(\text{NCS})_2]$ (tren, triaminotriethylamine) exhibits two CN stretching (2094 and 2060 cm^{-1}) and two CS stretching (818 and 745 cm^{-1}) bands. Based on this and other data, Raymond and Basolo²⁰⁶ suggested that this complex may contain both Cu—NCS and Cu—SCN linkages. Later, Jain and Lingafelter²⁰⁷ have shown from x-ray analysis that one NCS group is bonded to the metal through nitrogen and the other group is not coordinated with it. Bertini and Sabatini²⁰⁸ suggest that

may contain both types of coordination, since it exhibits two CN stretchings at 2120 and 2090 cm^{-1} .

The NCS group also forms a bridge between two metal atoms. The CN stretching frequency of a bridging group is generally higher than that of a terminal group. For example, $\text{HgCo}(\text{NCS})_4(\text{Co}-\text{NCS}-\text{Hg})$ absorbs at 2137 cm^{-1} , whereas $(\text{Et}_4\text{N})_2[\text{Co}(-\text{NCS})_4]$ absorbs at 2065 cm^{-1} . According to Chatt and Duncanson,²⁰⁹ the CN stretching frequencies of Pt(II) complexes are 2182–2150 cm^{-1} for the bridging and 2120–2100 cm^{-1} for the terminal NCS group. $[(\text{P}(n\text{-Pr})_3)_2\text{Pt}_2(\text{SCN})_2\text{Cl}_2]$ (compound I) exhibits one bridging CN stretching, whereas $[(\text{P}(n\text{-Pr})_3)_2\text{Pt}_2(\text{SCN})_4]$ (compound II) exhibits both bridging and terminal CN stretching bands. Thus the infrared spectra suggest the structure of each compound to be:

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3.A. Síntesis de $\text{Co}(\text{acac})_3$.

Add approximately 1.0 g of powdered $\text{Co}(\text{II})$ carbonate into a 50-mL Erlenmeyer flask and add 8 mL of acetylacetone. Drop in a small magnetic stir bar and cover with a small watch glass. Place the reaction vessel into a water bath on a stirring hot plate.

(Note: do not set the flask directly on the bottom of the water bath; instead, clamp it so that the flask is suspended within the water bath with its contents below the water line; refill the bath as needed as the water evaporates during heating).

Stir the reaction using the magnetic stirring bar and heat the reaction to about 80°C , maintaining it at this temperature for approximately 5 minutes. Slowly add 6 mL of 10% hydrogen peroxide (H_2O_2) drop-wise at a rate of 1-2 mL/min using a disposable pipette while continuing to stir the reaction mixture. **(CAUTION: 10% H_2O_2 will damage the surface of the skin! Wear gloves or be extra careful! The reaction will bubble vigorously.)** Reheat the mixture to 80°C after the addition of H_2O_2 is complete, then add an additional 6 mL of H_2O_2 in the same manner. After the second addition of H_2O_2 is complete, remove the water bath, place the flask directly on the hotplate surface and briefly bring the flask to boiling. Cool the mixture to room temperature and then place it in an ice bath for about 15 minutes or until crystallization appears to be complete. Filter the product using a Büchner funnel, washing remnants out of the flask with deionized water (2 x 5 mL). Wash with several small (2-3 mL) volumes of cold ethanol (make sure you chill it first or it will dissolve your compound!). Allow the compound to dry until the next lab period and then weigh the compound to get the yield of dry product.

3.B. Síntesis de $\text{Cr}(\text{acac})_3$.

Heat a water bath on a stirring hotplate. Add approximately 1.30 g of $\text{CrCl}_3 \cdot 6\text{H}_2\text{O}$ into a 50-mL Erlenmeyer flask. Drop in a small magnetic stirring bar, add 20.0 mL of distilled water and stir briefly. When the chromium compound has dissolved into solution, add 5.0 g of urea and 4.0 mL of acetylacetone. Clamp the flask in place in the water bath (see note in part A) and heat the water bath to a temperature greater than 90°C . Heat the mixture, uncovered and with vigorous stirring, and deep maroon crystals will begin to form. After one hour, allow the reaction to cool to room temperature, and then place it in an ice-bath to complete crystallization (Note: if you do not see solid after 1 h, continue to heat and stir until you do: sometimes it takes a bit longer!). Collect the product by vacuum filtration using a Büchner funnel. Wash out the flask and the crystals on the filter with plenty of deionized water (about 3 x 10 mL). Allow the product to dry until the next lab period and then weigh the dry final product to determine the yield of dry product.

3.C. Preparation of $\text{Mn}(\text{acac})_3$.

Add approximately 1.0 g of $\text{MnCl}_2 \cdot 4\text{H}_2\text{O}$ and 2.6 g of sodium acetate trihydrate into a 100-mL Erlenmeyer flask. Add a stirring bar, and dissolve the solids in a minimal amount of distilled water (add very small amounts at a time until the solids dissolve completely). When the solids are dissolved, add 4.0 mL of acetylacetone with continued stirring. Prepare another solution containing 200 mg of KMnO_4 in 10 mL of deionized water. Make sure you stir this for quite a while - KMnO_4 is difficult to dissolve and its color makes it hard to see if it has dissolved!- Add the KMnO_4 solution slowly and drop-wise, with stirring, to the above mixture. After addition is complete, a small portion of distilled water (1-2 mL) should be pipetted (using a disposable pipet) into the flask or beaker originally containing the KMnO_4 solution and then used to rinse the walls of the reaction flask to remove any permanganate that is sticking to the side. After stirring the mixture for about five minutes, add drop-wise a second 2.6 g portion of sodium acetate trihydrate dissolved in a minimal amount of water. Heat the resulting mixture directly on a hotplate for 10 minutes at about 70°C (put a thermometer directly in the reaction flask) and then set aside to cool to room temperature. Collect the dark

brown precipitate on a Büchner funnel. Use small portions of deionized water (2 x 5 mL) to wash any remaining crystals out of the flask. Allow the product to dry until the next lab period and then weigh the dry final product to determine the yield of dry product.

3.D. Preparación de Fe(acac)₃.

Add approximately 1.3 g of Fe(NO₃)₃•9H₂O into a 50-mL Erlenmeyer flask and dissolve in 7 mL of distilled water. Drop in a magnetic stir bar, stir the solution, and add a mixture of 1.0 mL of acetylacetone in 7 mL of methanol. Next, add a solution containing 1.3 g of sodium acetate trihydrate in 7 mL of distilled water, and briefly heat the mixture on a hotplate until some of the methanol has boiled away (you should see the volume of the reaction mixture go down). Cool to room temperature, and then place in an ice bath for about 15 minutes. Filter the red crystalline solid using a Büchner funnel, wash the flask and crystals, using 2 x 5 mL) of deionized water. Allow the product to dry until the next lab period and then weigh the dry final product to determine the yield of dry product.

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
(1) Complexes of Acetylacetonone

A number of β -diketones form metal chelate rings of type A:

Among them, acetylacetonone(acacH) is most common ($R_I = R_{III} = \text{CH}_3$ and $R_{II} = \text{H}$). According to x-ray analysis on $\text{Fe}(\text{acac})_3$,⁶⁰⁵ the chelate ring is planar and symmetrical (C_{2v} symmetry), and the two $\text{C}=\text{O}$ bonds are equivalent, as are the two $\text{C}=\text{C}$ bonds in the ring.

The infrared spectra of metal acetylacetonate complexes have been studied extensively. In 1956, Mecke and Funck⁶⁰⁶ assigned the infrared spectra of acetylacetonone (keto and enol forms) and its metal complexes on an empirical basis. In 1960, Nakamoto and Martell⁶⁰⁷ carried out an approximate normal coordinate analysis on $\text{Cu}(\text{acac})_2$. At the time of this earlier research, no spectral data were available below 400 cm^{-1} and no isotopically substituted acetylacetonate complexes were studied. Since then, a number of investigators^{608-610a} have extended their measurements to the far-infrared region. The spectra of metal complexes containing C^{13} - and O^{18} -labeled acetylacetonates have been obtained.^{611,612} A rigorous normal coordinate analysis considering all the atoms in the molecule has been made for $\text{M}(\text{acac})_2$ and $\text{M}(\text{acac})_3$ -type complexes.⁶¹³ Table III-49a lists the observed frequencies, band assignments, and M—O stretching force constants obtained by Mikami et al.⁶¹³ Figure III-33a shows the infrared spectra of six acetylacetonate complexes.

The nature of the 1577 and 1529 cm^{-1} bands of $\text{Cu}(\text{acac})_2$ has been a subject of controversy. Originally, Nakamoto and Martell⁶⁰⁷ assigned the former to a $\text{C}\cdots\text{C}$ stretching and the latter to a $\text{C}\cdots\text{O}$ stretching mode. Mikami et al.⁶¹³ also gave similar assignments although they showed that these two modes are coupled slightly with each other. However, the results of the C^{13} and O^{18} experiments^{611,612} suggest that the former is a pure $\text{C}\cdots\text{O}$ stretching and the latter is a pure $\text{C}\cdots\text{C}$ stretching. A further theoretical study is needed to solve this discrepancy.

The M—O stretching vibrations of acetylacetonate complexes are most important since they provide direct information about the strength of the M—O bonds. Thus far, they have been assigned from normal coordinate analyses^{607,613} without experimental evidence. As stated in Sec. III-9, the metal isotope technique³⁸⁹ provides a unique method to detect the metal-ligand vibrations. Nakamoto et al.^{614a}, therefore, applied this technique to acetylacetonate complexes of Fe(III), Cr(III), Pd(II), Cu(II), and Ni(II). Their results indicate that previous band assignments on the M—O stretching vibrations are essentially correct. As an example, Table III-49b gives the observed frequencies, isotopic shifts, and band assignments for $\text{Cr}(\text{acac})_3$. Figure III-33b shows an actual tracing of the infrared spectra of $\text{Cr}^{50}(\text{acac})_3$ and its Cr^{53} analog. Two bands at 463 and 358 cm^{-1} of the former must be assigned to the Cr—O stretching modes since they give large shifts relative to others. On the other hand, Pinchas et al.⁶¹¹ assigned the 592 cm^{-1} band to a pure Cr—O stretching since it gives the largest isotopic shift (19 cm^{-1}) by the $\text{O}^{16}\text{—O}^{18}$ substitution of $\text{Cr}(\text{acac})_3$. However, this band is shifted by only 0.7 cm^{-1} by the $\text{Cr}^{50}\text{—Cr}^{53}$ substitution, and cannot be assigned to a pure Cr—O stretching mode. It may be due to a ligand vibration in which the oxygen atom moves appreciably (e.g., $\text{C}\cdots\text{O}$ out-of-plane bending). It

TABLE III-49a. OBSERVED FREQUENCIES AND BAND ASSIGNMENTS OF ACETYLACETONATO COMPLEXES⁶¹³ (CM⁻¹)

Cu(acac) ₂	Pd(acac) ₂	Fe(acac) ₃	Predominant Mode
3072	3070	3062	$\nu(\text{CH})$
2987	2990	2895	$\nu(\text{CH}_3)$
2969	2965	2965	
2920	2920	2920	
1577	1569	1570	
1552	1549		
1529	1524	1525	$\nu(\text{C}\cdots\text{O}) + \nu(\text{C}\cdots\text{C})$
1461	(1425)	1445	$\delta(\text{CH}) + \nu(\text{C}\cdots\text{C})$
1413	1394	1425	$\delta_d(\text{CH}_3)$
1353	1358	1385	$\delta_s(\text{CH}_3)$
		1360	
1274	1272	1274	$\nu(\text{C}-\text{CH}_3) + \nu(\text{C}\cdots\text{C})$
1189	1199	1188	$\delta(\text{CH}) + \nu(\text{C}-\text{CH}_3)$
1019	1022	1022	$\rho_r(\text{CH}_3)$
936	937	930	$\nu(\text{C}\cdots\text{C}) + \nu(\text{C}\cdots\text{O})$
780	786	801	$\pi(\text{CH})$
	779	780	
		771	
684	700	671 ^a	$\nu(\text{C}-\text{CH}_3) + \text{ring def.}$ $+ \nu(\text{M}-\text{O})$
		667	
653	678	656	$\pi\left(\text{CH}_3-\text{C}\begin{array}{l} \diagup \text{C} \\ \diagdown \text{O} \end{array}\right)$
612	661	559 ^a	ring def. + $\nu(\text{M}-\text{O})$
		548	
451	463	433	$\nu(\text{M}-\text{O}) + \nu(\text{C}-\text{CH}_3)$
431	441	415	ring def.
		408	
291	294	298	$\nu(\text{M}-\text{O})$
1.45	1.85	1.30	K(M—O), mdyn/A (UBF)

^a Pure ring deformation.

should be noted that an isotopic substitution of the α -atom of the ligand (atom bonded directly to a metal) shifts the metal-ligand as well as ligand vibrations involving the motion of the α -atom. Thus, the metal-ligand vibrations cannot be uniquely assigned by this method.

Behnke and Nakamoto⁶¹⁵ prepared three deuterio analogs of $\text{K}[\text{Pt}(\text{acac})\text{Cl}_2]$, and calculated a set of UBF force constants to fit the observed frequencies of four isotopically substituted species. Table III-50 lists the observed

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You will prepare one member of the $Y_xGd_{3-x}Fe_5O_{12}$ ($0 \leq x \leq 3$) family, that is, you will be assigned a particular value of x . A description of the synthesis of $Y_{1.5}Gd_{1.5}Fe_5O_{12}$ via a mixed metal hydroxide precursor follows. Other compositions in the series can be prepared by varying the volumes of $Gd(NO_3)_3$ and $Y(NO_3)_3$ used, in such a way that the amount of $Gd(NO_3)_3$ plus the amount of $Y(NO_3)_3$ totals 6 mL.

Place 10 mL of 0.3 M $FeCl_3$ in a beaker. Add 3 mL of 0.3 M $Gd(NO_3)_3$ solution and 3 mL of 0.3 M $Y(NO_3)_3$ solution. Add 5-10 mL of 2 M NaOH dropwise to the metal ion solution in order to precipitate a reddish-brown solid. Decant the solution, and wash the remaining solid repeatedly with water until the wash is no longer basic (test with pH paper). Filter and dry the solid in a drying oven (120 °C) overnight.

Press the powder into 1/2-inch-diameter pellets using a standard pellet press, and fire in a furnace at 900 °C for 18-24 hours. A visible color change from reddish-brown to olive green after firing indicates that a reaction has taken place.

Determine the strength of the attraction of your pellet of garnet to a strong magnet at room temperature, dry-ice temperature, and liquid-nitrogen temperature. Discuss your procedure with other students in the class and agree on a method that will allow comparison of the magnetic effects as a function of the garnet composition.

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Rare Earth Iron Garnets

Their Synthesis and Magnetic Properties

Margret J. Geselbracht, Ann M. Cappellari, and Arthur B. Ellis¹
University of Wisconsin–Madison, Madison, WI 53706

Maria A. Rzeznik
University of California Berkeley, Berkeley, CA 94720

Brian J. Johnson
St. John's University, St. Joseph, MN 56374

The interaction of solids with magnetic fields is a topic that is discussed rarely in introductory chemistry courses. This occurs in spite of the proliferation of magnetic materials in today's society, including magnetic recording media and data storage. We believe that a presentation of the magnetic properties of solids fits well into the general chemistry course, because it provides a vehicle for discussing the existence of electrons, the electronic configurations of ions, and the role of solid state structure. Our experiment involving the preparation and study of a series of magnetic garnets with the composition $Y_xGd_{3-x}Fe_5O_{12}$ ($0 \leq x \leq 3$) that illustrates these topics follows. In addition, this experiment provides the opportunity to discuss the different types of magnetism in solids, to perform some chemistry using rare earth elements, and to learn about solid solutions and the tunability of properties by varying the composition of a solid solution. Furthermore, it can be performed in such a way that each student has a unique compound to prepare. The class data can then be compiled to obtain a picture of the magnetic properties of the entire range of composition of $Y_xGd_{3-x}Fe_5O_{12}$.

The majority of molecules and solids that we encounter are nonmagnetic. All of the electrons in these materials are paired, and the molecule or solid is said to be diamagnetic. The simplest type of magnetism, paramagnetism, often is demonstrated in introductory chemistry courses by observing the behavior of liquid oxygen in a magnetic field. Due to the presence of unpaired electrons in O_2 , the liquid is

attracted to a strong magnet, and if a large horseshoe magnet is used, the liquid oxygen appears to hover, trapped between the poles of the magnet (*1*).

Many solids exhibit interesting magnetic behavior as well, often more complex than simple paramagnetism, due to the cooperative effect of many electrons in the solid acting in concert. Ferromagnetism is a property of solids that is technologically important and exploited in the use of permanent magnets, magnetic recording media, and

Addendum

Description of the March 1994 Cover

The cover of the March 1994 issue needs further amplification to be useful for teachers. It shows a sample of zinc rod on the left-hand side of the page. The bottom portion of the rod was etched by human handling (this portion of the rod was exposed as the rod was passed around classes over a period of several semesters; the top portion of the rod had a protective cover) and reveals the boundaries between different grains. As discussed in the accompanying paper (Geselbracht, M. J.; Ellis, A. B.; Penn, R. L.; Lisensky, G. C.; and Stone, D. S. *J. Chem. Educ.* **1994**, *71*, 254) the bubble raft on the right-hand side of the cover (also shown schematically in Figs. 2 and 3 in the paper) provides a model for the grains, which are close-packed arrays of atoms having different relative orientations.

¹Author to whom correspondence should be addressed.


Figure 1. X-ray powder diffraction pattern of a sample of $Y_{1.5}Gd_{1.5}Fe_5O_{12}$ prepared via a mixed metal hydroxide precursor. Diffraction peaks that are marked with an asterisk correspond to the impurity $RFeO_3$ ($R = Y$ or Gd).

transformers. We believe that a discussion of the chemistry of magnetic solids in introductory chemistry courses not only illustrates the effect of unpaired electrons, but also introduces the concept that electrons in a solid can communicate with one another, leading to properties that are vastly different from those stemming from isolated, localized electrons.

What do we mean by the chemistry of magnetic solids? The magnetic properties of a solid depend on composition (the presence of magnetic ions), on the underlying crystal structure (the geometrical arrangement of the ions in three dimensions), and on temperature. In order to eliminate the influence of structure and to isolate the effects of composition and temperature on magnetic properties, ideally, we would like to change the composition of a magnetic solid without disrupting the crystal structure. In other words, we would like to study the magnetic properties of a solid solution.

In this paper, we will describe a general synthesis for compositions in the solid solution series $Y_xGd_{3-x}Fe_5O_{12}$ ($0 \leq x \leq 3$), and a simple demonstration that illustrates the differing magnetic properties of these materials. This experiment serves several purposes. Students gain experience in synthetic techniques commonly used to prepare many solid state oxides that are important in materials chemistry today. If X-ray powder diffraction is available, the X-ray diffraction patterns can be analyzed to discuss the structural characterization of solid solutions and illustrate the evolution of the lattice parameter with composition. Finally, these materials display some unique magnetic properties, that open the discussion to cooperative effects in solids and the applications of magnetic materials to technology.

Experimental

Synthesis

Compounds in the $Y_xGd_{3-x}Fe_5O_{12}$ ($0 \leq x \leq 3$) solid solution series were synthesized via a mixed metal hydroxide precursor (2). A 1 M stock solution containing yttrium nitrate was prepared by dissolving $Y(NO_3)_3 \cdot 5H_2O$ (Aldrich, 99.9%) in deionized water. Likewise, 1 M stock solutions were prepared from $Gd(NO_3)_3 \cdot 6H_2O$ (Aldrich, 99.9%) and $FeCl_3 \cdot 6H_2O$ (Alfa, ACS grade). A description of the synthesis of $Y_{1.5}Gd_{1.5}Fe_5O_{12}$ follows. Other compositions in the series can be prepared by varying the volumes of $Gd(NO_3)_3$ and $Y(NO_3)_3$ used.

A 3-mL aliquot of 1 M $Gd(NO_3)_3$ solution and a 3-mL aliquot of 1 M $Y(NO_3)_3$ solution were added to 10 mL of 1 M $FeCl_3$ in a beaker. Excess 6 M NaOH (5–10 mL) was then added dropwise to the metal ion solution in order to precipitate a reddish-brown solid. After precipitation was complete, the solution was decanted, and the remaining solid was washed repeatedly with water until the supernatant solution no longer tested basic using standard pH paper. Then the reddish-brown solid was filtered and dried in a drying oven (120 °C) overnight. The resultant powder was pressed into 1/2-in. diameter pellets using a standard pellet press and fired in a furnace at 900 °C for 18–24 h. A visible color change from reddish-brown to olive green after firing indicated that a reaction had taken place.

Preparing the material in the form of a pellet is necessary to carry out the magnetic tests described below. If a pellet press is not available, we have found that the loose powder can be packed into an empty gelatin capsule (available at a pharmacy counter) before carrying out the magnetic tests. Note that if a press is not available, the loose powder may require longer or repeated firings in order to complete the reaction.

X-Ray Powder Diffraction (optional)


Some instructors may wish to have their students verify the identity of the olive-green solid by X-ray powder diffraction. X-ray powder diffraction patterns of the samples were recorded from 10–100° 2θ on a Nicolet X-ray powder diffractometer. The powder diffraction pattern of the sample with the nominal composition $Y_{1.5}Gd_{1.5}Fe_5O_{12}$ is shown in Figure 1. The major product of the reaction was the garnet phase, $R_2Fe_5O_{12}$ ($R =$ rare earth element). The peaks indicated with an asterisk in Figure 1 are due to $RFeO_3$ ($R =$ rare earth element), which is a common impurity in samples of $R_2Fe_5O_{12}$. Grinding the olive-green solid, repelletizing, and firing a second time at 900 °C (or higher temperatures) for 24 h led to a noticeable decrease in the intensities of the peaks of the impurity phase in the diffraction pattern. It should be noted that the presence of a small amount of the $RFeO_3$ phase in the samples does not adversely affect the investigation of the magnetic properties in this experiment.

To obtain the lattice parameters for the samples, silicon was added as an internal standard prior to collecting the diffraction data, and the line positions were corrected for systematic errors by applying a linear correction, derived from the Si data, to the observed angles. Lattice parameters were refined from 15–18 reflections using standard least-squares refinement techniques. Using line positions that have not been corrected led to the same overall trend in lattice parameter with composition, although there was more scatter in the data (see Fig. 4).

Studies of Magnetic Properties

To demonstrate the magnetic properties of these materials we have used a qualitative test of the strength of the magnetism: the strength of the attraction of a pellet of garnet to a cow magnet (Edmund Scientific). Pellets that exhibit strong magnetism can be picked up by the cow magnet, while pellets that are only weakly magnetic cannot be picked up by the cow magnet. The temperature dependence of the magnetic properties can be demonstrated by studying the magnetic behavior of the pellets of the different compositions at room temperature, and after cooling the pellets in dry ice or liquid nitrogen. For more data points, students can prepare and use other low temperature baths (3).

More detailed investigations of the magnetism were carried out at the University of California, Berkeley with a Quantum Design SQUID magnetometer. The magnetiza-

Metal ion location within one octant of the unit cell:


Figure 2. Portion of the unit cell of the garnet structure of $\text{Gd}_3\text{Fe}_5\text{O}_{12}$ is shown at the top. All of the metal ions in one-eighth of the unit cell are drawn along with the oxygen ions that complete the octahedral or tetrahedral coordination around the iron atoms. The rest of the atoms have been deleted for clarity. At the bottom, horizontal sections perpendicular to the z direction are shown to indicate the precise positions of the iron and gadolinium ions within this octant of the unit cell. Note that in this figure, the $z = 0$ section corresponds to the top face of the cube and z increases in the downward direction.

tion was measured as a function of temperature in an applied magnetic field of 40 kGauss.

Suggested Use in Lab

Due to the need to fire these samples in a high temperature furnace for 18–24 h, we would suggest that the instructor allow two to three laboratory periods for this experiment, perhaps running concurrently with another experiment. One period is necessary for the preparation of the hydroxide precursor, and depending on the availability of furnaces, the firing can be done in the week between two meetings. In the second laboratory period, students can ascertain whether a second firing is necessary to complete the reaction. The magnetic studies then can be carried out in a second or third laboratory period. If the instructor includes X-ray powder diffraction in this experiment, the time must be adjusted accordingly.

The students can be given individual projects by asking them to synthesize a particular member of the $\text{Y}_x\text{Gd}_{3-x}\text{Fe}_5\text{O}_{12}$ family of compounds. This might be done by giving each student the formula (or x value) for their target compound and informing them that 1 M solutions of $\text{Y}(\text{NO}_3)_3$, $\text{Gd}(\text{NO}_3)_3$, and FeCl_3 are available in the lab and that everyone should use 10 mL of 1M FeCl_3 so that ap-



Figure 3. The three polyhedral coordination environments that are present in the garnet structure. In $\text{Gd}_3\text{Fe}_5\text{O}_{12}$, gadolinium ions are in the dodecahedral sites, and iron atoms are in both the octahedral and the tetrahedral sites.

proximately the same amount of material is prepared by each student.

After the synthesis, the students can perform the magnetic studies on their own compounds. Students may be asked to design their own tests of magnetic strength to get a more quantitative measurement of the relative magnetization of two different samples at a given temperature. This might be done by measuring the distance between the magnet and the sample when an interaction is first noticed. Finally, the results of the class can be tabulated to get an overall picture of the magnetic properties of these compounds.

Results and Discussion

The magnetic properties of rare earth iron garnets have been of great interest for the past two decades. Rare earth iron garnets, $\text{R}_3\text{Fe}_5\text{O}_{12}$, can be prepared with all of the rare earth ions except $\text{R} = \text{La}, \text{Ce}, \text{Pr}, \text{or Nd}$. Detailed reviews of the crystal chemistry and the physical properties of rare earth iron garnets abound in the literature (4).

Structure

Gadolinium iron garnet (GIG) and yttrium iron garnet (YIG) are members of the garnet structural family of complex oxides. The general formula for a garnet is $\{\text{C}_3\}[\text{A}_2](\text{D}_3)\text{O}_{12}$ where the C cations occupy dodecahedral sites, the A cations occupy octahedral sites, and the D cat-

Table 1. Variation in Cubic Lattice Parameter, a , with Composition for $\text{Y}_x\text{Gd}_{3-x}\text{Fe}_5\text{O}_{12}$

Composition	Lattice Parameter (\AA)
$\text{Y}_3\text{Fe}_5\text{O}_{12}$	12.370 (1)
$\text{Y}_{2.5}\text{Gd}_{0.5}\text{Fe}_5\text{O}_{12}$	12.382 (2)
$\text{Y}_2\text{Gd}_1\text{Fe}_5\text{O}_{12}$	12.402 (2)
$\text{Y}_{1.5}\text{Gd}_{1.5}\text{Fe}_5\text{O}_{12}$	12.423 (2)
$\text{Y}_1\text{Gd}_2\text{Fe}_5\text{O}_{12}$	12.437 (2)
$\text{Y}_{0.5}\text{Gd}_{2.5}\text{Fe}_5\text{O}_{12}$	12.450 (1)
$\text{Gd}_3\text{Fe}_5\text{O}_{12}$	12.468 (2)


Figure 4. The variation in the cubic lattice parameter, a , with composition in $Y_xGd_{3-x}Fe_5O_{12}$. The corrected data were obtained by adding silicon as an internal standard prior to collecting the diffraction data. The line positions of $Y_xGd_{3-x}Fe_5O_{12}$ were corrected for systematic errors by applying a linear correction, derived from the Si data, to the observed angles.

ions occupy tetrahedral sites in the crystal structure. The unit cell of the garnet structure has cubic symmetry and contains eight formula units for a total of 160 atoms (5). A portion of the unit cell of $Gd_3Fe_5O_{12}$ is drawn in Figure 2 with only the cations and selected oxygens shown for clarity. The three different polyhedral coordination environments for the cations in the garnet structure are drawn in Figure 3. In the case of $Gd_3Fe_5O_{12}$, Gd^{3+} occupies dodecahedral sites and Fe^{3+} occupies both octahedral and tetrahedral sites in the structure. However, a wide variety of cations in different valence states can reside in the cation sites. The primary consideration for site occupancy is ionic size. This leads to many compositions that form the garnet structure.

Because of the similarity in ionic radii of the rare earth ions, many rare earth iron garnets will form solid solution phases. In the case of yttrium (ionic radius = 0.900 Å) and gadolinium (ionic radius = 0.938 Å), the complete solid solution is stable:

Substitution of yttrium for gadolinium occurs on the dodecahedral site. Because Gd^{3+} is slightly larger than Y^{3+} , we would expect that the lattice parameter, a (the length of the side of the cubic unit cell), would increase with increasing gadolinium content. This trend is clearly seen in the data presented in Table 1 and graphically depicted in Figure 4.²

Many students have difficulty understanding the difference between a single-phase solid solution compound, such as $Y_{1.5}Gd_{1.5}Fe_5O_{12}$, where there is a 50% probability that the atom on the rare earth site is Y or Gd, and a physical two-phase mixture containing equal amounts of $Y_3Fe_5O_{12}$ and $Gd_3Fe_5O_{12}$. The X-ray powder patterns illustrate beautifully the difference between these two cases. Because the



Figure 5. A portion of the X-ray powder diffraction pattern of a) a two-phase mixture of $Y_3Fe_5O_{12}$ + $Gd_3Fe_5O_{12}$ and b) a single-phase solid solution compound $Y_{1.5}Gd_{1.5}Fe_5O_{12}$.

lattice parameters of $Y_3Fe_5O_{12}$ and $Gd_3Fe_5O_{12}$ differ, the positions of the diffraction peaks also should differ. Thus, in a two-phase mixture, formed for this purpose by physically mixing equal amounts of $Y_3Fe_5O_{12}$ and $Gd_3Fe_5O_{12}$, the diffraction pattern consists of two distinct patterns superimposed as shown in Figure 5a. However, the solid solution compound, $Y_{1.5}Gd_{1.5}Fe_5O_{12}$, should have a single lattice parameter that is intermediate between the pure yttrium-containing material and the pure gadolinium-containing material. The diffraction pattern of the solid solution compound should consist of one set of peaks that are midway between the two sets for the pure materials. This is shown in Figure 5b. The point should be made to students that a two-phase mixture, such as that described in Figure 5a, could be physically separated into the two components provided that one had "tweezers" small enough to pick out the particles of pure $Y_3Fe_5O_{12}$ and leave the pure $Gd_3Fe_5O_{12}$ behind. However, in practice, this physical separation is impossible to carry out in most cases, so we must design the syntheses to lead to a single-phase product.

Magnetic Properties

Electrons in a solid can behave cooperatively, leading to different properties from those observed in the gas or liquid states. An excellent example of cooperative behavior in solids is ferromagnetism.³ The difference between simple paramagnetism and ferromagnetism is shown in Figure 6. Unpaired electrons or spins are depicted in this figure by

²The shifts in the positions of the lines in the X-ray powder diffraction pattern are not clearly visible to the eye. However, by accurately determining the position of the lines and refining the lattice parameter, the shifts in the lattice parameters are significant.

³Theoretically, ferromagnetism is feasible in a liquid, although a liquid ferromagnet has not, to our knowledge, been reported. Ferrofluids have been made that are colloidal suspensions of tiny solid ferromagnetic particles (~100 Å in diameter) suspended in a liquid medium.

arrows. There are unpaired spins in a simple paramagnet (Fig. 6a), but in the absence of a magnetic field, the spins are oriented randomly relative to one another due to thermal motion. In a ferromagnet, the unpaired electrons communicate with one another and align in large regions known as magnetic domains⁴ (Fig. 6b). The sizes of the domains vary with the material but are typically on the order of hundreds of angstroms. In the absence of a magnetic field, the composite spins of individual domains are ordered randomly so that the net magnetization of a macroscopic piece of ferromagnet is zero. The application of a magnetic field (Fig. 6c) aligns all of the spins in the direction of the magnetic field, leading to an overall net magnetization. The application of a magnetic field to a paramagnet also would align the spins in the direction of the magnetic field, but the overall magnetization is less than in the case of a ferromagnet, because there is no communication or magnetic ordering between the spins.

The observed temperature-dependent magnetic behavior of a ferromagnet is derived from the interplay of thermal energy, which randomizes the spins, and the energy stabilization gained by aligning the spins in the ferromagnetic state. Ferromagnetism is an example of a solid-state phase transition with a characteristic transition temperature. Above an ordering temperature, known as the Curie

temperature (T_c), thermal energy is sufficient to break the alignment of the spins, and the material exhibits simple paramagnetism. But below T_c , the aligning forces overcome thermal randomization, the spins are locked in alignment, and the material becomes ferromagnetic.

Another type of magnetic ordering that can occur in solids is antiferromagnetism. In this type of magnetic behavior, the spins align opposite to each other below the ordering temperature⁵ to lessen the overall magnetism.

Ferrimagnetism, the type of magnetic ordering that occurs in the garnets, differs from ferromagnetism in that ions with different local moments are present in ferrimagnets. These unequal spins align antiparallel (antiferromagnetically) below the Curie temperature (T_c), but because the spins have different magnitudes, the cancellation is not complete. A simple schematic representing the differences between the magnetic ordering in ferromagnetism, antiferromagnetism, and ferrimagnetism is shown in Figure 7.


The magnetic properties of the rare earth iron garnets can be understood by considering that each type of coordination site in the crystal structure leads to a separate magnetic sublattice. Thus, for $Gd_3Fe_5O_{12}$, there are three magnetic sublattices comprised of the unpaired spins of the

⁴The application of magnetic properties in technology relies on the existence of domain walls and the ability to control the movement of the domain walls. For example, in permanent magnets, the domain walls are not moved easily. However, in the promising technology of magnetic bubble domain memory storage, the ability to move the domain walls underlies the writing capabilities in the read/write mode.

⁵The ordering temperature of an antiferromagnet is called the Neel temperature, T_N .


A. Paramagnet

Spins are randomized by thermal energy.

B. Ferromagnet

Spins are ordered in magnetic domains.

C. Ferromagnet

Spins are aligned with an applied magnetic field.


Figure 6. (a) In the absence of an applied magnetic field, the unpaired electrons (or spins) of a paramagnet are randomly oriented due to thermal energy. In contrast, the spins of a ferromagnet (b) spontaneously align with each other within regions known as magnetic domains. (c) In an applied magnetic field, the individual magnetic domains all align with the field, resulting in a large, overall magnetization.


Figure 7. The magnetic ordering of linear arrays of spins (unpaired electrons) in a ferromagnet, antiferromagnet, and a ferrimagnet. Ferrimagnetism is the type of magnetic ordering that occurs in the rare earth iron garnets.


Figure 8. The temperature dependence of the magnetization per gram for $Y_3Fe_5O_{12}$.

different cations: the gadolinium sublattice, the octahedral iron sublattice, and the tetrahedral iron sublattice. The total magnetic moment can be represented by eq 1:

$$M(T) = M_d(T) - M_a(T) - M_c(T) \quad (1)$$

where $M_d(T)$ represents the temperature-dependent magnetization of the tetrahedral sublattice, $M_a(T)$ represents the temperature-dependent magnetization of the octahedral sublattice, and $M_c(T)$ represents the temperature-dependent magnetization of the dodecahedral sublattice.⁶ Each term, such as $M_d(T)$, represents a sum of the magnetic moments at that particular crystallographic site (i.e., There are three tetrahedral sites in each formula unit.). The negative signs in eq 1 result from the fact that the spins in the octahedral sublattice (M_a) and the spins in the dodecahedral sublattice (M_c) align parallel to each other but antiparallel to the spins in the tetrahedral sublattice (M_d). The different magnetic sublattices can have different temperature dependencies for their magnetization that may lead to $M(T) = 0$ at a specific temperature below T_c . This temperature is the compensation temperature, T_{comp} .

For the material $Y_3Fe_5O_{12}$, the Fe^{3+} ions (with five unpaired electrons) in the octahedral holes order antiferromagnetically (spins aligned antiparallel) with the Fe^{3+} ions (also with five unpaired electrons) in the tetrahedral holes. However, because there are three tetrahedral sites and two octahedral sites in the garnet formula, there is a net magnetic moment of five unpaired electrons. There is

no magnetic contribution from the closed shell yttrium ion. Thus, $Y_3Fe_5O_{12}$ is strongly magnetic at all temperatures. This behavior can be demonstrated by the fact that the pellet of $Y_3Fe_5O_{12}$ can be picked up by the cow magnet at room temperature, dry ice temperature, and liquid nitrogen temperature. A plot of the magnetization of $Y_3Fe_5O_{12}$ versus temperature, shown in Figure 8, agrees with this qualitative behavior.

At the other extreme in the solid solution, $Gd_3Fe_5O_{12}$ has the same net five unpaired electrons from the two iron sublattices. But, in addition, each Gd^{3+} ion (there are three in the garnet formula) has seven unpaired electrons. The electrons in the gadolinium sublattice order antiferromagnetically (antiparallel) to the net five unpaired electrons of the tetrahedral iron. This suggests that the magnetic moment of $Gd_3Fe_5O_{12}$ should be $16 (3 \times 7 - 5)$ unpaired electrons. This is true if the moment is measured at very low temperatures. However, the gadolinium sublattice and the iron sublattice thermally randomize at different temperature dependencies. The result is a compensation temperature, T_{comp} , where the net magnetization is zero. In $Gd_3Fe_5O_{12}$, the compensation temperature is just below

⁶The identity of the subscripts is derived from the standard notation for the garnet formula, $C_3A_2D_3O_{12}$, where C represents the dodecahedral site, A represents the octahedral site, and D represents the tetrahedral site.


Figure 9. The temperature dependence of the magnetization per gram for $Gd_3Fe_5O_{12}$. The inset contains an enlargement of the region from 260 to 300 K in order to identify the compensation temperature ($T_{comp} = 289$ K).

Y₂GdFe₅O₁₂


Figure 10. The temperature dependence of the magnetization per gram for Y₂Gd₁Fe₅O₁₂. The inset contains an enlargement of the region from 50 to 90 K in order to identify the compensation temperature ($T_{\text{comp}} = 69$ K).

room temperature. This can be verified qualitatively by the observation that the Gd₃Fe₅O₁₂ pellet is not attracted to the cow magnet at room temperature. However, if this pellet is cooled in dry ice or liquid nitrogen, for example, the pellet becomes strongly magnetic and can be picked up with the cow magnet.⁷ The corresponding plot of magnetization versus temperature for this material is shown in Figure 9. Note that at the compensation temperature (289 K) and at room temperature, the magnetization is very small,⁸ but it increases dramatically as the temperature is lowered. Also note that at low temperatures, the magnetization of Gd₃Fe₅O₁₂ is roughly three times larger than the magnetization of Y₃Fe₅O₁₂, as is predicted from the described simple calculation of the number of unpaired electrons.

⁷This unusual magnetic behavior suggests an interesting experiment: If a magnet is secured at the proper distance above a cooling bath, upon cooling, the pellet will automatically jump up to the magnet. As the pellet warms, it should fall back into the cooling bath, where upon cooling again, the pellet will repeat this motion cycle. This behavior should continue indefinitely, provided that the cooling bath is replenished periodically and provided that the pellet does not fall out of the range of the magnet.

⁸In this particular measurement, the magnetization is not exactly zero, probably due to the presence of a small amount of GdFeO₃ in the sample.

The effect of combining both yttrium and gadolinium in the garnet structure in the solid solution compounds is that the compensation temperature becomes a strong function of composition (6). This is obvious from the magneti-

Table 2. Variation in Compensation Temperature, T_{comp} , with Composition for Y_xGd_{3-x}Fe₅O₁₂

Composition	T_{comp} (K)	Reference
Y ₃ Fe ₅ O ₁₂	— ^a	This work
Y _{2.7} Gd _{0.3} Fe ₅ O ₁₂	—	6
	(only measured above 73 K)	
Y _{2.1} Gd _{0.9} Fe ₅ O ₁₂	—	6
	(only measured above 73 K)	
Y ₂ Gd ₁ Fe ₅ O ₁₂	69	This work
Y _{1.8} Gd _{1.2} Fe ₅ O ₁₂	93	6
Y _{1.2} Gd _{1.8} Fe ₅ O ₁₂	158	6
Y ₁ Gd ₂ Fe ₅ O ₁₂	175	This work
Y _{0.6} Gd _{2.4} Fe ₅ O ₁₂	218	6
Gd ₃ Fe ₅ O ₁₂	289	This work

^aNo compensation temperature was observed down to 5 K.

Y₁Gd₂Fe₅O₁₂


Figure 11. The temperature dependence of the magnetization per gram for $Y_1Gd_2Fe_5O_{12}$. The inset contains an enlargement of the region from 130 to 220 K in order to identify the compensation temperature ($T_{comp} = 175$ K).

zation behavior plotted in Figures 10 and 11 for the two intermediate compositions $Y_2Gd_1Fe_5O_{12}$ and $Y_1Gd_2Fe_5O_{12}$. From the data, it is clear that the compensation point for $Y_2Gd_1Fe_5O_{12}$ is ~70 K (Fig. 10) and for $Y_1Gd_2Fe_5O_{12}$ is ~175 K (Fig. 11). Again, these can be verified qualitatively with the cow magnet. If $Y_2Gd_1Fe_5O_{12}$ is cooled in liquid nitrogen, the pellet cannot be picked up with the cow magnet. Likewise, if $Y_1Gd_2Fe_5O_{12}$ is cooled in dry ice, the pellet only weakly responds to the cow magnet. This information is summarized in Table 2, which lists the compensation temperatures for various compositions in the solid solution series, $Y_xGd_{3-x}Fe_5O_{12}$.

Conclusions

This laboratory presents an opportunity to introduce solid-state chemistry into the introductory chemistry course. The synthesis and qualitative magnetic measurements described above are not difficult, yet they provide insight into preparative chemistry and the magnetic properties of solids. In addition, the experiment is such that the students have their "own" individual compound to prepare in order to contribute to the overall understanding of the

system. We believe that this team approach generates a more positive laboratory experience.

Acknowledgment

The authors would like to thank Frank DiSalvo at Cornell University for providing the initial suggestion for this experiment, Debbie Stafstien for early experimental investigations, and the Dreyfus Foundation and the National Science Foundation (grant USE-9150484) for financial support.

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Place 500 mg of WO₃ into a 250 mL beaker. Carefully pour about 50 mL dilute hydrochloric acid (3 molar) onto the WO₃. Does anything happen? Now very carefully drop some zinc filings (*ca.* 3 g) into the acid. After all reaction has ceased filter off the solid product. After washing, place it on a filter paper and observe what happens when you leave it exposed to the air. Write balanced equations for both the reduction of WO₃ and its subsequent re-oxidation.

Make pellets of WO₃ and H_xWO₃ in a pellet press and measure the electrical resistance using an ohmmeter. Is there any difference?

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The Conversion of Triphenylphosphine to the Oxide: A Simple Experiment to Illustrate Catalytic Oxotransfer

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Oxygen atom transfer reactions involving dioxomolybdenum(VI) complexes are of current interest because of their relevance in some important biological processes (1). Most synthetic complexes that mimic oxotransferases contain sterically demanding sulfur ligands because two of the criteria for relevant models are the sulfur environment and hindrance to irreversible μ -oxo dimer formation. However, some molybdenum complexes lacking these requisites show a remarkable ability to act as oxotransfer agents (2).

Here we describe the facile synthesis of a simple, very active, catalyst ($\text{MoO}_2\text{Cl}_2(\text{DMSO})_2$, (DMSO = dimethylsulfoxide)) and its application in the oxidation of triphenylphosphine (PPh_3) to triphenylphosphine oxide (OPPh_3) by DMSO.

This oxidation of PPh_3 provides a good example of the catalytic oxotransfer process because PPh_3 is inexpensive, resistant to aerial oxidation, widely used in oxotransfer studies, and has a low volatility. Conducting the experiment in an aprotic medium with simple reagents minimizes the number of species involved in the overall process, which facilitates the interpretation of the results at the time that dramatically evidences the role of the catalyst (cleavage of strong multiple X-O bonds).

Preparing the Catalyst

The procedure presented here is related to that recently reported for the multigram scale synthesis of $\text{MoO}_2\text{Cl}_2(\text{DMSO})_2$ from aqueous hydrochloric solutions of MoO_3 (3) although the use of $(\text{NH}_4)_6\text{Mo}_7\text{O}_{24}\cdot 4\text{H}_2\text{O}$, in place of MoO_3 , saves time and gives a product that is not contaminated appreciably with NH_4Cl . A typical run is described:

In a 5-mL vial are placed 0.25 g (0.2 mmole) of $(\text{NH}_4)_6\text{Mo}_7\text{O}_{24}\cdot 4\text{H}_2\text{O}$, 1 g H_2O and 2 g (20 mmole) 37% HCl . The mixture is heated for 10 min at -100°C (e.g., boiling water bath) with stirring, and the resulting solution is cooled to room temperature. Then, 0.5 g (6.4 mmole) DMSO is added with stirring and the mixture is cooled again. The white precipitate is filtered, washed with acetone (2×1 mL) and peroxide-free diethyl ether (2×1 mL), and dried in air. 0.37 g (~75% yield) of $\text{MoO}_2\text{Cl}_2(\text{DMSO})_2$ is obtained as a white microcrystalline powder.

The compound melts at $172 \pm 1^\circ\text{C}$ and is very stable at ambient conditions. The IR spectrum is characterized by two strong bands at 921 and 892 cm^{-1} corresponding to the stretching vibrations of the *cis*- MoO_2 arrangement (3, 4).

Preparing OPPh_3 via Catalytic Oxotransfer

In a medium-size test tube (e.g., 100×15 mm) is placed 50 mg (0.14 mmole) of $\text{MoO}_2\text{Cl}_2(\text{DMSO})_2$, 500 mg (1.9 mmole) PPh_3 and 2.5 g (32 mmole) DMSO . The mixture is heated at -100°C for 15 min in a well-ventilated hood because the foul-smelling dimethylsulfide (DMS) is released. The resulting colorless solution is poured into a 25-mL vessel containing 12


mL of ~ 1 M NaOH solution (this solution can be prepared by dissolving 0.5 g of NaOH in 12 mL of water). The precipitate is filtered, washed with water (3×2 mL), and dried in an oven at -120°C for 30 min. One obtains 0.49 g (92% yield) of white OPPh_3 . This product melts at $157 \pm 1^\circ\text{C}$ and displays an IR spectrum characteristic of OPPh_3 (e.g., a strong $\nu(\text{P}=\text{O})$ band at 1182 cm^{-1}) (4).

Discussion

The above procedures can be explained, in a simplified manner, as follows.

With regard to the preparation of the catalyst:

1. Ammonium heptamolybdate reacts with warm hydrochloric acid to form NH_4Cl , H_2O , and *cis*, *trans*, *cis*- $\text{MoO}_2\text{Cl}_2(\text{H}_2\text{O})_2$.


This species has been detected in solution (5), isolated as $\text{MoO}_2\text{Cl}_2(\text{H}_2\text{O})_2 \cdot \text{H}_2\text{O} \cdot (\text{C}_2\text{H}_5)_4\text{NCl}$ (6) and, more recently, as the outer-sphere diglyme adduct $[\text{MoO}_2\text{Cl}_2(\text{H}_2\text{O})_2] \cdot (\text{diglyme})_2$ (7). Molybdenum trioxide, previously dissolved in the minimum amount of $\sim 12\%$ NH_3 , can be used in place of ammonium molybdate in this preparation.

2. DMSO is able to displace H_2O from the coordination sphere of many metallic centers producing a more hydrophobic environment around the metal. Thus, the $\text{HCl-NH}_4\text{Cl-MoO}_2\text{Cl}_2(\text{H}_2\text{O})_2$ aqueous solution reacts with excess DMSO , and the resulting *cis*, *trans*, *cis*- $\text{MoO}_2\text{Cl}_2(\text{DMSO})_2$ precipitates while NH_4Cl remains dissolved.

With regard to the preparation of OPPh_3 :

1. PPh_3 reacts with DMSO in the presence of $\text{Mo}^{\text{VI}}\text{O}_2\text{Cl}_2(\text{DMSO})_2$ to give OPPh_3 and DMS as represented in the catalytic cycle:

2. The lack of color in the resulting DMSO solution, which contains mainly OPPh_3 , the catalyst, and some DMS , suggests that Mo(V) or Mo(IV) species are not present in significant amounts (8).

3. The DMSO solution is poured into sufficient aqueous NaOH so that OPPh_3 precipitates as the $\text{MoO}_2\text{Cl}_2(\text{DMSO})_2$ is converted into soluble sodium chloride and molybdate, thus avoiding contamination of the OPPh_3 .

Final Considerations

To highlight the role of the molybdenum complex it is advisable to conduct a blank experiment. To minimize

the microscale laboratory

waste, the procedure for the preparation of OPPh_3 can be modified so that a solution of PPh_3 in DMSO is prepared and heated while the catalyst is prepared. Then, a sample of ~ 0.5 mL is taken off and treated with water to precipitate the unreacted PPh_3 , which may be further characterized (mp 80°C). Finally, the catalyst is added, and the solution is worked up as described.

These experiments can be performed in a three-hour laboratory period, although this catalyzed oxotransfer reaction also proceeds at room temperature over a longer period. If necessary, significant time can be saved by washing the wet OPPh_3 with small portions of ethanol (e.g., 2×1 mL) and diethyl ether (e.g., 2×1 mL), followed by drying at room temperature, although this reduces the yield considerably.

Acknowledgment

To the editor for helpful comments and to DGICYT (Grant PS90-0134) for financial support.

A198 Journal of Chemical Education

Reprinted from Journal of Chemical Education, Vol. 72, Page A7-A8, January 1995.

Copyright 1995 by the Division of Chemical Education of the American Chemical Society, and reprinted by permission of the copyright owner.

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$\text{MoO}_2(\text{acac})_2$ (acac = 2,4-pentanedionate) is one of the most important starting materials for preparing other dioxomolybdenum compounds of current interest (1) because:

- acac can be displaced from the coordination sphere of molybdenum by numerous ligands
- acacH is volatile
- acacH is soluble in many organic solvents and can be removed from the product
- $\text{MoO}_2(\text{acac})_2$ works well in fundamental IR and ^1H NMR experiments

The synthesis described here is a modification of that recently reported (2). The main alterations are (in addition to the scale and molar ratio of reactants):

- starting from MoO_3 powder, which is cheaper than the conventional $(\text{NH}_4)_6\text{Mo}_7\text{O}_{24} \cdot 4\text{H}_2\text{O}$.
- using the common 24% (13 M) instead of 15% (8 M) NH_3 solution.
- adding the reagents sequentially by weighing, which facilitates the work up.
- ending the washings with diethyl ether so that drying the product is more rapid.

A Typical Run


In a 5-mL conical vial (on an electronic balance) is placed 0.25 g (1.7 mmole) of MoO_3 powder, 0.40 g H_2O and 0.40 g (5.5 mmole) 24% NH_3 (Pasteur pipet). The mixture is stirred with a glass rod for about one min to obtain a clear colorless solution (exothermic reaction). Then, 0.75 g (7.5 mmole) of 2,4-pentanedione is added, and the mixture is stirred; however, two phases remain. Finally, 1 g (10 mmole) of 63% HNO_3 is added dropwise (exothermic reaction). The mixture is cooled to room temperature while stirring and the solid formed is ground until a yellow pow-

Literature Cited

1. (a) Holm, R. H. *Chem. Rev.* 1967, 67, 1401, and references therein. (b) Holm, R. H. *Covul. Chem. Rev.* 1990, 100, 163, and references therein. (c) Enemark, F. H.; Young, C. G. *Adv. Inorg. Chem.* 1994, 40, 1, and references therein.
2. (a) Roberts, S. A.; Young, C. G.; Kipke, C. A.; Cleland, W. E., Jr.; Yamanouchi, K.; Carducci, M. D.; Enemark, F. H. *Inorg. Chem.* 1990, 29, 3650. (b) Arzumanyan, H.; López, R.; Agrifoglio, G. *Inorg. Chem.* 1994, 33, 3177. (c) Arnáiz, F. J.; Aguado, R.; Martínez de Harbaya, J. M. *Polyhedron* 1994, 13, 3515.
3. Arnáiz, F. J. *Inorg. Synth.*, in press.
4. Butcher, R. J.; Oms, H. P.; MacLagan, R. G. A. R.; Kipton, H.; Powell, J.; Wilkins, C. J.; Hian, Y. S. *J. Chem. Soc., Dalton Trans.* 1978, 1223.
5. Coalington, J. M.; Taylor, M. J. *J. Chem. Soc., Dalton Trans.* 1990, 1, 41.
6. Taylor, M. J.; Jirong, W.; Rickard, C. E. F. *Polyhedron* 1993, 12, 1433.
7. Arnáiz, F. J.; Aguado, R.; Sanz-Aparicio, J.; Martínez-Ripoll, M. *Polyhedron* 1994, 13, 2745.
8. Greenwood, N. N.; Earnshaw, A. *Chemistry of the Elements*; Pergamon Press: Oxford, 1984.

der results. (Powdering the initially lumpy precipitate improves the purity of the product and requires less than 5 min.) The product is filtered through a Hirsch funnel, washed consecutively with distilled water (4×2 mL), ethanol (2×2 mL) and peroxide-free diethyl ether (2×2 mL), and dried in air. Washings also may be done by decanting. The powder settles sufficiently in less than 1 min so that more than 75% of the liquid can be removed. Using a Pasteur pipet is advisable in this case to reduce product losses because $\text{MoO}_2(\text{acac})_2$ forms a supernatant skin every time water is added.

The entire process requires about 30 min. About 0.19 g of a light yellow finely powdered product melting at $179 \pm 1^\circ\text{C}$ (the reported value) is obtained. The IR spectrum of the compound shows two strong bands at 900 and 930 cm^{-1} corresponding to the stretching vibrations characteristic of the *cis*- MoO_2 arrangement. Singlets at 2.12 and 2.14 δ (CD_2Cl_2 , room temperature), assignable to the two sets of inequivalent methyl groups (see drawing), are observable in the ^1H NMR spectrum.

Acknowledgement


I am grateful to the editor for helpful comments.

Literature Cited

1. See, e. g.: (a) Srivastava, A.; Ma, Y.; Pankavatselvan, R.; Dinges, W.; Nicholas, K. M. *J. Chem. Soc., Chem. Commun.* 1992, 853. (b) Bustos, C.; Manzur, C.; Carrillo, D.; Robert, F.; Gouzerh, P. *Inorg. Chem.* 1994, 33, 1427.
2. Chakravorti, M. C.; Bandyopadhyay, D. *Inorg. Synth.* 1992, 29, 129.

(2) Complexes of Sulfoxides and Related Compounds

Cotton et al.⁴¹⁰ studied the infrared spectra of sulfoxide complexes to see whether coordination occurs through oxygen or sulfur. The electronic structure of sulfoxides may be represented by a resonance hybrid of the structures:

If coordination occurs through oxygen, contribution of structure II will decrease and result in a decrease of the SO stretching frequency. If coordination occurs through sulfur, contribution of structure I will decrease and may result in an increase of the SO stretching frequency.

TABLE III-34. SHIFTS OF SO STRETCHING BANDS IN DPSO AND DMSO COMPLEXES (CM⁻¹)⁴¹²

Metal	DPSO Complex		DMSO Complex
	$\nu(\text{SO})$	Shift	Shift
Ca(II)	1012-1035	0 - (-23)	-
Mg(II)	1012	-23	-
Mn(II)	983-991	-45	-41
Zn(II)	987-988	-47	-
Fe(II)	987	-48	-
Ni(II)	979-982	-55	-45
Co(II)	978-980	-56	-51
Cu(II)	1012, 948	-23, -87	-58
Al(III)	942	-93	-
Fe(III)	931	-104	-

The far-infrared spectra of DMSO complexes have been studied by Johnson and Walton,⁴¹³ who assigned the M—O stretching bands at 500–400 cm⁻¹. The S-bonded DMSO complexes (Pt(II) and Pd(II)) also exhibit a strong band in the same region. Hence it is difficult to differentiate the O-

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

Their actual goal was the synthesis of fulvalene:


Structures which were originally proposed:

1952 E. O. Fischer: “**Double-cone structure**”

based on:

- X-ray structural analysis (compare the electron density contour lines)
- diamagnetism
- chemical behavior

1952 G. Wilkinson, R. B. Woodward: “**Sandwich structure**”

based on:

- IR spectroscopy
- diamagnetism
- dipole moment = 0

Woodward: The cyclopentadienyl rings of $\text{Fe}(\text{C}_5\text{H}_5)_2$ are amenable to electrophilic substitutions. This similarity to the aromatic behavior of benzene led to the trivial name **ferrocene**, later being extended to the designation **metallocenes** for compounds $\text{M}(\eta^5\text{-C}_5\text{H}_5)_2$ in general.


15.4.3.1 Binary Cyclopentadienyl Transition-Metal Complexes

PREPARATION


1. Metal Salt + Cyclopentadienyl Reagent

Dicyclopentadiene has first to be “cracked” in a retro-Diels-Alder reaction to give monomeric C_5H_6 . Cyclopentadiene is a weak acid, $\text{p}K_a \approx 15$; it can be deprotonated by strong bases or by alkali metals. Cyclopentadienylsodium (NaCp) is the most common reagent for the introduction of cyclopentadienyl ligands.

Ferrocene, ruthenocene and osmocene are susceptible to **electrophilic substitution**. Compared to benzene, ferrocene reacts 3×10^6 times faster.

The electrophile should not be an oxidizing agent, as substitution would then be suppressed by oxidation to the ferricinium ion $[FeCp_2]^+$ which is inert to attack by electrophiles. That electrophilic substitution of ferrocene does not involve direct participation of the metal has been established by the intramolecular acylation of two isomeric ferrocene carboxylic acids. It was found that the *exo*-isomer, in which the acylium ion cannot directly interact with the metal, cyclizes faster than the *endo*-isomer (Rosenblum, 1966):

This result suggests that electrophilic attack occurs from the *exo*-direction.

Harold T. McKone
St. Joseph College
West Hartford, CT 06117

Acylation of Ferrocene

Effect of temperature on reactivity as measured by reverse phase high performance liquid chromatography

The Friedel-Crafts Reaction, a widely used synthetic tool, is an integral part of most organic laboratory courses. Formally, benzene was the preferred substrate for this electrophilic substitution reaction. However, recent reports confirming benzene's carcinogenicity necessitate using a safer aromatic compound. Ferrocene (I) is a reasonable alternative since it is relatively non-toxic, stable, more reactive than benzene, inexpensive in the quantities needed per student, and affords an interesting array of acylation products (1-2). In fact, the Friedel-Crafts acylation of ferrocene with product identification by thin layer and open column chromatography is currently a popular undergraduate experiment (3-10). Haworth and Liu recently reported a high performance liquid chromatographic (lc) separation of acylated ferrocenes using "normal-phase" conditions of silica packing and an ether: methanol mobile phase (11-12). The increasing use of lc in separating components in organic reaction mixtures encouraged us to investigate incorporating this lc separation of ferrocenes into our undergraduate laboratory. However, "normal-phase" lc has several disadvantages over the more versatile "reverse-phase" chromatography. The former requires relatively long column equilibration times and generally necessitates using non-aqueous mobile phases. We have developed a reverse-phase separation of these ferrocenes using a C₁₈ bonded to silica column (μ BONDAPAK/C₁₈-Waters Associates) and a far less bothersome mobile phase of methanol:water. A complete introduction to reverse phase lc may be found in references (13) and (14)

Instrumentation

A Waters ALC/GPC-204 Liquid Chromatographic with a Model U6K Injector and a Houston Omni-Scribe Model 5211-1 Recorder were used. The prepacked analytical column consisted of C₁₈ bonded to 10 μ m silica particles (Waters Associates, Inc., Milford, MA 01757), and the mobile phase was a mixture of methanol:water in a ratio of 3:1. The flow rate was 0.9 ml/min and the absorbance detector was at 254 nm.

Experimental Conditions

Since ferrocene acylates readily, the reaction may be investigated at a variety of temperatures from 45°-100°C. A styrofoam cup serves as an excellent individual constant temperature bath. The results described here were obtained at 100°C with each student assigned a different reaction time.

Procedure

Place 0.2 g (~0.001 mole) of ferrocene, 1.0 ml of acetic anhydride and 0.1 ml of 85% phosphoric acid into a small test tube, cork, shake, and immediately place in a constant temperature bath. At designated intervals (15 sec, 30 sec, etc.) quench the reaction by pouring onto about 5 g of crushed ice. Allow the red-brown mixture to stand (with occasional stirring) for 15 min. Filter, wash the solid with 5 ml of cold


Figure 1. Acylation of ferrocene.


Figure 2. Liquid Chromatogram of the Separation of Standards of Ferrocene (I), Acetylferrocene (II), and 1,1-Diacetylferrocene (III).

Table 1. *k'* and *R_f* Values for Ferrocene and Some Acylated Derivatives

Compound	<i>k'</i>	<i>R_f</i>
Ferrocene (I)	2.14	0.9
Acetylferrocene (II)	1.00	0.6
1,1-Diacetylferrocene (III)	0.64	0.3

Table 2. LC Analysis of Products of Acylation of Ferrocene at 100°C

Reaction Time	Product Composition			
	Ferrocene	Acetyl-Ferrocene	Diacetyl-Ferrocene	Other
0.1 min	(+)	trace	(-)	(-)
0.5 min	(+)	(+)	(-)	(-)
1.5 min	trace	(+)	trace	(-)
10 min	(-)	(+)	(+)	trace
25 min	(-)	(+)	(+)	(+)

water, and air-dry. Dissolve a small amount of the dry solid in 5 ml of methanol, filter¹, and separate by tlc and lc.

Results and Discussion

Ferrocene acylates readily, yielding a variety of products as is indicated in Figure 1. A typical chromatogram, depicting the separation of ferrocene (I) and two major acylation products (acetylferrocene, II, and 1,1-diacetylferrocene, III) is shown in Figure 2. Under reverse phase conditions, the most polar compound, III, elutes first while the relatively non-polar ferrocene (I), which is adsorbed more strongly on the non-polar packing, elutes last. The separation is complete in less than 10 min allowing many injections during a single laboratory session. Table 1 summarizes the results of the lc analysis of reactions quenched at different times and Table 2 lists *k'* and *R_f* values for

¹To extend column life, all samples should be filtered prior to lc separation.

²*k'*, the capacity factor, may be defined in terms of the difference between the time (*t*₁) required for a component to travel the length of the column and the time (*t*₀) required for solvent (or other non-retained molecules) to traverse the column, divided by *t*₀. See Figure 2.

ferrocene and its mono and diacylated products². Thin layer chromatograms were run on silica using a developing solution of toluene:methanol in proportions of 20:1.

Conclusions

Monitoring the ferrocene acylation reaction by reverse phase lc is an excellent method of integrating this important separations technique into the undergraduate organic laboratory and of providing students with hands-on experience with lc. Because of the short reaction time, simple product work-up, and rapid lc separation, many students have the opportunity to complete this experiment in one laboratory session.

Acknowledgment

The author wishes to thank the National Science Foundation (Undergraduate Research Participation Grants #GY5892 and GY7241) and Ms. Gerry Granda, Ms. Helen Grzybowski, and Ms. Diane Sommerer for their work in tlc

analysis. Thanks are also extended to Waters Associates, Milford, MA, for providing materials and financial support.

Literature Cited

- (1) Rosenblum, M., "Chemistry of the Iron Group Metallocenes: Ferrocene, Ruthenocene and Osmocene," John Wiley and Sons, Inc., New York, 1965.
- (2) Dauben, W. G., "Organic Reactions," 17, John Wiley and Sons, Inc., New York, 1969.
- (3) Clark, A. K. and Bell, C. E., "Laboratory Experience in Organic Chemistry," Charles E. Merrill Publishing Co., Columbus, Ohio, 1972, p. 55.
- (4) Fieser, L. F., and Williamson, K. L., "Organic Experiments," 4th Ed., D.C. Heath and Co., Lexington, MA, 1975, p. 169.
- (5) Adams, R., *et al.*, "Laboratory Experiments in Organic Chemistry," 7th Ed. Macmillan Publishing Co., New York, 1979, p. 462.
- (6) Wade, L. G., *J. CHEM. EDUC.*, **55**, 208 (1978).
- (7) Bozak, R. E., *J. CHEM. EDUC.*, **44**, 78 (1969).
- (8) Herz, J. E., *J. CHEM. EDUC.*, **43**, 599 (1966).
- (9) Bohan, J. M., *et al.*, *J. CHEM. EDUC.*, **50**, 367 (1973).
- (10) Gilbert, J. C., and Monti, S. A., *J. CHEM. EDUC.*, **50**, 369 (1973).
- (11) Haworth, D. T., and Liu, T. J., *J. CHEM. EDUC.*, **53**, 730 (1976).
- (12) Haworth, D. T., and Liu, T. J., *J. Chromatogr.*, **14**, 519 (1976).
- (13) Colin, H., and Guiochon, G., *J. Chromatogr.*, **141**, 289 (1977).
- (14) Gloor, R., and Johnson, E. L., *J. Chromatogr.*, **15**, 419 (1977).

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All reagents and solvents are used as supplied except 1,2-diaminobenzene (*o*-phenylenediamine) that should be purified by a single recrystallization from EtOH in the presence of decolorizing charcoal or by sublimation. Dry the *n*-butanol over 4 Å molecular sieves.

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Place 2.00 g of Ni(OAc)₂·4H₂O (8.03 mmol), 1.73 g of *o*-phenylenediamine (1.60 mmol), 35 mL of **dry** butan-1-ol, 1.7 mL of 2,4-pentanedione (~1.66 mmole), and a stirring bar in a 100-mL round-bottom flask. Adapt a reflux condenser to the flask and set it in a heating mantle. Bring the mixture to a brisk reflux and stir for 1.5 hours, noting any color changes. Then, concentrate the solution (by slow distillation of the BuOH at ordinary pressure) to approx. 20 mL and reflux it 35 minutes more.

Remove the flask from the heat source and allow it to cool until it is just warm to the touch. Add 30 mL of methanol and cool the mixture in an ice-salt bath for at least 15 min to precipitate the nonpolar macrocyclic complex. Filter the solid and wash it with methanol (about 2 x 10 mL; the washings should be colorless to pale green). The dark purple microcrystalline macrocycle is generally obtained in a 37-40% yield. Retain a small sample of product for IR analysis. Use the remaining product to synthesize the free macrocycle as described below.

Preparación de [H₄(Macrocycle)][NiCl₄]

Prepare a suspensión of 1.0 g of the nickel complex in 30 mL of ethanol in a 100-mL round-bottom flask. Bubble HCl gas through the suspension at a moderate flow and stir it from time to time. (See note at the end of the experimental part). **Caution: the**

mixture gets quite warm. Once a copious turquoise precipitate forms, filter it off and wash it with ethanol (2×5 mL) and ether to obtain about ~95% yield of product.

Aislamiento de la sal del macrociclo (Mac) $[H_4(\text{Mac})][BF_4]_2$

Dissolve the turquoise $[H_4(\text{Mac})][NiCl_4]$ salt in 10 mL water. Occasionally white, water insoluble, impurities are present at this stage and can be removed by filtration. Add 1.0 g of solid NH_4BF_4 to the filtrate and swirl until a large amount of white precipitate forms. Filter the $[H_4(\text{Mac})][BF_4]_2$ in a sintered glass filter and wash with water until the product is pale green. (In this step the $(NH_4)_2NiCl_4$ byproduct needs to be washed away to prevent nickel from re-coordinating into the macrocycle upon basification). Transfer the sticky solid to a 25 mL Erlenmeyer flask, rinsing the sintered glass filter with methanol several times. Enough methanol should be present in the flask to cover the salt. Add triethylamine dropwise to the suspension and swirl. Filter the free yellow macrocycle, wash with 5 mL of methanol, and air dry. Record the IR spectrum.

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
Charge a 10-mL round-bottom flask with $Cu(OAc)_2 \cdot H_2O$ (0.1 g, 0.6 mmol), a stirring bar, and acetonitrile (4 mL). Attach a condenser and reflux the suspension until all the $Cu(OAc)_2 \cdot H_2O$ dissolves. Remove the heater and, with a micropipete add slowly a mixture of the macrocycle (0.1 g, 0.56 mmol), acetonitrile (2 mL), and triethylamine (0.20 mL). Reflux the resulting dark-green solution for about 10 min, then cool to room temperature, and finally keep it in an ice bath for 20 min. Filter the dark green microcrystalline solid and wash it with a 2 mL of acetonitrile and then 2 mL of ether.

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2. Synthesis of cobalt complex of encapsulating ligand L₁

Procedure

Dissolve 2.45 g of $[\text{Co}(\text{en})_3]\text{Cl}_3$ and 1.2 g of Na_2CO_3 in 25 cm³ water and add 18 cm³ of 40% aqueous formaldehyde (with CAUTION!) and 2.5 cm³ of nitromethane. Allow the mixture to stand at room temperature for 4 hours and add 5 cm³ ethanol to induce crystallization. Filter the resulting solid off using suction and wash it with ethanol (10 cm³) and diethyl ether (10 cm³). Dry the solid in air and store it.

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Referencias

E.C. Constable, *Coordination chemistry of macrocyclic compounds*, Oxford Chemistry Primers, 1999

K. Nakamoto, *Infrared and raman spectra of inorganic and coordination compounds*, 4th edition, John Wiley and Sons: New York, 1986.

A. Gordon, *The chemist's companion: A handbook of practical data, techniques, and references*, John Wiley and Sons: New York, 1972.

Facile Synthesis of a Macrobicyclic Hexamine Cobalt(III) Complex Based on Tris(Ethylenediamine)Cobalt(III)

An Advanced Undergraduate Inorganic Laboratory Experiment

Jack MacB. Harrowfield

The University of Western Australia, Nedlands 6009, Western Australia

Geoffrey A. Lawrance

The University of Newcastle, New South Wales 2308, Australia

Alan M. Sargeson

The Australian National University, Canberra 2601, Australia

The synthesis and resolution of tris(ethylenediamine)-cobalt(III) is a classical undergraduate experiment in coordination chemistry, which had been described in THIS JOURNAL and elsewhere (1-4). Chelation of cobalt(III) by multidentate ligands has been extended in recent years to include many examples of macromonocyclic ligands (5), and recently examples of saturated macrobicyclic ligands (6). The latter complexes totally encapsulate the metal ion, and are prepared by an organic template reaction about the inert Co(en)_3^{3+} precursor (en = ethylenediamine = $\text{H}_2\text{NCH}_2\text{CH}_2\text{NH}_2$), e.g.





Figure 1. Visible absorption spectra of (A) $[\text{Co(dinosar)}]\text{Br}_3$ and (B) $[\text{Co(en)}_3]\text{Br}_3$.

The encapsulation of the cobalt ion in the macrobicyclic ligand renders it substitution inert, even in the normally labile cobalt(II) state. The synthesis is a good example of reactions of coordinated ligands, and it proceeds in a stepwise fashion, most likely as illustrated in the reaction scheme below. The constraint and ordering of reactants about the metal ion and the small ring closures involved make the template synthesis entropically favorable.

This type of synthesis introduces a new class of coordination compounds as direct descendants from classical complexes. The synthesis of an aza-capped analog of $\text{Co}(\text{dinosar})^{3+}$ has been described (7), but that synthesis involves an isolation procedure less facile, and less suitable for subsequent experiments, than is the case of $\text{Co}(\text{dinosar})^{3+}$ (dinosar = 1,8-dinitro-3,6,10,13,16,19-hexaazabicyclo-(6.6.6)icosane). The described synthesis from $\text{Co}(\text{en})_3^{3+}$ can be completed within 2 to 3 h, and is inexpensive and safe. Various aspects of visible, infrared and nuclear magnetic resonance spectroscopy, and of electrochemistry and optical activity can be incorporated into the program as desired.

Synthesis of the Macrobicyclic Complex

$\text{Co}(\text{en})_3\text{Br}_3$ (2.5 g) and anhydrous Na_2CO_3 (1.5 g) are dissolved in water (90 mL) in a 200-mL beaker, and nitromethane (5 mL) and 37% aqueous formaldehyde (40 mL) added. (Caution: Reaction should be performed in a fume hood.) After stirring briefly to mix reagents, the solution is placed (not stirred) in a waterbath and maintained at $35^\circ\text{--}40^\circ\text{C}$ for approximately 1 to $1\frac{1}{2}$ h. The color rapidly darkens from yellow to brown, and after ~ 1 h an appreciable amount of yellow crystals form. The solution is stirred briefly and occasionally with a glass rod to help induce crystallization. The solution is removed from the waterbath after the $1\text{--}1\frac{1}{2}$ h period, and cooled immediately to $\sim 5^\circ\text{C}$ in an ice bath. The crystals are collected on a sintered frit, under suction. The filtrate is discarded (fume hood sink). If crystallization upon cooling is not substantial, addition of up to an equal volume of ethanol to induce precipitation may be employed but is normally not required.

The collected solid is dissolved on the frit using portions of hot dilute acetic acid (a total of ~ 30 mL of ~ 1 M). The hot yellow filtrate is treated immediately with 15 mL of concentrated HBr. A copious fine yellow crystalline precipitate forms. After cooling to room temperature, this is collected on a frit under suction, washed with cold 1:1 methanol:water, methanol, then ether, and dried at the pump and, if required, in a vacuum desiccator (yield 1.95g, 55%).

Characterization of the Macrobicyclic Complex

Electronic Spectrum: The $\text{Co}(\text{dinosar})^{3+}$ complex has absorption maxima in the visible region at 475 nm (ϵ 148 $\text{M}^{-1}\text{cm}^{-1}$) and 342 nm (ϵ 133), typical of an octahedral saturated CoN_6^{3+} chromophore, and similar to the $\text{Co}(\text{en})_3^{3+}$ precursor (466 nm (ϵ 75), 338 nm (ϵ 68)), Figure 1.

Infrared Spectra: The "fingerprint" pattern of the macrobicyclic ligand is quite different from that of the $\text{Co}(\text{en})_3\text{Br}_3$ precursor. The characteristically strong absorbance of $-\text{NO}_2$ can be seen at $\sim 1560\text{ cm}^{-1}$ in a KBr disc spectrum of the macrobicyclic complex.

NMR Spectra: The ^{13}C NMR spectrum (JEOL FX-60 spectrometer, D_2O , shifts relative to dioxan, +ve downfield) illustrates the essentially D_3 symmetry in the macrobicyclic complex. Three signals are observed at -11.3 ppm (due to $-\text{CH}_2-$ groups of the en ring residues), -15.1 ppm (due to $-\text{CH}_2-$ groups of the capping units), and a resonance at $+21.4$

ppm (due to the tertiary $\geq\text{C}-\text{NO}_2$ of the caps). For $\text{Co}(\text{en})_3^{3+}$, a single resonance occurs at -21.8 ppm. The high symmetry of the macrobicyclic is also evident from the ^1H NMR spectrum (200 MHz spectrometer, D_2O , N-deuterated, shifts relative to DSS); mainly two overlapping quartets are observed, associated with AB (or AA'BB') patterns of the cap and en methylene groups, Figure 2.

Electrochemistry: In water, a reversible 1 e^- couple at -0.20 V (versus SCE, mercury electrode, 0.1 M NaClO_4 electrolyte) is observed for $\text{Co}(\text{dinosar})^{3+}$ in both dc (at DME) and cyclic (at HMDE) voltammetry. Further, an irreversible 8 e^- reduction occurs in neutral solution at -0.7 V, due to reduction of both $-\text{NO}_2$ capping groups to $-\text{NHOH}$ (4 e^- each). Cyclic voltammetry at a hanging mercury drop electrode clearly defines the reversible nature of the first ($\text{Co}(\text{III})/\text{Co}(\text{II})$) reduction, and the irreversible nature of the second ($-\text{NO}_2 \rightarrow -\text{NHOH}$) reduction. By comparison, $\text{Co}(\text{en})_3^{3+}$ exhibits only a single wave at -0.41 V, and this reduction is not fully reversible. Detailed electrochemical studies of a range of such macrobicyclics have appeared (8).

Further Experiments: Reduction of the $\text{Co}(\text{dinosar})^{3+}$ complex chemically with Zn/HCl and air oxidation subsequently yields the $(\text{Co}(\text{diansarH}_2)\text{Cl})_5$ complex from HCl solution with fully reduced ($-\text{NH}_3^+$) groups on the apical caps (9). The absence of the $-\text{NO}_2$ groups can be readily established by IR spectroscopy and electrochemistry. The original synthesis can be performed with optically active $\text{Co}(\text{en})_3^{3+}$, whence the product is optically active. Reduction with Zn/HCl and subsequent reoxidation with air (the cobalt(II) salts are air sensitive but can be isolated and their properties studied) yields a complex which retains its chirality. In fact, reduction/reoxidation of the $\text{Co}(\text{diansarH}_2)^{5+}$ ($E_{1/2} -0.18$ V) or $\text{Co}(\text{diansar})^{3+}$ ($E_{1/2} -0.54$ V versus SCE) product can be performed at room temperature with total retention of ligation and chirality; the same experiment performed with chiral $\text{Co}(\text{en})_3^{3+}$ leads to complete loss of chirality. This is a simple test of the unusual stability of the cobalt(II) macrobicyclics which arises as a consequence of the encapsulation of the metal ion within the ligand cavity.


Figure 2. ^1H NMR spectrum in D_2O (N-deuterated) of $\text{Co}(\text{dinosar})^{3+}$.

Literature Cited

- Pass, G., and Sutcliffe, H., "Practical Inorganic Chemistry," 2nd ed., Chapman and Hall, London, 1974.
- Angelici, R. J., "Synthesis and Technique in Inorganic Chemistry," W. B. Saunders Co., Philadelphia, 1969.
- Bell, C. J., "Synthesis and Physical Studies of Inorganic Compounds," Pergamon, Oxford, 1972.

- Krause, R. A., and Magargle, E. A., *J. Chem. Educ.*, **53**, 667 (1976).
- Curtis, N. F., *Coord. Chem. Rev.*, **3**, 3 (1968); Lindoy, L. F., and Busch, D. H., *Preparative Inorg. React.*, **6**, 1 (1971).
- Sargeson, A. M., *Chem. Brit.*, **15**, 23 (1979).
- Harrowfield, J. MacB., Herlt, A. J., and Sargeson, A. M., *Inorg. Synth.*, **20**, 85 (1980).
- Bond, A. M., Lawrance, G. A., Lay, P. A., and Sargeson, A. M., *Inorg. Chem.*, **22**, 2010 (1983).
- Geve, R. J., Hambley, T. W., Harrowfield, J. MacB., Sargeson, A. M., and Snow, M. R., *Inorg. Chem.*, **23**, *J. Amer. Chem. Soc.*, **106**, 5478 (1984).

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