Handylaser sprint

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Índice

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
GB Operating Instructions

Use of the charger

Charger/discharger for nickel cadmium (NiCd) and nickel metal hydride (NiMH) battery packs of 1-10 cells (1.2-12.0 V) with a capacity of 800-7,200 mAh.

Features

- For worldwide use thanks to a switch mode power supply (100-240 V AC) and exchangeable primary plug set
- Microprocessor controlled charging
- Testing cycle at the beginning of the charging in order to identify the number of cells and to identify and report defective battery packs
- Short circuit detection and electronic protection against reversed polarity
- Battery pack charge status at the beginning of the charging is not important
- Microprocessor controlled charge status monitoring during the whole charging time
- Includes safety features such as voltage gradient supervision and -delta V switch off, plus a safety timer
- Optional discharging of the battery pack before use by pressing the discharge button. Automatic switching over to fast charging after discharge
- LED status indication
- Automatic switching over to trickle charge once fully charged

Indicators

LED red „Power“ (1): Steady light shows that the charger is ready for use. It lights once the charger is plugged in.

LED red „Charge“ (2): Steady light indicates the fast charging process after connecting the battery pack.

LED green „Ready“ (3): Steady light indicates the battery pack is fully charged. After approx 2 minutes it switches to flashing, which indicates trickle charging.

Once the battery pack is connected, the testing cycle is initiated and the green LED flashes for around 1 minute.

If it carries on flashing and the red LED „Charge“ does not light, it indicates that the battery pack is not connected properly and maybe has an incorrect polarity.

LED yellow „Discharge“ (4): Steady light indicates the discharging process after pressing the yellow button. The LED „Ready“ flashes for around 1 minute indicating the testing cycle.

Controls and Accessories

Discharge button (5): Press this button for about two seconds to start the optional discharging process of the battery pack.

Attention

Charge only nickel/cadmium (NiCd) or nickel/metal-hydride (NiMH) battery packs. Danger of explosion if other types of batteries are inserted.

Safety Instructions

Do not attempt to open the charger
Keep the charger in a dry place (indoor use only).
In order to avoid the risk of fire and/or electric shock, the charger must be

protected against high humidity and water.

Do not plug in the charger if there are any signs of damage to the housing, mains pins cables or connectors. In case of a defect please return to an authorised service centre.

Keep the charger out of reach of children.

If the safety instructions are not followed, this may lead to damage to the charger or batteries or even to serious injury to the user.

Operation

Connect the charger to the mains: With the exchangeable primary plug set and the electronic power supply (100-240 V AC) the charger can be used worldwide. To change the primary plug, unlock the mechanism on the back of the unit towards the arrow. Attach the right primary plug to the unit until it is clicks in place. Once the charger is connected to the mains, the power indicator lights and the charger is ready for use.

Connect the battery pack to the charger: Usually the battery pack is brought into contact via the output plugs which are enclosed. When connecting the plug to the secondary cable, pay attention to the correct polarity (see figure 6). The red LED „Charge“ lights and indicates the charging process. During the testing cycle the green LED „Ready“ flashes also but goes out after around 1 minute when the test phase is over. After termination of the fast charge the charger switches automatically to trickle charge. The red LED „Charge“ is off and the green LED „Ready“ is constantly on for approx. 2 minutes. After around 2 min. the indicator changes to green flashing. The battery pack can be removed at this time or left connected to keep topped up until use.

If the green LED „Ready“ flashes immediately after connecting the battery pack and the red LED „Charge“ lights up sporadically, this means that the polarity (+/-) at the charging cable is not set correctly. Please reverse the polarity (see figure 6). If both LEDs are still flashing after changing the polarity and the steady light „Ready“ does not light permanently, the battery pack is defective and cannot be recharged at all. In this case the pack has to be replaced.

The optional discharging procedure can be started by pressing the button for discharging (5) for about 2 seconds. The yellow LED „Discharge“ (4) lights and indicates the discharging procedure. During the first minute the green LED „Ready“ (3) flashes but should turn off after the test phase. After discharging, which can in some cases last for several hours, the charger automatically switches over to fast charging.

Environment


Rechargeable batteries are not to be disposed in domestic waste. Return used batteries to your dealer or to a battery recycling collection point.

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

Electromagnetic compatibility – Warning instructions

Precautionary measures

There are special precautionary measures for medical electric devices as far as electromagnetic compatibility is concerned. This device may be used only for the purpose described in the manual, set up and put into operation according to the notes for electromagnetic compatibility.

Influence of mobile and portable communication equipment

High-frequency energy radiation from mobile communication equipment may influence on the operation of the medical electric device. It is not allowed to use those devices (for instance mobile phones and GSM telephones) close to the medical electric device.

Guidelines and manufacturer's declaration – Electromagnetic radiation

The laser device is designed for the operation only with original accessories in ambient conditions indicated below. Client and user have to make sure that the device is operated in such ambient conditions.

Radiated disturbance measurements	Conformity	Electromagnetic ambient conditions – guideline
HF-radiation to CICPR 11	Group 1	The laser device uses HF-energy only for its internal function. Therefore HF radiation is very low. The influence on neighbouring electric devices is unlikely...
HF-radiation to CICPR 11	Class B	The laser device may be used in all kinds of buildings, including residential buildings and those directly connected to mains for the public supply which are also servicing buildings used for residential purposes.
Radiation of harmonic oscillation according to IEC 61000-3-2	Conformar	
Radiation of voltage variations/ Flicker effect according to IEC 61000-3-3	Not applicable	

Guidelines and manufacturer's declaration – Electromagnetic immunity

Immunity tests	IEC 60601-Test level	Conformity level	Electromagnetic ambient conditions – guidelines
Discharge of static electricity (ESD) according to IEC 61000-4-2	± 6 kV Contact discharge ± 8 kV Air discharge	conforms	The floor should be made from wood or concrete or covered with ceramic tiles. If the floor is covered with synthetic material the relative humidity should be minimum 30%
Fast transient electric disturbance variable / Bursts according to IEC 61000-4-4	± 2kV for network leads ± 1 kV for input and output leads	conforms	The quality of the supply voltage should be typical for business or hospital ambient conditions.
Impulse voltage (Surges) according to IEC 61000-4-5	± 1 kV opposed mode voltage ± 2 kV common mode voltage	conforms	The quality of the supply voltage should be typical for business or hospital ambient conditions.
Voltage dips, short interruptions and voltage variations according to IEC 61000-4-11	< 5% U_T (> 95% dip of U_T) for 1/2 period 40% U_T (60% dip of U_T) for 5 periods 70% U_T (30% dip of U_T) for 25 periods < 5% U_T (> 95% dip of U_T) for 5 s	conforms	The quality of the supply voltage should be typical for business or hospital ambient conditions. If the user wants to continue operation in case of an interrupted energy supply it is recommended to feed the laser device from a system without interruption or a battery.
Magnetic field at supply frequency (50/60 Hz) according to IEC 61000-4-8	3 A/m	conforms	Magnetic fields at mains frequency should have the typical values for business and hospital ambient conditions.

Note: U_T - alternating voltage mains before the application of the test level.

Guidelines and manufacturer's declaration – Electromagnetic immunity – for the laser device without life supporting function

Immunity tests	IEC 60601–test level	Conformity level	Electromagnetic ambient conditions – guidelines
HF-disturbance variables (transmission) according to IEC 61000-4-3	3 Veff 150 kHz up to 80 MHz	conforms	Portable and mobile radio sets must not be operated in a distance to the laser device and the leads less than the recommended protective distance what was calculated according to the equation for the transmission frequency. Recommended protective distance: $d = 1,2 \sqrt{P}$ $d = 1,2 \sqrt{P}$ for 80 MHz up to 800 MHz $d = 2,3 \sqrt{P}$ for 800 MHz up to 2,5 GHz with P as power rating of the transmitter in Watt (W) according to data from the manufacturer of the transmitter and as recommended protective distance in meter (m). According to a local ^a investigation the field intensity for all frequencies of immobile radio transmitters should be less than the conformity level ^b . Interference is possible in the surroundings of devices carrying the following signs. 
HF-disturbance variables (radiation) according to IEC 61000-4-3	3 V/m 80 MHz up to 2,5 GHz	conforms	

Note 1 At 80 MHz and 800 MHz the higher frequency range is applied.

Note 2 These guidelines are not applicable in all cases. The propagation of electromagnetic quantities is influenced by the absorption and reflection of buildings, items and human beings.

^a It is not possible to predict theoretically the field strength of immobile transmitters like base stations of mobile phones and mobile land radio sets, amateur radio stations, AM and FM television and radio broadcasting transmitters. In order to determine electromagnetic ambient conditions with regard to immobile transmitters a study of the location should be taken into consideration. If the field strength measured in the place where the laser device is operated exceeds the above mentioned conformity level the laser device should be watched to proof its function in accordance to the requirements. Unusual performance data may require additional measures, for instance change of orientation or different location of the laser device.

^b Above the frequency range from 150 kHz to 80 MHz the field strength should be less than 3 V/m.

Recommended protective distance between portable and mobile HF-telecommunication devices and the laser device without life supporting function

Recommended protective distance between portable and mobile HF-telecommunication devices and the laser device			
The laser device is designed for the operation in electromagnetic ambient conditions where HF-disturbance variables are controlled. Client or user of the laser device can help to avoid electromagnetic disturbances by keeping the minimum distance between portable and mobile HF-telecommunication devices (transmitters) and the laser device as indicated below, dependent on the power output of the communication device.			
Power rating of the transmitter W	Protective distance dependent on the transmitting frequency m		
	150 kHz up to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz up to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz up to 2,5 GHz $d = 2,3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
If the maximum power rating for certain transmitters is not given in the above table the recommended protective distance d in metres (m) can be determined by using the equation indicated in the respective column, where P is the maximum power rating of the transmitter in Watt (W) according to the manufacturer of the transmitter.			
Note 1 At 80 MHz and 800 MHz the higher frequency range is applied			
Note 2 These guidelines are not applicable in all cases. The propagation of electromagnetic quantities is influenced by the absorption and reflection of buildings, items and human beings.			

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REIMERS & JANSSEN GmbH