Dear EvenCare G2 Owner,

Thank you for choosing the EvenCare G2 Blood Glucose Monitoring System to help you in monitoring your diabetes. We designed the EvenCare G2 Blood Glucose Monitoring System to be accurate and easy to use. The added convenience of the voice feature can be used as a verbal aid for using the meter and hearing test results by visually impaired users. This manual contains all of the information needed to use and to maintain your new blood glucose meter. Please read it carefully before use.

If you need assistance, please call Medline's Technical Service Center at 1-866-274-2131 Monday through Friday between 8:00 a.m. and 5:00 p.m.(Central Time). Contact your Healthcare Provider at all other times.

Intended Use

The EvenCare G2 Blood Glucose Monitoring System is intended for the quantitative measurement of glucose in fresh capillary whole blood samples drawn from the fingertips, forearm, or palm. Testing is done outside the body (In Vitro diagnostic use). The meter includes voice functionality to assist visually impaired users. It is intended for both lay use by people with diabetes and in a clinical setting by healthcare professionals, as an aid to monitoring levels in Diabetes Mellitus. It is not intended for the diagnosis or screening of diabetes or for neonatal use.

IMPORTANT:

- Alternate Site Testing (AST) with the EvenCare G2 Blood Glucose Monitoring System can be used only during steady-state blood glucose conditions.
- Dehydration Severe dehydration may lead to inaccurate blood glucose test results. If you suspect you are severely dehydrated, contact your healthcare professional immediately.
- Hematocrit range A hematocrit range that is higher than 55% or lower than 30% can cause inaccurate blood glucose test results.

Table of Contents

Understanding Your New Blood Glucose System	
The EvenCareG2 Meter	6
Meter Display Screen	8
The EvenCareG2 Test Strip	10
EvenCareG2 Glucose Control Solution (Recommended)	11
Setting Up Your New System	
The Set Up Screen	12
Setting the Clock	13
Setting the Alarms	16
Voice and Beeper Function	18
Language Selection	19
Voice Function Table	20
Control Solution Testing	
Performing a Control Solution Test	25
Control Solution Trouble Shooting	29
Testing Your Blood Glucose	
Preparing Your Lancing Device	30
Important Information on Alternate Site Testing (AST)	32
Performing Blood Glucose Test	34
Understanding Your Blood Glucose Test Results	38
Viewing Your Meter's Memory	39
Caring for the Meter	
Changing the Battery	42
Cleaning the Meter	43

Ξ	venCareG2 Table of Contents		
	Storage and Precautions	43	
	Solving Problems	44	
	Product Warranty	47	
	Specifications	48	

Understanding Your New Blood Glucose System

The EvenCare G2 Blood Glucose Monitoring System contains the following items if it arrives as a meter only system:

- EvenCareG2 Meter
- User's Guide
- Carrying Case
- Log Book

The EvenCare G2 Blood Glucose Monitoring System contains the following items if is arrives as a Starter Kit:

- EvenCareG2 Meter
- Lancing Device
- 10 Lancets
- 10 Test Strips
- User's Guide
- Carrying Case
- Log Book

The EvenCareG2 Meter

Test Strip Port Insert test strip (see page 34) to turn on the meter.

Meter Display

Shows your blood glucose test results, memory values, averages, and other messages.

Enter (O) and $\blacktriangle \bigtriangledown$ buttons Press and hold Enter (O) for 2 seconds to turn on or off the Meter, or press Enter (O) to go into setup mode when not testing (see page 12)

Press ▲ or ▼ to enter control mode with test strip inserted (see page 26), during meter setting (see page 12), or to navigate the stored test values and averages (see page 39)

This function is not yet available.

Meter Display Screen

Displays the date, time, am/pm (in 24-hours format, see page 12), or type of day averages (7, 14, or 30, see page 39)

BEEPER Set sound volume (see page 18).

TEMPERATURE Appears when it is either too hot or too cold to test (outside the ranges of 50°F~104°F).

-CLOCK Indicates the meter is ready to set up time and date (see page 13).

VOICE Indicates the voice feature is on (see page 18).

-UNITS Unit of measurement for your blood glucose is ma/dL.

BLOOD DROP and ARROW

Indicate the meter is ready for testing (see page 34).

The EvenCareG2 Blood Glucose Test Strip

The EvenCare G2 Test Strips are glucose specific, biosensor-based test strips that require very little blood sample to quickly test glucose in fresh capillary whole blood. The test results are plasma-calibrated for easy comparison to lab results, however results may be 10% to 15% higher than results obtained from whole-blood referenced blood glucose assays. The test strip is also equipped with under-fill detection to alert you when there is not enough blood to perform a test, so you can be assured that every reading you get is an accurate and meaningful result.

IMPORTANT:

- Immediately close the vial cap of the test strip bottle tightly after each use. Keep the unused test strips in its original bottle.
- · DO NOT leave any test strips outside the bottle while not in use.
- DO NOT reuse test strips. Test strips are for single use only.
- Carefully discard used test strips and lancets in proper waste containers.

EvenCare G2 Glucose Control Solution (Recommended)

Call your local supplier or pharmacy to order control solution. If they do not stock it, please contact us at 1-866-274-2131 between 8:00 am and 5:00 pm (central time), Monday through Friday.

EvenCare G2 Glucose Control Solution is the standard glucose concentration solution required to perform a Glucose Control Solution Test (See Control Solution Testing). It is available at 2 levels. EvenCare G2 LOW Glucose Control Solution helps you to validate the performance of your system at a low range blood glucose reading, while HIGH solution helps you to validate your system performance at a high range blood glucose reading. Run the control test to make sure the test strips and the meter are working together properly, to practice testing procedure, and when using a new lot of EvenCare G2 Blood Glucose Test Strips.

Setting Up Your New System

The Setup Screen

The Setup Screen of the EvenCare G2 Meter lets you take advantage of the advance functions of this meter. Simply use $\blacktriangle \lor$ buttons to navigate to the desired function icon, and press the Enter (\diamondsuit) button to enter the function when the icon is flashing. Navigate to End at the display screen to exit the Setup Screen. The Meter will also announce each function when the voice feature is on.

memPress

Enter (ϕ) to enter Memory Recall to retrieve stored results. Your meter can memorize up to 300 memories with date and time. (see page 39)

✓ ④ ………Press Enter (也) to enter PC Link Mode for data transfer onto a PC*

*This function is not available yet.

- ⁽∭)Press Enter (⁽⊕) to enter Time and Date Setup (see page 13).
- ⁽¹⁾Press Enter (⁽Φ) to enter Alarm Setup (see page 16). Your meter can have as many as 3 alarm settings.

IIPress Enter (Φ) to enter Volume Setup (see page 18).

ඳිාPress Enter (也) to enter Language Selection (see page 19).

EndExit Setup Screen and return to testing.

Setting the Clock

You will need: $\sqrt{}$ Your EvenCare G2 Meter

NOTE

If date and time on the EvenCare G2 Meter has not been set, test or control result will not be stored in memory.

Step 1: Press and hold Enter (${\mbox{$\varPhi$}}$) to turn on the meter.

Step 2:

Step 3-Set the Year

The current year will flash. If the voice function is on, the meter will say, "Set Year." press \blacktriangle or \forall to select the correct year. Press Enter (0) to confirm your choice and advance to set the month.

EvenCareG2 Setting Up Your New System

Step 4 - Set the Month:

The current month will flash. If the voice function is on, the meter will say, "Set Month." press \blacktriangle or \blacktriangledown to select the correct month. Press Enter (ϑ) to confirm your choice and advance to set the date.

Step 5 - Set the Date:

The current date will flash. Press $\blacktriangle \text{ or } \forall$ to select the correct day. If the voice function is on, the meter will say, "Set Day." Press Enter (\circlearrowright) to confirm your choice and advance to set time format.

Step 6 - Set Time Format to 12-hour or 24-hour: The current time format will flash. Press \blacktriangle or \blacktriangledown to select for either 12-hour format or 24- hour format. Press Enter (0) to confirm your choice and advance to set hour.

Step 7 - Set the Hour:

The current hour will flash. If the voice function is on, the meter will say, "Set Hour." press \blacktriangle or \lor to select the correct hour. Press Enter (\oint) to confirm your choice and advance to set minutes.

Step 8 - Set Minute:

The current minutes will flash. If the voice function is on, the meter will say, "Set Minutes." Press \blacktriangle or \blacksquare to select the correct minutes. Press Enter (ϑ) to confirm your choice and finish setting the clock.

NOTE:

- Anytime during set up, inserting a EvenCare G2 Test Strip will begin testing. Any changes made so far are stored.
- · You may need to reset time after changing batteries.

After date and time has been set, the meter will announce the date and time if the voice feature is on. Begin testing by inserting an EvenCare G2 Blood Glucose Test Strip (see page 30, Testing Your Blood Glucose), or press and hold Enter (0) to turn off the meter, or see the next section on how to set up the alarm feature on your meter.

EvenCareG2 Setting Up Your New System

Setting the Alarms

You can set up to three alarms on the EvenCare G2 Meter. You must set the clock before you can set the alarms. When an alarm sounds, the meter will turn on, and a melody will sound for 30 seconds. Pressing Enter (0) or inserting a test strip will silence the alarm. You can set up one or all three alarms.

You will need: $\sqrt{}$ Your EvenCare G2 Meter

Step 1: Press and hold Enter (υ) to turn on the meter.

Step 2: Press Enter (Φ) to enter setup screen.

Step 3: Press \blacktriangle or \lor until 0 starts to flash. If the voice function is on, the meter will say, "Set testing Alarm." Press Enter (0) to enter alarm setting.

Step 4: The display screen will indicate alarm 1 status (the default is OFF). You can turn it on or off by pressing \blacktriangle or \blacktriangledown

Setting Up Your New System EvenCareG2

Step 5:

Step 6:

The current minutes will flash. Press $\blacktriangle or \forall$ to select the correct minutes. Press Enter (\diamond) to confirm your choice and advance to set alarm 2.

Step 7:

Set alarm 2 and 3 following steps 4 through 6. After you have finished setting up alarm 3, press Enter (ϕ) to exit alarm setup.

NOTE:

- Alarm will not sound during testing.
- Alarm will follow the time format (12-hour or 24-hour) you choose when setting up the clock.
- · Alarm settings will not be erased when changing batteries.

Begin testing by inserting a EvenCare G2 Blood Glucose Test Strip (see page 30, Testing Your Blood Glucose), or press and hold Enter (\odot) to turn off the meter, or see the next section on how to set up the voice feature on your meter.

Voice and Beeper Function

Your meter has a voice function that can be used as a verbal aid for using the meter and hearing test results. The default setting for the voice function is on, and can be switched off to conserve battery power. When the voice function is off, the meter will beep to alert you to important messages.

You will need: $\sqrt{}$ Your EvenCare G2 Meter

Setting Volume

Step 1: Press and hold Enter (0) to turn on the meter.

Step 2:

Press Enter (ϕ) to go into setup screen. Press \blacktriangle or \forall until \triangleleft starts to flash. Press Enter (ϕ) to enter volume setting.

Step3:

The display screen will show the current sound volume (default is 2).

Press ▲ or ▼ to adjust the sound volume from 1 (lowest) to 5 (highest). A melody will sound

according to the sound volume you have selected. Setting to 0 will turn off the voice function. The meter will instead beep to alert you to important messages. Once you have determined the suitable volume, press Enter (Φ) again to exit volume setting.

Language Selection

Step 1:

Press and hold Enter (ϕ) to turn on the meter.

Step 2: Press Enter ($\boldsymbol{\vartheta}$) to enter setup screen.

Step 3: Press ▲or▼ until இ starts to flash. The meter will announce "Language Selection."

Step 4: The default language (1 on the display screen) is announced. Press (\bullet) to change the language. Press \blacktriangle or \triangledown to select another language. Press (ϕ) when you have finished selecting the language and exit language selection. English and Spanish are the two options.

NOTE:

- Any time during set up you can insert an EvenCare G2 Test Strip and begin testing. Any changes made so far are stored.
- Your choice of language and volume settings will not be erased when changing batteries.

Begin testing by inserting an EvenCare G2 Blood Glucose Test Strip (see page 30, Testing Your Blood Glucose) or press and hold Enter (Φ) to turn off the meter.

Voice Function Table

The voice function can be used as an aid for you to hear the test results and use the meter. The following table lists WHEN the meter will speak and WHAT the meter will say.

Testing blood and control solutions (see page 25 , Control Solution Testing, and page 30, Testing your blood glucose)

WHEN the meter will speak	WHAT the meter will say
Turned on.	Greeting Tone will sound.
Prompt user to insert test strip.	"Please insert test strip"
Test strip is inserted and the meter is ready to test (with flashing arrow and blood drop).	"Blood glucose test. Please apply blood to the test strip."
Blood Glucose test is completed and the result is shown on screen.	"Your blood glucose reading is (number) milligram per deciliter."
Blood Glucose test result is BELOW the measurement range of 20 mg/dL.	"Your blood glucose reading is below 20 mg/dL."
Blood Glucose test result is ABOVE the measurement range of 600 mg/dL.	"Your blood glucose reading is above 600 mg/dL."

WHEN the meter will speak	WHAT the meter will say
Temperature is out of range (be- tween 50°F~104°F).	"Temperature is out of range"
When the meter is turned off.	Melody will sound.

Test Strip is inserted and press \blacktriangle or \blacktriangledown (CTL showing next to the test strip icon) .

"L1 control solution test. Please apply control solution onto the test strip."

"Blood glucose test. Please apply blood onto the test strip."

"L2 control solution test. Please apply control solution."

Setting Meter Functions via the Setup Screen

(see page 12, setting up your new system)

Pressing Enter (\underline{o}) to enter setup screen (after the meter has been turned on)

WHEN the meter will speak	WHAT the meter will say
During date and time setup, setting up the Minutes	"Set Minutes"
After finishing setting up date and time	"The date and time is (year) (month) (date) (hour) (minutes)"
Setting alarm 1, 2, or 3	"Set Alarm (number)"
When sending data to PC *This function is not available yet.	"Sending data to PC"
Finish sending data to PC *This function is not available yet.	"Data transmission complete"

Other warning messages (see page 44, Solving Problems)

WHEN the meter will speak	WHAT the meter will say
Meter date and time has not been set	"Please set date and time"
Meter failure	"Error 1 Meter Failure"
Used test strip	"Error 3 used test strip"
Not enough blood or control solution	"Error 4 not enough sample"

Control Solution Testing

The purpose of the control solution testing is to make sure the Even-Care G2 Meter and the EvenCare G2 Test Strip are working properly. You should perform control solution testing when:

- · Using the meter for the first time
- · Using a new bottle of EvenCare G2 Blood Glucose Test Strips
- · You left the test strip bottle cap open for awhile
- · You dropped the meter
- · You suspect your meter and test strips are not working properly
- The blood glucose test results do not reflect how you feel
- · You want to practice the testing procedure

IMPORTANT:

- Use only the EvenCare G2 Glucose Control Solutions with the EvenCare G2 Test Strips. Other brands of control solution will produce inaccurate results.
- Always check the expiration date. DO NOT use control solutions if they are expired.
- DO NOT FREEZE. Store the control solutions at room temperature.

You will need:

- √ EvenCare G2 Control Solution (Low and High)
- √ Your EvenCare G2 Meter
- √ A new EvenCare G2 Blood Glucose Test Strip

Step 1:

Take out an EvenCare G2 Test Strip from the test strip bottle and close the bottle immediately. Insert the test strip to turn on the meter.

IMPORTANT:

- Check the expiration date printed on the test strip bottle and its package. Do not use expired test strips.
- · Use each test strip immediately after removing it from the bottle.
- Do not use wet or damaged test strips.
- Keep away from direct sunlight and heat. Store the test strip bottle in a dry, cool place.
- Record the "date opened" on the bottle label. Discard the bottle and any remaining test strip after 6 months from date of opening.
- Test under the temperature range of 50°F to 104°F. Your meter will not work outside these ranges.
- · Must insert test strip properly.
- Must see a flashing blood drop icon if the test strip has been fully inserted.

Step 2:

Wait until the flashing blood drop and arrow icons appear on the meter display screen. The meter will also announce, "Please apply blood onto the test strip" if you have the voice feature turned on. Press \blacktriangle or \blacksquare to enter L1(LOW) control solution testing.

Step 3:

ctl icon will appear next to the test strip icon and L1 will appear on the meter display screen. The meter will also announce "L1 control solution test, please apply control solution onto the test strip." If you decide not to do a control solution test, Press ▲ or ▼until ctl disappears from the meter display.

Step 4:

Squeeze a drop of control solution onto a clean, dry, nonabsorbent surface. Do not apply control solution to the test strip directly from the bottle. Replace the bottle cap on the control solution bottle immediately after use.

DO NOT touch the Sampling End of the test strip with the bottle of control solution. See illustration of Sampling End below.

🕬 EvenCare.G2. 📖

Gently touch the drop of the control solution to the tip of the test strip. The meter will beep indicating it has received the control solution.

Step 6:

The screen will start to count down. After 6 seconds, the control test result will appear on the meter's display screen.

Step 7:

Compare the reading on the screen to the Low range printed on the test strip bottle or its package. If the reading falls outside the Low range printed on the test strip bottle or its package, see Control Solution Trouble Shooting, on page 29.

NOTE:

- Control solution testing results will be stored into the meter memory and indicated by ct icon.
- Control solution testing results will not be used for calculating averages.

Sampling End

EvenCareG2 Control Solution Testing

Step 8:

Removed the used test strip. The meter will automatically turn off. Discard the used test strip.

IMPORTANT:

Do not reuse test strips.

Follow steps 3 through 8 to do High control solution testing and compare the result with the High range printed on the test strip bottle.

Control Solution Trouble Shooting

If your control solution testing is out of range (too high or too low), it may be caused by the following:

	Possible Causes		What you can do
•	Wrong brand of Control Solution being used	•	Make sure you are using EvenCare G2 Control Solution.
•	Control solution not at room temperature	•	Make sure the testing environ- ment is between 50°F to 104°F.
•	Expired or contaminated control solution or damaged test strip.	•	Check the expiration and open bottle date on both the control solution and test strips. Repeat testing using a new test strip. If the result is still out of range, use a new bottle of control solution and retest.
•	Meter malfunction	•	If the problem persists, please call Medline's Technical Service Center at 1-866-274-2131 Monday through Friday between 8:00 a.m. and 5:00 p.m.(Central Time). Contact your Healthcare Provider at all other times.

Testing Your Blood Glucose

You will need:

- √ Your EvenCare G2 Meter
- $\sqrt{}$ A new EvenCare G2 Blood Glucose Test Strip
- $\sqrt{}$ Lancing device with a clean, unused lancet (or safety lancets)
- $\sqrt{}$ Clear cap for Alternate Site Testing (AST) on palm and forearm

Before you begin, make sure:

- To set up your meter properly and run a control test. See Setting Up Your New System on page 12, and Control Solution Testing on page 25 for details.
- Wash your hands and the testing site thoroughly with soap and warm water, and dry well.
- You are testing in an area between the temperature ranges of 50°F~104°F. Your meter will not test outside of this range and will display thermometer icon with verbal indication if the voice feature is turned on. Move the meter into an area that is between 50°F~104°F, and let it sit for 10 to 15 minutes before testing again.

Preparing Your Lancing Device (for Home Testing)

Step 1:

Unscrew the adjustable cap of the lancing device and insert the lancet by pushing down firmly until it is fully seated.

Step 2:

Twist the protective cover off the lancet. Twist the adjustable cover clockwise back onto the lancing device.

NOTE:

- For fingertip testing, use the opaque adjustable cap.
- For alternate site testing on palm and forearm, use the clear adjustable cap.

Step 3:

Adjust the puncture setting on the adjustable cap for the puncture depth level (1 is the lightest and 5 is the deepest).

Step 4:

With one hand holding the adjustable cap, pull back the lancing device with your other hand until you hear a click sound. The lancing device is now ready for blood sampling.

IMPORTANT:

Use a new sterile lancet every time you test to avoid cross contamination. If alcohol wipes are used to cleanse the fingers, make sure the area is dry before the blood sample is obtained.

Important Information on Alternate Site Testing (AST)

The EvenCare G2 Blood Glucose Monitoring System can test for blood glucose from areas other then your fingertip such as palm and forearm (alternate site testing, or AST). Alternate site testing can be less painful then fingertip testing, but because of the physiological differences between your fingertip and palm and forearm¹; alternate site test results maybe significantly different

then results from fingertip testing under certain conditions.

You should consult with your doctor or healthcare professional before using alternate site testing.

DO AST ONLY in the following intervals:

- In a pre-meal or fasting state (more than 2 hours since the last meal)
- 2 hours or more after taking insulin.
- · 2 hours or more after exercise.

Alternate Site Testing SHOULD NOT be used when:

- You are Hypoglycemic unaware (not able to tell if you have low blood sugar)
- · Within 2 hours of a meal, exercise, or medication
- · You will be operating machinery or driving a car.
- During times of stress.
- You are sick
- · You think your blood glucose is low.
- Your AST results do not match the way you feel.
- · You are testing for hyperglycemia.
- · Your routine glucose results are often fluctuating.

NOTE:

- Rubbing forearm until warm before lancing will help minimize the difference with fingertip test results. The palm (at the base of the thumb) does not require rubbing before testing.
- If results from alternate site testing do not agree with how you feel, use fingertip testing instead.
- Jungheim K., Koschinsky T.: "Risky delay of hypoglycemia detection by glucose monitoring at the arm." Diabetes Care 2001:24(7):1303-04.

Performing a Blood Glucose Test

Step 1:

Take out an EvenCare G2 Test Strip from the test strip bottle and close the bottle immediately. Insert the test strip to turn on the meter.

IMPORTANT:

- Check the expiration date printed on the test strip bottle and its package. Do not use expired test strips.
- · Use each test strip immediately after removing it from the bottle.
- · Do not use wet or damaged test strips.
- Keep away from direct sunlight and heat. Store the test strip bottle in a dry, cool place.
- Record the "date opened" on the bottle label. Discard the bottle and any remaining test strip after 6 months from date of opening.
- Must see a flashing blood drop icon if the test strip has been inserted properly.

Step 2:

Wait for a flashing blood drop and an arrow icon to appearon the meter display screen. The meter will also announce, "Please apply blood onto the test strip" if you have the voice feature turned on.

Step 3: Obtain blood sample. For Fingertip Testing:

Hold the lancing device (use the opaque adjustable cap or safety lancet) against the side of your fingertip and press the release button to create a puncture.

TIP:

- Gently massage your hand and finger toward the puncture site to form a drop of blood. Do not "milk," or squeeze around the puncture site.
- Lance the side of your fingertip to avoid soreness. To avoid calluses, choose a different lancing site each time.

EvenCareG2 Testing Your Blood Glucose

For Testing on Palm and Forearm:

Hold the lancing device (use the clear adjustable cap) against the test

site (base of palm or forearm). Avoid a test site that has veins, moles, excessive hair, or bone.

Press the lancing device firmly against the test site and press the release button, do not lift up the lancing device.

For forearm test sites, apply and release pressure ("pumping") for a few seconds while keeping the lancing device in constant contact with the skin. Palm testing (at the base of the thumb) does not require pumping. When the blood is approximately () size, lift the lancing device straight up without smearing the blood.

NOTE:

- Rubbing forearm until warm before lancing will help minimize the difference with fingertip test results. The palm (at the base of the thumb) does not require rubbing before testing.
- If results from alternate site testing do not agree with how you feel, use fingertip testing instead.
Step 4:

Gently bring the test strip and touch the drop of blood at a slight angle. The test strip acts like a straw to pull the blood in. Keep the test strip in the blood drop until the meter beeps to indicate the test strip has enough blood to test.

Step 5:

The screen will start to count down. After 6 seconds, your glucose testing result will appear on the meter display screen and will be announced verbally if you have the voice feature turned on. Your test result will be automatically stored into the meter memory.

CAUTION:

 If you see "HI" or "LO" displayed, your blood glucose level may be above 600 mg/dL or below 20 mg/dL. Test again using fingertip testing, DO NOT test on forearm. If you still receive the same result, call your physician or healthcare professional immediately.

Step 6:

Properly discard the used test strip. The meter will automatically shut off.

Step 7:

Remove the used lancet from your lancing device according to instructions and discard into proper waste container.

IMPORTANT:

 Used lancets and strips are biohazardous materials and can transmit bloodborne diseases. Dispose of them according to local government regulation to avoid injury or contamination.

Understanding Your Blood Glucose Test Results

Blood glucose values will vary depending on food intake, medication, health, stress, and exercise. The ideal range for adults with and without diabetes should be¹:

- · 70 to 100 mg/dL before meals, and
- · 70 to 140 mg/dL after meals

It is important to consult with your physician or healthcare professional to determine an appropriate target range for you.

What to do if you get a high or low reading

If the Meter displays results that are HI or LO, or you get a result that is more than 250 mg/dL or below 50 mg/dL AND you feel ill:

• Treat your diabetes according to the instruction from your doctor and/or consult your healthcare provider.

If the Meter displays results that are "HI" or "LO," or you get a result that is more than 250 mg/dL or below 50 mg/dL AND you DO NOT feel ill:

- Test your meter with a control solution, refer to page 25, Control Solution Testing.
- · Test again using fingertip.

If you still get a high or low reading, contact your health care professional.

NOTE:

- Inaccurate results may occur in severely hypotensive individuals or patients in shock.
- Inaccurate low results may occur for individuals experiencing a hyperglycemic-hyperosmolar state, with or without ketosis.
- · Critically ill patients should not be tested with blood glucose meters.

Reference:

1. American Diabetes Association. Standards of Medical Care in Diabetes—2008. Diabetes Care. 2008;31:S12–S54.

Viewing Your Meter's Memory

Your EvenCare G2 Meter automatically stores up to 300 test and control results with date and time. The meter also provides 7, 14, and 30 day averaging to help track your blood glucose trend. You must set up the date and time on your meter properly to use the memory and day averaging function, see Setting Up Your New System on page 12. The meter will not memorize test or control results if the date and time are not set.

Viewing Stored Readings and 7, 14, 30 Days Average from-Memory

You will need: $\sqrt{}$ Your EvenCare G2 Meter

Step 1: Press and hold Enter () to turn on the meter.

Step 2: Press Enter () to enter setup screen.

Step 3:

Press ▲or▼ until **MEM** starts to flash. If the voice function is on, the meter will say "memory recall." Press Enter () to enter.

MEM will appear at the upper left corner of the display screen together with your most recent testing result. Press \blacktriangle to see the next result, and the next, and so on. Holding \blacktriangle will flash through the results.

END will appear on the display screen after the last recorded test result if less then 300 readings are stored. You can press ▼ to go back viewing from the oldest to the most recent readings.

Step 4.

When you reach the most recent result, continue to press ▼ to start viewing day averages. The meter will display the 7 day's average, press ▼ again to see 14 day's average, and press ▼ again to see 30 day's average.

The number of results included in the averages is shown right next to the icon. Control solution results are not included in the day average calculations.

NOTE:

- Control testing results are flagged by ctl on the bottom left corner of the display screen in addition to **MEM**.
- When the meter memory is blank, the display screen will show nil.
- When there are no day average data available, the display screen will show 4 dashes.
- If the meter memory is full, adding a new test result will cause the oldest one to be deleted.
- Memorized test results and day averages will not be erased when changing batteries.
- You can begin testing anytime by inserting a test strip.

EvenCareG2 Caring for the Meter

Caring for the Meter

Caring for your EvenCare G2 Meter is easy. Follow these simple guidelines to keep your meter working properly.

Changing the Battery

Batteries need to be changed when the meter display.

icon appears on the

Materials you will need:

- ✓ Two Alkaline AAA Batteries
- √ Your EvenCare G2 Meter

Step 1.

Press firmly on the battery cover and slide in the direction of the arrow.

Step 2. Insert new batteries . Slide the battery cover back into the meter firmly.

NOTE:

- Meter will prompt you to reconfirm the date and time when turned on. See page 13 for Setting the Clock.
- The date, stored results, alarm settings, volume setting, and language choice stay in memory when the batteries are changed.

Cleaning the Meter

- If the meter gets dirty, use a moist (NOT WET) lint-free cloth dampened with a mild detergent.
- Do not get water inside the meter or hold it under running water.
- · Do not use glass or household cleaners on the meter.
- · Do not try to clean the test strip holder.

Storage and Precautions

- Handle the Meter with care- severe shock, such as dropping the meter, could damage the electronics.
- The meter and test strips are designed to be used within the temperature ranges between 50°F to 104°F.
- Avoid leaving the meter in extremely hot or cold place, such as near a heat source or in an extremely hot or cold car.
- Do not store or use the meter or test strips in high humidity levels, such as bathroom or kitchen.
- Always close bottle cap immediately after removing a test strip and make sure it's closed tightly.
- Do not take the meter apart as doing so will void the warranty. If there are technical problems or questions, please call Medline's Technical Service Center at 1-866-274-2131 Monday through Friday between 8:00 a.m. and 5:00 p.m.(Central Time). Contact your Healthcare Provider at all other times.

Solving Problems

This section details the significant display screen messages and error codes you will encounter when using your EvenCare G2 meter and EvenCare G2 test strips.

What	What the "Đ roo rSee Meter Failure"	What it Timeterrestarys electronic or test strip is damaged	 What Should Remove/leabat/eny Brod turn on the meter again. Remove the test strip and insert a new test strip again. If the problem persists, please call 1-866-274-2131 Monday through Friday, 8:00 a.m 5:00 p.m.(CST). Contact your Healthcare Provider at all other times.
Err3	"Error 3 used test strip"	Used or damaged test strip	 Repeat the test with a new test strip. Remember to add sample only when you see the prompt on the display screen or the meter verbally prompts you to add sample.
Erry	"Error 4 not enough sample"	Not enough sample on the test strip to start	• Remove the Test Strip and repeat the test with a new test strip. See Testing Your Blood Glu- cose on page 30.

EvenCareG2 Table of Contents

E N	"Low Battery" Vhat itWhat Sh /leansYou Do	Batteries are low ould You See	Change the battery ac- cording to page 42, Changing the Battery.
	"Tempera- ture is out of range"	Temperature out of range.	Move the meter into an area that is between $50^{\circ}F$ ~104°F, and allow 10 to 15 minutes for it to reach the new temperature.
H 12-18 12-12m	"Your blood glucose reading is above 600 mg/dL."	Test result is higher than 600 mg/dL	Wash and dry your hands and repeat the test on your fingertip with a new test strip. If the result is still "HI," contact your physician or healthcare professional immediately.
L Ø 12-18 12:12m	"Your blood glucose reading is below 20 mg/dL."	Test result is lower than 20 mg/dL	Wash and dry your hands and repeat the test on your fingertip with a new test strip. If the result is still "LO," contact your physician or healthcare professional immediately.

mem	Not Applicable What Vhat Should You See	There are no memorized results in the meter What the meter says	 Check if the date and time on your meter is set up. See Setting Up Your New System on page What it Start testing your blood glucose on page 30.

Product Warranty

The EvenCare G2 Blood Glucose Meter is warranted to be free of defects in workmanship and materials for a period of five (5) years from the date of purchase. Medline's liability for warranty claims is limited to repair or replacement, and in no event shall Medline be liable for indirect or consequential damages, or for any loss arising from misuse, improper use, abuse, neglect, unauthorized repair or modification.

This warranty is void and of no force and effect in the event of product misuse, improper use, abuse, neglect, unauthorized repair or modification. This warranty specifically excludes EvenCare G2 Test Strips and EvenCare G2 Control Solution.

THIS WARRANTY IS THE SOLE AND EXCLUSIVE WARRANTY TO THE EXCLUSION OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, ALL OF WHICH ARE WAIVED AND DISCLAIMED.

All warranty claims must be directed to Medline's Technical Service Center at 1-866-274-2131. This warranty is extended only to the original purchaser of the meter.

Specifications

Test Strips:	EvenCare G2 Blood Glucose Test Strips
Test Range:	20 - 600 mg/dL
Calibration:	Plasma
Test Time:	6 Seconds
Sample Size:	Larger than 1uL
Hematocrit Range:	30-55%
Display Type:	Large LCD screen
Memory:	300 blood test results with date and time
Result Averaging:	7, 14, and 30 days averaging
Dimension:	93L x 58W x 20.5H (mm)
Weight:	79g
Battery:	2 Alkaline AAA Battery
Battery Life:	1000 tests of continuous use or 1 year
Automatic Power-off:	After 2 minutes of nonuse
Operating Temperature:	50°F to 104°F.
Relative Humidity:	Less than 85%
Storage Condition:	Meter at 39°F ~ 122°F

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