ENDO-PASS

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ENDO-PASS

Endodontic Penta-Aggregate Silicate Sealer

MULTIMINERAL PORTLAND CEMENT FOR DENTAL PURPOSES

Product for dental use
Keep away from children
Date of the specifications: 03/2015

AVAILABLE PACKAGES:

1 blister of 12 capsules of 900 mg/cad. (COD. EP 1)

CHARACTERISTICS

Endo-PASS is made of a mixture of mineral components suitably mixed and processed in order to obtain a cement paste easily workable and quick setting times. Endo-PASS must be mixed with bidistilled water.

COMPOSITION

Endo-PASS is constituted by a mixture of calcium silicates (tricalciumsilicate, dicalciumsilicate, tricalciumaluminate) and phyllosilicates (smectite and hydrotalcite), zirconium dioxide and barium sulfate.

INDICATIONS

- Permanent root canal fillings in combination with gutta-percha
- Endodontic surgery (retrograde filling)
- Perforation
- Internal and external resorption
- Apexification and Apexogenesis
- Retrograde fillings
- Pulp capping (direct and indirect)
- Repair of decayed tissue deep (sandwich technique with adhesives and composites)

CONTRAINDICATIONS

Hypersensitivity to the components. There are no known adverse reactions to the product in clinical use.

WAY OF USE

Endo-PASS must be handled with care. It is good to respect the doses of cement and bidistilled water required by the manufacturer and respect the time of mixing and processing (ratio cementitious powder / bidistilled water preferably 2: 1, modified according to the needs of the intervention). The cement must be tooled incorporating the product in bidistilled water and making sure to incorporate all the particles in order to obtain a creamy consistency.

In the presence of a non-permeating area particularly wet, as, for example, in the case of retrograde filling, the amount of bidistilled water may be lower.

Endo-PASS may be used for all situations where it is necessary a hardening in presence of moisture. For this reason, it is advisable to maintain a moist field. The cements consist of mineral aggregates (type MTA), require a curing time of about 12- 24 hours, resulting in at least two operational sessions. The new cement EndoPASS, hardens in 7-10 minutes.

The powder is white and makes it suitable for cosmetic sectors. The first step is the mixing of the powder, which must be done in neutral pH environment, in the absence of ions, therefore with bidistilled water pure.

Anesthetics should be avoided, because they tend to usually strangers or starches in basic solutions; should be avoided physiological solution of NaCl, for the presence of chlorine ions and sodium; should be avoided water in the dentist armchair for the presence of hypochlorous acid.

The bidistilled water temperature influences the rate of hardening. The powder-liquid mixing occurs in a sheet of plexiglass or glass plate. To mix use a spatula for cement.

Prepare the powder on one end of the floor of mixing and two or three drops of liquid on the other end. It starts mixing incorporating the cement little by little in the liquid, taking care to mix and simultaneously press the mixture. Once you reach the required viscosity (preferably a pasty appearance at low flow), the cement can be applied in the operative area.

SPECIFICATIONS

Endo-PASS is a mineral material completely biocompatible thanks to its bioactivity promotes the repair of the dentin.

In fact, when used as a material for direct and indirect pulp capping, it induces a reaction that leads to the formation of dentin repair that favors the maintenance of the vitality of the pulp.

Endo-PASS is easy to use and requires no procedure of preparation of the surfaces, it does not need any type of adhesive and it shows excellent sealing properties in order to avoid the risk of bacterial penetration and thus ensuring the absence of post-operative sensitivity and relapse.

The mechanical properties of Endo-PASS greatly exceed those of mineral cements based on aggregate of Portland cements, making Endo-PASS particularly suitable as endodontic cement

in combination with warm gutta percha and for endodontic repair (perforations, internal or external resorptions, apexifications, retrograde filling). The eventual extrusion of the concrete beyond the apex promotes the healing bone, or, in the presence of bone defects, apical, allows the rapid formation of new bone tissue.

Endo-PASS can also be used in the coronal part of damaged teeth, as material for the direct and indirect pulp capping, for the repair of carious lesions deep in pulpectomy.

Endo-PASS can be coupled to any technical adhesive / composite in order to finalize the definitive reconstructions even in the same session, making them extremely durable and aesthetic. The mixing of the Endo-PASS leads to the formation of a paste of creamy appearance with hardening reaction that starts immediately after the dust comes in contact with water or other biological fluids (saliva, blood, cleivocoral fluid and dentin).

This process (hydration reaction) initially produces calcium hydroxide and other silicates thus resulting in an alkaline environment (about pH = 13) strongly antibacterial. Hydration of the various silicates in the formulation depends on the principal physical, chemical and mechanical.

To influence the development of these properties, in addition to the ratio of tricalcium silicate and dicalcium silicate, is the presence of zirconium oxide and phyllosilicates of cationic and anionic surfactants, which encourage greater strength of the cement due to a more efficient distribution of moisture.

It also presents a refinement of particle size which is particularly cared in the present preparation. In fact, a more concrete end has a greater specific surface area and therefore greater speed of hydration.

Most speed cementation is also demonstrated by NMR analysis. The presence of barium sulphate in the formulation determines the necessary radio-opacity without presenting toxicity.

WARNINGS

Smaller quantities can be disposed of with household waste, while any overdue packaging must be disposed of in accordance with the provisions of law.

CONSERVATION

Keep the product in and/or dry boxes and fresh.

Via Torino, 765 - 21020 Mercallo (Va) - Italy
Phone +39 0331 969270 - Fax +39 0331 969271
www.deitalia.it

ENDO-PASS

Endodontic Penta-Aggregate Silicate Sealer

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