Estimado cliente

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Dear Customer

Thank you for having chosen a quality product from Dentaurum.

In order to use this product at its best for you and your patients, it is important to study and follow these directions for use carefully.

The written instructions cannot cover all eventualities during operation. For this reason our Hotline is available to answer any other questions and ideas that may arise.

Due to a constant development we recommend, even when you use the same products frequently, to study the enclosed latest issue of the instructions for use or refer to our website at [www.dentaurum.de](#).

1. Manufacturer

Dentaurum GmbH & Co. KG
Turnstraße 31 · 75228 Ispringen · Germany

2. General information

The light-curing retainer-adhesive ConTec LCR is a flowable, radiopaque microhybride composite. It is used with the Primer ConTec Prime. Due to its low viscosity ConTec LCR provides very good wetting and adaptation to the tooth substance. The product shows high durability and color stability. ConTec LCR contains about 64 % inorganic filler material (highly disperse silicon dioxide, micro fine glass 0.7 µm).

3. Application

To fasten retainers.

3.1. Preparation of enamel surfaces

Clean enamel areas to be etched with a fluoride-free polishing paste and water. Then dry thoroughly with oil-free air.

Etch enamel with the etching gel ConTec Etch for a minimum of 15 seconds Prolong the etching time by 10 seconds if enamel surfaces have been fluoridated.

Rinse enamel thoroughly with water for about 10 seconds and dry with air. The etched and dried enamel surfaces should have a chalky-white appearance. Create a dry work field. If etched enamel is contaminated with saliva, repeat the etching (10 seconds).

3.2. Primer application

Dispense a small amount of ConTec Prime on a mixing pad. With a brush, apply a thin, even layer to the etched enamel, distribute with a weak stream of air and polymerize for 20 seconds with a light source.

3.3. Application

Flow caps: Remove protective cap of the flow cap. Insert the caps into the opening of the dispenser (ConTector), click it into place under slight pressure and pull it to the front. Bend the metal cannula according to the indication into the desired direction. Apply ConTec LCR directly by pressing on the dispenser slowly and evenly.

Use flow caps for one patient only!

Apply ConTec LCR to the ends of the preformed retainer wire, which will then be positioned in the patient’s mouth. The material is then to be polymerized with a light source for 40 seconds at the end of the wire for pre-fixing. Apply the adhesive between the wire and etched tooth surface of the other teeth. Polymerize the adhesive for 40 seconds with a light source.

Position the light source close to the adhesive to polymerize. If the distance is more than 5 mm the polymerization will be less effective. This can cause discoloration and insufficient adhesion to tooth surface, which in turn can loosen the retainer.

After polymerization, the removal of access material must be done immediately under cooling, (such as, extra fine diamond discs, or polishing discs).

3.4. Attention

- Do not use cleaning agents or disinfectants to clean the flow caps.
- ConTec LCR has to be at room temperature before use.
- Apply ConTec LCR with even but not too strong pressure.
- To light-cure this material conventional polymerization devices are suited. The light output should be a minimum of 500 mW/cm² on halogen polymerization devices and 300 mW/cm² on LED devices.

4. Information about composition / ingredients

The composition /ingredients are included in the Safety Data Sheet. → [www.dentaurum.de](#)

5. Safety instructions

- The adhesive and the primer contain Bis-GMA, urethan dimethacrylates, TEGDMA, hydroxyethyl-methacrylate, BHT, and amines. Do not use in case of known hypersensitivity to ingredients (allergies).
- No known side effects. Sensitisation of hypersensitive persons cannot be excluded.
- When treating with cationic mouth-rinses, plaque indicator solution or chlorhexidin, discoloration may occur.

6. Warnings

Avoid long-term and repeated contact of the unpolymerized material with the skin and the mucosa. On contact these parts must be cleaned thoroughly. Upon contact with the eyes immediately rinse with water and consult an ophthalmologist.

7. Storage

Store the containers tightly closed, at room temperature. Do not use after expiry date.

8. Availability

ConTec LCR Set	REF 163-530-00
Content	
10 x ConTec LCR flow caps adhesive 0.25 g each	
1 x syringe etching gel ConTec Etch 5 ml + 2 curved application cartridges	
1 x ConTec Prime Primer 5 ml in bottle	
50 x brushes	
1 x brush holder	
2 x mixing trays	
1 x instructions for use	


Individual products

Content		
25 x ConTec LCR flow caps adhesive 0.25 g each		REF 163-531-00
1 x syringe etching gel ConTec Etch 5 ml + 2 curved application cartridges		REF 163-512-00
1 x ConTec Prime Primer 5 ml in bottle		REF 163-513-00
50 x brushes		REF 163-417-00
5 x brush holders		REF 163-419-00
2 x mixing pads		REF 163-416-00
The product REF 163-530-00 does not contain:		
1 x ConTector, applicator for flow caps		REF 163-532-00

9. With reference to our quality standards

Dentaurum ensures a faultless quality of the products manufactured by us. These recommendations are based upon our own experiences. The user himself is responsible for the processing of the products. Responsibility for failures cannot be taken, as we have no influence on the processing on site.

10. Explanation of symbols used on the label

	Order number
	Batch description
	CE Mark MP Class I
	CE Mark MP Class II
 → www.dentaurum.de	Instructions for use available in the Internet
	Codablock barcode
	Matrix code
	F/Flammable
	Xi/Irritant
	C/Corrosive
	Prescribed storage temperature
	Use by... (Best before date)

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