RANGE OF MOTION ADJUSTMENTS

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ENGLISH

The OActive 2 comes with standard 0° extension stops. Additional flexion and extension stops may be installed if they are prescribed by your physician or needed for optimal performance or fit of the brace. These should be inserted by the technician who fits your brace.

Optional 5°, 10°, 15°, 20°, and 25° extension stops and optional 45°, 60°, 75° and 90° flexion stops are included with each brace.

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CALF SUSPENSION WRAP APPLICATION INSTRUCTIONS

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ENGLISH

A suspension wrap may be provided with your brace if necessary. This wrap is designed to aid suspension and brace position by wrapping more surface area and applying even compression to the soft tissue of the calf.

STEP 1.

Attach the hook Velcro® strip sewn into the edge opposite the finger pull to the inside two straps of the lower half of the brace (1).

STEP 2.

Pull the wrap around the outside of the leg and brace (2).

STEP 3.

Continue the wrap all around the brace and leg and attach it to itself with the finger pull end (3).

NOTE: Wrap should be positioned so that the top edge is ½ inch above the strap right below the knee (Strap #1, Fig. 4).

CARE AND MAINTENANCE

Wrap should be removed and hand washed as needed using cold water and mild detergent. Rinse thoroughly and air-dry (do not machine-dry).

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ENGLISH

INDICATIONS:

- Mild to severe unicompartmental knee osteoarthritis
- Medial or lateral compartment

WARNINGS:

If you experience pain, swelling, sensation changes, or unusual reactions while using this product, contact VQ OrthoCare's Patient Care department at 800.452.7993 or consult a physician.

FOR SINGLE PATIENT USE ONLY.

NOTICE: While every effort has been made using state-of-the-art techniques to obtain maximum compatibility of function, strength, durability and comfort, this device is only one element in the overall treatment program administered by a medical professional. There is no guarantee that injury will be prevented through use of this device.

CAUTION: Applicable federal, state, and/or local regulations may restrict this device to sale by or on the order of a practitioner licensed by law of the State in which he/she practices to use or order the use of this device.

WARRANTY: VQ OrthoCare warranties all its products from the original date of purchase against defects in materials and workmanship. Normal wear and tear during use of a product is not considered a defect. Contact VQ OrthoCare at **949.261.3000** for specific product warranty information.

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
Subject to one or more patents listed on
www.vqorthocare.com/patents

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










MANUFACTURED AT
1390 Decision St, Suite A
Vista, CA 92081, USA

© 2015 VisionQuest Industries Inc.
VQO151830REVB

OACTIVE 2 APPLICATION INSTRUCTIONS

STEP 1.

- Open the brace by lifting up on the buckles (1a) and unhooking the buckles from the frame (1b).
- Loosen the straps by peeling up the strap tab (1c) and re-applying at the end of the strap (1d).

STEP 2.

- While sitting at edge of chair, bend knee slightly and apply brace to lateral side (outside) of leg.
- Align center of hinge with top of kneecap and position slightly behind midline of leg (2).

Important note for initial fitting: When positioning the **OActive 2 off-the-shelf brace**, it should first be adjusted to neutral. This is achieved by setting the correction angle (see step 5) so the cuffs and hinge rest on leg without gapping or pressure. Medial or lateral compartment correction will be adjusted in step 5. The **OActive 2 Custom** has correction pre-set and **should not** be set to neutral.

STEP 3.

Attach strap buckles in sequence indicated on strap ends (3a) by hooking the buckle on the frame (3b) and snapping the buckle down (3c). Be sure to securely snap buckle into place. Hold brace at hinge to prevent rotation while tightening straps (3d).

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- Fasten strap 1 snugly just above calf muscle. Strap 1 is important in maintaining brace position.
- Fasten strap 2 around middle of calf.
- Fasten strap 3 around lower thigh.
- Fasten strap 4 around upper thigh.

All straps should be comfortably tightened to properly secure brace and provide relief. Do not overtighten thigh straps.

STEP 4.

- Once straps are secure, center of hinge should be aligned with top of kneecap and positioned slightly behind midline of leg.
- To ensure proper fit for activity, walk a short distance. It is normal for brace to drop or settle into position. Make strap and position adjustments as necessary.
- The straps are now set for proper fit and do not need to be loosened and tightened for each application. The brace can now be put on and taken off using just the buckles (4).

STEP 5.

Warning: Initial adjustment of the OActive 2 brace should be done by a trained bracing professional in accordance with a Doctor's prescription indicating the compartment to unload. Subsequent adjustments, though rarely needed, should also be made by a trained bracing professional (5a, 5b).

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For initial fitting of OActive 2 off-the-shelf (for OActive 2 Custom skip to section 5.1): With the brace in neutral (from step 4), bend the knee to 90°. With the provided patient adjustment tool, turn the Q-hinge screw 1/8 turn clockwise (important) to increase lateral compartment relief or 1/8 turn counterclockwise to increase medial compartment relief. This is a good starting point for initial use.

STEP 5.1 FOR INITIAL FITTING AND ALL ADJUSTMENTS

- Stand and walk a short distance to ensure comfort and proper correction after each adjustment or upon initial fitting of OActive 2. If there is any discomfort in the knee joint or leg, turn the angle adjustment back slightly.
- If more relief is required, turn Q-hinge screw no more than 1/8 turn. Walk after the adjustment and check as above.
- Over-adjusting the Q-hinge may cause discomfort, rotation or skin pressure issues. Adjust conservatively.

“WEAR TIME”

Gradually increase brace wear time to allow your body to adjust to the brace. Start with 30 minutes on day one and gradually increase. Pain caused by osteoarthritis may take time to subside. It is important that you routinely wear the brace during weight bearing activities (e.g. when you're on your feet) as it will reduce joint loading and allow you to begin feeling relief. It may take a number of weeks for relief to become significant.

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CARE AND MAINTENANCE

Salt water, sand, dirt and other elements can cause damage. Thoroughly rinse brace with fresh water after exposure to these conditions. Towel or air-dry. Liners should be removed and hand-washed as needed using cold water and mild detergent. Rinse thoroughly and air-dry (do not machine dry.) Lubricate hinges with dry silicone spray as needed.

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