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009095EX-20150422

Multisil–Mask hart

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Instructions for use

Multisil–Mask hard

Before using the product, please read these istructions for use carefully.

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

1. Description of the product

Multisil–Mask hard is a self-curing pink resin in the double-mix cartridge for simple fabrication of a gingival mask on a working model.

2. Indication

Multisil–Mask hard serves for the integration (reconstruction) of the gingiva for over-impressions in implantology and dental prosthetics. Multisil–Mask hard is used to fabricate a removable, torsion-resistant gingival mask made of special resin. Fixation on the model is achieved using the Varío-Stud-Snap vks-oc system.

3. Contraindication

Based on the current state of knowledge there are no known contraindications.

4. Hazard and safety information

Signal word:	
Attention	
Hazard statements:	
H319	Causes serious eye irritation.
H412	Harmful to aquatic life with long lasting effects.
Precautionary statements:	
Prevention:	
P264	Wash thoroughly after handling.
P273	Avoid release to the environment.
P280	Wear protective gloves / eye protection.
Response:	
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337+P313	If eye irritation persists: Get medical advice / attention.
Disposal:	
P501	Dispose of contents / container in accordance with local / regional / national / international regulations.
Product contains:	dibenzoyl peroxide. May produce an allergic reaction.

5. Storage and durability

Multisil–Mask hard should be stored at temperatures from 15 to 28 °C. Protect against sunlight and heat. Please observe the expiration date printed on the label.

6. Processing

6.1 Preparation of the cartridge

The cartridge is inserted correctly into the dispensing device with the notch facing downward. The cap is turned by 90° to remove it and the respective mixing cannula is attached. It must be ensured that the individual components will not contaminate each other. If new cartridges are used, first discard a strip of a length of 1 cm since due to manufacturing-technical reasons there may be slight variations in the quantity of material filled into the cartridge halves. After the fabrication of the gingival mask, remove the mixing cannula and place the cap on the cartridge again.

6.2 Preparation of the impression

The impression is disinfected, blood and saliva residues are removed and then the impression is dried. The implant impression abutments are screwed on the laboratory analog, placed in the impression and the correct fit is checked. The Implant-compensating cones (REF 430 0739 2) are pushed over the respective laboratory implant (analog) in the impression to avoid any undercut.

6.3 Injecting Multisil–Mask

Multisil–Mask hard can be processed at room temperature for approx. 7–8 minutes. Multisil–Mask hard is filled into the impression by exerting pressure onto the lever of the dispensing device. It must be ensured that the material is applied without the formation of bubbles. Multisil–Mask hard must not be applied into undercuts. In order to fix the gingival mask in the model later on, vks-oc plastic matrices (diameter 1.7 mm) are pressed into the gingival mask as long as it is still soft. Once the application has been completed, Multisil–Mask should harden under atmospheric pressure. The surface can be smoothed by vibrating slightly. Polymerization is completed after 10 minutes at 23°C and the gingival mask can be removed from the impression; if required the impression abutment can be screwed on to remove the impression.

6.4 Finishing the gingival maks

After removing from the impression, the dispersion layer is removed and the outer contour is ground conically to allow removal from the plaster model later on. After grinding, the gingival mask and the laboratory implants (analog) are placed back on the model. The corresponding metal transfer patrices 1.7 are inserted into the plastic matrices vks-oc 1.7. Then the impression is filled with plaster. Before removing the gingival mask, the compensating cones must be heated with hot water to ensure they do not form any resistance when lifting the gingival mask; then the cones are boiled out. Finally, transparent the varnish (REF 540 0100 6) should be applied to the gingival mask to obtain a smooth surface.

7. Technical data

Processing time incl. mixing time: approx. 7 minutes.
Time until removal from the mould: approx. 20 minutes.

8. Symbols

8.1 Symbols according to DIN EN ISO 15223-1

	Reference number	
	Lot number	
	Can be used until YYYY-MM	
	Temperature limits	Lower limit: 15 °C Upper limit: 28 °C
	Manufacturer	
	Attention!	

8.2 Symbols according to GHS regulation

	GHS02, flame
--	--------------

8.3 Symbole nach NFPA (USA)

	HEALTH 1 FIRE 1 REACTIVITY 0	Hazard diamond
--	------------------------------------	----------------

9. Other indications

The product may only be used by dentists, dental technicians and adequately trained dental staff. Only original tools and components may be used for processing. The lot numbers of all components should be recorded to ensure traceability and to handle complaints. These instructions for use are based on state-of-the-art methods and equipment and our own experience. The product may only be used for the indication described under item 2. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on further processing. Any occurring claims for damages may only be made up to the value of our products.

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