**DISTRIBUTEURS / DISTRIBUTORS :**

FR	Laboratoires Septodont 58, rue du Pont de Créteil 94107 Saint-Maur-des-Fossés Cedex
GB	SEPTODONT Ltd. Units R & S - Orchard Business Centre St Barnabas Close Allington, Maidstone Kent ME 16 0JZ
ES	LABORATORIOS MAB DENTAL, S.A. C/Córcega 116 SP - 08029 - BARCELONA
DE	SEPTODONT GmbH Felix - Wankel - Str. 9 53859 Niederkassel - Mondorf
IT	GIOVANNI OGNA & FIGLI SpA Via Figini 41 20053 MUGGIO (MI)
BE	Distributeur/Verdeler BENELUX SEPTODONT NV-SA Grondwetlaan 87 - Av. de la Constitution 87 Brussel 1083 Bruxelles

PL	SEPTODONT-POLSKA SP. Z.O.O. Taneczna 18 A 02-829 Warszawa
SE	BIGMAN AB Järnvägsgatan 5 SE-172 75 Sundbyberg
US	SEPTODONT Louisville, CO, 80027
CA	SEPTODONT Cambridge, ON N1R 6X3, CANADA
AU	SPONSOR : Specialities SEPTODONT PTY Ltd. 14 Cliffbrook Crescent, Emu Plains, NSW 2750

**Ultra
Safety Plus**

PART B

**Ultra
Safety Plus XL**

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EN

INSTRUCTIONS FOR USE

Part B (handle) must be assembled with Part A (needle). This handle is designed to be used with the Sterile Single Use Ultra Safety Plus Device (Part A) for maximum safety to the practitioner and staff against the risk of needle stick injury and cross contamination.

FEATURES:

- Ergonomic design and ease of use.
- Locking on the barrel of the Ultra Safety Plus Syringe device Part A.
- Enables aspiration.
- Can be used with 1.7, 1.8, 2.2 ml cartridges.

HOW TO USE:

Step 1 : Insert the cartridge into the Ultra Safety Plus Syringe barrel.

Step 2 : Grip the handle plunger and put the thumb behind the finger holder. Introduce the handle tip into the barrel of the Part A, behind the cartridge.

Step 3: Now slide the sheath protecting the needle backwards, towards the handle until it CLICKS (the click is made as the sheath hits the handle and locks the unit together).

Step 4 : Perform the injection.

Self aspiration: At the base of the cartridge barrel of the injection device you will see there is a small protuberance. It appears as a plastic tip holding the centred back end of the needle that penetrates the diaphragm of the cartridge when inserted. At the beginning of the injection, the diaphragm presses against the protuberance, and aspiration occurs as the diaphragm moves back away from the protuberance when released.

Aspiration : The pliable O-Ring situated at the end of the handle piston creates a vacuum when the practitioner pulls back the piston.

The cartridge stopper follows the piston providing active aspiration. It is best observed when a minimum of 0.25 ml to 0.35 ml of solution has been previously expelled from the cartridge (providing space).

Step 5 : After you have finished the injection, remove your fingers from the O-Ring Handle and grip the finger holder, slide the protective sheath down the barrel of the syringe until you reach the first holding position, then grip the barrel of the syringe and pull the protective sheath down further to the final locking position. Now the injection device (Part A + Part B) can be disposed of safely.

To remove the handle, hold the barrel with one hand and using the other hand, place a finger in the ring of the handle and pull backwards until the plunger is fully retracted. After you have fully retracted the plunger, peel off the handle in one movement. Dispose of the Syringe into the sharps container. The handle can be disposed of in the clinical waste bag or in any other container according to the regulations in force in the country involved.

WARNINGS :

- 1)** This is designed to be a single use product. Any reuse constitutes a risk to the safety of the user and the patient.
 - Do not resterilise.
 - Do not use for more than one patient.
 - Use a new handle with any new sterile injection device (Part A).
- 2)** Prior to use, please check the presence of the pliable o-ring at the end of the plunger.
- 3)** When assembling the handle with Part A, be sure that it is perfectly locked (audible CLICK).

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