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Instructions for use

Cast-on prosthetic elements

Before using the product, please read these instructions for use carefully.

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

1. Description of the product

The cast-on prosthetic elements are used for prosthetic restorations made from SKY classic and blueSKY implants and ensure the accurate fabrication of individual suprastructures.

2. Indication / Application areas

Individual abutment solutions can be achieved using the cast-on SKY abutments (REF SKY-PV00). The SKY bar abutment HL (REF SKY-SH11) of the modular system is available for the fabrication of bars.

3. Contraindication

The SKY implant system is not suitable for standard contraindications in implant dentistry such as patients with cardi-ovascular problems, bruxism, osteoporosis, heavy smokers, alcoholics, diabetes and metabolic diseases and insufficient primary stability of the implants. The product should not be used for patients with a suspected allergy until allergy tests have been completed to confirm that the patient is not allergic to the materials contained in the product.

4. Storage and durability

The product is supplied in the non-sterile condition and dry and dust-free storage in the original packaging must be ensured.

5. Processing / Use

Generally, we recommend thorough surgical planning and selection of the prosthetic components depending of the situation and type of restoration.

The retention screws used in the laboratory are not suitable for the clinical use, since they are subject to stress in the laboratory; as a result long-term and reliable fit can no longer be guaranteed. bredent medical offers special laboratory screws to be used in the dental laboratory.

All components must be disinfected prior to the insertion (see item 8).

SKY abutment, cast-on (REF SKY-PV00):

- Place the cast-on prosthetic element in the model and shorten the modelling sleeve in accordance with the anatomical situation but no more than down to the metal base. Wax must be used to attach the modelling sleeve to the metal base and in order to ensure adequate fit of the modelling sleeved during investing.
- Use aluminium oxide (110 µ) to roughen the modelling sleeve and model the abutment in accordance with the anatomical situation. The model can be adapted to the gingival course in cases of individually formed gingiva.
- Use investment material, which enables careful divesting.
- When casting on, adhere to the instructions of the manufacturer of the cast-on alloy. Please observe the technical data (see item 6).
- When divesting, use only glass beads (abrasive material) and exert minimal pressure; the torx must not be sand-blasted.

Attention: Mechanical processing of the basal surface and the collar of the metal base must be avoided!

6. If horizontal screw retention is intended, threads should be tapped into the coping.

7. Fabricate the prosthetic restoration using the standard dental procedures and techniques.

When selecting the cast-on material it must be verified whether the material is suitable for the fabrication of implant suprastructures. The manufacturers' instructions must be observed.

6. Technical data

REF	Product	Material	platform
	SKY abutment, cast-on/ SKY bar abutment HL		
SKY-PV00, SKY-SH11		<u>Material - metal base:</u> Au 60 %, PD 20 %, Pt 19 %, Ir 1 % <u>Coefficient of thermal expansion</u> 25-100° C / 25-600° C CTE 11,9 - 12,2 10 ⁻⁶ K ⁻¹ <u>Melting range:</u> 1.400-1490° C <u>Material - modelling sleeve:</u> PMMA	regular
	SKY Prosthetics accessories		
SKYLP522	SKY laboratory screw 2.2 grey	Titanium grade 4, grey	

7. Symbols

-  Lot number
-  Observe instructions for use
-  Single use
-  Non-sterile
-  Dry storage

8. Additional information

The product may only be used by dentists, dental technicians and adequately trained dental staff.

Only original tools and components may be used for processing.

The lot numbers of all components should be recorded to ensure traceability and to handle complaints.

Prior to the integration into the patient's mouth, all components must be disinfected. Only standard disinfectants for dental instruments may be used, e.g. Dürr ID 212 forte. The manufacturer's instructions on the use of the disinfectant must be observed.

After disinfecting, disinfectant residues must be thoroughly removed from the components (see the instructions of the manufacturer of the disinfectant; for example under running water - for ID212 forte).

During patient recall appointments, proper fit of the prosthetic restoration or retention element needs to be checked since loose screws may cause fatigue fracture and fragments may be aspirated or swallowed.

At least one recall appointment must be made each year.

All products labelled  are intended for single use only. However, if this product is used several times, the following risks may arise: cross contamination, malfunction, inaccurate fit, etc.

These instructions for use are based on state-of-the-art methods and equipment and our own experience. The product may only be used for the indication described under item 2. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on further processing. Any occurring claims for damages may only be made up to the value of our products.

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