INSTRUCTIONS FOR USE

DESCRIPTION

MACRO-LOCK POST™ ILLUSION™ X-RO® is a serrated taper post, length 17,5 mm.

MACRO-LOCK POST™ ILLUSION™ X-RO® is a light yellow translucent fiber post with a colored resin matrix which becomes translucent in the mouth. The post color will re-appear on command under cold air, water spray if re-treatment is required.

The radiopacity level of the MACRO-LOCK ILLUSION™ X-RO® is 340 % Al equ. (ISO 4049)

	• 1	• 2	• 3	• 4	• 5	• 6
Ø apical tip	0.80	0.80	1.00	1.00	1.30	1.30
Ø Post head	1.35	1.47	1.67	1.83	2.05	2.20

INDICATIONS FOR USE

In case of insufficient residual tooth structure (<4mm), the post is needed to support and secure the coronal restoration.

CONTRAINDICATIONS

Insufficient residual dentin: at least 1.5 mm of tooth structure around the preparation (ferrule effect) is required.

PRECAUTIONS

The post must be cleaned with alcohol before insertion in the canal. Avoid touching the posts with your fingers after cleaning. Shortening the post should be done outside of the mouth. The use of rubber dam is recommended. However, sometimes the rubber dam may prevent the temperature elevation required for the color disappearance until the dam is removed. The crown must cover the preparation for at least 1.5 mm of healthy dentin in order to get the ferrule effect.

ADVERSE REACTIONS

None known to date.

STEP-BY-STEP INSTRUCTIONS

1. Root canal preparation.
2. Select the correct MACRO-LOCK POST™ ILLUSION™ X-RO®. Determine the post size according to the anatomical situation of the tooth using the radiograph and the calibrating paper.
3. Determine the length of the coronal build-up, which will account for 1/3 of the final post length. Remove the root filling with the starter Drill (rotation speed 1000-2000 rpm) as deep as necessary for the post to be inserted to 2/3 of the

root length. At least 4 mm of root canal filling should remain in the apical region. In case of curved canal, this calculated length is not suitable.

4. Shape the canal with the Finishing Drill corresponding to the selected MACRO-LOCK POST™ ILLUSION™ X-RO® size, in increasing size sequence (rotation speed 1000-2000 rpm).
5. Try in the post, to check that the post is properly seated in the canal. Shorten the post to its final length with a diamond disc outside of the mouth. **Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post.**
6. Clean post with alcohol.
7. Apply etchant in the post space and to the exposed dentin for 15 seconds*. Rinse for 10 seconds. Remove the excess water with paper points but leave surface moist. Do not put the etchant in contact with the gingiva.
8. Apply 2 coats* of the primer with a thin brush (e.g. Composi-brush RTD) in the post space. **Remove excess using paper points and gently air dry all surfaces.** Enamel and dentin surfaces should have a uniform, glossy appearance. If not, repeat application. Light-cure the primer* inside the post space for 10 to 20 seconds.
9. Apply a single coat of primer to the post. Gently air dry 5 seconds and light-cure the post for 10 to 20 seconds (if needed) outside the mouth*.
10. Mix the cement* and apply it in post space and to the post using a Lentulo spiral or a syringe tip. Seat the post immediately. Remove excess with appropriate instruments. If dual cure cement, light-cure for 40 to 60 seconds, applying gentle pressure with the tip of the curing light probe.

11. Apply 2 coats of bonding agent* to the exposed post, cement and coronal areas. Remove excess and air-dry with jet of air. Light-cure for 20 seconds (if needed).
12. Directly model the build-up using a composite core build-up material*. R.T.D. recommends that the coronal end of the post be covered by core composite. In cases where a prosthetic crown is not indicated, a minimum of 1 mm of composite restorative should cover the post.

OTHER RECOMMENDATIONS

MACRO-LOCK POST™ ILLUSION™ X-RO® can be disinfected with medical alcohol.

They can neither be sterilised nor disinfected in thermo disinfectant. **The use of autoclave is not recommended.**

Do not store in direct sun light or artificial light MACRO-LOCK POST™ ILLUSION™ X-RO® Starter and Finishing Drills: the drills are supplied non-sterile. Disinfect and sterilize drills before each use.

Disinfect drills with disinfecting or cleaning agents adapted for rotary instruments and follow the manufacturer's instructions for use.

Do not disinfect drills in thermo disinfectant.

Control drills visually before each use: change drills when they are damaged, corroded or dull.

Use a new drill after 12 - 15 use.

SAFETY AND LIABILITY

RTD declines any liability or compensation for possible damage due to:

- Use of instruments which do not belong to the system and which may impair the function.
- Failure to observe the Instructions For Use The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the instructions for use. For Dental use only.

* Refer to the manufacturer's instructions for use.

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