PROTEINAS / PROTEINS / PROTÉINES

Componente Component Composant	Método Method Méthode	Valor Value Valeur	Rango Range Intervalle	Unidades Units Unités
Albumina Albumin Albumine	Método BCG / BCG method / Méthode BCG	3.09 30.9	2.63 – 3.55 26.3 – 35.5	g/dl g/l
Fructosamina / Fructosamine	Método NBT / NBT method / Méthode NBT	426	362 – 490	µmol/l
Proteínas Totales Total Proteins Protéines Totales	Biuret	4.83 48.3	3.86 – 5.80 38.6 – 58.0	g/dl g/l

OTROS METABOLITOS / OTHER METABOLITES / AUTRES MÉTABOLITES

Ácido Úrico Acide Urique Uric Acid	Uricase POD / Uricase POD	9.56 569	7.65 – 11.5 455 – 684	mg/dl µmol/L
Bilirrubina Total Bilirubine totale Total Bilirubin	Jendrassik – Gróf Manual / Manuelle Jendrassik – Gróf Méthode Automático / Automatic / Automatique	4.44 75.9 4.54 77.6	3.54 – 5.34 60.5 – 91.3 3.63 – 5.45 62.1 – 93.2	mg/dl µmol/L mg/dl µmol/L
Bilirrubina Directa Bilirubine Directe Direct Bilirubin	Jendrassik – Gróf Méthode Automático / Automatic / Automatique	2.17 37.1	1.74 – 2.60 29.7 – 44.5	mg/dl µmol/L
Creatinina Creatinine Créatinine	Jaffé cinético Kinetic Jaffé Jaffé cinétique	4.02 355	3.29 – 4.74 291 – 419	mg/dl µmol/L
Glucosa / Glucose	GOD-POD	278 15.4	223 – 333 12.4 – 18.5	mg/dl µmol/L
Urea / Urée	Ureasa / Urease / Uréase GLDH Ureasa / Urease / Uréase Berthelot	110 18.3 119 19.8	90 – 131 15.0 – 21.8 96 – 141 16.0 – 23.5	mg/dl µmol/L mg/dl µmol/L

SERISCANN® ANORMAL

SUERO CONTROL ANORMAL / ABNORMAL CONTROL SERUM / SÉRUM CONTRÔLE ANORMAL

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Product description

The control serum is prepared from human serum with human and nonhuman enzymes and nonprotein constituents added. Bacteriostatic agents have been added. The value and expected range for each constituent are derived from interlaboratory data. The mean of several determinations may not duplicate the value printed on the insert but should fall within the expected range. To determine the accuracy and precision of a certain analytical method, it is advisable to run a Normal as well as an Abnormal control serum samples.

Caution

Human serum was used in the manufacture of this product. Each donor unit was tested with licensed reagents and found negative for HBsAg and nonreactive for the HIV antibody. Because of no test method can offer complete assurance that products derived from human blood will not transmit infection agents, it is recommended that this product be handled with the same precautions used for patient specimens.

Procedure

Add exactly 5 ml. of deionised water to each vial. Allow to sit 5 to 10 minutes. Swirl the contents until homogeneous. This control serum should be handled as if a sample from human origin. Analytical values should fall in between the assay range.

Storage and stability

The product remains stable, when stored at 2 - 8° C., until the expiration date stated on the label. Once rehydrated, the constituents are stable for ten days, when stored at 2 - 8° C., except for Acid and Alkaline phosphatases only stable for 48 h. and GOT, GPT only stable for 4 days at 2 - 8° C. Since Bilirubin is a light-sensitive metabolite it is advisable, for the sake of better storage conditions, to keep the control in the dark.

Configuration

Kit 10 x 5 ml. (Ref. 99 93 29). Contents: 10 Freeze-dried vials of abnormal control serum.

1 x 5 ml. (Ref. 99 46 85). Contents: 1 Freeze-dried vial of abnormal control serum.

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