CE 0483

Gebrauchsanweisung
Instructions for use
Mode d'emploi
Modo de empleo
Modalità d'uso

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Instruction for use Self-ligating brackets – discovery® sl 2.0

Dear Customer

Thank you for choosing this high-quality product from Dentaurum.

In order to achieve the best results with this product, it is important to carefully study and follow these instructions for use.

The written instructions cannot cover all possible situations that can occur during use. For this reason, our hotline is available to answer any other questions or concerns that may arise.

Due to constant developments, we recommend, even when using the same products regularly, that you study the enclosed current instructions for use or refer to our website at www.dentaurum.de.

1. Manufacturer

Dentaurum GmbH & Co. KG · Turnstraße 31 · 75228 Ispringen · Germany

2. General description

discovery® sl 2.0 are the smallest self-ligating brackets worldwide. Developed on the basis of the well-proven discovery® brackets discovery® sl 2.0 unites the sum of many years of experience with the modern standard of orthodontic treatment.

Per CAD (computer aided design) the brackets are ideally adapted to the tooth shape and offer ideal handling. The lasered FDI-marking on the base and the lasered cross marking on the lid help to easily identify and position the brackets. The patented laser-structured base ensures the optimal bond and guarantees easy debonding.

3. Working stages

3.1. Patient consultation

Inform the patient about the possible consequences and risks involved in treatment using brackets, in particular during debonding.

3.2. Tooth conditioning

Check each tooth thoroughly for any existing damage to the enamel to avoid being "surprised" by damage to the tooth during debonding. Adhere to the bracket adhesive manufacturer's instructions for conditioning the teeth.

3.3. Bracket adhesive selection

discovery® sl 2.0 brackets can be bonded using all commercially available bracket adhesives. We recommend our ConTec SE adhesive (REF 163-550-00).

Adhere to the manufacturer's instructions for the bracket adhesive as well as the Adhesive System Table in our orthodontic catalog.

3.4. Removal and handling

Grip the discovery® sl 2.0 brackets using clamping tweezers (e.g. REF 025-277-00) mesially-distally, to remove them from the packaging and handle according to the following steps. The base of the discovery® sl 2.0 brackets should not be touched with your fingers or contaminated in any way, as this could otherwise detrimentally affect the bond.

3.5. Application of the bracket adhesive

Prepare the bracket adhesive according to the manufacturer's instructions. Apply a small amount of adhesive to the center of the bracket base, but enough to cover the entire bracket base. This guarantees an optimum bond.

3.6. Placement and alignment of the bracket

After applying the adhesive, place the discovery® sl 2.0 bracket in the usual way onto the tooth.

Optimally align the discovery® sl 2.0 bracket with the angulation aid of the clamping tweezers and press it onto the tooth so that the surplus adhesive exudes from the sides of the bracket base.

Use the laser markings on the top of the brackets (Fig. 1) as a positioning guide. The vertical line marks the bracket center along the relevant long axis of the tooth crown (angulation). The horizontal line marks the middle of the slot (archwire level). **Please note that the slot does not lie in the center of the bracket, and take this into consideration when aligning the bracket!**

It is essential to avoid contact with the antagonist due to the risk of abrasion and/or premature bracket loss. If necessary a palatal bite raise can be used e.g. Bite-Bumper™ (REF 705-000-00).

It is essential to ensure that no adhesive material enters the bracket fastening mechanism on the gingival and occlusal side of the bracket, as this could restrict or destroy the functionality of the lid.

Completely remove the exuded excess adhesive before it has cured. This is particularly important to ensure easy debonding at the end of treatment.

3.7. How the fastening mechanism works / securing the arch wire

The slot of the discovery® sl 2.0 brackets is closed off with a lid. All lids, whether they are designed for the upper or lower jaw, open from the occlusal side towards gingival (Fig. 2). This is carried out using only the probe-shaped end of the sl instrument (REF 023-277-00) (Fig. 3).

The lid must be opened in order to insert the arch wire. Place the pointed probe-shaded end of the sl instrument from the labial/buccal side into the notch in the lid between the occlusal bracket wings (Fig. 4a and 4b) and flip the lid open, using slight leverage towards vestibular-lingual (Fig. 5). The lid will spring open.

The discovery® sl 2.0 brackets include a mechanism (CapSafe®), which enables the lid to remain open while the arch wire is being inserted. In order to activate CapSafe®, open the lid and press it back in a gingival direction with your finger or the sl instrument until there is a feeling of resistance (Fig. 6). The lid is now secured in the open position (Fig. 7).

Place the arch wire into the slot (Fig. 8). If necessary, use the forked end of the sl instrument to aid in this step. To close the lid, use your finger or the sl instrument to press the lid from the gingival to the occlusal position (Fig. 9). The lid should close with a distinct click, which guarantees the lid is now secure (Fig. 10).

In order to prevent lip and cheek irritations during arch wire change, we recommend to use lip retractors (REF 075-600-00 or 075-601-00).

3.8. Bracket debonding

The excellent bonding characteristics of discovery® sl 2.0 brackets means the process of debonding requires special attention: When debonding discovery® sl 2.0 brackets it is advisable to use a Weingart plier (e.g. REF 003-120-00) or in the molar region, an angled bracket removing plier (REF 004-347-00). Grasp the bracket mesio-distally on the tie wings, and twist toward the left or the right without pulling away from the tooth.

4. Note for doctor and patient

Please consider that teeth must never be brushed with the bracket lids open e.g. for arch wire change! This leads to torsion of the delicate lids, which distracts the function of the closing mechanism so that the lid cannot be optimally closed any more.

5. Information about the composition

Information about the composition is included in the materials list, see catalog or www.dentaurum.de.

6. Safety information

- The product should not be used if there is a known allergy to one or more of the ingredients.
- As the material 1.4542 used in the self-ligating discovery® sl 2.0 brackets is one of the slightly magnetizable metals, we recommend to take the brackets off before the magnetic resonance imaging examination.

7. Storage and shelf-life

Store discovery® sl 2.0 brackets in closed packaging.

8. Information for single-use products

discovery® sl 2.0 brackets are intended for single use only. Reconditioning of discovery® sl 2.0 brackets that have already been used (recycling) and reuse on patients is not permitted. Due to the mechanical wear during treatment and debonding neither dimensional accuracy nor bonding characteristics of the brackets can be guaranteed after a recycling process.

9. Additional tips

- Technical information such as torque, angulation, size etc. can be found on the relevant product label or in the catalog information.
- There are no particular precautions to be taken when disposing of discovery® sl 2.0 brackets after treatment.

10. Delivery program

For the complete range of discovery® sl 2.0 products please see our orthodontics catalog.

11. With reference to our quality standards

Dentaurum ensures a faultless quality of the products manufactured by us. These recommendations are based upon our own experiences. The user himself is responsible for the processing of the products. Responsibility for failures cannot be taken, as we have no influence on the processing on site.

12. Explanation of symbols used on the label

REF	Order number
LOT	Batch description
CE 0483	CE Mark MP Class 2a
	Instructions for use available in the internet
	Not for reuse
	Codablock barcode
	Matrix code

Date of information: 01/12
Subject to modifications

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Dentaaurum GmbH & Co. KG · Turnstraße 31 · 75228 Ispingen · Germania