Anclasil onewash

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E

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Editing date: 13.09.2010

Instructions for use

Precision impression silicone, addition curing, low consistency – light bodied – fast setting

Fig. 1

1. Tray preparation

Select a tray type suited for desired application (standard-size or custom, in metal or hard plastics). To avoid teeth contact with the tray, place stops from wax or light-curing resin into the tray at non-relevant areas of the impression. For optimal adhesion of impression material, we recommend applying a thin layer of an adhesive to tray (Fig. 1).

Let dry for approx. 1 – 2 minutes.

Fig. 2

2. Mixing and dispensing

Insert cartridge into Automix2-System mixing gun (Fig. 2). Remove and dispose of closure cap. For exact flow control extrude slightly material until uniform flow from both orifices is achieved (Fig. 3). Attach mixing cannula by aligning notched side to cartridge and turning clockwise until fastened (Fig. 4). Dispenser is now ready to mix silicone in any amount needed. After each use leave the used mixing cannula in place to seal the cartridge.

Fig. 3

3. Using the sandwich technique

First apply the mixed **Anclasil oneputty** onto the impression tray, giving the teeth area a distinctly trough-like shape (Fig. 5). Then apply a layer of **Anclasil onewash** into the shaped depression (Fig. 6). For applications in the groove area use intra-oral-tip or disposable syringe (Fig. 7). Insert the filled tray into patient's mouth and position it within 3 – 4 seconds, under light pressure. Then seat the impression without pressure, until the product has set.

Fig. 4

4. Disinfecting

After removing tray from the mouth, rinse impression under running, lukewarm water. If desired, impression can then be disinfected by immersing into 2% glutaraldehyde for 15 minutes.


Fig. 5

5. Pouring the model

We recommend waiting at least 15 minutes before pouring the model. The impression will retain its dimensional stability for up to 14 days. All class III and IV dental plasters and standard modelling resins can be used.


Fig. 6

6. Cleaning the tray


After setting, material can be removed with a blunt instrument. Remaining film of the adhesive is dissolved by immersion in acetone or similar cleaning agents (be sure to use only in well-ventilated rooms). Tray can then be cleaned and disinfected in the usual manner.

Fig. 7

important working hints

- Substances such as retraction liquids can adversely affect the setting reaction and should therefore be thoroughly rinsed off before taking the impression.
- Do not combine with condensation curing silicones.
- Contact with latex examination gloves or other latex material may impair the setting reaction of the product.
- Be sure to remove all traces of impression material from patient's mouth.
- Cured impression materials are chemically inert – spots on clothing should be avoided.
- We recommend use of **Anclasil oneputty** as tray material in the double-mix technique.

Indications for use:

- Impressions with sandwich and double-mix technique
- Crown- / bridge work
- Inlays, onlays, partial crowns
- Dowel crowns, veneers

Technical data:

DIN EN ISO 4823 - Type 3

- Mixed volume:** 50 ml cartridge
- Mixing ratio:** 1:1
- Product colour:** base: violet catalyst: purple
- Mixing time:** omitted (Automix2-System)
- Working time:** approx. 1 min. 30 sec.*
- Setting time:** approx. 3 min. 30 sec.*
- Time in the mouth:** approx. 2 min.
- Strain in compression:** < 6,5 %
- Recovery from deformation:** > 99,8 %
- Linear dimensional change:** < 0,2 %
- Application:** At 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5% rel. humidity.

*from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5% rel. humidity. Increased temperatures accelerate, decreased temperatures retard these times.

Ordering information:

Anclasil oneputty A-040005
1 x 250 ml jar of putty
1 x 250 ml jar of catalyst
2 dispenser spoons

Anclasil onewash A-040006
2 x 50 ml cartridges
12 mixing cannulas
6 intraoral ends

Mixing cannulas A-040007
Bag with 50 mixing cannulas

Intraoral ends A-040008
Bag with 50 intraoral ends

Ancladén

Ancladén, S.L. | Rubinstein, 2 Entlo 1ª | 08022 Barcelona
93 417 29 84 | ancladen@ancladen.com | www.ancladen.com