onetime perfect foilbag

GB
Instructions for use

**Precision impression silicone, addition curing,
high consistency – fast setting**

1. Preparing the impression tray

Select the appropriate impression tray (custom or stock tray made from metal or plastic resin) for the impression technique. Insert wax or light-curing resin stops in non-functional areas of the impression to prevent excessive pressure on the impression material. For optimum adhesion we recommend coating all trays with a thin layer of **Detaseal® Adhesive rapid** (Fig. 1). Allow to dry for approx. 1 - 2 min.

Fig. 1

2. Operating the cartridge holder

Before using the foil bag for the first time, remove the white security pin on the activator cap of the large **onetime perfect** foil bag by twisting it to the right (Fig. 2). Insert both **onetime perfect** foil bags with the activator caps attached into the **DETAX cartridge holder**. Ensure that the notches in the cartridge holder and activator cap are correctly aligned (Fig. 3). Press the activator caps firmly onto the cartridge holder (Fig. 4 + 5); the **onetime perfect** foil bags are automatically pierced by a pin.

Fig. 2

3. Mixing and measuring

Insert the loaded **DETAX cartridge holder** into the dispensing unit. Start the unit and, if new foil bags are being used, dispense the material until it is homogeneous. Discard the initial inhomogeneous material. Then attach the dynamic mixer (Fig. 6) and secure the mixer with the locking plate on the cartridge holder. Fill the tray with the required amount by holding the tray at an angle and applying a little pressure. When loading the tray, keep the mixer in the material (Fig. 7). Use the filled mixer as a seal for the cartridge holder. Before reuse, loosen the locking plate to remove the mixer and check that the material can flow freely through the outlets in the activator caps. Attach a new mixer and ensure that it is in the correct position. Close the locking plate and proceed as usual, discarding the material dispensed initially. To remove the cartridge from the mixing unit, take hold of the top of the closed locking plate and remove the complete cartridge (Fig. 8).

Fig. 3

4. Procedure with the sandwich technique

Load the mixed **onetime perfect** onto the impression tray and cover it with a layer of **onetime perfect wash** (Fig. 9). After intraoral application of **onetime perfect wash** (Fig. 10), the loaded tray is positioned in the patient's mouth using light pressure for approx. 3 - 4 sec. Hold the impression in position without applying pressure until the **onetime perfect** has set.

Fig. 4

5. Disinfecting

After removing the impression from the mouth, rinse it under tepid running water. The impression can then be disinfected for 15 min. in 2% glutaraldehyde.

Fig. 5

6. Fabricating the model

The impression should be left for 15 min. before pouring. Dimensional stability is guaranteed for 14 days. We recommend using Class III and IV dental stone and standard model resins as model materials.

7. Cleaning the tray

Remove any set material mechanically with a blunt instrument. Place the tray in acetone or a similar solvent to soften the **Detaseal® Adhesive rapid**. Use only in well-ventilated areas. Then clean and disinfect the trays in the usual manner.

Fig. 6

Important working hints

- Ensure that solutions used prior to taking the impression (e.g. retraction solutions) have been completely removed, as they can affect the setting reaction.
- Do not use in conjunction with condensation-cured silicones.
- Latex gloves and surfaces contaminated with latex can affect the setting of **onetime perfect**.
- Remove any residual impression material from the oral cavity.
- Cured impression materials are resistant to chemicals – avoid stains on clothing.
- With the multiple mix technique we recommend using **onetime perfect wash** as a light bodied material for one-step impression taking.
- We recommend always using DETAX dynamic mixers to ensure optimum mixing.
- Adhere to the mixing unit manufacturer's instructions for use.
- Any large undercuts or wide interdental spaces should be blocked out in the usual way before taking an impression.
- When reloading the cartridge holder with new foil bags, always use base and catalyst with the same batch no.

Fig. 7

Fig. 8

Fig. 9


Fig. 10

Areas of application:

- Tray material for the sandwich and one-step putty-wash techniques
- Crown and bridge work
- Inlays, onlays, partial crowns
- Posts, veneers, abutments

Technical data:

DIN EN ISO 4823 - Type 1

- **Mixed volume:** 300 ml base 62 ml catalyst
- **Mixing ratio:** 5:1
- **Colour:** Base: light green Catalyst: white
- **Working time:** Approx. 1 min. 30 sec.*
- **Setting time:** Approx. 3 min. 30 sec.*
- **Intraoral setting time:** Approx. 2 min.
- **Strain in compression:** < 3.0 %
- **Recovery from deformation:** > 99,8 %
- **Linear dimensional change:** < 0.2 %
- **Application:** At 23°C ± 2°C / 73°F ± 4°F, 50 ± 5% relative humidity

* From beginning of mixing at 23°C ± 2°C / 73°F ± 4°F, 50 ± 5% relative humidity. Increased temperatures accelerate, decreased temperatures retard these times.

Ordering information:

onetime perfect foilbag

Basic set 02063
1 x 300 ml base
1 x 62 ml catalyst
1 cartridge holder
10 dynamic mixers

Refill kit 02068
2 x 300 ml base
2 x 62 ml catalyst

onetime perfect wash

Economy pack 02152
8 x 50 ml cartridges
48 mixing tips, pink
48 intraoral tips, white

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