FORMULATION

The unique resin based matrix contains BisEMA, UDMA, and the new Dimer Dicarbamate Dimethacrylate, DDCDMA. The filler system contains a radiopaque glass and ceramic filler ranging in size from nanosize to 13 microns (58wt%) with less than 1wt% catalysts, inhibitors and pigments.

DESCRIPTION

N'Durance® Dimer Core is a radiopaque dual-cure resin based dental restorative material containing Sodium Fluoride. The material is a creamy two paste system that is easily premixed and dispensed through a mixing tip and produces a tooth colored composite restoration. The light-curing cure occurs within 40 seconds and the chemical cure response will allow 2 minutes of working time, but will set completely in about 4 minutes even where the light does not reach. There are two shades of N'Durance® Dimer Core, "Natural" and "Bleach White". It contains a new, nanosize ytterbium trifluoride. The unique combination of monomers can achieve a high conversion.

INDICATIONS

N'Durance® Dimer Core is indicated for direct build-up of the crowns of teeth that have been badly broken down and are in need of an indirect restoration such as a crown in order to regain function. The "Natural" shade is especially indicated as a universal base coloring under esthetic full porcelain and composite crowns, onlays and inlays where the underlying tooth structure and core restoration will show through the esthetic veneer. The "Bleach White" shade is used where higher contrast with tooth structure is desired or where the final esthetics demands a white core. It can be used on vital or de-vital endodontically treated teeth. It may be used with pins and posts and should be bonded to the tooth with a dentin bonding agent.

DIRECTIONS FOR USE

1. The tooth should be cleaned by scaling and prophylaxis with flour of pumice.
2. Isolate the teeth. Use of a rubber dam is highly recommended.
3. Follow usual procedures for tooth or cavity preparation. Be certain the prepared cavity is clean and free of all debris.
4. Base the preparation as needed. Biodentine® or calcium hydroxide cement may be used for covering small pulpal exposures. Any additional dentinal coverage desired can be achieved by using a glass ionomer lining cement or other type of compatible cement. Eugenol containing cements must be avoided.
5. Pins and posts can be used to enhance the retention of the core build-up composite. Professional judgment should be used in the decision on when and where to place such additional retention mechanisms and on what kind of pin or post to use.
6. Apply a matrix band or removable shell crown around the preparation in the tooth as needed for application directly to the tooth.
7. Use a dentin/enamel bonding system to seal this composite to tooth structure. Alternative bonding systems can be used to bond this composite to restorative materials already found in or on the tooth.

8. Place the mixing/dispensing tip on the dual barreled syringe.
 - a) Remove the end cap carefully from the end of the dual barreled syringe by turning 90° counterclockwise and pulling the cap off. Do NOT discard the cap.
 - b) Line up the notch on the mixing nozzle with the corresponding notch on the syringe, push to place and turn 90° clockwise.
 - c) Attach the dispensing tip by pushing to place until it snaps into position. Dispense a tiny first increment of material onto a pad, to remove any initial inconsistencies in mixing, before dispensing directly into the tooth preparation.
 - d) After dispensing the N'Durance® Dimer Core, remove and discard the mixing nozzle from the syringe, and carefully re-align and replace the cap, being careful not to cross-contaminate the two exposed ends of the barrels.
9. Place composite in the cavity assuring good adaptation to walls of the preparation and to any pins and posts. Complete placement and freehand contouring within 2 minutes. The material will flow to place but can be contoured with little to no slumping with a dental hand instrument. A non-metallic instrument would help to avoid leaving any gray marks on any set composite material.
10. Light cure for 40 seconds with a conventional quartz-tungsten-halogen or LED dental curing light, minimum radiance of 400 mW/cm². Cutting the crown preparation can begin after the light-curing period. Complete setting of the material not reached by light will occur in 4 minutes.
11. Remove flash with a sharp instrument. Adjust N'Durance® Dimer Core to proper margins, and contours. Finish crown preparation using diamonds or carbide burs.

WARNINGS

- Dentists and assistants should wear gloves and protective eye wear. Patients also should wear eye protection.
- Uncured monomers in this resin material may cause an allergic reaction or skin irritation in some individuals. Do not use in individuals with allergies to these materials.
- Any non-sterilizable item that is handled in the delivery of the dental service should be disinfected by standard dental office hygiene procedures.
- Do not allow saliva or water contamination of material during the placement and setting procedure.
- Protect the composite material from direct light exposure from the operator light when removing the nozzle and replacing the cap.
- Do not use after expiration date.

STORAGE

- The system is designed to be stored between 15°C – 24°C / 59°F – 75°F.
- Do not store materials in proximity of eugenol-containing products.
- Do not store the materials in intense light, or under wet conditions.

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Nano-Dimer Conversion Technology®

Composite de reconstitution coronaradiculaire radio-opaque dual-cure

Made in USA
 Manufacturer/ Fabricant/
 Fabricante:
 SEPTODONT
 416S Taylor avenue,
 Louisville, CO80027, USA

01/12

C 05 14 273 00 00

N'Durance® Dimer Core

Nano-Dimer Conversion Technology®

Dual-Cure Radiopaque Core
 Build-Up Composite

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