F

## PROPRIÉTÉS

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UK

## PROPERTIES

Although careful mechanical root canal preparation is essential, it must be completed by chemical treatment.

One of the application techniques for this product consists in introducing into the canal a relatively large amount of antiseptic with an alkaline reaction, viz. sodium hypochlorite which, in addition to its bactericidal activity, has the capability to disintegrate necrotic tissues and to clear out the debris stuck to the canal walls.

- Irrigation should not be irritant to the periapex. That is the reason why the recommended sodium hypochlorite solution is diluted.
- So as to be effective, you should use a relatively large amount of liquid (10 to 20 ml) in order to remove all the debris from the canal walls and flush them away.
- To avoid injecting sodium hypochlorite into periapical tissues, it is recommended to use an appropriate needle so that the injected liquid may not, in any case, penetrate into the periapex.

The recommended needle has no opening at the tip: the vent is on the side of the gauge. The advantage of side-vented needles is the following: unlike irrigation needles with a distal opening, side-vented needles prevent the liquid under pressure from irritating the periodontal ligament, causing very painful acute periodontitis. The vent diameter is relatively small so that the pressure is sufficient to wash up the dental and pulpal debris that remained in the root canal.

## INSTRUCTIONS FOR USE

Suck the required amount of sodium hypochlorite out of the bottle. Then attach the needle to the syringe and carry out irrigation.

As with any other endodontic instrument, the needle should be carefully inserted into the canal so as not to injure the periapex during irrigation. The distance between the syringe tip and the apex should be at least 1.5 mm since the irrigation solution has a low surface tension and the exerted pressure is low. Before inserting the needle into the canal, it can be bent to conform to the canal shape. A perforated disk can also be slipped onto the needle so as to gauge the needle insertion depth.

## WARNINGS

KIT PARCAN instruments must not be heat-sterilized.

## PRESENTATION

The kit contains:

- 1 x 10 ml syringe
- 10 side-vented, blunted needles (Ø 0.50 x 40 mm)

Refill: 25 side-vented, blunted needles (Ø 0.50 x 40 mm)

*For professional use only in the practice of dentistry.*

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