F

#### COMPOSITION

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#### COMPOSITION

Paste A: butyleneglycol salicylate, zinc oxide, calcium phosphate, excipients.

Paste B: calcium hydroxide, zinc oxide, excipients.

#### INDICATIONS

Cavity liner.

#### PROPERTIES

Self-curing calcium hydroxide preparation indicated for direct or indirect pulp capping and cavity isolation, under filling materials, cements, composites and other materials.

This product acts as a protective barrier between the dentin and the pulp on the one hand, and acid generating materials on the other hand, and promotes the formation of secondary dentin in contact with the pulp.

It does not impede the polymerization of acrylic and composite restorations. The product obtained by mixing pastes A and B is radiopaque and hardens rapidly.

The composition gives the product optimal spreading properties.

#### INSTRUCTIONS FOR USE

Place an equal volume of paste A and B on the mixing pad. Mix rapidly (10-15 seconds) until a homogeneous mixture is obtained (even colour).

Immediately introduce the so-obtained mixture into the dry cavity using a ball-end instrument. The setting time in the mouth is about 2 minutes.

#### WARNINGS AND PRECAUTIONS FOR USE

- Do not change the ratio of paste A and B: the 1:1 mixture is designed for a mixing time of 15 seconds and a setting time of about 2 minutes.
- The setting time depends on the hygrometric conditions and the temperature. It is shorter in the mouth because of the high temperature and moisture content. Conversely, the setting time will be longer at low temperatures and moisture levels.
- Avoid touching the edges of the cavity or introducing an excessive amount of material. Any excess must be removed.
- This product cannot be used as a temporary filling material, and does not interfere with composites and other restorative materials.
- Avoid interchanging caps on the paste tubes.
- Any temperature rise is likely to cause deterioration of the products.

Xi – Irritant.

N – Dangerous for the environment.

Contains calcium hydroxide (Paste B / catalyst).

Toxic to aquatic organisms, may cause long-term adverse effects in the environment.

Irritating to eyes, respiratory system and skin.

Risk of serious damage to eyes.

In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

Wear suitable protective clothing, gloves and eye/face protection. In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible). Avoid release to the environment. Refer to special instructions/safety data sheets. This material and its container must be disposed of as hazardous waste.

#### STORAGE

Do not store above 25°C.

#### PRESENTATION

Kit containing:

- 1 tube of paste A: 13 g

- 1 tube of paste B: 11 g

- 1 mixing pad

*For professional dental use only*

#### Distributor:

SEPTODONT Ltd.

Units R & S - Orchard Business Center  
St Barnabas Close Allington, Maidstone  
Kent ME16 0JZ - UK

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