# DEI® EXPERIENCE

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# DEI® EXPERIENCE

THE NEW COMPOSITE SYSTEM FOR UNIVERSAL RESTORATIONS AND PREPARATIONS OF BRIDGES, CROWNS, INLAYS, FACETS, AND CHARACTERIZATION OF MOBILE PROSTHESES

Product for dental use  
Keep away from children  
Date of the specifications: 11/2010

## Characteristics

The planning of DEI® experience issues from the consideration that the natural tooth owns a complex chromatic characterization, due either to the shape and surface texture, or to the superimposition of different structures carrying various translucencies and opacities that must be “copied” for the purpose of restoring a correct visual perception of natural and aesthetic conditions.

The simple technique of stratification used with DEI® experience permits the achievement of ever got before aesthetic results.

Peculiarity of the product is the presence of three masses having different opacities, that enable the achievement of natural, highly aesthetic restorations, using an easy and quick technique.

The presence of the mass “enamel” with its special grade of transparency, available in 6 tints, gives the tooth the natural translucence brought about by passing of light, that is reflected on the inner mass of “dentin”, available in the 16 Vita® colours.

The supply of 6 cervical tints more opaque than the enamels and dentins, to be placed in the neck area and used in presence of reduced depths for easier making of the colour, increases the realistic chromatic effect of the natural tooth.

## Composition

Resins based on Bis GMA, UDMA, TEDMA, EBPA-DMA (22,5 %); vitreous fillers, pyrogenic silica (77%), catalysts, stabilizers, pigments and fluorescent agents.

## Technical data

Resistance to compression: 350 MPa  
Transverse resistance (flexion): 130 MPa  
Transverse resistance modulus: 7.400 MPa  
Water absorption: 25 µg/mm<sup>2</sup>  
Colour stability: optimum  
Polishing: high  
Radiopaqueness: 200% AL  
Cervical masses translucency: 20-35%  
Dentin masses translucency: 34-46%  
Enamel and incisal masses translucency: 53-62%  
Particle size: 0.02 - 1 micron

## Way of use

### Dentist® surgeries

Proceed as usual preparing the cavity, protecting the pulp, etching and applying the adhesive system. Use of the adhesive system Extra Bond (not included in the kit) is recommended. Apply the composite material stratifying the various masses. Every layer must be cured for 40 seconds for depths not thicker than about 1.5 mm. Photo-polymerize with curing lamps with light performance of at least 500 mW/cm<sup>2</sup>. If it is impossible to put the light source close enough to the material, it is necessary to lengthen the exposure time.

Halogen lamps are preferable because offering wider light spectrum compared to LED lamps, and consequently are efficient on any photosensible material.

The observance of photo-polymerization times and depths are fundamental for achieving complete conversion of the materials, and ensuring resistance and duration in the time. The composites harden also very few seconds, but are not polymerized. In this phase saving time could cause failure of the restoration.

Maximum Depth of polymerization:

for enamel mass: 3.5 mm  
for cervical mass: 1.7 mm  
for dentin mass: 2.5 mm  
Work time at full light: 60 sec.

## For laboratories

Experience proceeding on metal:

- Prepare a substructure identical to that used for metal/ceramic. Because of the highest adhesion got among the various components: metal / Extra Bond / opaque / DEI® experience composite, it is possible the use of micro-retentions positioned as rings.  
- Sandblast the substructure with aluminium dioxide, 90 microns grains.  
- Let fade carefully or clean with ethylic alcohol for 5 minutes under ultrasounds.

Since both the Extra Bond adhesive system and the experience opaque are self-and photo-curable, when the elements to treat are many, we recommend the repetition of the listed below procedures, rather than to increase the doses.

- Application of Extra Bond: mix for 15 sec. 1 drop of Extra Bond with 1 drop of Extra Catalyst and apply evenly with a dry small brush on the substructure. After 3 minutes proceed applying the opaque.

## Utilization of the opaque

- Deliver a drop of liquid for opaque on a suitable plate and add, mixing, the powder for opaque up to the wanted consistency. If you wish to increase the opaque hardness add one drop of Seal Coat Fast for each drop of opaque liquid in the mixing.  
- Considering the very high covering capacity of the experience opaque, the consistency must not be too much thick; apply the opaque with a clean small brush.

- Put the so prepared substructure for 3-8 minutes in the light oven (3 minut light strobo-plasma; 8 minut light alogenous) in a photo-polymerizer oven; afterwards, proceed with the application of the composite.

All these procedures (primer-opaque-composite) should be made uninterruptedly; when it is impossible, follow the listed below instructions:

- Reactivation of primer: carry out a second application of Extra Bond primer.

To get perfect adhesion between already cured and finished composite and composite not yet cured (for example in case of repairs): roughen a little the already cured composite, clean it with ethylic alcohol and clean brush; apply on the so prepared surface a thin layer of DEI® experience Seal Coat Fast; don't cure and add new composite.

The photo-polymerization times indicated refer to the experience Lamp Strobo oven and are valid also for the dental color xs kulzer. For other typologies of light and ovens, the times can considerably change.

## Finishing technique

Before the photo-polymerization of last layer of composite, apply DEI® Surface Hardener and cure for 40 seconds, wash and dry. The hardening of the composite shall be complete. Using this product avoids the necessity of eliminating the superficial layer of inhibition with burs or other products, with shorter work times.

Should increases of appearance and resistance of the reconstruction be wanted, apply DEI® experience Seal Coat Fast. Cure for 60 seconds and let polymerize for other 20 seconds. Take away the superficial layer of DEI® experience Seal Coat Fast using abrasive paste and small brush made with goat hair, small brush of cotton and Polish Paste.

After finishing, an extremely hard, resistant and shiny surface is obtained. Further, this technique considerably decreases the build-up of tartar and bacterial plaque on the restored tooth.

## Finishing in laboratory

After the photo-polymerization of DEI® experience Seal Coat Fast (in suitable light oven), it is necessary to spread with goat hair some clean and not recycled aluminium dioxide kneaded with clean water. The micro-engine must run at LOW NUMBER of revolutions. The ideal granulometry of the dioxide must be in the range of 50 - 90 microns. After elimination of DEI® experience Seal Coat Fast (the dry surface must not sparkle) apply the cotton brush with Polish Paste.

## Warnings

In restorations with composite don't use products containing eugenol. Keep refrigerated at constant temperature of 4°C, or in environment with fresh and constant temperature lower than 25°C. Take the pack out from the refrigerator some hours before use. Don't expose to sources of light or heat. After use seal again the syringe with its suitable cap. To release the pressure of the material, turn the piston halfway anticlockwise. Avoid contact with the skin, as in some cases a hypersensitivity may occur: in these cases don't use. For the utilization of the accessory products, read always the instructions for use of the specific product.

Attention: the liquid for opaque contains substances that must be handled with care, and is inflammable. Don't smoke during use and don't keep close to sources of heat; avoid electrostatic charges. Keep in airy environment and with cap well sealed. In case of contact with skin, wash with plenty of water and soap; in case of contact with eyes, rinse with plenty of water and consult a specialist. The peculiar characteristics of DEI® experience allow modifications of the viscosity of the product either heating the masses for making it more fluid, or cooling down them for making it more compact.

## Conservation and validity

If correctly kept at temperature lower than 25°C, the expiration date of the product is 3 years from the date of production. Correctly kept far from damp and at temperature lower than 25°C, the liquid and powder opaquer have an expiration date of 2 years.

\* trade mark registered by Vita Zahnfabrik

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