Light curing dental adhesive  
Lichthärtendes Dentaladhasiv  
Adhésif dentaire photopolymérisable  
Adesivo dentale fotoindurente  
Adhesivo dental fotopolimerizable

# Fiber-Bond

Directions for Use  
Gebrauchsanweisung  
Mode d'emploi  
Istruzioni per l'uso  
Instrucciones de uso

REV. 01.07.2013

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**1) Description of the medical product Fiber-Bond**

Light curing dental adhesive.

**2) Composition**

Bis-GMA, TEGDMA, Dodecane diol dimethacrylate, HPMA, Butane diol dimethacrylate.

**3) Indications**

Dental adhesive (bonding agent) for Fiber-Splint, Fiber-Splint ML and Fiber-Splint Ortho when using the enamel etching technique. Wetting agent when using Fiber-Splint, Fiber-Splint ML, Fiber-Splint Lab Ret and Fiber-Splint Lab Bridge for repairing and/or strengthening dental acrylic.

**4) Side-effects**

Although, to date, Fiber-Bond has not been observed to have a toxic effect on the pulp, irritation of the pulp cannot be ruled out if Fiber-Bond is placed in direct contact with the dentine.

**5) Interactions**

Preparations containing eugenol inhibit polymerization.

**6) Warning**

As Fiber-Bond is blue-sensitive, polymerization is triggered off by daylight or the operatory light (blue filters slow this process down considerably).

Light curing materials should therefore not be exposed to light unnecessarily (remove the material from the bottle prior to application and close it immediately afterwards) and the material should be placed and contoured rapidly.

**7) Application**

Soak the glass fiber strips (Fiber-Splint, Fiber-Splint ML, Fiber-Splint Ortho, Fiber-Splint Lab Ret, Fiber-Splint Lab Bridge) thoroughly with Fiber-Bond. The fabric must be fully soaked in liquid. The impregnated glass fiber strip is then cured with blue light. The hardening of the fiberglass ribbon (respectively of the Fiber-Bond) depends from the polymerisation light's performance. As a rule, a polymerisation of 20 seconds is sufficient for a strip length of 1 cm, however approx. 30-40 seconds are recommended to ensure evenness. This evaluation is based on a light intensity of approx. 800 mW/cm<sup>2</sup>, which is the usual intensity of lamps at a bandwidth of 400-500 nm (blue light). For polymerisation with LED light, a slightly higher power is recommended, as many LED units have one light of 430 nm only.

**8) Emergency measures, symptoms and antidotes**

Avoid contact with the skin or eyes. Should an accident occur, rinse immediately with copious amounts of water and, if necessary, consult a doctor.

**9) Shelf-life**

See indications on the packaging.

**10) Storage**

Refrigeration (+5°C/+8°C) is recommended, but do not exceed 25°C.

**11) Presentation and packaging**

3.0 ml and 5.0 ml bottles - available separately and in assortments.

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