Fig. 1


Fig. 2


Fig. 3

Range of indications

Sealing in caries prophylaxis, conditioning the enamel surfaces with the etch technique in conjunction with **smartseal® & loc/F**.

1. Preparation of the application syringe

Remove the sealing cap of the application syringe, attach one of the blue application needles and turn the needle to lock it in position. If necessary, bend the needle a little to facilitate application.

CAUTION! Before using for the first time, dispense a small amount of **smartseal® etch** from the application syringe to prevent any etch gel spurting out due to air blows in the application syringe.

2. Application

Apply **smartseal® etch** to the tooth enamel only, avoiding contact with the dentine and edges of the tooth. Allow the etching gel to take effect for a minimum of 15 - 30 sec. (Fig. 1). We recommend increasing the reaction time to 90 - 120 sec. for deciduous teeth which have been treated with fluoride lacquer. At the end of the reaction time the etching gel is rinsed off with running water for approx. 20 seconds (Fig. 2). The etched area is then dried with an air jet uncontaminated with water or oil (Fig. 3). The etched enamel should then have a milky, calcium white colour. Close the **smartseal® etch** properly after use.

Important working hints

- Apply using the blue application needles. These slender needles ensure safe and precise application.
- If the etched surface is contaminated, e.g. with saliva, it should be re-etched, then rinsed again with running water and dried as described above.
- Avoid **smartseal® etch** coming into contact with the mucosal tissue, eyes or skin. On accidental contact, rinse thoroughly with plenty of running water.
- Before using for the first time, dispense a small amount of **smartseal® etch** from the application syringe to prevent any etch gel spurting out due to air blows in the application syringe.
- Adhere to the manufacturer's instructions for use when using **smartseal® etch** in conjunction with adhesives, bonders and other commonly used materials.
- In case of following application of **smartseal® & loc/F** the dentine surface should not be dried finally, it must remain moist.

Technical data:

- **Volume:** 2 ml
- **Colour:** transparent blue
- **Dosage:** direct application using a syringe and application needles or application brush
- **Reaction time:** 15 - 30 sec.
- **Application:** at 23°C ± 2°C / 73°F ± 4°F, 50% ± 5% relative humidity
- **Storage:**

Ordering information:

smartseal® etch

Standard packing 02326

1 x 2 ml smartseal® etch application syringe
10 application needles, blue

Application needles 02458

25 needles, blue

smartseal® & loc

Combi Pack 02612

1 x 1,5 g smartseal® & loc application syringe
1 x 1,5 g smartseal® & loc F application syringe
12 application needles, black, long

smartseal® & loc

transparent

Standard packing 02573

2 x 1,5 g smartseal® & loc application syringe
12 application needles, black, long

smartseal® & loc F

natural opaque

Standard packing 02581

2 x 1,5 g smartseal® & loc F application syringe
12 application needles, black, long

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