DEI® RAINBOW EASYTEMP

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Barium glass silica, ca-
talyist, stabilizer, pigments,
polyfunctional acrylates
and polyfunctional me-
thacrylates, unsaturated
esteri, malonyl urea deri-
vatives.

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DEI® RAINBOW EASYTEMP

DUAL CURE COMPOSITE FOR TEMPORARY BRIDGES AND CROWNS

Available packages

One 75 g cartridge in the colours:

A2 (COD. EASY 2)

A3 (COD. EASY 10)

A3,5 (COD. EASY 3)

50 mixing cannulas (COD. EASY 4)

Dispenser gun (COD. EASY 5)

Product for dental use

Keep away from children

Date of the specifications: 04/2011

Characteristics

DEI® Rainbow Easytemp is a self mixing, dual curing material consisting of two components. It is indicated for production of temporary crowns and bridges.

It offers very important advantages:

1. It swiftly and easily meets the requirements of aesthetic and functional temporary items.
2. The patient feels gratified by the prompt operation.
3. Working time is only few minutes (only 1-2 minute in the mouth).
4. After extraction from oral cavity it maintains during 5-6 minutes

an elastic consistency useful for easier removal and elimination of eventual surplus.

5. As the duration of work is determined by self-mixing and capacity of photopolymerization, it is very easy to use.
6. It doesn't heat in the hardening phase (38°C).
7. Highly resistant to abrasion and flexure.
8. Top protection of the preparation.
9. Best polishing.
10. In case, repairable using the same material as a glue, after treatment of the two parts to connect with ethylic alcohol and DEI® Seal Coat Fast (not cured).

Composition

Barium glass silica, catalyst, stabilizer, pigments, poly-functional acrylates and polyfunctional methacrylates, unsaturated esters, malonyl-urea derivatives.

Technical data

Resistance to compression: 290 MPa

Resistance to flexure: 80 MPa

Way of use

1. Take the impression either before the preparation of the stumps, or before the extraction; eliminate the interdental details from the impression.
2. Control that the outlets on the cartridge are free from hardened material; eventually, remove it using a probe. The material must be extruded slowly and continuously.
3. Insert the cannula on the cartridge

rotating fully. The cannula shall remain inserted on the cartridge up to the next use, acting as a top.

4. After discard of the first 2 mm, distribute the product inside the impression; place the impression inside the oral cavity within 30" and leave it there 1-2 minutes.
5. Extract the provisional that will be elastic for other 3-4 minutes and proceed to elimination of excesses and proximal undercuts. Wanting immediate hardening, photopolymerize every element for 60 sec with halogen light.
6. Take out the surface inhibiting layer with ethylic alcohol or acetone.
7. Finish using burs or rubber tops.
8. To get excellent aesthetic performance, apply DEI® Seal Coat Fast (not included in the package) and photopolymerize. Take out the surface part using abrasive paste and goat hair brush.
9. Fix with DEI® Easycem Hard.

Cautions

In case of individual sensitivity to product's components don't use.

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Composición

Barium glass silica, catalyst, stabilizer, pigments, polyfunctional acrylates and polyfunctional methacrylates, unsaturated esters, malonyl urea derivatives.

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