Ball-bearing turbines
Turbine a cuscinetti a sfere
Turbine mit Kugellagern
Turbines à roulements à billes
Turbinas de cojinetes de bolas
Turbinas de rolamentos de bolas
Γωνιακή αερότορ με ρουλεμάν

Operating instructions
Istruzioni per l'uso
Gebrauchsanleitung
Mode d'emploi
Modo de empleo
Instruções de uso
Οδηγίες χρήσεως

TKD

TeKne Dental

REF 186.30
REF 186.40

CE
0051

Fig. 1


Fig. 2

Fig. 3

Fig. 4

ACCESSORIES AND SPARE PARTS / ACCESSORI E RICAMBI / ZUBEHÖR UND ER-SATZTEILE / ACCESSOIRES ET RÉCHANGES / ACCESORIOS Y RECÁMBIOS / ACESSÓRIOS E PEÇAS SOBRESSELENTES / ΕΞΑΡΤΗΜΑΤΑ ΚΑΙ ΑΝΤΑΛΛΑΚΤΙΚΑ

REF 186.90	Rotor / Rotore / Rotor / Rotor / Rotor / Rotor / Ρότορας
REF 120.00	SO2000 Service Oil / Lubrificante / Schmiermittel / Lubrifiant / Lubricante / Lubrificante / Λάδι λιπανσης
REF 120.01	Nozzle / Beccuccio / Kanüle / Buse / Boquilla / Bocal / Ρύγχος
REF 129.40	Air pressure gauge for 4-hole connection / Manometro per attacco 4 vie / Manometer für 4-Wege-Anschluss / Manomètre pour raccordement 4 voies / Manómetro para enlace 4 orificios / Manómetro de 4 vías / Μανόμετρο για κορδόνι με 4 οπές

DESCRIPTION

High-torque ball-bearing turbine for professional dental use, compliant to ISO 7785-1 standard.
 Weight: 59 g.
 Maximum speed of rotation: 355000 rpm.
 Maximum noise level: 66 dBA (@2,3 bar).
 Instrument classification: Type 2 according to ISO 7785-1.

This medical device meets the requirements of the European Directive 93/42 EEC (Class IIa).

CONNECTION

The THALIA® turbine with 3-hole connection meets the requirements of Standard ISO 9168 Tipo A.
 The THALIA® turbine with 4-hole connection meets the requirements of Standard ISO 9168 Tipo B.
 The turbines must be directly connected to the supply hose.

AIR PRESSURE

Recommended air pressure: 2,1 ÷ 2,3 bar (30,5 ÷ 33,5 psi).
 Air consumption: 46 NI/min (@2,3 bar).
 Check air pressure by means of a gauge (REF 129.40) connected between the turbine and the supply hose (Fig. 2).
 Supplied air must be perfectly dry and well-filtered, kept at constant pressure.

SPRAY

One-jet type.
 Water pressure: 0,8 ÷ 2,5 bar (11,6 ÷ 36,3 psi). Water must be well-filtered.
 Spray air pressure: 1,0 ÷ 4,0 bar (14,5 ÷ 58,0 psi).
 To clean the spray duct, use the supplied wire cleaner, as shown on Fig. 4
 The turbine is equipped with a non-return valve for water.

BUR TYPE

Shank diameter: Ø1,60 mm, Type 3 according to ISO 1797.
 Maximum length: 25 mm, code 4 or 5 or 6 (6th digit) according to ISO 6360-1.
 Minimum fitting length of shank: 11 mm.

BUR CHANGING

The turbine is equipped with a sophisticated push-button chuck: press the button with the thumb and subsequently insert or extract the bur (Fig. 1).

⚠ **When fitting on a bur make sure it is fully inserted.**
 Before inserting a bur, make sure the bur shank is clean.
 The turbine must not be started without a bur inserted into the chuck.

⚠ To avoid overheating of the push-button, which could cause burns, the button should not be pressed while the turbine is running. Internal tissues (tongue, cheeks, etc..) must then be protected by a retractor or a dental mirror.

LUBRICATION BY OIL

It is possible to lubricate instrument by means of the SO2000 service oil (REF 120.00) carrying out the following operations:
 1) Disconnect turbine from the supply hose;
 2) Insert the special nozzle (REF 120.01) onto the lubricant can;
 3) Holding turbine by means of a cloth, fully insert the nozzle into the turbine (Fig. 3). Lubricate by pushing once only for one second.

In order for the instrument you have purchased to be helpful in your profession for a very long time, carefully repeat the described operations **before every sterilization cycle and, anyway, at least TWICE A DAY.**

STERILIZATION

⚠ **The turbine can be sterilized in autoclave at 135 °C / 2 bar (275°F / 29 psi).**
 Before sterilization, remove residues and carefully clean the instrument.
 The device is supplied "not sterile".

CLEANING AND DISINFECTION

Clean with a cloth soaked in a disinfectant containing:
 - Glutaraldehyde up to 3%
 - O-phenylphenol
 - Alcohols based on ethanol
 Products not recommended: disinfectants containing benzalkonium chloride or acetone.
 The instrument cannot be immersed in any kind of liquid.

SERVICE

In case of irregular noise, head vibrations, bur swinging or any other malfunction of the instrument, and for any overhaul and repair work, please contact your usual supplier, an authorized service centre, or directly to TKD.

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