Eigenschappen van de zool:

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
Mode of action:

Pressure relief

The 4 twin chambers are designed to support optimally the natural rolling action of the foot.

When walking

When walking the air cushions first absorb the weight. The air flows to the forefoot to support the vault of the foot. The heel gets in full contact with the floor.

When rolling

When rolling the foot the air flows backwards through the twin chambers and lifts the heel, thus the foot arch is optimally supported.

Structure of the insole:

- Altogether 4 air chambers placed in pairs
- Wear-resistant textile surface
- Color: beige, 5 sizes
- Machine washable at 30°C

You will feel the advantages:

- Walking becomes softer, the cushioning takes pressure off the foot, knee and hip joints
- Good massage effect thanks to the air exchange from one air chamber to the next (additional reflex zone massage)
- Support of the natural rolling action of the foot by a conducted air exchange between the twin chambers
- Reduces the tiredness of your feet
- Even pressure on the insole of your foot
- Effective alleviation of pressure (no more pressure points)
- Air cushions have an insulation effect at low temperatures

Instructions for use:

The air cushion insoles are marked “L” for left and “R” for right side. Place the insoles correspondingly into your shoes so that the air chambers face downwards. In sports shoes or hiking shoes place the insoles under the removable shaped insoles of your shoes. Arch supports in shoes have to be removed in order to achieve correct positioning of the insoles.

Indications:

Shock absorption.

Contraindications:

Sensory loss and blood flow disorders in the affected body part, skin disorders in the treated part of the body, Diabetes mellitus

Washing instructions:

Please take the air cushion insoles out of your shoes to air them whenever you used them. Wet insoles must never be dried on the heater. That could lead to damages of the material. From time to time you may wipe clean the air cushion insoles with a wet cloth or wash them in the washing machine at a temperature of 30 °C.

Important advice: The heel caps of the shoes should be high enough to avoid slipping out of the heels. Avoid direct contact of the air cushion insole with injured skin! Little grains of sand, stones and other dirt have to be removed from the shoes as otherwise they might damage the material. Do not cut the air cushion insoles.

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