# Prophylux

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# Prophylux

Light-curing composite used preventively for filling  
pits and fissures of premolar and molar grooves

## Composition

Bis-GMA, triethylene glycol dimethacrylate, glyceryl propoxytriacrylate,  
sodium fluoride, excipient.

## Indications

Filling pits and fissures of molar and premolar grooves.

## Instructions for use

### **Tooth preparation :**

- Prepare the tooth with a non-fluoride, oil-free paste. Wash with an oil-free air/water spray.
- Isolate the tooth with cotton rolls or rubber dam.
- Blow-dry completely with an oil/moisture-free air source.

### **Etching :**

- Apply the etchant to the surface of the tooth with a brush. Do not rub the enamel.
- Allow the etchant to remain for 30 seconds.
- Wash area completely for 20 seconds. Saliva must not touch the etched area.
- Blow-dry with oil/moisture-free air.
- Etched surfaces of the tooth should have a dull, white frosted appearance. If not, repeat the process allowing etchant to remain on the tooth for an additional 15 seconds. Several repetitions of this etching procedure may be needed especially on primary teeth or teeth high in fluoride content.

### **Application of the composite :**

- Place 2 drops of PROPHYLUX on the groove to be filled.
- Using a brush, apply the sealant to each etched dry occlusal surface, taking care to avoid air entrapment.
- Expose all coated areas to a light-curing source for 40 seconds.

### **Finishing :**

- After curing, wipe the occlusal surfaces with a wad of cotton to remove the non-polymerized layer from the surface of the composite.
- Check for complete coverage with an explorer.
- Incomplete coverage requires re-etching and reapplication of PROPHYLUX.

## Shelf-life and storage

- This system is designed to be used at room temperature (21-24°C). Shelf-life at room temperature is two years for unopened syringes and one year for opened syringes when caps are replaced immediately after dispensing.
- An opened syringe can be refrigerated to 4°C to extend shelf life. If stored in refrigerator, the syringe should be taken out to equilibrate to room temperature prior to use.
- Do not store PROPHYLUX in the proximity of eugenol-containing products.
- Do not expose PROPHYLUX to intense light or to high temperature.

## Presentation

2 syringes of 3.5 g composite,  
20 tips.

*For professional use only in the practice of dentistry.*