EFO

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◆ Content

- 6 Description
- 6 Indication
- 6 Contraindication
- 6 Function
- 7 Measurement/Size
- 7 Application
- 8 Important information
- 9 Material
- 9 Product care/Cleaning/
Maintenance
- 9 Storage/Disposal

◆ Description

EFO is a bandage for elastic guidance and stabilisation of the hip joint. The leg is guided in its movements. Positions that might damage the joint are resisted by the bandage and thus avoided by the patient.

The degree of stabilisation is individually defined by the continuously adjustable circular closure and the lateral coil springs.

The EFO is supplied in 4 sizes, for the left or right side only, and for both sides.

- Elastic guidance of the leg movements
- Stabilisation continuously adjustable
- Invisible underneath the patient's clothes

◆ Indication

- Stimulus situations with arthrosis/
arthritis
- Postsurgical

◆ Contraindication

- A lack of patient compliance

◆ Function

- The hip is stabilised by an even, circular compression of the soft parts

◆ Measurement/Size

EFO						
REF left side only	REF right side only	REF both sides	Size	Height/cm	Hip circumference/cm	Thigh circumference/cm
40 367 101	40 367 201	40 367 301	1	48	60-80	40-55
40 367 102	40 367 202	40 367 302	2	48	70-100	50-65
40 367 103	40 367 203	40 367 303	3	48	90-110	55-70
40 367 104	40 367 204	40 367 304	4	48	100-130	60-75

◆ Application

Measurements required:

Hip circumference/cm (circumference 1)

Thigh circumference/cm (circumference 2)

1. The bandage must be applied without crinkles and without the edges pressing into the skin. The lateral pressure plate with the coils lies to the side of the hip joint.
2. The degree of compression is correct when it has a stabilizing effect on the hip joint without restraining the patient's movements more than necessary.
3. In rare cases, individual adjustment to the patient may be required.
4. The degree of stabilisation can also be regulated by removing some of the coil springs. Patients who spend most of their time in a sitting position may not need the anterior coil spring.

◆ Important information

- Before using this device an individual adaptation of the orthosis and a familiarisation in its use by an experienced orthopaedic technician according to the medical instructions is essential.
 - Incorrect selection, adaptation, application and/or failure to monitor the correct fit of this device can cause health damages.
 - Never wear the orthosis on bruised skin.
 - Check orthosis for damage (e.g. loose connections, fasteners, belts, etc) each time it is put on.
 - Should any performance malfunctions (e.g. apparent reduction in or lack of corrective action) be suspected, immediately call at an orthopaedist or an orthopaedic technician.
 - The orthosis must be checked for function, correct fit and possible damage by an orthopaedist or orthopaedic technician at regular intervals. The intervals between checks are to be individually determined for each patient.
 - Should the product have been subject to disproportionately strong forces, it must be checked for possible damages by an orthopaedic technician before further use.
 - Sharp edges or the like, which may occur when machining the product, must be deburred in order to avoid skin injuries.
 - Avoid excessive heat accumulation. Wear cotton clothes under the orthosis if necessary.
- Due to the required individual adaptation (productdependent forming and/or machining) of the orthosis it is not suitable for reuse by other persons!

◆ Material

- Cotton
- Spandex
- Rayon
- Polyamide compounds
- Foam Latex

◆ Product care/Cleaning/ Maintenance

- Please close the Velcro closures before washing
- Hand wash
- Do not iron
- Do not tumble dry
- Do not bleach
- Do not dry-clean

◆ Storage/Disposal

Product does not contain any environmentally hazardous substances. It can be disposed of with household refuse.

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Wilhelm Julius Teufel GmbH