tomas®-uprighting spring

CE 0483

Gebrauchsanweisung
Instructions for use
Mode d'emploi
Modo de empleo
Modalità d'uso

D
DENTAURUM

Turnstraße 31 · 75228 Ispringen · Germany · Telefon +49 72 31 / 803-550
Fax +49 72 31 / 803-295 · www.dentaurum.de · E-Mail: info@dentaurum.de

Germany 02/10BR1-3
Printed by Dentaurum
989-651-00


Fig. 1:
Aufrichten ohne
Intrusion und Extrusion
Uprighting without
intrusion or extrusion
Redressement sans
intrusion ni extrusion
Enderezamiento sin
intrusión ni extrusión
Uprighting senza intru-
sione ed estrusione


Fig. 2:
Aufrichten ohne
Intrusion und Extrusion
Uprighting without
intrusion or extrusion
Redressement sans
intrusion ni extrusion
Enderezamiento sin
intrusión ni extrusión
Uprighting senza intru-
sione ed estrusione


Fig 3a:
Aufrichten mit
Intrusion, Variante 1
Uprighting with
intrusion, version 1
Redressement avec
intrusion, variante 1
Enderezamiento con
intrusión, variante 1
Uprighting con intru-
sione, variante 1


Fig. 3b:
Aufrichten mit
Intrusion, Variante 2
Uprighting with
intrusion, version 2
Redressement avec
intrusion, variante 2
Enderezamiento con
intrusión, variante 2
Uprighting con intru-
sione, variante 2


Fig. 4a:
Aufrichten mit
Extrusion, Variante 1
Uprighting with
extrusion, version 1
Redressement avec
extrusion, variante 1
Enderezamiento con
extrusión, variante 1
Uprighting con estru-
sione, variante 1


Fig. 4b:
Aufrichten mit
Extrusion Variante 2
Uprighting with
extrusion, version 2
Redressement avec
extrusion, variante 2
Enderezamiento con
extrusión, variante 2
Uprighting con estru-
sione, variante 2


Fig. 5:
Aufrichten mit mesia-
lem Wurzel torque
Uprighting with mesial
root torque
Redressement avec tor-
que radiculaire mérial
Enderezamiento con
torque mesial de la raíz
Uprighting con torque-
mesiale della radice

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Stand: 02/08
Hotline: +49 72 31 / 803-550

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Instructions for use for the tomas®-uprighting spring REF 302-009-00

Dear Customer

Thank you for having chosen a quality product from Dentaurum.

In order to use this product at its best for you and your patients, it is important to study and follow these directions for use carefully.

The written instructions cannot cover all eventualities during operation. For this reason our Hotline is available to answer any other questions and ideas that may arise.

Due to a constant development we recommend, even when you use the same products frequently, to study the enclosed latest issue of the instructions for use or refer to our website at www.dentaurum.de.

1. Manufacturer

Dentaurum GmbH & Co. KG · Turnstraße 31 · 75228 Ispringen · Germany

2. With reference to our quality standards

Dentaurum ensures a faultless quality of the products manufactured by us. These recommendations are based upon our own experiences. The user himself is responsible for the processing of the products. Responsibility for failures cannot be taken, as we have no influence on the processing on site.

3. Brief description

The tomas®-uprighting spring is for uprighting teeth. The uprighting action can be combined with intrusion or extrusion, depending on how the spring is set. The tomas®-uprighting spring consists of a super-elastic nickel titanium spring (dimensions 0.42 mm x 0.62 mm; 16 x 24) and a steel wire (dimensions 0.46 mm x 0.64 mm; 18 x 25). Both wire units are connected with a clamping tube. The nickel titanium spring is pushed into the clamping tube which is then crimped together to fix the wires in the desired position.

If the tooth is tilted at approx. 40° the uprighting moment is between 10 and 20 Nmm. Depending on how the Alpha curve is formed, the tomas®-uprighting spring has a large plateau in the range of 8 to 15 Nmm and creates an intrusive force of approx. 0.5 to 1.0 N.

The tomas®-uprighting spring is for single use only. Recycling the tomas®-uprighting spring or repeated use on a patient is not permitted.

4. Application

For uprighting lower molars, it is advisable to place a tomas®-pin or tomas®-pin SD between the first and second lower premolars. The pin head must be in close proximity with the attached gingiva.

The tooth being uprighted can be provided with a buccal tube for the adhesive technique or conventionally with a band and welded buccal tube. Depending on the structure of the tooth, the fit of the buccal tube for the adhesive technique, the strength of the adhesive bond and the power of the transferring forces, a buccal tube may detach from the tooth.

Step 1: Fitting the tomas®-uprighting spring

The effectiveness of the nickel titanium spring can be influenced by sliding it in the clamping tube. If the nickel titanium spring is slid as far as possible into the tube (refer to Fig. 1, distance C₀), the maximum alpha and beta moments are attained. The same applies to the intruding and extruding forces. Drawing the nickel titanium spring from the tube lowers the alpha and beta moments as well as the intruding and extruding forces (C < C₀). Once the desired position has been set, the tube is crimped with pliers to fix the nickel titanium spring in place.

Slide the nickel titanium spring up to its first curve into the slot of the buccal tube. Hold the tomas®-uprighting spring parallel to the occlusal plane. Drop a perpendicular from the vertical slot in the tomas®-pin or tomas®-pin SD to the steel wire and mark the wire accordingly (Fig. 1). Depending on the objective of the treatment (refer to step 2), the steel wire is then curved at the mark and shortened as required.

Step 2: Setting the desired treatment objective

The tomas®-uprighting spring can be used for uprighting teeth and intruding, extruding or distalizing them simultaneously. Teeth can also be uprighted by applying mesial root torque only.

Distance A₀ is the horizontal distance between the vertical slot in the tomas®-pin and the mesial edge of the buccal tube.

Distance B₀ is the vertical distance between the steel wire and the occlusal edge of the slot in the pin. Angle alpha is located occlusodistally from the intersection between the perpendicular of the pin slot and the steel wire. In the initial condition $\alpha = 90^\circ$ (Fig. 1). The individual values can be varied to produce different results.

Uprighting (Fig. 2)

When uprighting teeth, the above mentioned values should remain:

$$\alpha = 90^\circ$$

$$B = B_0$$

$$A = A_0$$

Uprighting, with intrusion

When uprighting and intruding teeth, the following parameters should be modified:

Version 1	Version 2
$\alpha > 90^\circ$ (max. 135°)	$\alpha = 90^\circ$

$$B < B_0$$

$$A = A_0$$

$$A = A_0$$

(Fig. 3a) (Fig. 3b)

Uprighting, with extrusion

When uprighting and extruding teeth, the following parameters should be modified:

Version 1	Version 2
$\alpha < 90^\circ$ (min. 45°)	$\alpha = 90^\circ$

$$B > B_0$$

$$A = A_0$$

$$A = A_0$$

(Fig. 4a) (Fig. 4b)

Uprighting teeth, with mesial root torque (Fig. 5)

When uprighting teeth, the original values should remain:

$$\alpha = 90^\circ$$

$$B = B_0$$

$$A = A_0$$

If an additional ligature is attached, the tooth can be uprighted with mesial root torque.

Step 3: Fitting the tomas®-uprighting spring

The nickel titanium spring must be fixed securely in the clamping tube. To facilitate care by the patient and depending on the setting, it is advisable to remove the nickel titanium wire overlapping the steel wire.

Slide the nickel titanium spring into the buccal tube. Place the angled steel wire in the slot of the tomas®-pin or tomas®-pin SD and secure it in position with light-curing adhesive such as ConTec LC (REF 163-510-00). Reactivation during the treatment is usually not required.

Step 4: Removing the tomas®-uprighting spring

As soon as the desired treatment objective has been attained, the tomas®-uprighting spring is removed by carefully squeezing the adhesive on the head of the tomas®-pin or tomas®-pin SD with Weingart forceps (REF 003-120-00) until the adhesive releases. The steel wire is then removed from the slot and the NiTi spring pulled out of the slot in the buccal tube.

If the tomas®-pin or tomas®-pin SD is not needed for further retention purposes, it can also be removed. Please refer to the instructions for use of the tomas®-pin and tomas®-pin SD (REF 989-631-20).

5. Explanation of symbols

	Non-sterile
	Batch description
	Not for reuse
	Refer to instructions for use
	Order number

To date: 02/08
Hotline: +49 72 31 / 803-550

Subject to modifications