Compex®

www.compex.info

Compex®

duofit

User and applications manual

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**Compex Médical SA
Switzerland
www.compex.info**

English

Directions for use	10
Table with placements of electrodes and stimulation positions	19
Programs and specific applications	21

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CONTENTS

I. WARNINGS	4
1. Contra-indications	4
2. Safety measures	4
II. PRESENTATION	7
1. Reception of equipment and accessories	7
2. Presentation of the device	7
3. Guarantee	8
4. Maintenance	8
5. Storage and transportation conditions	8
6. Use conditions	9
7. Elimination	9
8. Standards	9
9. Patent	9
10. Normalized symbols	9
11. Technical information	10
III. DIRECTIONS FOR USE AND INSTRUCTIONS	10
1. Inserting and replacing the batteries	10
2. Connections	10
3. Preliminary settings: language, light intensity and sound volume	11
4. Selecting your motivation	11
5. Selecting a program	12
6. Personalizing a program	12
7. During the stimulation session	13
8. Electricity consumption	14
9. Problems and solutions	15
IV. HOW DOES ELECTROSTIMULATION WORK ?	15
1. Motor nerve stimulation (EMS)	16
2. Stimulation of the sensitive nerves	16
V. USAGE PRINCIPLES	17
1. Placement of electrodes	17
2. Stimulation position	18
3. Setting the stimulation energy level	18
4. Progression of the levels	18
VI. TABLE WITH PLACEMENTS OF ELECTRODES AND STIMULATION POSITIONS	19
VII. PROGRAMS AND SPECIFIC APPLICATIONS	21
Your motivation: Fitness	22
Your motivation: Wellness	26

I. WARNINGS

1. Contra-indications

Important contra-indications

- Cardiac stimulator (pacemaker)
- Epilepsy
- Pregnancy (do not position in the abdominal area)
- Serious arterial circulatory problems in the lower limbs
- Abdominal or inguinal hernia

Precautions when using the Compex

- After trauma or a recent operation (less than 6 months)
- Muscular atrophy
- Persistent pains
- Need for muscular rehabilitation

Osteosynthesis equipment

The presence of osteosynthesis equipment (metallic equipment in contact with the bone: pins, screws, plates, prostheses, etc.) is not a contra-indication for the use of Compex programs. The electrical current of the Compex is specially designed to have no harmful effect with regard to osteosynthesis equipment.

In all of the above cases, you are advised to:

- Never use the Compex for prolonged periods without medical advice;

- Consult your doctor if you are in any doubt whatsoever;
- Read this manual carefully, in particular chapter VII, which provides information concerning the effects and indications of each stimulation program.

2. Safety measures

What you should not do with the Compex

- Do not use the Compex in water or in a humid atmosphere (sauna, hydrotherapy, etc.).
- Never carry out an initial stimulation session on a person who is standing. The first five minutes of stimulation must always be performed on a person who is sitting or lying down. In rare instances, very highly strung people may have a vagal reaction. This is of psychological origin and is connected with a fear of the muscle stimulation as well as surprise at seeing one of their muscles contract without having intentionally contracted it themselves. This vagal reaction produces a feeling of weakness with a fainting tendency, slowing down of the heart rate and a reduction in arterial blood pressure. If this does occur, all that is required is to stop the stimulation and lie down with the legs raised for the time it takes for the feeling of weakness to disappear (5 to 10 minutes).
- Never allow movement resulting from muscular contraction during a stimulation session. You should always stimulate isometrically; this means that the extremities of the limb in which a muscle

is being stimulated must be firmly fixed, so as to prevent the movement that results from contraction.

- Do not use the device within 1.5 metres from a short-wave device, a micro-wave device or high frequency (H.F.) surgical equipment to avoid the risk of causing skin irritations or burns under the electrodes. If you are in any doubt about using the Compex in close proximity to any other medical device, please consult its manufacturer or your doctor.

- Only use the electrode cables supplied by Compex.

- Do not disconnect the stimulator electrode cables during the session without first switching off the stimulator.

- Never connect stimulation cables to an external electric power supply because of risk of electric shock.

- Never use the Compex if it is damaged (case, cables, etc.) or if the battery compartment is open. There is a risk of electric discharge.

- Keep the Compex and its accessories out of the reach of children.

- Do not allow any foreign items (soil, water, metal, etc.) to enter into the Compex and the battery compartment.

- Sudden temperature variations can create droplets in the device due to condensation. Use the device only once it reached ambient temperature.

- Do not use the Compex while driving or operating machinery

- Do not use the stimulator in the mountains at an altitude higher than 3.000 metres.

Where you should never apply electrodes...

- To the head.
- Counter-laterally: Do not use the two poles of the same channel (one electrode connected to the positive/red pole and one electrode connected to the negative/black pole) on either side of the body centreline (eg: right arm and left arm).
- To or in the vicinity of skin lesions of any kind (wounds, inflammations, burns, irritations, eczema, etc.).

Precautions when using the electrodes

- Only use the electrodes supplied by Compex. Other electrodes may have electrical properties that are unsuitable for the Compex stimulator.

- Do not remove or move the electrodes during a stimulation session without first switching off the instrument.

- Do not plunge the electrodes into water.

- Do not apply a solvent of any kind to the electrodes.

- Before applying the electrodes, it is recommended to wash and degrease the skin, and then dry it.

- Apply the whole surface of the electrodes to the skin.

■ For very important reasons of hygiene, each user must have his own electrode set. Do not use the same electrodes on different people.

■ Do not use a set of electrodes for more than fifteen sessions, as the quality of the contact between the electrode and the skin, which is important for the comfort and effectiveness of stimulation, progressively degrades.

■ In certain people with very sensitive skin, redness may be observed under the electrodes after a stimulation session. Generally, this redness is completely harmless and disappears after 10 to 20 minutes. However, avoid starting a stimulation session on the same area until the redness is no longer evident.

Precautions when using the batteries

■ Do not open the battery compartment without first switching off the stimulator.

■ Use only alkaline batteries 1.5 V, size AAA, LR03. Never use batteries with other characteristics.

■ All the four batteries must be replaced at the same time.

■ Batteries must be positioned according to the diagram inside the battery compartment.

■ Do not allow any foreign items (object, water, sand, dust, etc.) to enter into the battery compartment.

■ Always close the battery compartment before switching on the device.

■ Remove immediately the batteries if the Compex “beeps” continuously, in case of abnormal heating or smell, or if smoke comes from the device.

■ Never leave the batteries inside your device if you do not intend to use it for a long period of time. Batteries may leak and damage your unit. If the batteries show signs of a leak, avoid contact with skin, eyes and clothing. In case of contact, wash surface with water and contact your doctor.

II. PRESENTATION

1. Reception of equipment and accessories

Your stimulator is supplied with:

- The stimulator
- 4 batteries
- One set of electrode cables (blue, yellow)
- Bags containing the electrodes (small (5 x 5 cm) and large (5 x 10 cm) electrodes)
- A user and applications manual
- A poster with the pictures for the placements of electrodes and the stimulation positions
- Belt clip

2. Presentation of the device

A – “On/Off” switch

B – Sockets for the 2 electrode cables

C – Electrode cables
channel 1 = blue channel 2 = yellow


D – Selection key

E – Validation key

F – “+”/“-” keys of the 2 stimulation channels

G – Battery compartment

H – Belt clip socket

3. Guarantee

The Compex stimulators are contractually covered with a guarantee of 2 (two) years.

The Compex guarantee comes into effect on the date of purchase of the device.

The Compex guarantee applies to the stimulator (goods and work) and does not cover cables and electrodes. It covers all the defects resulting from a problem of quality of the material or from a defective manufacture. The guarantee does not apply if the device was damaged further to a shock, an accident, a false operation, an insufficient protection against the humidity, a dumping or a repair not made by our after-sales services.

 *The guarantee is valid only on presentation of proof of purchase.*

Legal rights are not affected by this guarantee.

4. Maintenance

To clean your unit, use a soft duster and an alcohol-based cleaning product, which does not contain any solvents.

The user must not attempt any repairs to the device or any of its accessories. Never dismantle the Compex containing high-voltage parts because of risk of electric discharge.

Compex Médical SA declines all responsibility for damages and consequences resulting from any attempt to open, modify or repair the device or any of its

components by a person or a service centre not officially approved by Compex Médical SA.

Compex stimulators do not require calibration or verification of performance parameters. The characteristics are systematically verified and validated for each device manufactured. These characteristics are stable and do not vary when used under normal conditions.

If your device seems not to function as expected, regardless of the situation, contact an official Compex service centre for assistance.

Medical and health professionals must refer to local legislation for information related to maintenance. Normally, these laws require verification of certain criteria at regular intervals.

5. Storage and transportation conditions

The storage and transportation conditions must not exceed the following figures:

Storage and transportation temperature from -20°C (-4°F) to 45°C (113°F)

Max. relative humidity: 75%

Atmospheric pressure: from 700 hPa to 1060 hPa

 *Never leave the batteries inside your device if you do not intend to use it for a long period of time.*

6. Use conditions

Temperature of use: from 0°C (33°F) to 40°C (104°F)

Max. relative humidity: from 30% to 75%

Atmospheric pressure: from 700 hPa to 1060 hPa

 *Do not use in an explosion risk area.*

7. Elimination

The Directive 2002/96/CEE (WEEE) has a first priority, the prevention of waste electrical and electronic equipment, and in addition, the re-use, recycling and other forms of recovery of such wastes so as to reduce the disposal of waste.

The pictogram wheeled bin barred means that the equipment can not be thrown with the household refuse, but that it makes the object of a selective collection. The equipment has to be given to a suitable collection point for the treatment.

By this way, you contribute to the safeguarding of the natural resources and the human health protection.

Batteries must be disposed of in accordance with your country's national laws governing the disposal of such items.

 *Batteries must never be disposed of in a fire. There is a risk of explosion.*

8. Standards

The Compex is based directly on medical technology.

To guarantee your safety, the design, manufacturing and distribution of Compex are in conformity with the requirements of the European Directive 93/42/CEE.

The device is in conformity with the standard for general safety rules for electromedical devices IEC 60601-1, the standard for electromagnetic compatibility IEC 60601-1-2 and the standard for special safety rules for nerve and muscle stimulators IEC 60601-2-10.

Current international standards (IEC 60601-2-10 AM1 2001) require that a warning be given concerning the application of electrodes to the thorax (increased risk of cardiac fibrillation).

Directive 2002/96/CEE Waste electrical and electronic equipment (WEEE).

9. Patent

Electronic protection system: patent pending

10. Normalized symbols

 **IMPORTANT:** Under some conditions, the effective figure for the stimulation pulses can exceed 10 mA and 10 V. The information given in this handbook must be strictly observed.

 The Compex is a class II device unit with its own internal electric power, with type BF applied sections.

 Waste electrical and electronic equipment (WEEE).

 The "On/Off" switch is a multi-function key:

Functions Symbol N° (according to CEI 878)

On/Off (2 positions, stable) 01-03

Waiting or preparation for part of the unit 01-06

Stop (switching off) 01-10

11. Technical information**General points****P**ower supply: 4 alkaline batteries 1.5 V, size AAA, LR03**C**ase: Plastic**Neurostimulation****A**ll electrical specifications are given for a charge between 500 and 1000 ohms per channel.

Channels: 2 independent and individually adjustable channels electrically insulated from each other and from ground

Impulse shape: Constant rectangular current compensated to eliminate any direct current component to avoid any residual skin polarisation

Maximum current of an impulse: 80 milliamperes



Intensity increase interval:
– Manual stimulation intensity adjustment: 0-389 (energy units)
– Minimum interval: 0.5 mA

Impulse duration: From 60 to 350 microseconds

Maximum quantity of electricity per impulse: 56 microcoulombs (2 x 28 µC compensated)

Typical rise time for an impulse: 3 microseconds (between 20% and 80% of maximum current)

Impulse frequency: From 1 to 150 Hertz

*ted to the nerves by self-adhesive electrodes.**The choice of the size, the connection and the correct positioning of the electrodes are essential factors for ensuring effective and comfortable stimulation. So they require particular care.**To do all this – and for information on the recommended stimulation positions – refer to the pictures and pictograms in the poster, and to the “Table with placements of electrodes and stimulation positions” in chapter VI. The specific applications contained in chapter VII also provide useful relevant information.***Connecting the electrodes and cables****T**he electrode cables are connected to the stimulator through sockets on the base of the device.**T**wo cables can be connected simultaneously to the two channels of the device.**F**or easier use and better identification of the two channels, we advise you to follow the colour of the electrode cables and the stimulator sockets.**3. Preliminary settings: language, light intensity and sound volume****B**efore using the unit for the first time, you should select the working language of the device which is displayed on the options screen. Afterwards, for the greatest comfort, Compex offers you a number of setting options (operating language selection, display contrast setting and sound volume setting).**T**o do this, display the options screen by pressing the “On/Off” switch on the left of the Compex and hold it down for a few seconds.**Fig. 3** 3a 3b 3c**3a** To select the language of your choice, press the selection key.**3b** To set the contrast of the display, press the “+”/“-” key for channel 1.**3c** To set the sound volume, press the “+”/“-” key for channel 2.**3d** To confirm the parameters selected, press the “On/Off” switch.**4. Selecting your motivation****Fig. 4** 4a 4b**T**o switch on your stimulator, briefly press the “On/Off” switch on the left of the Compex.**B**efore choosing a programme, you should select the desired motivation.**4a** To select the desired motivation, press the selection key or channel 2, until the desired motivation is displayed in a frame.**4b** To confirm your choice, press the validation key.**III. DIRECTIONS FOR USE AND INSTRUCTIONS**

 Before using the unit for the first time, you are strongly advised to take careful note of the counter-indications and safety measures detailed at the beginning of this manual (chapter I: “Warnings”), as this powerful equipment is neither a toy nor a gadget!

1. Inserting and replacing the batteries**B**efore using the unit for the first time and to make it work, you must insert the four batteries delivered with your device into the battery compartment on the back of your device.

 Use only alkaline batteries 1.5 V, size AAA, LR03. Never use batteries with other characteristics.

2. Connections

 The electrical impulses generated by the Compex are transmitted

5. Selecting a program

 To choose a program, it is particularly useful to consult chapter VII of this manual ("Programs and specific applications").

Regardless of the motivation you selected in the previous step, a list containing a select number of programs appears on the screen.

5a To select the desired program, press the "+"/"-" key for channel 2.

5b Press the "On/Off" key to return to the previous screen.

5c Press the validation key which, depending on the program selected, displays the symbol **OK** or **START**


a) **START** = the stimulation session starts immediately;

b) **OK** = a parameter setting screen displays.

6. Personalizing a program

 The program personalisation screen is not available for all programs!

6a Some programs require manual selection of the muscular group that you want to stimulate. To select your desired group, press the selection key.

6b The warming-up sequence suggested for some programs is active by default (an animated rising convection symbol above the heater). If you prefer to skip this sequence, press the "+"/"-" key for channel 1.


6c Some programs allow you to adjust the level of work. To do this, press the "+"/"-" key for channel 2, until the desired level of work is displayed.

6d Total programme time in minutes.

6e By pressing the "On/Off" key, you can return to the previous screen.

6f To confirm your selection and start the program, press the validation key.


 Current international standards require that a warning be given concerning the application of electrodes to the thorax (increased risk of cardiac fibrillation).

7. During the stimulation session

Starting a program

When the program starts, the Compex prompts you to increase the stimulation energy, a key factor for the effectiveness of the stimulation.

7a The Compex "beeps" and the symbols of the two channels flash, changing from "+" to "000". You must increase the stimulation by pressing the "+" of the keys of the relevant channels, until the desired setting is reached.

7b Remaining time in minutes and seconds.


7c The animated rising convection symbol represented the warming-up sequence.

7d Session progression scale.

7e By pressing the "On/Off" key, you can return to the previous screen.

Program progression

Stimulation does really start when the stimulation energy has been increased. The examples reproduced below are sufficient to understand the general rules.

7f Remaining time (in minutes and seconds) before the end of the program.

7g A small animated symbol indicating the motivation to which the program belongs.

7h Session progression scale: The progression scale is composed of three distinct parts: the program consists of three stimulation sequences, i.e., warming up (light hatching/left), the actual work (black area/centre) and relaxation at the end of the program (heavier hatching/right). The progression scale can be also entirely hatched: the program then consists of a single stimulation sequence.

The small cursor located just below the scale indicates precisely the state of progression of the program.

7i Contraction and active rest time bar: This bar is displayed only during the work sequence. It indicates the progression of the contraction (black part) and active rest (hatched part) times.


7j Bar charts:

The muscular contraction phase is represented by the black bar charts. The active rest phase is represented by the hatched bar charts.


■  The stimulation energy during the active rest phase is automatically set at 50% of the stimulation energy during the muscular contraction.

7k You can interrupt the program momentarily (**STOP** = "PAUSE") by pressing the "On/Off" key. To restart the session, simply press the "+"/"-" key for channel 1 or 2.

■  After an interruption, stimulation restarts at 80% of the stimulation energy set before the **STOP**.

7l During the session, the Compex "beeps" and some  symbols start to flash above the active channels (over the energy indications): the Compex is prompting you to increase the stimulation energy level. If you cannot support an increase of the energy, just ignore this message.

End of program


7m At the end of the session, you hear a musical signal and a small flag is displayed.

8. Electricity consumption


 Use only alkaline batteries 1.5 V, size AAA, LR03. Never leave the batteries inside your device if you do not intend to use it for a long period of time. Batteries may leak and damage your unit. If the batteries show signs of a leak, avoid contact with skin, eyes and clothing. In case of contact, wash surface with water and contact your doctor.

The Compex runs on four alkaline batteries 1.5 V. Their working life depends on the programs and stimulation energy used.

The symbol of a small battery indicates the charge level of the batteries.


8a The small battery is empty and it is blinking; batteries are dead. You must replace them

8b The **START** symbol normally displayed above the validation key has disappeared (): the batteries are completely dead. It is no longer possible to use the device. Change the batteries immediately

9. Problems and solutions

Electrode fault


9a The Compex "beeps" and alternatively displays the symbol of a couple of electrodes and an arrow pointing to the channel where a problem has been detected. This message may mean that:

- 1) there is no electrode connected to this channel;
- 2) the electrodes are old, worn out and/or the contact is poor: try using new electrodes;
- 3) the electrode cable is defective: try it using another channel. If cable shows a default again, please replace.

Device fault

In the case a fault happens during use, the following screen is displayed:

If this message appears, switch off the device and remove the batteries. Please contact the hotline, noting this error code.

If, for any reason, you experience any other trouble with your device, please contact the hotline mentioned and approved by Compex Médical SA.

IV. HOW DOES IT WORK ?

The principle of electrostimulation is to stimulate nerve fibres by means of electrical impulses transmitted by electrodes.

The electrical impulses generated by Compex stimulators are high-quality impulses – ensuring safety, comfort and efficiency – that stimulate different types of nerve fibres:

1. the motor nerves to stimulate a muscular response. The quantity of work and the benefits obtained depend on the stimulation parameters. This is known as electromuscular stimulation (EMS).


2. certain types of sensitive nerve fibres to obtain analgesic or pain-relieving effects.

1. Motor nerve stimulation (EMS)

In voluntary activity, the order for muscular work comes from the brain, which sends a command to the nerve fibres in the form of an electrical signal. This signal is then transmitted to the muscular fibres, which contract.

The principle of electrostimulation accurately reproduces the process observed during a voluntary contraction. The stimulator sends an electrical current impulse to the nerve fibres, exciting them. This excitation is then transmitted to the muscular fibres causing a basic mechanical response (= muscular twitch). The latter constitutes the basic requirement for muscular contraction.

This muscular response is completely identical to muscular work controlled by the brain. In other words, the muscle cannot distinguish whether the command comes from the brain or from the stimulator.

The parameters of the Compex programs (number of impulses per second, contraction time, rest time, total program time)

subject the muscles to different types of work, according to muscular fibres. In fact, different types of muscular fibres may be distinguished according to their respective contraction speed: slow, intermediate and fast fibres. With a good knowledge of human physiology and a perfect mastery of the stimulation parameters of the various programs, the muscular work can be directed very precisely towards the desired goal (muscular reinforcement, increased blood flow, firming up, etc.).

2. Stimulation of the sensitive nerves

The electrical impulses can also excite the sensitive nerve fibres to obtain an analgesic or pain-relieving effect.

The stimulation of the tactile sensitive nerve fibres blocks the transmission of pain by the nervous system. The stimulation of another type of sensitive fibres creates an increase in the production of endorphins and, therefore, a reduction of pain.

With pain relief programs, electrostimulation can be used to treat localized sharp or chronic pains as well as muscular pains.

Do not use the pain relief programs for a long period without medical advice.

Benefits of electrostimulation

Electrostimulation is a very effective way to make your muscles work:

- with significant improvement of different muscular qualities,
- without cardio-vascular or mental fatigue,
- with limited stress on the joints and tendons.

Electrostimulation thus allows a greater quantity of work by the muscles compared with voluntary activity.

To be effective, this work must involve the greatest possible number of mus-

cular fibres. The number of fibres working depends on the stimulation energy. It will therefore be necessary to use the maximum tolerable energy. The user controls this aspect of stimulation. The higher the stimulation energy, the greater the number of muscular fibres that are working and, therefore, the more significant the progress achieved.

To maximize your results, Compex recommends that you complement your electrostimulation sessions with other efforts, like:

- regular exercise,
- proper and healthy nutrition,
- balanced lifestyle.

V. USAGE PRINCIPLES

The usage principles presented in this section should be considered general rules. For all programs, it is recommended that you read carefully the information and advice on use presented in chapter VII of this manual ("Programs and specific applications").

1. Placement of electrodes

It is recommended that you comply with the suggested placements. To do this, refer to the pictures and pictograms shown on the poster and to the "Table with placements of electrodes and stimulation positions" in chapter VI.

A stimulation cable consists of two poles:

a positive pole (+) = red connection

a negative pole (-) = black connection


A different electrode must be connected to each of the two poles.

Note: It is possible and normal to have an electrode arrangement that leads one electrode connection free from a cable.

For all muscle electrostimulation programs, i.e., programs involving muscle contractions, it is important to place the positive polarity electrode on the muscle motor point.

Always use the size of electrodes shown in the pictures. Unless you have other specific medical instructions, always follow the placement directions in the pictures. If necessary, find the best position by slightly moving the positive polarity electrode, so as to obtain

the best muscle contraction or the position that seems most comfortable.

2. Stimulation positions

To determine the stimulation position to adopt depending on the placement of the electrodes and the chosen program, refer to the pictures and pictograms on the poster and to the "Table with placements of electrodes and stimulation positions" in chapter VI.

The stimulation position depends on the muscular group to be stimulated and the program chosen.

For the programs involving powerful muscular contractions, the muscle should always be stimulated in an isometric fashion. In this way, you provide maximum resistance to the movement and prevent any shortening of the muscle during the contraction, which could create cramp pains and serious stiffness after the session.

For the other types of programs (for example, the Wellness programs and the **Active recovery** program), which do not involve powerful muscular contractions, position yourself as comfortably as possible.

3. Setting the stimulation energy level

In a stimulated muscle, the number of recruited fibres depends on the stimulation energy. For the programs involving powerful muscular contractions, it is therefore absolutely necessary to use maximum stimulation energy, always at the limit of what you can tolerate, so as to involve the greatest possible number of fibres.

The "Table with placements of electrodes and stimulation positions" in chapter VI provides information on the energy to use, depending on the selected program.

4. Progression for the levels

In general, it is not advisable to go through the different levels quickly with the intention of reaching level 3 as fast as possible. In fact, the different levels correspond to progress with electrostimulation.

The simplest and most usual procedure is to start with level 1 and raise the level when changing to a new stimulation cycle.

At the end of a cycle, you may either start a new cycle at the next level up or do some maintenance at the rate of 1 session a week at the last level reached.

VI. TABLE WITH PLACEMENTS OF ELECTRODES AND STIMULATION POSITIONS

The following table will help you to determine which electrode placement and associated stimulation position to adopt, depending on the chosen program.

It should be read together with the poster:

- the numbers on the table refer to the electrode placement pictures in the poster;
- the letters refer to the pictograms of the corresponding body positions.

The table also indicates the stimulation energy to use, according to the program.

Example:

Muscle building 13 I

FITNESS	Energy	01	02	03	04	05	06	07	08	09	10	11	12	13	14	15	16	
Toning		0	A	A	A	-	IL	-	-	-	-	-	-	EGL	I	B	B	B
Firming		0	A	A	A	-	IL	-	-	-	-	-	-	EGL	I	B	B	B
Abs		0	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	I
Buttocks		0	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	L
Muscle building		0	A	A	A	-	IL	-	-	-	-	-	-	EGL	I	B	B	B

The diagram shows a Compex device with a screen displaying a person sitting. A line connects the device to a poster titled "Compex" which shows various body parts with electrode placement numbers (1-16) and pictograms (I, L, B, E). A specific example is highlighted: "Muscle building 13 I". This points to the "Muscle building" row in the table above, which shows energy levels 0, A, A, A, -, IL, -, -, -, -, EGL, I, B, B, B. The "13" points to the 13th energy level (EGL) and the "I" points to the pictogram of a person lying down.

Table with placements of electrodes and stimulation positions

	Energy	01.	02.	03.	04.	05.	06.	07.	08.	09.	10.	11.	12.	13.	14.	15.	16.	17.	18.	19.	20.	21.	22.	23.	24.	25.	26.	27.	28.	29.	30.	31.	32.	33.	34			
FITNESS																																						
Toning	●	A	A	A	IL				E	G	L	I	B	B	B	B	H	C	C	D	D	D	D	D	F	F												
Firming	●	A	A	A	IL				E	G	L	I	B	B	B	B	H	C	C	D	D	D	D	D	F	F												
Abs	●											I																										
Buttocks	●											L																										
Muscle building	●	A	A	A	IL				E	G	L	I	B	B	B	B	H	C	C	D	D	D	D	F	F													
Jogging	●											G																										
Active recovery	*																																					
Power	●	A	A	A	IL				E	G	L	I	B	B	B	B	H	C	C	D	D	D	D	F	F													
Combined workout	●	A	A	A	IL				E	G	L	I	B	B	B	B	H	C	C	D	D	D	D	F	F													
WELLNESS																																						
Relaxing massage	*																																					
Heavy legs	*											JL																										
Pain management TENS	☼																																					
Chronic pain	*																																					
Cramp prevention	*																																					

● = Maximum tolerable energy
 * = Pronounced twitches
 ☼ = Tickling sensation

◇ = Comfortable position
 . = Non applicable

Letters in the table refer to pictograms for the body positions in the poster

Cover the painful area (eg.: 26, 28, 31, 32, 33, 34); Comfortable position

Positive electrode on the painful point (eg.: 14, 15, 16); Comfortable position

www.compex.info

You will find additional applications on our Internet site (www.compex.info). These protocols will help you to determine the best procedure to follow depending on your needs (choice of the program, muscular group, duration, placement of the electrodes, body position).

VII. PROGRAMS AND SPECIFIC APPLICATIONS

The applications that follow are given by way of example. They provide a better understanding of how electrostimulation sessions can be combined with voluntary activity.

You can of course choose a muscle group different from the one indicated, depending on the activities you practice or the body area you wish to stimulate. The specific applications provide you information on the placement of electrodes and the stimulation position to adopt. This information is represented by numbers (placement of electrodes) and letters (stimulation position). They refer to the "Table with placements of elec-

trodes and stimulation positions" in chapter VI and to the poster with the pictures for the placements of electrodes and the pictograms for the body positions. The table also indicates the stimulation energy to use, according to the program. You must select in the device the muscle group indicated on the figurine. When no figurine appears (→), it is not necessary to choose the muscle group: the selection is automatic.

Your motivation: Fitness

Programs

Toning

The use of this program will help you to tone up your muscles and prepare them before the more intensive work of firming. You will quickly recover a tonic, firm and well shaped body.

Firming

This program is very useful if you want to increase firmness and restore the elastic quality of the supporting tissues. It also compensates for insufficient muscular activity by specifically putting to work some flabby areas.

Abs

Specifically adapted to the abdominal muscles, this program will enable you to reduce your waist line, while improving the muscular support of the abdominal wall.

Buttocks

This program subjects the buttocks muscles to the most appropriate work for improving your figure. Its hyper-tensor effect is the ideal solution for restoring or improving the tone of the buttocks region.

Muscle building

This program will help you not only to obtain firm and well shaped muscles, but also to develop in a balanced manner your musculature in general (strength, volume, tone)

Jogging

This program is useful as preparation and support for your voluntary jogging sessions. It increases your comfort level during this type of effort and significantly improves your endurance.

Active recovery

Through a strong increase of the blood flow and an accelerated elimination of muscular contraction waste, this program stimulates and activates your muscular recovery. It has a relaxing and soothing effect after particularly sustained muscular work or strenuous physical activity.

Power

This program allows you to improve your muscular strength with a slight increase in volume.

Combined workout

By subjecting your muscles to different working routines, this program has been designed to maintain your overall physical condition or to reactivate your muscles after a period of rest and before a more intense and targeted physical or sports activity.

Specific applications

■ *At the end of a cycle, you may either start a new cycle at the next level up or do some maintenance at the rate of 1 session a week at the last level reached*

You want to firm your thighs

The **Toning** and **Firming** programs impose perfectly adapted and progressive work on the muscles. This very intensive muscular activity (hundreds of seconds of sustained contractions) first tones then firms your muscles.

In this example, the person wants to work on her thighs. It is of course possible to stimulate other muscles as well. Also, the same stimulation routine can be applied to several muscular groups during the same cycle.

Cycle duration: 8 weeks, 4x/week

weeks 1-4: 4 x **Toning** 11G
👤/week

weeks 5-8: 4 x **Firming** 11G
👤/week

You want to reduce your waist and tone up your abdominal belt

The **Abs** program is especially appropriate for improving the qualities of support and firmness of the abdominal belt muscles. The various work routines, associated with the tonic massage performed during the rest phases, guarantee beneficial results on the waist line.

Cycle duration: 3 weeks, 4 x/week

Program: **Abs** 13 I –

Your objective is to improve the firmness and shape of your buttocks

The considerable amount of progressive and specific work of the buttocks muscles produced by the **Buttocks** program effectively increases the firmness of these muscles. The shape of the buttocks improves progressively, with beneficial effects on the whole figure.

Cycle duration: 3 weeks, 4 x/week

Program: **Buttocks** 12A, B or L –

You want to develop harmonious arm muscles with a moderate increase in muscle volume

Most voluntary physical activities, such as jogging and cycling, require the arm muscles to do only limited work. That is why it is particularly important to compensate this under-utilisation with electrostimulation sessions.

The **Muscle building** program allows to impose a lot of specific work to the muscles of the upper part of your body and thus ensures the harmonious development of your arms, with firm and well shaped muscles.

Cycle duration: 5 weeks, 4 x/week

Eq for 1 week

Mon: 45'-1h of voluntary physical activity (jogging, swimming, cycling, fitness activity, etc.), then 1 x **Muscle building** 21D, 22D or 23D 👤

Tue: Rest

Wed: 1 x **Muscle building** 21D, 22D or 23D 

Thu: 45'-1h of voluntary physical activity (jogging, swimming, cycling, fitness activity, etc.)

Fri: 1 x **Muscle building** 21D, 22D or 23D 

Sat: Rest

Sun: 1 x **Muscle building** 21D, 22D or 23D 

You wish to progress in endurance during your jogging sessions

By using the **Jogging** program, you will feel more comfortable when running and will be capable of doing long runs more easily.

Moreover, this program is very useful for you if you jog regularly for recreation and cannot do a training session: you can usefully replace the voluntary training with a stimulation session.

Cycle duration: 8 weeks, 2-3 x/week

Eg for 1 week

Mon: 30'-1h of jogging or physical activity in a fitness centre

Tue: 1 x **Jogging** 11G –

Wed: 45'-1h of physical activity in a fitness centre

Thu: 1 x **Jogging** 11G –

Fri: Rest

Sat: 1h of voluntary physical activity (jogging, swimming, cycling, etc.), then 1 x **Jogging** 11G – (optional)

Sun: Rest

You want to recover better and faster after an effort

A session with the **Active recovery** program, which accelerates the drainage of toxins, will help you to relax your muscles and to restore your muscular balance faster, whatever type of activity you practice.

The program is usually used between 15 minutes and 3 hours after a physical effort. If the result is insufficient, the session can be repeated immediately a second time. It is also possible to repeat a session the following day, if aches and/or spasms persist.

An example of how to use this program is suggested under the application “You want to use a varied muscle training approach”.

You want to develop and strengthen the large pectoral muscles

The development and strengthening of the pectoral muscles can be useful for sports activities requiring efficient upper limb and shoulder strength. The **Power** program is the perfect solution for working on these muscles to improve their contraction strength. Using this program, you can also achieve a harmonious development of muscular volume.

Cycle duration: 4 weeks, 3 x/week


Program: **Power** 20 C 

You want to use a varied muscle training approach

Through these eight sequences that automatically follow each other, the **Combined workout** program imposes various kinds of stimulation during the same session. It does not emphasise any particular kind of muscular performance, but rather allow you to reach the right level in all kinds of muscular work (endurance, resistance, strength).

Cycle duration: 6 weeks, 2-3 x/week


Eg for 1 week

Mon: 1 x **Combined workout** 11G , then 1 x **Active recovery** 11G 

Tue: Normal training

Wed: 1 x **Combined workout** 11G , then 1 x **Active recovery** 11G 

Thu: Normal training

Fri: 1 x **Combined workout** 11G , then 1 x **Active recovery** 11G  (optional)

Sat: Normal training or other outdoor physical activity

Sun: Rest

Your motivation: Wellness

Programs

Relaxing massage

This program is indicated for reducing unpleasant muscular tensions and creating a relaxing sensation of well-being, in particular after a tiring day. It induces a state of very pleasant optimal relaxation.

Heavy legs

This program is intended to eliminate the sensation of weight and discomfort in the legs that is felt in particular circumstances (following prolonged standing positions, excessively high ambient temperature, hormonal imbalance related to the menstrual cycle, etc.). It quickly produces a sensation of lightness and suppleness.

Pain management TENS

All painful phenomena and, in particular, pains due to movements-related pathological conditions, can be effectively eased using this program that blocks the transmission of pain by the nervous system.

Chronic pain

Thanks to its analgesic action, this program provides rapid relief due to chronic muscular pains (back, lower back pains, etc.). It reactivates the circulation, oxygenates and loosens up tense muscles.

Cramp prevention

By restoring the balance of the blood flow, this specific program prevents the onset of cramps during the night or after physical effort.

Specific applications

You want to prevent uncomfortable muscular tension in the back of the neck

Remaining seated for long periods, associated with repetitive arm movements (as is frequently the case in front of a computer screen), often provokes increased tension in the back of the neck muscles. Any other condition of stress may also create situations leading to a state of excessive muscular tension, which is often responsible for painful or uncomfortable sensations.

The in-depth effects produced by the **Relaxing massage** program ease such painful sensations, with a particularly relaxing effect.

Cycle duration: To be used on the back of the neck muscles, whenever an occasional sensation of painful muscular tension is felt; to be repeated, if necessary, in cases of particularly strong muscular tension

Program: **Relaxing massage** 16◇

You want to prevent the feeling of heaviness and discomfort in the legs

The feeling of heaviness in the legs is due to a temporary deficiency in the venous return. The insufficient supply of oxygen to the tissue, and in particular the muscles, that results from this situation, causes the unpleasant feeling in the legs.

Use of the **Heavy legs** program accelerates the return of blood to the veins and has an important relaxing effect on the sore muscles.

Cycle duration: To be used whenever an occasional feeling of heaviness in the legs occurs

Program: **Heavy legs** 5 J –

You have hand and wrist pains (eg: pains following prolonged use of a computer mouse)

Remaining for a long time in certain inappropriate working positions and the micro-traumas caused by the repetitive use of a computer mouse can lead to a wide range of painful phenomena.

Hand and wrist pains are the most frequent and can be treated very effectively using the **Pain management TENS** program.

Cycle duration: 1 week, min. 1 x/day, then adapt in relation to painful episodes

According to requirements, the **Pain management TENS** program can be repeated several times during the same day

Program: **Pain management TENS** 34◇ –


*The use of **Pain management TENS** and **Chronic pain** programs should not be prolonged without medical advice.*

You have muscle tensions in the back of the neck, middle of the back or lower back region

Pains in the back area or the back of the neck are due to a stiffness of the muscles, with accumulation of toxins. The pain can then become continuous or appear after only a few minutes spent working in an unfavourable position (for example, working in front of a computer screen).

These chronic muscle pains can be treated effectively with the **Chronic pain** program. In addition to its pain-reducing effect, it reactivates the circulation and relaxes the stiff muscles.

Cycle duration: 4 week, 2 x/day, with a 10 minute break between the 2 sessions

Program: **Chronic pain** 16  (back of the neck), 15  (middle of the back) or 14  (lower back region)

You want to prevent cramps in the calf muscles

Many people suffer from cramp in the calf muscles, which can appear spontaneously during rest at night or as a result of prolonged muscular effort. The phenomenon of cramp may in part be caused by a disequilibrium of blood flow in the muscles.

To improve the blood flow and prevent the onset of cramp, the Complex offers a specific stimulation program that will allow you to obtain positive results.

Cycle duration: 5 weeks, 1 x/day, at the end of the day or at evening

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
The diagram illustrates the workflow: a user selects a program (e.g., 'Spiertraining 13 I') on the Compex device. This selection points to a table of energy levels for various programs. The table indicates the energy level for 'Spiertraining' is '13'. This level points to a poster showing electrode placement positions (numbered 1-16) and body postures (lettered I, II, III, IV). The device is shown connected to a person in posture 'I', with arrows indicating the placement of the two electrodes on the back.

FITNESS	Energie	01	02	03	04	05	06	07	08	09	10	11	12	13	14	15	16
Tonus		0	A	A	A	IL	---	E	GL	I	B	B	B	B	B	B	B
Versterking		0	A	A	A	IL	---	E	GL	I	B	B	B	B	B	B	B
Buik		0	---	---	---	---	---	---	---	L	---	---	---	---	---	---	---
Billen		0	A	A	A	IL	---	E	GL	I	B	B	B	B	B	B	B

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