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- One-Step Cleaner and Disinfectant
- Reduces the Risk of Infection and Cross-Contamination

KEEP OUT OF REACH OF CHILDREN.
CAUTION: LIQUID CAUSES MODERATE EYE IRRITATION.
SEE BACK LABEL FOR ADDITIONAL PRECAUTIONARY STATEMENTS.

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ACTIVE INGREDIENT:
 Sodium Hypochlorite.....0.55%
 OTHER INGREDIENTS:.....99.45%
 TOTAL:.....100.00%

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150 PREMOISTENED WIPES
6" x 5" (15.2 cm x 12.7 cm)
NET WT 1 LB 11 OZ (765 g)

R0924030

Artwork shown at approximately 110%. PDF prints at approximately 100% from acrobat reader or when "Fit to printable area" settings are selected.

PATHOGEN/PATÓGENO	CONTACT TIME/TIEMPO DE CONTACTO
<i>Clostridium difficile</i> / <i>Clostridium difficile</i>	3 minutes/3 minutos
Bacteria/Bacteria	30 seconds/30 segundos
†Viruses/†Virus	1 minute/1 minuto
‡Bloodborne Pathogens/ ‡Patógenos de Transmisión Sanguínea	1 minute/1 minuto
TB/TB	3 minutes/3 minutos
Parvoviruses/Parvovirus	3 minutes/3 minutos
Fungi/Fungo	3 minutes/3 minutos

ORGANISMS:	
Bacteria:	<ul style="list-style-type: none"> • <i>Acinetobacter baumannii</i>* • <i>Bordetella pertussis</i>* • <i>Campylobacter jejuni</i>* • Carbapenem resistant <i>Klebsiella pneumoniae</i> (CRKP)* • <i>Clostridium difficile</i> spores***§ • Community Acquired Methicillin resistant <i>Staphylococcus aureus</i> (NARSA NRS123)* • Community Acquired Methicillin resistant <i>Staphylococcus aureus</i> (NARSA NRS384)* • <i>Enterobacter aerogenes</i>* • <i>Enterobacter cloacae</i>* • Vancomycin resistant <i>Enterococcus faecalis</i> (VRE)* • ESBL (Extended Spectrum Beta Lactamase) producing <i>Klebsiella pneumoniae</i> (ESBL producing <i>K. pneumoniae</i>)* • ESBL (Extended Spectrum Beta Lactamase) producing <i>Escherichia coli</i> (ESBL producing <i>E. coli</i>)* • <i>Escherichia coli</i>* • <i>Escherichia coli</i> New Delhi Metallo-Beta Lactamase-1 (NDM-1 <i>E. coli</i>)* • <i>Escherichia coli</i> O157:H7* • <i>Klebsiella pneumoniae</i>* • <i>Klebsiella pneumoniae</i> New Delhi
<ul style="list-style-type: none"> • Metallo-Beta Lactamase-1 (NDM-1 <i>K. pneumoniae</i>)* • <i>Legionella pneumophila</i>* • Linezolid resistant • <i>Staphylococcus aureus</i> (LRSA)* • <i>Listeria monocytogenes</i>* • Methicillin resistant <i>Staphylococcus aureus</i>* • Multi-drug resistant <i>Enterococcus faecium</i> (MDR <i>E. faecium</i>)* • <i>Proteus mirabilis</i>* • <i>Pseudomonas aeruginosa</i>* • <i>Salmonella enterica</i>* • <i>Serratia marcescens</i>* • <i>Shigella dysenteriae</i>* • <i>Staphylococcus aureus</i>* • <i>Streptococcus pneumoniae</i>* • <i>Streptococcus pyogenes</i>* • Vancomycin intermediate resistant <i>Staphylococcus aureus</i> (VISA)* • Vancomycin resistant <i>Staphylococcus aureus</i> (VRSA)* 	<ul style="list-style-type: none"> • †Human Hepatitis C** • †Influenza A** • †Norovirus** • †Poliovirus type 1** • †Respiratory syncytial virus** • †Rhinovirus type 37** • †Rotavirus**
Bloodborne Pathogens:	<ul style="list-style-type: none"> • ††HIV type 1* • †††Human Hepatitis B** • †††Human Hepatitis C**
Parvoviruses:	<ul style="list-style-type: none"> • Canine parvovirus*** • Feline parvovirus***
Fungi:	<ul style="list-style-type: none"> • <i>Trichophyton mentagrophytes</i>*** • <i>Candida albicans</i>***
	<ul style="list-style-type: none"> * 30-second contact time ** 1-minute contact time *** 3-minute contact time †† Kills Avian Influenza A virus on precleaned hard, nonporous surfaces. § Follow special instructions for cleaning prior to disinfection.
Mycobacterium:	<ul style="list-style-type: none"> • <i>Mycobacterium bovis</i> (BCG) (TB)**
Viruses:	<ul style="list-style-type: none"> • †Adenovirus type 2** • †Avian Influenza A**, †† • †Herpes simplex virus type 2** • †HIV type 1* • †Human Coronavirus** • †Human Hepatitis A** • †Human Hepatitis B**

Use Sites:

- Ambulances
- BP monitors
- Carts
- Critical care units
- Dialysis clinics
- Emergency rooms
- Examination rooms
- Glucometers
- Headboards
- Hospitals
- Intensive care units
- IV stands
- Nurse-call device
- Nurses' stations
- Nursing homes
- ORs
- Patient restrooms
- Patient rooms
- Radiology rooms
- Recovery rooms

Medical Surfaces:

- Cabinets
 - Computer keyboards
 - Doorknobs
 - External ultrasound transducers
 - Gurneys
 - Hard, nonporous medical surfaces
 - High touch surfaces
 - IV poles
 - Stretchers
 - Support bars
 - Telephones
 - Mobile phones
 - Other telecommunications equipment
- Surface Materials:**
- Chrome
 - Enamel
 - Glass
 - Glazed ceramic

- Glazed porcelain
 - Laminated surfaces
 - Plastic
 - Porcelain enamel
 - Stainless steel
 - Synthetic marble
 - Vinyl
- Avoid Contact with:**
- Brass
 - Chipped enamel
 - Clear plastic
 - Clothes
 - Copper
 - Fabric
 - Natural marble
 - Painted surfaces
 - Paper surfaces
 - Rubber
 - Sealed granite
 - Silver
 - Unfinished wood
 - Wood

DIRECTIONS FOR USE: It is a violation of Federal law to use this product in a manner inconsistent with its labeling. **Not for cleaning or sanitizing skin. Do not use as a diaper wipe or for personal cleansing. For surfaces that may come in contact with food, a potable water rinse is required.** Avoid contact with clothes, fabric, wood, natural rubber, and painted and paper surfaces. Do not use on silver or chipped enamel. **APPLY CAREFULLY** — This product contains bleach. Do not use on glassware, utensils or dishes. Always refer to manufacturer's care instructions before using on equipment or devices.

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DISPENSER DIRECTIONS: Remove lid and discard seal. Insert corner of center sheet through flap opening on lid. DO NOT push finger through opening. Replace lid and pull wipe at an angle. The next sheet pops up automatically. When finished, close lid flap to retain moisture.

DISINFECTION / PARA DESINFECTAR

TO CLEAN AND DISINFECT AND DEODORIZE HARD, NONPOROUS SURFACES: Wipe hard, nonporous surface to be disinfected. Use enough wipes for treated surface to remain visibly wet for the contact time listed on label. Let air dry. Gross filth must be removed prior to disinfecting. If streaking is observed, wipe with a clean, damp cloth or paper towel after appropriate contact time has expired. **Killing *Clostridium difficile*:** Preclean hard, nonporous surfaces by removing gross filth. Wipe surface to be disinfected. Use enough wipes for treated surface to remain visibly wet for 3 minutes. Let air dry.

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Special Instructions for Cleaning Prior to Disinfection against *Clostridium difficile* spores: **Personal Protection:** Wear appropriate barrier protection such as gloves, gowns, masks, and eye covering. **Cleaning Procedure:** Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with this product. Cleaning is to include vigorous wiping and/or scrubbing, until all visible soil is removed. Special attention is needed for high-touch surfaces. Surfaces in patient rooms are to be cleaned in an appropriate manner, such as from right to left or left to right on horizontal surfaces, and top to bottom on vertical surfaces, to minimize spreading of the spores. Restrooms are to be cleaned last. Do not reuse soiled cloths. **Infectious Materials Disposal:** Materials used in the cleaning process that may contain feces/wastes are to be disposed of immediately in accordance with local regulations for infectious materials disposal. **Special Instructions for Use Against HIV-1, HBV and HCV:** This product kills HIV-1, HBV and HCV on precleaned environmental surfaces/objects previously soiled with blood/body fluids in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus (HIV-1) (associated with AIDS), Human Hepatitis B Virus (HBV) and Human Hepatitis C Virus (HCV). **Special instructions for using this product to clean and decontaminate against HIV-1, HBV and HCV on surfaces/objects soiled with blood/body fluids:** **Personal Protection:** When handling items soiled with blood or body fluids, use disposable impervious gloves, gowns, masks and eye coverings. **Cleaning Procedure:** Blood and other body fluids must be thoroughly cleaned from surfaces and other objects before applying this product. **Contact Time:** Allow surfaces to remain wet for 1 minute, let air dry. For all other organisms, see directions for contact times. **Disposal of Infectious Materials:** Use disposable impervious gloves, gowns, masks and eye coverings. Blood and other body fluids must be autoclaved and disposed of according to local regulations for infectious waste disposal. This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream, or normally sterile areas of the body, or (2) contacts intact mucous membranes, but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

PRECAUTIONARY STATEMENTS: HAZARDS TO HUMANS AND DOMESTIC ANIMALS.

CAUTION: Liquid causes moderate eye irritation. Do not get in eyes or on clothing. Avoid contact with clothing. Wear protective eyewear. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. For sensitive skin or prolonged use, wear gloves. **FIRST AID: IF IN EYES:** Hold eye open and rinse slowly and gently with water for 15–20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing. Call a poison control center or doctor for treatment advice. Have the product container or label with you when calling a poison control center or doctor or going for treatment. **QUESTIONS?** Call (800) 234-7700. **PHYSICAL AND CHEMICAL HAZARDS:** This product contains bleach. Do not use this product with other chemicals, such as ammonia, toilet bowl cleaners, rust removers or acid, as this releases hazardous gases.

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STORAGE AND DISPOSAL: Store this product in a cool dry area, to avoid deterioration. Dispose of wipe in trash after use. Do not flush. Nonrefillable container. Do not reuse or refill empty container. Recycle empty container or discard in trash. CONTAINS NO PHOSPHORUS.

QUESTIONS OR COMMENTS? Call (800) 234-7700. For MSDS information, please go to www.cloroxhealthcare.com.

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Mfd. for Clorox Professional Products Company, 1221 Broadway, Oakland, CA 94612.

© 2012 Made in the U.S.A. Clorox Healthcare is a trademark of Clorox Professional Products Company.

U.S. Patents 7,008,600 and 7,070,737. EPA Reg. No. 67619-12. EPA Est. No. 8251-OH-1 (Q6).

Beginning of batch code indicates Est. No.

R0924030

NI-18491 0 00000 00000 0