SCALPELS, DISKS, EXPANDERS, BONE TAPS, DRILLS, SINUS TREPHINES, MOTORIZED WRENCHES, SCREWDRIVERS AND DRILL EXTENSION.

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Rx only

SURGICAL MATERIAL
SCALPELS, DISKS, EXPANDERS, BONE TAPS, DRILLS, SINUS TREPHINES, MOTORIZED WRENCHES, SCREWDRIVERS AND DRILL EXTENSION.

1. GENERAL COMMENTS
KLOCKNER® IMPLANT SYSTEM supplies its instrumentation in a non-sterile manner and, therefore, the material must be washed and sterilized before use (See [1.2.]). This procedure will be cyclically renewed before all surgery.
During surgery, it is wise to submerge contaminated instruments that have been used for their specific purpose in a detergent solution in order to avoid residue incrustation.

1.1. INSTRUCTIONS FOR USE, MAINTENANCE AND SAFETY OF KLOCKNER® IMPLANT SYSTEM PRODUCTS DEVELOPED BY SOADCO, S.L.
1.1.1. USE
- It is very important to check the assembly with the receptor pieces.
- In the case of motor-connected instruments, check that they are well secured by making the instrument turn before use in surgery.
- All instruments can be used repeatedly, although the degree of use must be decided by the professional, depending on their wear and their state of preservation.
- In the case of drills and trephines, read the instructions regarding this.
- Surgical instruments can be used repeatedly, but the professional must decide of the degree of use. Damaged instruments should be replaced by new instruments to avoid damaging the receiving parts and ensure proper functionality.

1.1.2. MAINTENANCE: CLEANING AND STERILIZING SURGICAL INSTRUMENTS
a. To disinfect and clean the material, closely follow the manufacturer's instructions on their use, contact time and properties of the disinfection and cleaning substances for certain types of instruments. Surgical instruments are cleaned using an aldehyde-free detergent solution with active quaternary ammonium and guanidine compound ingredients [2 % dilution, contact time 1 hour] and with a soft nylon brush, in a dish washer or in an ultrasonic solution.
b. Check that there are no bone fragments, sections or hematoma remains. If using drills with internal irrigation, is very important to note how clean the irrigation port of the drill as well as the supplementary material.
c. Pour ethanol onto the instruments to remove any remains of soap and minerals. Avoid rusting material.
d. Dry the instruments properly, using compressed air. Point the jet of air at the internal cooling conduits and concavities and make sure that drying is thorough.
e. Place the instruments on the surgical tray and wrap the material in paper or bags suitable for use in autoclave [see the manufacturer's comments]. All instrumentation can be sterilized through the use of an autoclave in accordance with the EN ISO 17665-1 standard on sterilisation: i.e. 4 mins at 134 °C]. Do not remove until the process is fully completed. On the day the surgery is to be performed, place the box in the autoclave and sterilise the entire set. Use the autoclave drying option if there is one.
f. Check the instruments are in good working before each surgical operation.

1.1.3. SAFETY
a. We do not recommend using dry heat sterilisers. You should also avoid temperatures greater than 180 °C and the use of chemical sterilisers.
b. The KLOCKNER® IMPLANT SYSTEM instrumentation must only be used by qualified professionals in the field of surgery, including the diagnostic, planning and surgical technique.
c. Improper use of the instruments can cause premature wear of the same and expose the patient to unnecessary physical danger.
d. Avoid damage caused by the overheating of instruments during use.
e. It is advisable to wear safety goggles during the reaming process.
f. We recommend using gloves during the cleaning procedures of contaminated instruments.
g. The instruments must not come into contact with H₂O₂ [hydrogen peroxide] for long periods of time].

2. SCALPELS
Circular scalpels designed for the technique of fitting transmuscular implants, for dentistry biopsies and follow-up surgery. They minimize surgical trauma by avoiding suture. They are used mechanically with the counter-angle at high revs [1.000 r.p.m.].
The depth marks mean that the thickness of existing soft tissue can be assessed.
[See fig. 1]
The scalpels can be used repeatedly, although the amount of use must be decided by a professional, depending on the wear and conservation thereof [are recommended no more than 12 uses].

3. DISCS
The KLOCKNER® IMPLANT SYSTEM represents an essential tool in the crest expansion techniques and the lateral sinus elevations. Their cutting capacity with a minimum thickness allows the volume to be retained.
The initial disc with the smallest diameter and a cutting depth of 1.3 mm commences the osteotomy. Then, if required, the depth of the cut is increased, with the second disc measuring 2.8 mm and the third 4.8 mm.
We recommend using KLOCKNER® IMPLANT SYSTEM discs making a crestal incision and, if required, two transverse incisions to relieve tension in the mesial and distal vestibular area.
We recommend working at high revolutions [5.000 / 10.000 r.p.m.] with plenty of rinsing.
All instruments can be used repeatedly, although the times of use must be decided by a professional, depending on the wear and conservation thereof [are recommended no more than 12 uses].

4. EXPANDERS
The KLOCKNER® IMPLANT SYSTEM expanders share the peri-implant socket, minimizing bone loss by replacing the use of drills. They enable progressive expansion of the osseous crest, with sequenced work possible. They act as a bone compactor, increasing primary stability of the implant in D3-D4 type bones.

4.1. INSTRUCTIONS FOR THE USE, MAINTENANCE AND SAFETY OF EXPANDERS DEVELOPED BY SOADCO, S.L.
In the case of the peri-implant bone being cancellous, type D3 - D4, use of the last drill should be replaced with the expander corresponding to the diameter of the implant to be fitted. [See the surgical procedure].
Work with expanders is recommended where existing bone is narrow, starting drilling with initiation reamers and then widening in line with the expanders sequence or the professional's criteria, working with the drilling sequence and replacing the last drill with the expanders corresponding to the diameter of the implant to be fitted.
In all cases, the expanders must be used after drilling to the required distance [minimum drill diameter 1.8 mm].
The threaded expanders can be used with the hand wrench, the socket wrench and/or the torque wrench. [See the surgical system brochure].

4.2. IDENTIFICATION MARKINGS
The color code links the diameter of each expander to the corresponding implant diameter. [See fig. 5]
The expanders have 2 mm depth marks made by laser. The first mark made by laser is located in the area between 6 mm and 8 mm and the second between 12 mm and 14 mm. To complete it there is a small laser mark at 10 mm. [See fig. 3]
These markings perceptibly help swift assessment of the expander insertion depth.

5. BONE TAPS
The bone taps have been designed for use in oral implantology based on the KLOCKNER ® IMPLANT SYSTEM for the best results and for easier insertion of implants into D1 - D2 bones.

5.1. INSTRUCTIONS FOR THE USE, MAINTENANCE AND SAFETY OF BONE TAPS DEVELOPED BY SOADCO, S.L.
The threaded bone taps can be used with the hand wrench, socket wrench and motorized wrench, respecting the indicated surgical sequence 15 - 30 rpm, at least for the first 30 rpm. [See the surgical procedure].
In the case of high bone resistance, use the hand wrench with tightening and loosening movements until the socket is fully worked.
Do not force tightening excessively, as instrument resistance is not unlimited (maximum safety stress 70 Ncm).
Excessive pressure can accelerate wear of the bone tap and lead to bone overheating.

5.2. IDENTIFICATION MARKINGS
The color code links the diameter of each one with the corresponding implant diameter. [See fig. 5].
The bone taps have 2 mm depth marks made by laser. The first mark made by laser is located in the area between 6 mm and 8 mm and the second between 12 mm and 14 mm. To complete it there is a small laser mark at 10 mm. [See fig. 3]
These markings perceptibly help swift assessment of the bone tap insertion depth.
All the initiators measure absolute distances from the tip of the instrument to the corresponding marking.
[See the surgical system brochure].

6. DRILLS
The S.F.K. reamers have been designed to be used in dental implants using the KLOCKNER® IMPLANT SYSTEM designed for each model of KLOCKNER® IMPLANT SYSTEM for the purpose of obtaining the best results with regard to the cutting, cooling and visibility of the markings.

6.1. INSTRUCTIONS FOR THE USE, MAINTENANCE AND SAFETY OF DRILLS DEVELOPED BY SOADCO, S.L. [S.F.K.]
Bear in mind that while the bone is able to accept surgical trauma it cannot withstand heat trauma. We therefore recommend the use of plenty of cooling during the drilling process, applying gentle alternating pressure to avoid bone necrosis and subsequent osseous attachment failure.

WARNING: IN THE CASE OF NOT USING COOLING WE RECOMMEND NOT EXCEEDING 50 R.P.M.

Also take into account that the leverage forces could cause the reamer to break.
All the reamers can be used repeatedly, although the degree of use must be decided by the professional, depending on the wear of same and its state of preservation [we recommend no more than 5 times].

6.2. IDENTIFICATION MARKINGS
The S.F.K. system drills have a double mark [laser and mechanical].
This marking perceptibly helps swift assessment of the depth of the drill insertion.
The first mechanical mark is at 6 mm and the second at 8 mm. There is a dark, continuous laser coloring between both distances. This process is repeated at the third mechanical mark at 12 mm and the fourth mechanical mark that indicates 14 mm.
These are completed with two small laser markings and 10 mm and 16 mm respectively. The long reamers have a mechanical marking at 18 mm.
All the reamers measure absolute distances from the tip of the instrument to the corresponding mark. [See fig. 4]
The new reamers have a laser mark on their mandrel the reference and diameter.
Each reamer has a color code, which links it to the corresponding implant diameter. [See fig. 5]
* Except the countersunk and the lanceolate reamers.

6.3. REAMER DESCRIPTION
The reamers KLOCKNER® IMPLANT SYSTEM are internal irrigation, double diameter reamers that adapt to the cone design of the implant core. Both the initiation reamer and the pilot reamer are helicoidal, whereas the finishing reamers are double diameter straight reamers.

6.3.1. SK2-NK2, S3M, S4-S6
The normal drilling sequence is designed for manual implant fitting and for immediate Loads. It is also ideal for D3 - D4 type bones.

6.3.1.1. EASY DRILLS SEQUENCE
These reamers have been specially designed for motor positioning of the NK2, SK2, S4 and S6 implants [except Ø 5.5 mm implants].
The EASY DRILLS reamers replace the normal sequence finishing reamers.

6.3.2. ESSENTIAL REAMERS
The ESSENTIAL drilling sequence is designed for motorized implant fitting. Ideal for implant insertion in D1 - D2 type bones.

6.3.2.1. C.I. REAMERS
These reamers have been designed for manual implant fitting and for immediate Loads. It is also ideal for D3 - D4 type bones.

6.3.3. DRILLS FOR GUIDED SURGERY
These drills are specially designed for use in guided surgeries, as the intermediate ring prevents the drill inserting further into the bed when in contact with the ferrule or cap.

6.3.4. DRILLS WITH STOPPER
These reamers are external irrigation. It must employ an intermittent technique for drilling and milling. You must ensure a good external irrigation throughout the operation.
To use the drill stops alveolar crest should be uniform, otherwise you will not be able to achieve the desired safety depth.
Drills are ready to receive and hold the stop during surgery. Make sure you use the right cap for each drill. There is differentiation by color and laser marking.
All drill caps can be used repeatedly, although the times of use must be decided by a professional, depending on the wear and conservation thereof [are recommended no more than 50 uses].
* For more information see the stops systematic "Drill stops."

7. SINUS TREPHINES
Trephines are recommended for use in sinus lifts using the PADROS technique. They measure 8 and 10 mm in diameter and are internally cooled. They have three polyoxethylene stops that allow an appropriate cut depth of 0.8 mm, 1.6 mm and 2.8 mm to be selected, avoiding the risk of tearing the sinus membrane.
The trephines must be used to check the stops to be checked to ensure the cutting depth.
Trephines must be used at low revolutions [500 - 500 r.p.m.] with substantial internal irrigation.
The trephines can be used repeatedly, although the times of use must be decided by a professional, depending on the wear and conservation [are recommended no more than 12 uses].

7.1. INSTRUCTIONS FOR THE USE, MAINTENANCE AND SAFETY OF TREPHINES DEVELOPED BY SOADCO, S.L.
We recommend approaching the sinus antrum from the front with the help of the trephine guard [Ref. Q.T03], which protects the surrounding mucous membrane and minimises the vibration of the trephine when it comes into contact with the bone.
They must be used at low speed [500 -500 r.p.m.] especially at the start.
It is important to use the stops [20-1] with plenty of internal cooling, which reduces the risk of tearing the sinus membrane. All trephines can be used repeatedly, although the degree of use must be decided by the professional, depending on the wear of the trephine and its state of preservation. [for use no more than 12 times].
Excessive pressure can accelerate deterioration of the stops at the same time as leading to overheating.

8. MOTORIZED WRENCHES
Motorized wrenches can be used with bone taps and for filling implants by motor.
KLOCKNER® IMPLANT SYSTEM motorized wrenches should be used at a very low speed of 15 - 30 r.p.m., ensuring that the torque does not exceed 30 Ncm. [see instructions for the counter-angle reducer and surgical motor].

9. SCREWDRIVERS
These are used to screw in the different components requiring counter-angle lighting.
KLOCKNER ® IMPLANT SYSTEM screwdrivers should be used at a very low speed of 15 - 30 r.p.m., ensuring that the torque does not exceed 30 Ncm. [SEE PROSTHETIC SYSTEMS FOR FURTHER INFORMATION].

9.1. INSTRUCTIONS FOR THE USE, MAINTENANCE AND SAFETY OF MOTORIZED SCREWDRIVERS DEVELOPED BY SOADCO, S.L.
Ensure the correct adjustment for better preservation of the screwdriver and threadable parts.

10. DRILL EXTENSION
The drill extension is indicated for use when there are problems with the length of the neighbouring teeth and the instrument cannot be inserted into the necessary place.

11. WARNINGS
The use of component brands other than KLOCKNER® IMPLANT SYSTEM, exempts SOADCO, S.L. from any kind of compensation and/or guarantee.
Product specifications alone do not guarantee correct product usage, and, thus, prior consultation, either with professionals experienced in the system or through the demonstrations and courses periodically given by KLOCKNER® IMPLANT SYSTEM, is advisable. Professional users are therefore obliged to ensure that there is up to date on the current status of the technology and its applications.
Users shall be exclusively responsible for the handling and use of KLOCKNER® IMPLANT SYSTEM products, as these activities are beyond the control of SOADCO, S.L. Consequently, the company shall not be liable for any damage arising from the improper use of its products.
We provide all users with details as to the instructions for the handling and use of our products, which can be consulted in our leaflets and prosthetic and surgical systems in www.klocknerimplantssystem.com or requested from our corresponding national representative.
The use and handling of KLOCKNER® IMPLANT SYSTEM products are solely and exclusively subject to user responsibility. SOADCO, S.L. is exempt from any civil liability for damages arising in this area that are outside its area of control. In compliance with EC and FDA standards, it will not be responsible for the refunding of any of its products where these have been previously manipulated.
Magnetic resonance environment: The KLOCKNER® IMPLANT SYSTEM products have not been evaluated for safety and compatibility nor have they been tested for heating and migration in the magnetic resonance environment.

SOADCO, S.L. Shall only recognize the use of its products by accredited professionals.
WAT SYSTEM and KLOCKNER® IMPLANT SYSTEM products are on sale in every country.
PRECAUTION: US Law restricts the sale of these products to dentists or specialists.

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