AESTHETIC UNIVERSAL Bonder

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Instructions for use

Mode d'emploi

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 CARMEN®  Triceram®

CE 0483

ESPRI  ENT

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Contents

	Page
Availability	8
Properties	9
Application	10
General instructions	12
Firing chart	13

Availability

Order No. **292-200-15** Set (syringe 4g, brush, liquid, small bowl)

Order No. **292-200-10** Syringe 4 g

Order No. **260-905-00** Brush No. 4

Order No. **299-200-10** Liquid

Properties

AESTHETIC UNIVERSAL bonder paste consists of pure gold and organic components with ceramic particle deposits. The ceramic particles produce a rougher surface texture after firing. Following several firings, these particles fuse with the bonding ceramic to form a stronger mechanical bond as well as a chemical bond.

Special properties of the AESTHETIC UNIVERSAL bonder:

- stronger bond between the metal and ceramics (e.g. from 42 MPa to over 55 MPa with titanium)
- the warm hue of the AESTHETIC UNIVERSAL bonder enhances the optical shade effect of the bonding ceramic
- optimum shade in the marginal area of the crown.

Our technical advisors will be pleased to help you with any questions regarding the use of our products.

Hotline tel. no. + 49 7231/803-440

Information and instructions for using Esprident ceramic systems can be found on the Internet at **www.esprident.com**

AESTHETIC UNIVERSAL Bonder is a class IIa medical product and certified with CE 0483.

Application

Note: When using AESTHETIC UNIVERSAL bonder on a titanium framework, Triceram® titanium bonder should no longer be applied!

1. Prepare the metal framework according to the bonding ceramic or alloy manufacturer's directions (Triceram® instructions for use, page 9 / CARMEN®, page 30). Sandblast the metal surface and clean thoroughly. Warm the metal framework gently using the heat from an open ceramic furnace.

Storage Information:

During long periods of storage the gold content inside the syringe can separate from the paste. For this reason we recommend turning the syringe from time to time. Should a dispersion of the mixture be visible, please proceed as follows:

- Remove the lid
- Withdraw the plunger as far as possible and stand the syringe upright with the plunger at the base.
- The paste will sink toward the base.
- Push the plunger back to its original position until the paste can be seen in the syringe head again. Repeat this procedure two or three times.
- If the AESTHETIC UNIVERSAL Bonder becomes difficult to spread it is possible to add a **small** amount of AESTHETIC UNIVERSAL Bonder Liquid, order no.

299-200-10. Place the paste as required in a small bowl and mix with a **small** amount of liquid.

Important: Use the thinner **sparingly!** Once in use, the object must be dried for a longer period of time, otherwise danger of bubbling.

2. Apply the AESTHETIC UNIVERSAL Bonder in a thin layer to the metal framework using a clean brush (Esprident brush no. 4, order no. 260-905-00). Use gentle pressure and light rotation of the brush on the surface to cover the metal framework uniformly.
3. Apply a thin covering layer of AESTHETIC UNIVERSAL bonder. Apply AESTHETIC UNIVERSAL bonder uniformly to ensure that no air becomes trapped.
4. After application, dry the covered surface for at least 5 mins in front of an opened ceramic furnace at the standby temperature of 400°C / 752°F. The shiny paste will turn matt. The object must not exceed 150°C / 302°F. If necessary, a second layer of AESTHETIC UNIVERSAL Bonder may be applied on top of the dried bonder. Again dry thoroughly.

Important: It is essential to dry the AESTHETIC UNIVERSAL Bonder thoroughly before firing. Otherwise the remaining bonder liquid within the superficially dried paste will break free later during ceramic firing. These inclusions can cause bubbles to occur. During drying, the objects should reach a temperature of between 100 – 120°C / 212 – 248°F irrespective of the furnace program, but never exceed 150°C / 302°F. If the open furnace does not generate enough heat, raise the standby temperature. The pre-heating furnace can also

be used as an alternative (increase heat slowly to 130°C / 266°F and hold for 5 minutes) or with a hot air blower (set at 120°C / 248°F for 5 minutes).

5. After firing, an even distribution of ceramic particles should be visible on the metal surface (rough surface). If necessary, apply a second layer.
6. Then apply the opaquer according to the manufacturer's instructions. Only one firing is generally required.

General instructions

- Always keep the brush clean.
- Use alcohol to clean paste off the brush.
- Replace the caps on the brush and AESTHETIC UNIVERSAL bonder syringe after use.
- If the consistency of the AESTHETIC UNIVERSAL bonder is too thick, it can be immersed in a container of warm water to restore its spreading consistency.

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Esprident GmbH · aesthetic dental products