## Nutriline TwinFlo

### INSTRUCTIONS FOR USE Code 1252.

The catheter is made of radio-opaque PUR. It has its own fixation wing with integral distal and proximal extensions, complete with Luer-Locking hubs. The catheter is introduced through a split cannula introducer.

The catheter is used for mid to long term infusion of TPN or intravenous medications.

This catheter is suitable for insertion into peripheral veins. Recommended insertion sites include the basilica vein, the median cubital at the elbow, the external jugular vein, the long saphenous veins of the ankle, and, in neonates, the temporal or femoral veins.

Strict aseptic technique must be observed. In order to avoid air embolism the catheter lumens must be flushed before it is introduced.

Do not stretch the catheter or subject it to pressure above 21,75 psi (1,5 bar, 1.140mmHg). Do not use syringes smaller than 10 cc, as these can generate very high pressures. It is possible to generate 4 or 5 times the maximum safety pressure, with any size of hand held syringe. Subjecting the catheter to pressure above 21,75 psi can result in catheter rupture and embolism. Small syringes generate higher pressures than larger ones.

Assess the patient and determine the site for introduction of the catheter.

Determine the desired position for the placement of the catheter tip and note the length of catheter to be introduced.

Wash hands thoroughly, don gloves, and rinse gloves free of starch with sterile water.

**Carefully clean the insertion site and position drapes. The puncture can only be made after the disinfectant has been completely dried on the skin. Flush both catheter lumens with saline.**

**Be aware that organic solvents such as alcohol or acetone may interact with catheter material and weaken it.**

Perform venipuncture with the split sheath. Venipuncture is confirmed by free flow of blood in the flash back chamber.

Holding the wings, advance the cannula and remove the needle (placing a gloved finger over the cannula opening to prevent excessive bleeding and possible air aspiration).

Pass the catheter through the split cannula using non-toothed forceps.

If difficulty is experienced advancing the catheter, gentle flushing through one of the lumens may assist in opening venous valves, which may be obstructing advancement.

Whilst maintaining a forward pressure on the catheter, remove the cannula gently from the vein.

Press down on the catheter with the middle finger of one hand in front of the introducer tip and with the thumb of the same hand directly behind the wings. With the other hand, press the cannula wings with the thumb and the forefinger. This will activate the opening of the split cannula.

To complete the removal of the cannula, maintain your hold on the wings, slide the cannula along the catheter whilst simultaneously removing it up with a gentle lateral pull (making sure you do not dislodge the catheter during this process).

Anchor the catheter, using the fixation wings. If a vein in the arm is used, a small loop should be left in the catheter to prevent kinking of the catheter if the patient flexes or bends the arm.

#### **IMPORTANT: ARRANGE FOR A CHECK X-RAY TO CONFIRM CORRECT PLACEMENT PRIOR TO STARTING ANY COURSE OF TREATMENT VIA THE CATHETER.**

If the catheter is not used immediately, or its use is temporarily discontinued, both lumens should be heparin locked to prevent clotting or blockage.

#### **NURSING CARE:**

The catheter and insertion site should be inspected daily for any indications of a developing complication. Signs and symptoms at the insertion site include warmth, redness, swelling, pain or tenderness and discharges. If any of these are present, notify the appropriate clinician and record in the nursing records. The catheter's general condition should also be noted.

Catheter cleansing and dressing policy should be in accordance with the policy of the individual facility.

Prior to catheter removal, the catheter site should be cleaned carefully to remove all sticky residues.

The catheter should be removed by gentle sustained traction close to the exit site.

#### **IMPORTANT CAUTIONS:**

Federal law restricts this device to sale by or on the order of a physician.

Incompatible drug delivery: Some studies suggest that offset exit holes at the distal end of the catheter would be beneficial when delivery incompatible drugs. The small size of the lumens of the Nutriline TwinFlo dual lumen Nutriline TwinFlo prevents the offsetting of the holes. Therefore care should be taken to ensuring that precipitation does not occur when infusing incompatible drugs; this can be done by ensuring accurate catheter tip placement in the superior vena cava.

Do not over stretch the catheter as it may rupture, and rebound into the insertion site, causing a catheter embolism.

The lumens are not suitable for blood aspiration.

Do not use this product for diagnosis, control, or treatment of a defect of the heart or the central venous system.

#### **PRECAUTIONS:**

Always confirm correct placement of the catheter tip by x-ray.

Do not use small syringes as these can generate very high pressures.

The catheter must not be advanced into the heart (right atrium). The location of the catheter within the heart may cause cardiac tamponade, or cardiac arrhythmias. It is advisable to make checks of the catheter tip position at regular intervals throughout the usage of the catheter.

Avoid any contact of the catheter tubing to alcohol containing disinfectants. This may irreversibly damage the catheter.

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