**PPS® Flow+****SAFETY HUBER NEEDLE****Instructions for use****I/ DESCRIPTION**

PPS® Flow+ are curved safety Huber needles mounted with one connection line. PPS® Flow+ curved safety Huber needles are available:  
 • in various lengths and diameters  
 • Some PPS® Flow+ curved safety Huber needles are available with a closed BD Q-SyTe™ Luer access system (pre-mounted to the Luer lock of the tubular):  
 • at the distal tip of the tubing for catalog numbers with no "Y" connector,  
 • or only at the lateral site of injection for catalog numbers with a "Y" connector.  
 They are provided with:  
 • Luer lock connectors on which are screwed  
 - closure caps (tend of the main line) which must not be used as sealing caps  
 - injection sites with a silicone septum on the "Y" connectors  
 - and clamps.  
 Tubing is DEHP (di-2 ethylhexyl-phthalate) free. PPS® Flow+ are packaged in 12 unit-cartons.

**II/ INDICATIONS**  
 PPS® Flow+ safety Huber needles are indicated for the administration or withdrawal of fluids through implantable ports.

PPS® Flow+ are designed to prevent blood-borne pathogen exposures caused by accidental needle punctures. POLYPER® Safe does not protect against other routes of blood-borne pathogen transmission.

Some PPS® Flow+ items have a closed Luer BD Q-SyTe™ access system at the tip of the line. The safety Huber needle remains always closed thanks to this device sealed on its tubing. A universal system, BD Q-SyTe™ constitutes a barrier against infections and gases, thus improving the access to intravenous routes. Its design ensures that the catheter remains aseptic after a syringe or infusion line removal.

Already pre-mounted on the Huber needles, the closed BD Q-SyTe™ system is instrumental in care procedure simplifications, and good practice compliance.

Single hand removal of the safety Huber needle PPS® Flow+ automatically and simultaneously generates a positive pressure (saline injection) in the implantable port, thus limiting the risk of blood clot obstruction of the distal tip of the catheter.

**III/ CONTRA-INDICATIONS**

On the skin, present if:  
 - Erythema, exanthem  
 - Ulcerous pain  
 - Oozing, fluid discharge, fistula, abscess, any sign or suspicion of localized infection  
 Any sign or suspicion of the reservoir rolling over  
 - Difficulty or impossibility to locate the septum of the implantable port

**IV/ POSSIBLE COMPLICATIONS**

The most frequently described complications, linked with the use of Huber needles are:  
 - During the puncture or chemotherapy treatment  
 - Infection  
 - Instability, turning over, expulsion with or without extravasation  
 - Low flow, or needle obstruction  
 - Cutaneous fistula due to repeated punctures at the same cutaneous area  
 - Distal bending of the needle (abutting against the bottom of the reservoir, or crushing due to a compression under the dressing).

**V/ WARNING AND PRECAUTIONS**

This product does not provide protection against leakage of contaminated fluids.  
 • Use of the device or needle or a needle/reconnection attempt may cause an accidental puncture with a contaminated needle.  
 • Sterilized with ethylene oxide. **DO NOT RESTERILIZE**.

**Do not use the device if the packaging is damaged, opened, or after the expiry date.**  
 Tighten the injection Y site cap before each use.

Failure to appropriately attach the Luer connector to the device after removing the injection cap can result in air embolism.

Inflammatory puncture slightly move up from the last puncture area in order to avoid sludge or fibrin formation between the septum and the patient skin.

**Pain during the infusion: immediately stop the perfusion and seek medical advice and supervision.**

**Check that the needle length suits the port and patient subcutaneous thickness:** If it is too short, the needle or the implantable port may be damaged during puncture; if it is too short, the septum of the implantable port may not be entirely punctured by the needle and the medication might be administered in the surrounding tissues and/or the needle could be obstructed.

Do not press too hard the tip of the needle on the bottom of the reservoir (a potential damage of the septum or a bend of the needle when the needle is removed).

If PPS® Flow+ are used during the infusion, the port (if plunger to a vertical position and totally locked in the port until a click is heard), the needle tip could come out of the plunger, potentially leading to an accidental puncture with a contaminated needle.

Light-sensitive solutions, the tubing must be protected from light.

PPS® Flow+ should not be used for power injection of contrast media.

**B/ WARNINGS AND PRECAUTIONS CONCERNING THE LUER BD Q-SYTEM™ CLOSED ACCESS SYSTEM RELATING TO CERTAIN PPS® FLOW+ SAFETY HUBER NEEDLES**

The device is intended for use **ONLY** with:  
 - Intravascular lines  
 - any type of Luer-slip and Luer-lock male connectors compliant with ISO standards

In order to use correctly this product, health care providers should be properly trained and properly trained to the use of BD Q-SyTe™. Its use should be preceded by the application of a facility protocol.

• Maintaining the integrity of the line by following the aseptic techniques established in the facility protocol.

• The BD Q-SyTe™ device should be disinfected with an appropriate antiseptic, such as 70% isopropanol, prior to each use. Do not use acetone.

In order to minimize the blood reflux in the intravenous catheter, clamp the tubing before connecting.

• Changing the BD Q-SyTe™ according to the facility protocol, US CDC guidelines, or when product integrity has been compromised, is recommended. Do not activate (connection + disconnection) the BD Q-SyTe™ pre-split septum more than 100 times. The BD Q-SyTe™ device should not be used for sensitive, blind-ended, non-ISO compliant Luer connectors, or with ISO Luer connectors showing signs of wear.

If this warning is not followed, there could be leads and/or a poor performance of the device.

• Pain during the infusion: immediately stop the perfusion and seek medical advice and supervision.

• Check that the needle length suits the port and patient subcutaneous thickness: If it is too short, the needle or the implantable port may be damaged during puncture; if it is too short, the septum of the implantable port may not be entirely punctured by the needle and the medication might be administered in the surrounding tissues and/or the needle could be obstructed.

• Do not press too hard the tip of the needle on the bottom of the reservoir (a potential damage of the septum or a bend of the needle when the needle is removed).

If PPS® Flow+ are used during the infusion, the port (if plunger to a vertical position and totally locked in the port until a click is heard), the needle tip could come out of the plunger, potentially leading to an accidental puncture with a contaminated needle.

• Light-sensitive solutions, the tubing must be protected from light.

• PPSS® Flow+ should not be used for power injection of contrast media.

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