# Technovit® 2000 LC

## D Gebrauchsinformation Technovit 2000 LC System

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## GB User instruction Technovit 2000 LC system

### Application

Technovit 2000 LC is a light-curing resin for embedding of translucent ground sections within material testing. Technovit 2000 LC had been specifically developed for preparation of sensitive materials and micro components. The main application areas are: Semiconductor technique, microelectronics, medical equipment technique, optoelectronics and micro-system techniques.

Technovit 2000 LC varnish is a light-curing varnish to avoid the development of an inhibition layer on ground sections embedded in Technovit 2000 LC.

### Properties

The Technovit 2000 LC liquid is light-curing and is polymerized within the non-dangerous spectrum of the visible blue-light generated by the polymerization unit Technotray Power of Heraeus Kulzer. This allows an unlimited moving of a specimen within the embedding material (positioning or removal of air entrapings) during daylight or artificial illumination. If properly handled gmal can be sawed, milled, drilled, formed, ground and polished.

### Advantages

- dry and non-sticky surface of the ground section embedded in Technovit 2000 LC
- no loss of material (one-component resin)
- light-curing thus almost unlimited handling time
- optimum translucency
- exceptional low gap formation no air entrapments
- alcohol and acid resistant
- usable under vacuum, no gas evaporation (can be used under REM)
- odor free
- no UV-light but only blue light
- low polymerization temperature

### Handling

Surfaces exposed to Technovit 2000 LC must be clean and degreased for embedding. For this degreasing solvents such as light benzine or ester are recommended. Polyethylene moulds for embedding are recommended. Heraeus Kulzer embedding moulds are specifically suitable.

### Curing

- Polymerization is possible in block or in layers. Larger specimens should be embedded in layers. Polymerization takes place in the curing unit Technotray Power. This unit is specifically tuned to Technovit 2000 LC.
1. Pour about 5 mm high Technovit 2000 LC into embedding mould.
  2. Position specimen (if necessary fix specimen with Technovit 2000 LC fixing paste). Air bubbles must be avoided.
  3. Position embedding mould onto the glass disk of the drawer of the curing unit Technotray Power (please observe operating instructions of the Technotray Power unit).
  4. Irradiation time: 5 min.
  5. Fill up embedding mould with Technovit 2000 LC till 3 mm under the brim
  6. Final polymerization for 20 min takes place now.

The maximum polymerization temperature reached is 90 °C. If lower temperatures are desired curing in increments is recommended. Each layer must then be irradiated for approx. 3 min. Final polymerization time is always 20 min. Please observe polymerization times!

### Working with Technovit 2000 Inside Cure

In order to harden Technovit 2000 LC in porous samples, shady places or lighting-inaccessible places (internal threads, pipe sections, capacitors, etc.), Technovit 2000 LC is mixed with the Inside Cure additive in a mixture of 1000 ml to 40 ml. For this purpose pour the contents of the original packing of Technovit Inside Cure (40 ml) into an original packing of Technovit 2000 LC (1000 ml). Tightly seal the Technovit 2000 LC bottle and shake well.

As this mixture has a shorter shelf life than Technovit 2000 LC please place the "Inside Cure mixture" label accompanying the packing of the Technovit Inside Cure on to the Technovit 2000 LC bottle and note the life from the beginning of the mixture as the shelf life + 1 year.

Technovit 2000 LC can now be used as usual but it is not recommended to work in several layers or temperature-time programmes as then the polymerisation cannot take place in a regulated manner.

### Technovit 2000 LC varnish

Technovit 2000 LC varnish avoids the development of an inhibition layer on the surface of the ground section, embedded in Technovit 2000 LC.

### Application

1. Pour Technovit 2000 LC into the embedding mould leaving at least 3 mm space to the brim and polymerize for 10 min. minimum.
2. Take the drawer out of the Technotray Power polymerisation unit and fill up with a 2 – 3 mm thick layer of Technovit 2000 LC varnish.
3. The surface of the embedding has to be covered completely with the varnish.
4. Polymerize specimen for further 10 min.!

### Danger/Safety advices

Please follow the instructions on the product packaging and safety data sheets.  
www.kulzer-technik.com

### Delivery units

Technovit 2000 LC	– 1000 ml liquid
Technovit 2000 LC varnish	– 100 ml liquid
Technovit fixing paste	– 1 x 4 g
Technovit Inside Cure	– 1 x 40 ml liquid

® = registered trademark of Heraeus Kulzer GmbH, Wehrheim/Ts.

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® = marca registrada de Heraeus Kulzer GmbH, Wehrheim/Ts.