Gebrauchsanweisung | Instructions for use |
Mode d'emploi | Modo de empleo | Modalità d'uso

rematitan[®] Plus speed

Einbettmasse | investment material |
revêtement | evestimientto | rivestimiento

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Dear customer

Thank you for choosing a quality product from Dentaaurum.

It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit. Instructions for use cannot describe every eventuality and possible application. In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current instructions for use supplied with the product and stored in the internet at www.dentaaurum.com, even though you frequently use the same product.

1. Manufacturer

Dentaaurum GmbH & Co. KG | Turnstraße 31 | 75228 Ispringen | Germany

2. General description

rematitan® Plus speed is an investment material developed especially for titanium casting and covers a wide spectrum of indications. This product is a continued development of the well proven investment material rematitan® Plus and can be used for crown and bridgework as well as for metal denture casting. The speed characteristic reduces working time drastically. Only one mixing liquid is used for both crown and bridgework and for metal denture casting. Excellent results can be achieved by observing the following instructions for use. Users of Dentaaurum-rematitan®-Technology maintain the same style of waxing-up and attaching sprues. Users of other titanium casting systems must adapt this process within their usual system.

3. Safety Tips

When using the speed heating method, ensure that the furnace door remains closed continuously for at least 15 minutes after the casting ring has been inserted. (Danger of burning due to flames escaping)

Investments contain quartz. Avoid inhalation of dust as this may cause lung diseases (silicosis / lung cancer). Wearing of face mask type FFP2-EN 149:2001 is recommended. Open pouch with scissors and avoid dust development when pouring into mixing bowl. Rinse empty pouch with water before disposal.

Remove dust from working area only with a wet cloth.

To avoid dust when devesting, place ring into water until it is completely soaked after ring has reached room temperature.

When sandblasting, use suction with fine dust filter.

4. Availability

	rematitan® Plus Speed powder	80 x 250 g	20 kg	REF 107-620-00
	rematitan® Plus Speed mixing liquid	1 l		REF 107-625-00

5. Shelf life

Powder – 36 months in a dry storage area.
Mixing liquid – 18 months (sensitive to frost)

6. Material parameters of DIN EN ISO 9694 (Liquid 100%)

Initial setting time	approx. 5 – 6 min
Compressive strength	9 – 10 MPa
Fluidity	approx. 170 mm
Thermal expansion	approx. 0.65 %
Linear Setting Expansion (acc. to german standards DIN 13919 part 2)	approx. 0.05 %



7. Quality information








Dentaaurum ensures a faultless quality of the products manufactured by us. These recommendations are based upon our own experiences. The user himself is responsible for the processing of the products. Responsibility for failures cannot be taken, as we have no influence on the processing on site.

8. Explanation of symbols used on the label

Please refer to the label. Additional information can be found at www.dentaaurum.com (Explanation of the symbols used on the label REF 989-313-00).

C+B, Cast partial denture		Storage	20 °C  22 °C 	Maintain a storage temperature of 20 to 22°C/68 to 72°F! Cool the mixing liquid at higher temperatures! Do not expose the mixing liquid to frost! Caution with delivery in winter!
		Mixing ratio	250 g : 35 ml	500 g : 70 ml 750 g : 105 ml Abide by the mixing ratio exactly. Ensure of clean measuring beakers!
		Working time	4 – 5 min	Measured at working temperature (20 – 22°C/68 – 72°F). At higher temperatures cool the mixing liquid.
Crowns and bridges		Crowns and bridges	Observe the rematitan® casting system instructions for use! (Current instructions for use or instructions for dental use can be found on the internet at www.dentaurum.com)	When using the rematitan® casting system, follow the dental technical instructions regarding spruing and preparation techniques! Users of other casting systems should maintain their usual working methods!
		Wax tension reducing	Use Lubrofilm® or Lubrofilm® plus as a wax tension reducing	Dry at all costs! If no wax tension reducing is used, ensure of adequate vacuum mixing! Lubrofilm® plus (REF 112-053-00) is free of alcohol, superficial damage is therefore avoided.
		Casting system	rematitan® casting system: Casting ring base sizes 3 – 9 Metal casting rings or elastic casting rings	Metal rings: Size 3 casting ring require 1 mm thick, Kera-Vlies® (REF 127-250-00) Size 6 and 9 casting rings require 2 mm thick, Kera-Vlies® (REF 127-251-00) When using metal casting rings ensure that Kera-Vlies® protrudes over the ring edge. Scratch the surface of the investment material with a knife once it has set.
		Liquid concentration	Approx. 70 – 80% concentration for normal crowns and bridge work. For special uses such as telescopic work use up to 100%	When using the speed investment method, use a higher liquid concentration! (approx. 80 %) A higher liquid concentration produces wider fittings!
		Mixing	60 s mixing time under vacuum	Empty first liquid, then powder into the mixing bowl. Mix thoroughly by hand with a spatula. Ensure the mixing bowl is clean! No plaster residues! The mixing speed and shape of the mixing apparatus influence the mixture and its expansion. Ideal mixing speed: approx. 360 rpm.
Cast partial denture		Cast partial denture	Observe the rematitan® casting system instructions for use! (Current instructions for use or instructions for dental use can be found on the internet at www.dentaurum.com)	When using the rematitan casting system, follow the dental technical instructions regarding spruing and preparation techniques! Users of other casting systems should maintain their usual working methods.
		Duplication	Only use silicone duplicating material such as rema®-Sil or ecosil	
		Silicon tension reducing	 No	
		Liquid concentration for model and encasement	Approx. 70 – 80% liquid concentration for normal model casting techniques	When using the speed investment method, use a higher liquid concentration! (Approx. 80 %) A higher liquid concentration produces wider fittings!

Cast partial denture		Mixing	90 s mixing time under vacuum	Empty first liquid, then powder into the mixing bowl. Mix thoroughly by hand with a spatula. Ensure the mixing bowl is clean! No plaster residues! The mixing speed and shape of the mixing apparatus influence the mixture and its expansion. Ideal mixing speed: approx. 360 rpm.
		Setting time and model drying	 <p>Setting time: 40 min at 22 – 22 °C/68 – 72 °F Model drying: 30 min at 70 – 100°C/158 – 212°F</p>	Dry in a circulating air cupboard or a drying cupboard (70°C/158°F). Do not expose to any higher temperatures!
		Model hardening	 Not required	
		Casting system	rema®-Form casting rings and special funnel formers from the rematitan® casting system	Users of other casting systems may continue to use their accustomed casting rings.
Crowns + Bridges, Cast partial denture		Investing	Using gentle vibration	C+B: When using metal casting rings, fill to the edge of the ring liner. After the investment material has set, remove the upper surface with a knife or a trimmer. Cast partial denture: After the investment material has set, trim the entire ring base!
		Pressure investing	 No	
		Conventional heating		Setting time 40 min. Place the casting rings into the furnace with the opening at the base! Ideal heating speed 5 °C/min (41°F/min), holding time 250°C/60 min (482°F/60 min) End temperature: 950°C/1742°F, holding time 60 min Remove from the furnace and cool for 45 min at room temperature, then cast. Protect the casting rings left to stand over the weekend from drying out!
		Speed heating	 <p>Do not re-open the furnace for at least 15 min after having inserted a speed investment ring – danger of explosion!</p>	C+B: Standing time 15 min Cast partial denture: Standing time 30 min Then place immediately into the furnace at the end temperature 950°C/1742°F! Holding time 60 min Remove from the furnace and allow to cool for 45 min at room temperature and then cast.
		Cooling	 The ring must be quenched immediately after having been cast	
		Divesting	Sandblast using aluminium oxide, grain size 125 µm (C+B), and 125 – 250 µm (Cast partial denture)	Moisten the casting ring and then divest carefully. Do not strike with a hammer to divest!

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