Orthosin SP

Instructions

Orthosin SP is a self-curing acrylic resin for manufacturing removable orthodontic appliances.

Orthosin SP has been developed especially for the spray technique.

The advantages:

- immediately stabilizing - therefore quicker building up of the different layers of liquid/powder
- very well sprayable
- excellent brilliancy
- high transparency (no bubbles)
- no smearing
- colours available: pink-clear and clear.

Instructions:

1. Soak dry models for 5 min. in warm water. Isolate as usual with alginate isolation liquid (e.g. FuturaSep Plus of Schütz-Dental). Dry.

2. Put a little monomer to a plastic beaker.

3. Moisten the model surface with a brush and monomer and spray an equally thin coat of powder. Then add monomer until all dry spots of powder are soaked. Add more by the same procedure until desired thickness is achieved.

Finally smoothen the total surface with brush and monomer.

4. Polymerisation

Hardening is achieved in the polymeriser at 45°C warm water and 2-3 bar pressure.

Time: 10 minutes.

5. Finish and polish as usual.

Repairs:

Repairs are to be done with Orthosin SP. Before repairs are carried out, all fracture areas must be roughened, cleaned and wetted with monomer. Polymerize for 10 min. in a pressure polymerization unit at a pressure of 2-3 bar and a water temperature of 45°C.

Important:

Clean brushes after use. Never pour monomer back to the bottle.

Storage:

Store at room temperature (21°C). Please note date of expiry.

Hazard identification:

Avoid skin contact with the unpolymerized material and do not breathe in monomer vapours. Open containers: Do not use contents after date of expiry.

Side-effects:

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects - even if it is doubtful that the side-effect has been caused by our product - please kindly contact us.

Contra-indications/interactions:

If a patient has known allergies or hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. In such cases, we will supply the composition of our medical device upon request. The dentist should consider known interactions and crossreactions of the product with other materials already in the patient's mouth before using the product.

Note

Please supply the dentist with the above information, if this medical device is used to produce a special model. Please also consider the safety data sheets.

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