Suture Shuttle Needle**STERILE | R****Device Description**

The Smith & Nephew Suture Shuttle Needle is a single-use disposable device, designed for use in conjunction with the Smith & Nephew ELITE® PASS Premium Arthroscopic Suture Shuttle.


Figure 2

Figure 1

Indications for Use

The suture shuttle needle is indicated for passing suture through soft tissue for orthopedic repairs.

The suture shuttle needle is intended for use with any braided suture material equal in size to USP #0 through #2 suture. Use the suture shuttle needle with a 5.5 mm or larger cannula. The use of cannulae smaller than 5.5 mm may compromise the performance of or damage the instrument.

Contraindications

- Inufficient integrity of soft tissue to employ arthroscopic or limited-access suturing techniques.
- Large tears or lesions in soft tissue where arthroscopic or limited-access suturing techniques may provide poor durability of lesions.
- Soft tissue tears that cannot be properly accessed or effectively approached with the orientation of this instrument.
- Hard bone.

Warnings

Do not use if package is damaged. Do not use if the product sterilization barrier or its packaging is compromised.


- Contents are sterile unless package is opened or damaged. Do NOT RESTERILIZE. Discard any open, unused product. Do not use after the expiration date.
- The needle is single use only. DO NOT REUSE. Reuse can result in improper function or breakage. Broken needle parts may not be retrievable.
- It is the surgeon's responsibility to be familiar with the appropriate surgical techniques prior to use of this device.
- Read these instructions completely prior to use.
- While deploying the needle, do not open the grasper jaws as it may interfere with the performance of the instrument.
- Arthroscopic suturing requires specialized knowledge of knot tying and suture techniques. Knowledge of these and other appropriate techniques is an important consideration for successful utilization of this equipment.

Precautions

R U.S. Federal law restricts this device to sale by or on the order of a physician.

- Hazards associated with reuse of this device include, but are not limited to, patient infection and/or device malfunction.
- Prior to use, inspect the device to ensure it is not damaged. Do not use a damaged device.
- Careful attention must be paid to asepsis and avoidance of anatomical hazards.

Remove the disposable needle from the suture shuttle hand instrument and discard according to local and state regulations at the completion of the surgery (Figure 6).

**Warranty**

For single use only. This product is warranted to be free from defects in material and workmanship. Do not reuse.

For Further Information
If further information on this product is needed, contact Smith & Nephew Customer Service at +1 800 343 5717 in the U.S., or an authorized representative.

Instructions for Use

- To load the single use disposable needle into the instrument, insert the needle into the back of the instrument as shown in Figure 2. Place the tip of the straight portion of the needle shaft into the instrument bore. Retract the actuation lever and push the bent portion completely into the hole.

*Trademark of Smith & Nephew, registered U.S. Patent & Trademark Office.

©2004, 2010 Smith & Nephew, Inc. All rights reserved.

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*Varumärke som tillhör Smith & Nephew, och har registrerats vid U.S. Patent and Trademark Office (USA:s Patent och registreringsverk).
©2004, 2010 Smith & Nephew, Inc. Alle rättigheter förbeholdes.