addition protect XPress

GB
Instructions for use

Special silicone impression material for hearing protection, low-viscosity, thixotropic, quick-setting

Silicone ear impression material specially designed for fitting ear protection, fast and effective, e.g. with serially produced impressions. Ultra-short setting time and significantly less force required during application. Soft creamy consistency initially, good flow properties, does not run out of the ear. Easy to extract from mould and unproblematic when it comes to further processing.


Fig. 1


Fig. 2


Fig. 3


Fig. 4

Fig. 5

Fig. 6

Fig. 7


Fig. 8

1. Ear preparation

Before taking an impression, the auditory canal and the eardrum have to be examined (Fig. 1). If distinctive features are found (e.g. inflammations or perforation of the eardrum, tubes) impression must not be taken. The auditory canal has to be cleaned and depilated. In case of insufficient cleaning, remaining cerumen can cause retarded setting resp. inhibition of setting process. In order to achieve a microfine separating layer in addition to cleaning, the whole skin contacted area (auditory meatus, concha and marginal areas of the auricle) is moistened with **preclean** solution, using a wadding stick or a cotton pad. Conditioning of the ear with **preclean** solution offers the following advantages:

- easier removal of the otoplastics
- homogeneous, completely vulcanised surface, without any smearing layer (inhibition layer).

Before taking the impression, the eardrum must always be protected with an impression plug (e.g. blue secure) placed at the end of the external auditory canal in front of the eardrum (Fig. 2). Now **addition protect XPress** is injected directly into the prepared ear of the patient as described under point 2.

2. Mixing and dispensing

Place cartridge into Automix mixing gun (Fig. 3) and unscrew cap. Carefully extrude a small amount of material to ensure proper flow from both orifices (Fig. 4). Insert mixing cannula into guide grooves on cartridge and rotate anti-clockwise to a stop (Fig. 5). If necessary, mount also an Intra-Tip. Dispenser is now ready to mix silicone in any amount needed. Slowly apply the impression material into the auditory canal until the concha is filled completely (Fig. 6). Particularly light bodied **addition protect XPress** allows an almost pressure-free, steady flow into the ear. Keep the end of the cannula immersed on the surface of the material during dispensing. Leave mixing cannula on cartridge after use, thus sealing the cartridge. After curing remove the impression carefully from the ear (Fig. 7). Re-examine the ear (Fig. 8).

Important working hints

- Do not combine with condensation curing silicones.
- Cured impression materials are chemically inert – spots on clothing should be avoided.
- Be sure to remove all traces of impression material from the patient's ear.
- Latex gloves and latex contaminated surfaces as well as cerumen, cremes and resins may inhibit the setting reaction of **addition protect XPress** (we recommend standard gloves made of nitrile rubber or polyethylene).
- A careful preparation (point 1) as well as a thorough cleaning of the auditory canal is required before taking an impression.

Safety advice

The impression taking of the external auditory canal has to be executed by trained specialist staff only. The working instructions and precautionary measures have to be strictly observed. Non-observance might lead to irreparable damage of the ear or ear drum. DETAX is not liable for any damage caused by improper application of the impression material.

Further information: Silicone based materials are proven a million times. On condition of a proper application, undesired effects are not to be expected. However, reactions of the immune system like allergies, irritations, cannot be absolutely excluded. In case of doubt, we recommend to make an allergy test before the application of the material.

Indications for use:

- Ear impression for hearing protection
- Serially produced impressions

Technical data:

- **Mixing volume:** 50 ml (cartridge)
- **Mixing ratio:** 1:1
- **Colour code:** base: turquoise blue catalyst: light blue
- **Application time:** approx. 30 sec.
- **Time in the ear:** approx. 1:55 min.
- **Total set:** approx. 2:25 min. at body temperature
- **Final hardness:** approx. 38 Shore A
- **Strain in compression:** approx. 6 %
- **Recovery from deformation:** approx. 99.8 %
- **Linear dimensional change:** approx. 0.2 %
- **Application:** At 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5% rel. humidity
- **Storage:**

Ordering information:

addition protect XPress
Cartridge bag of 8 **03442**
8 x 50 ml cartridges base + catalyst

Cartridge box of 50 **03443**
50 x 50 ml cartridges base + catalyst

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