**ENGLISH INSTRUCTIONS FOR USE****I. INTRODUCTION**  
CLEARFIL PHOTO CORE PLT is a light-cure, radiopaque core build-up composite material with a special dispensing system for quick placement.**II. INDICATION**  
Core build-up of vital or non-vital tooth**III. CONTRAINDICATION**  
Patients with a history of hypersensitivity to methacrylate monomers**IV. INCOMPATIBILITY**  
Do not use eugenol containing materials for pulp protection or temporary sealing, since the eugenol could retard the bonding system curing process.**V. PRECAUTIONS**

- Avoid use of the product for patients with a history of hypersensitivity to methacrylate monomers.
- If any hypersensitivity occurs, such as dermatitis, discontinue use of the product and consult a physician.
- Wear gloves or take other appropriate protective measures to prevent the occurrence of hypersensitivity that may result from contact with methacrylate monomers.
- Use caution to prevent the product from coming in contact with the skin or getting into the eye. Before using the product, cover the patient's eyes with a towel to protect the patient's eyes from splashing material.
- If the product comes in contact with human body tissues, take the following actions:  
 - If the product gets in the eye:  
 Immediately wash the eye with copious amounts of water and consult a physician.  
 - If the product comes in contact with the skin:  
 Immediately wipe it off with a cotton pledgelet moistened with alcohol or gauze and wash with copious amounts of water.
- Use caution to prevent the patient from accidentally swallowing the product.
- Avoid looking directly at the dental curing unit when curing the product.
- For infection control reasons, PLT tips are for single use only.

**2. Handling and manipulation precautions**

- Do not use the product in conjunction with the other composite resins. Mixing materials may cause a change in physical properties, possibly a decrease, from the properties expected.
- When light curing the product, note the depth of cure in this instructions for use.
- The product polymerize under an operating light or natural light (sunlight from windows). Avoid operating light or natural light, and use within 5 minutes from dispensing.
- The emitting tip of the dental curing unit should be held as near and vertical to the resin surface as possible. If a large resin surface is to be light cured, it is advisable to divide the area into several sections and light-cure each section separately.
- Low intensity of light causes poor adhesion. Check the lamp for service life and the dental curing unit guide tip for contamination. It is advisable to check the dental curing light intensity using an appropriate light evaluating device at appropriate intervals.
- Squeeze the dispenser with a slow, steady pressure. Excessive force is not necessary.
- The use of the product is restricted to a licensed dental professional.

**3. Storage precautions**

- The product must be used by the expiration date indicated on the package.
- Keep away from extreme heat or direct sunlight.
- The product must be stored at 10–25°C/50–77°F when not in use.
- Do not refrigerate or store at low temperature. If the product is refrigerated or stored at a low temperature, it will become excessively soft.
- The product must be stored in proper places where only dental practitioners can access it.

**VI. COMPONENTS**

- Components  
Please see the outside of the package for contents and quantity.

**2) Ingredients**

- Principal ingredients:
  - Silanated silica filler
  - Silanated barium glass filler
  - Bisphenol A diglycidylmethacrylate (Bis-GMA)
  - Triethylenglycol dimethacrylate
  - dl-Camphorquinone

The total amount of inorganic filler is approx. 68vol%.  
The particle size of inorganic fillers ranges from 0.49 to 75 µm..

**VII. CLINICAL PROCEDURES****A. CORE BUILD-UP OF VITAL TOOTH****A-1. Cleaning tooth structure**  
Be sure the cavity is adequately cleaned. An adequately cleaned cavity assures maximum adhesive performance.**A-2. Moisture control**  
Avoid contamination of the treatment area from saliva or blood to produce optimal results. A rubber dam is recommended to keep the tooth clean and dry.**A-3. Cavity preparations**  
Remove any infected dentin and prepare the cavity in the usual manner.**A-4. Pulp protection**  
Any actual or near pulp exposure could be covered with a hard setting calcium hydroxide material. There is no need for cement lining or basing. Do not use eugenol materials for pulp protection.**A-5. Applying a matrix strip and wedges**  
Apply a matrix strip and wedges as needed.**A-6. Tooth surface treatment and bonding**  
Tooth surface treatment and bonding should be performed according to the Instructions for Use of the bonding system used (e.g. CLEARFIL TRI-S BOND or CLEARFIL SE BOND).**A-7. Placement and light curing of CLEARFIL PHOTO CORE PLT**

- Dispensing  
Place a PLT tip of the product into the dispenser barrel according to the Instructions for Use of the dispenser. The tip may be rotated to provide the proper angle for delivery into the cavity.

[CAUTION]  
Squeeze the dispenser with a slow, steady pressure. Excessive force is not necessary.  
Discard the tip after use and sterilize the dispenser according to the dispenser instructions for use. The product polymerize under an operating light or natural light (sunlight from windows). Avoid operating light or natural light, and use within 5 minutes from dispensing.

- Placement  
Incremental placement and light-curing each increment is strongly recommended especially in deep cavities.

- Curing  
Cure the resin with a dental curing unit. Hold the emitting tip as close to the resin as possible. The light curing time and depth of cure for each dental curing unit is used is shown below.

**Light curing time and depth of cure**

depth of cure	Dental curing unit	Light curing time			
		5 sec.	10 sec.	20 sec.	40 sec.
Conventional halogen*		—	—	3.5mm	4.0mm
Fast halogen*	4.0mm	4.5mm	—	—	—
Plasma arc*	5.0mm	6.0mm	—	—	—
LED*	—	—	4.0mm	4.5mm	—

*Dental curing unit		Wavelength range and light intensity			
Conventional halogen		Light intensity <sup>a</sup> of 300 - 550 mW/cm <sup>2</sup> in wavelength range from 400 - 515 nm			
Fast halogen		Light intensity <sup>a</sup> of more than 550 mW/cm <sup>2</sup> in wavelength range from 400 - 515 nm			
Plasma arc		Light intensity <sup>a</sup> of more than 2000 mW/cm <sup>2</sup> in wavelength range from 400 - 515 nm, and light intensity of more than 450 mW/cm <sup>2</sup> in wavelength range from 400 - 430 nm			
LED		Light intensity <sup>a</sup> of more than 300 mW/cm <sup>2</sup> in wavelength range from 400 - 515 nm			

- Peak of emission spectrum: 450 - 480 nm
- Evaluated according to ISO 10650-1.
- Wavelength distribution and light intensity values measured with a spectro-radiometer calibrated using an IEC or the NIST (National Institute of Standards and Technology) standard lamp

**A-8. Preparation of the abutment tooth**  
After curing the composite resin, prepare the abutment tooth.**B. CORE BUILD-UP OF NON-VITAL TOOTH****B-1. Cleaning tooth structure**  
Be sure the cavity is adequately cleaned. An adequately cleaned cavity assures maximum adhesive performance.**B-2. Moisture control**  
Avoid contamination of the treatment area from saliva or blood to produce optimal results. A rubber dam is recommended to keep the tooth clean and dry.**B-3. Cleaning cavity and preparing root canal**  
Remove the temporary sealing material from the cavity and filling material from the root canal. Using a Pico reamer, prepare and clean the root canal opening.**B-4. Try-in of the post**  
A post of the proper diameter is fitted to the prepared root canal and the length of the post adjusted.  
Sand blast the post as needed.**B-5. Either B-5a or B-5b**  
**B-5a. When PANAVIA F 2.0 is used.**

- Cement the post into the root canal according to the Instructions for Use of PANAVIA F 2.0.

[NOTA]  
Do not use CLEARFIL SE BOND when PANAVIA F 2.0 is used. Since PANAVIA F 2.0 plays the role as adhesive bonding agent, apply excess PANAVIA F 2.0 paste to the remaining tooth structure on which ED PRIMER II, a component of PANAVIA F 2.0 is applied and the surface of the post.

**B-5b. When another resin luting cement is used.**

- Cement the post into the root canal according to the Instructions for Use of the resin luting cement.

**B-6. Placement and light curing of CLEARFIL PHOTO CORE PLT**

- Dispensing  
Place a PLT tip of the product into the dispenser barrel according to the Instructions for Use of the dispenser. The tip may be rotated to provide the proper angle for delivery into the cavity.

[CAUTION]  
Squeeze the dispenser with a slow, steady pressure. Excessive force is not necessary.  
Discard the tip after use and sterilize the dispenser according to the instructions for use of the dispenser.

The product polymerize under an operating light or natural light (sunlight from windows). Avoid operating light or natural light, and use within 5 minutes from dispensing.

- Build-up  
Using appropriate instruments, build-up the resin to the appropriate anatomy being careful to avoid incorporation of air bubbles.

- Incremental placement and light-curing each increment is strongly recommended.

- Curing  
Refer to "Light curing time and depth of cure" in A-7.

**B-7. Preparation an abutment tooth**  
After curing the composite resin, prepare the abutment tooth.**[WARRANTY]**

Kuraray Noritake Dental Inc. will replace any product that is proved to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

[NOTE]  
"CLEARFIL", "CLEARFIL PHOTO CORE", "SE BOND", "CLEARFIL TRI-S BOND" and "PANAVIA" are trademarks of KURARAY CO.,LTD.

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"CLEARFIL", "CLEARFIL PHOTO CORE", "SE BOND", "CLEARFIL TRI-S BOND" e "PANAVIA" sono marchi di KURARAY CO.,LTD.