Drala® Phosphatzement

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For cementations:

Normal setting: 1,4 g of powder + 0,5 ml of liquid

The recommended by-weight ratio is 1,6 parts of powder to 1 part of liquid.

Fast setting: 1,6 g of powder + 0,5 ml of liquid

The recommended by-weight ratio is 1,9 parts of powder to 1 part of liquid.

For cavity lining:

Only normal setting: 1,9 g of powder + 0,5 ml of liquid

The recommended by-weight ratio is 2,2 parts of powder to 1 part of liquid.

Mixing

In order to achieve correct mixing of powder and liquid, a glass slab with a temperature of 18° - 21°C / 64° - 70°F and a plastic or stainless steel spatula should be used. Be careful to use only clean dry glass slab and spatula. The cement should not come into contact with any moisture! When mixing the cement, first a little of powder is added to the liquid and then gradually more until the correct consistency is reached. The powder is mixed in by using alternately circular and patting movements. Fresh amounts of liquid should never be added during mixing!

The cementing consistency is reached when, on lifting the spatula, the point so formed slowly sinks back again into the mixture.

The filling consistency is reached when, on lifting the spatula, the extended point folds over to form a small hook and doesn't sink back into the mixture.

Please note

- Increased plasticity or a longer setting period should not be obtained by using a thinner consistency but by cooling the mixing slab. A too thin mixture may impair the quality of the cement.
- Overdosing of powder or increased temperatures accelerate the setting time.

Application

The filling should be realized immediately after mixing. Be sure to avoid penetration of moisture!

Caution

The liquid contains phosphoric acid that may cause acid burns. Avoid contact with eyes, mucosa and skin. In case of accidental contact immediately flush affected body areas with generous amounts of water and soap. In case of contact with eyes immediately flush under running water for at least 15 minutes and seek immediately medical care.

Storage

Close recipients carefully after each use. Keep tightly closed at temperatures of 10 °C - 24 °C/50 °F - 75 °F.
Keep away from direct sunlight.

Indications

- Fixing cement for crowns, bridges, inlays.

Drala® Phosphate cement normal setting additionally:

- As a cavity lining / base under composite or silicate restorative materials, if a sufficient pulp protection is existing.

Technical

Data

DIN EN ISO 9917-1

Contra indications

- Drala® Phosphate cement is not allowed to be applied directly in close proximity to the unprotected pulp.

Attention

When working in close proximity to the pulp a calcium hydroxide lining should be applied to the cavity before applying the zincphosphate cement. Pulpal irritation may occur in rare cases with zincphosphate cements.

Mixing time

- 1½ minutes

Working time

- **For cementation:**
normal setting:
3 minutes*
fast setting :
2½ minutes*
- **For cavity lining:**
only normal setting:
2½ minutes*

Setting time

- **For cementation:**
normal setting :
5½ - 7 minutes**
fast setting :
3½ - 5 minutes**
- **For cavity lining:**
only normal setting:
3½ - 5 minutes**

* from beginning of mixing,
at 21 °C/70 °F, 50 ± 5 % rel.
humidity

** from end of mixing, at 37 °C/
99 °F, 50 ± 5 % rel. humidity

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GmbH & Co. KG