Collection liner for fluids suction canisters

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SUCTION CANISTER LINER

Warning

WARNING: DANGER INDICATION
Attention: Important indication

General information

- Read this manual carefully before using the FLOVAC® fluid suction system.
- AFTER UNPACKING, CHECK THAT THERE ARE NO PARTS MISSING OR DAMAGED.
- EACH TIME BEFORE USING THE DEVICE CARRY OUT THE OPERATIONS DESCRIBED IN THE CHAPTER “START-UP PROCEDURE”.
- The company will not accept any responsibility if these instructions are not observed.
- The medical device and its components or accessories do not contain parts in natural rubber latex.

Connections

- MAKE THE CONNECTIONS AND CHECK THE SEALS BETWEEN THE COMPONENTS AS DESCRIBED IN “INSTRUCTIONS FOR USE” SECTION. NON PERFORMANCE OF THESE CHECKS MAY COMPROMISE THE SAFETY AND FUNCTIONING OF THE DEVICE AND ANY CONNECTED ACCESSORIES (TUBES, PARTICLE CONTAINER, SUPPORT SYSTEMS, ETC).
- ANY SUCTION ADJUSTMENT DEVICE CONNECTED AND THE CONNECTION HOSES MUST COMPLY WITH EN ISO 10079-3 STANDARDS.
- INADVERTENT INVERSION OF THE CONNECTIONS COULD LEAD TO CONTAMINATION OF THE OPERATOR AND/OR OF THE VACUUM GENERATING SYSTEM.

Operation

- THE FLOVAC® SUCTION LINERS MUST ALWAYS BE USED WITH CARE AND ONLY BY PERSONNEL WHO ARE AWARE OF THE CONSEQUENCES OF THE ONGOING THERAPY.
- FULL CONTAINERS MUST BE HANDLED WITH GREAT CARE DURING TRANSFER TO THE DISPOSAL AREAS, FOLLOWING THE LOCAL PROCEDURES AND REGULATIONS.

The device has been designed and manufactured in conformity with the safety regulations established by standard EN ISO 10079-3.

Copyright[©]

Copyright® FLOW METER S.P.A. 2003-2010

All right reserved. None of the information contained in this publication must be used for purposes other than the original ones. No part of this publication may be reproduced without written permission. CA-MI S.r.l. devices are subjected to periodic revisions to maintain and improve the production standards, functioning and increase their reliability. For this reason, the contents of this publication may be modified at any time without prior warning. The device described in this publication is designed and manufactured by CA-MI S.r.l., Via Ugo La Malfa, 31 – 43010 PILASTRO (PR) – Italy.

Applications

The FLOVAC® suction liner is used to collect organic fluids by aspiration and has been designed and manufactured for "High flow and high vacuum" applications (EN ISO 10079-3).

It basically consists of a rigid supporting canister and of this disposable collection system, the latter formed by a liner bag hermetically welded to the lid.

The lid has connectors that make it possible to connect the suction liner with the vacuum line and with the patient. A hydrophobic shut off and antibacterial filter protects the device and/or the vacuum system. This filter also operates as an overflow valve, deactivating vacuum generation whenever the specified maximum filling limit is reached. The "PATIENT" port permits aspiration of fluids containing bone fragments once the right-angled reducing connector has been removed. A graduated scale printed on the external container makes it possible to monitor the volume of liquid collected. The size of the device to be used can be chosen on the basis of the real amount of fluid to be collected. Where a collection capacity greater than the maximum available is needed, the connector "TANDEM" can be employed. This permits connection of two or more equal devices to be joined in a stack, thus increasing the system's collection capacity.

The hermetic seal on the lid makes disposal particularly hygienic, simple and secure, for the patient and for the personnel dealing with the removal and disposal of the device.

Instructions for use

After connecting the "PATIENT" port to the suction tube and the "VACUUM" port to the vacuum line, the FLOVAC® can be used for the safe and aseptic collection of organic fluids.

Where specified, the patient connection tube and suction cannula are supplied with the container kit.

Start-up procedure

1. Insert the liner bag (A) and, pressing firmly all around the edge, apply the lid (B) to the reusable rigid container of the appropriate dimensions (C) (supplied separately), making sure that the system is completely sealed.
2. Close the "TANDEM" (D) connector using the supplied tap (E) and pressing it firmly.
3. Position the device in the appropriate ring support (F) with the sliding connection device (G), ensuring that the system is vertical and firmly connected.
4.  **Connect the vacuum feed source to the "VACUUM" port (H) using a tube with suitable size and characteristics, ensuring that the butterfly connection (I) is fully inserted in the cover and turned to "ON". when the canister is not in use, aspiration can be simply interrupted by turning the butterfly connector to "OFF".**
5.  **Connect a suitable patient tube (J) of appropriate dimensions and, if necessary, fitted with a cannula (K), to the "PATIENT" port (L) on the lid using, if necessary, the right-angled connector (M) supplied with the device. The "PATIENT" port permits aspiration of fluids containing bone fragments once the right-angled reducing connector has been removed. To prevent it from getting in the way, the tube may be easily attached to the support ring using the spring (P).**
6. When the right-angled connector is used (M), ensure that it is firmly held in port (L).

Start-up procedure

7. Before use, check all the seals and make sure that there are no leaks by turning on the suction and blocking the cannula (K) connected to the "PATIENT" port (L). There are no leaks in the system if the liner bag is fully extended, closely adhering to the wall of the rigid container and the lid is slightly buckled towards the inside of the canister.
8.  **DO NOT USE A SOURCE OF VACUUM OF LESS THAN -750 mbar (-570 mmHg) WITH THIS CANISTER.**
9. Turn on the suction and check regularly the level inside the canister. The vacuum check valve will interrupt suction if the fluid reaches the maximum level envisaged for the canister.
10.  **WHEN THE CHECK VALVE INTERVENES, IT IS NECESSARY TO DISCONNECT THE VACUUM WITHIN A MAXIMUM OF 5 MIN.**
11. Proceed as described in the instructions given in the "Instructions for disposal" section.

Use of several units in tandem

1. Fit the liner bags with the lids to all the canisters to be used, following the instructions given in point 1 of the previous paragraph.
2. Seal the "TANDEM" opening (D) on the last unit with the appropriate tap (E), pressing firmly.
3.  **Connect the vacuum feed source to the "VACUUM" port (H) of all units, using tubes with suitable size and characteristics and suitable line splitters (N) complete with exclusion taps (O).**
4.  **Connect the "TANDEM" port (D) of the first unit to the "PATIENT" port (L) of the second unit, using the right angled connector (M) supplied with the kit. Continue for all other canisters, using tubes (P) with suitable size and characteristics (e.g. TANDEM tube Ø8x11 L=0.38m – FLOW METER Ref. 970010120).**
5.  **Connect the patient tube (J), with suitable dimensions and characteristics and where necessary furnished with cannula (K), to "PATIENT" port (L) of the last canister, using if necessary the right angled connector (M) supplied with the kit.**
6. Carry out all the checks described in points 6 and 7 of the previous section.
7. Start up suction and proceed as described in points 8, 9, 10 and 11 of the previous section.

Use of the particle container

1. If necessary, a particle container can be fitted (FLOW METER ref. 000036100) to the FLOVAC device in order to hold the organic material necessary for laboratory analyses in a suitable containment bag.

Use of the particle container

2. After having fitted the bag with cover to the support container, following the instructions in the paragraph "Procedures for Use", remove the right-angled connector (M) from the port (L).
3. Insert the particle container (Z) in the adapter (R), pressing the cover firmly down around its entire edge, ensuring that the system is completely sealed. Connect the device to the bag with cover, inserting it firmly in the port (L).
4. **⚠ Connect the patient tube (J) with suitable dimensions and possibly equipped with a cannula (K) to the "PATIENT" port (S) of the cover, using if necessary the right angled connector (M) previously removed from the bag with cover. To avoid it causing an obstruction, the tube can be conveniently fixed to the support ring using the spring (P).**
5. If the right-angled connector (M) is used, ensure that it is firmly inserted in its port (S).
6. Perform checks as described in points 6 and 7 of the paragraph "Procedures for Use".
7. Initiate aspiration and proceed as described in points 8, 9, 10, and 11 of the paragraph "Procedures for Use".

Instructions for disposal

1. **⚠ THE DEVICE IS INTENDED FOR SINGLE PATIENT APPLICATION. AVOID RE.USE IN ORDER TO LIMIT THE RISK OF CROSS CONTAMINATION.**
2. **⚠ TURN OFF THE VACUUM AND REMOVE ALL THE TUBES CONNECTED TO THE LINER, GIVING PARTICULAR ATTENTION TO AVOID ACCIDENTAL CONTAMINATION. IF RIGHT-ANGLED CONNECTOR HAS BEEN USED, THIS MUST BE REMOVED.**
3. **⚠ IF THE PARTICLE CONTAINER IS USED, IT MUST BE REMOVED IMMEDIATELY AFTER THE CONTAINER'S DEACTIVATION TO AVOID DEGRADATION OF THE COLLECTED PARTICLES. THEN CONNECT THE APPROPRIATE STOPPERS (U – V) TO THE "PATIENT" CONNECTORS (S) AND (T), PRESSING HARD, BEING CAREFUL TO AVOID ACCIDENTAL CONTAMINATION. TRANSFER THE SUITABLY IDENTIFIED DEVICE TO THE LABORATORY, AND AFTER HAVING OPENED THE PARTICLE CONTAINER AND CUT THE INTERNAL CONTAINMENT BAG (Q), RESPECTING THE HOSPITAL'S APPLICABLE REGULATIONS, THE ORGANIC MATERIAL TO BE ANALYSED CAN BE REMOVED.**
4. **⚠ FIT THE APPROPRIATE PLUGS (E – P) TO THE "PATIENT" (L) AND "TANDEM" (D) PORTS, PRESSING THEM HOME FIRMLY, TAKING CARE TO AVOID ACCIDENTAL CONTAMINATION.**
5. **⚠ Turn the butterfly connector (I) to "OFF".**
6. **⚠ Remove the liner bag from the rigid container and transfer it to the waste disposal area, ensuring that all the openings are sealed, as well as any used tubes and cannulas, keeping in mind that the product is potentially infectious.**
7. The product must be disposed of in accordance with the current hospital regulations.

Instructions for disposal

7.  The reusable rigid canister must be cleaned and disinfected with water and neutral detergents and/or sterilised in an autoclave a 121 °C (relative pressure 1 bar) for 15 min. **DO NOT USE SOLVENTS OR ALCOHOL FOR CLEANING AND DISINFECTION: THE USE OF THESE PRODUCTS MAY DAMAGE THE UNIT.** The reusable canister mechanical resistance is guaranteed for 30 cleaning and sterilisation cycles under the specified conditions (EN ISO 10079-3). beyond this limit it may show signs of physical and mechanical damage to the plastic and its substitution is thus advised.

Technical features

Sales description	FLOVAC®
Reference code	31845 - 1.0 L
	31846 - 2.0 L
	31847 - 3.0 L
	000036030 - 1.0 L with 1.8m tube Ø 8x11
	000036031 - 2.0 L with 1.8m tube Ø 8x11
	000036032 - 3.0 L with 1.8m tube Ø 8x11
	000036050 - 1.0 L with 1.8m tube Ø 8x11 and vacuum stop
	000036051 - 2.0 L with 1.8m tubo Ø 8x11 and vacuum stop
	000036052 - 3.0 L with 1.8m tubo Ø 8x11 and vacuum stop
	31840 - Reusable 1.0 L canister
	31841 - Reusable 2.0 L canister
	31842 - Reusable 3.0 L canister
	970010120 - TANDEM tube
	000036100 - Particle container
"VACUUM" connector	Ø 9.0÷10.0 with ON-OFF switch
"PATIENT" connector	Ø 14.0÷15.5 (Ø 8.0÷9.2 with right-angled connector)
"TANDEM" connector	Ø 8.0÷9.2
"TANDEM" tube (optional)	Ø 8x11 L=0.38m
Particle container (optional)	"PATIENT" connector: Ø 14.0÷15.5 (Ø 8.0÷9.2 with right-angled connector) Outlet connector: Ø 14.0÷15.5 F.
Maximum vacuum value	-750 mbar (-570 mmHg)
Maximum flow value	42 LPM under the recommended conditions
Maximum graduation interval	50 ml
Gradation precision	±50 ml
Patient tube dimensions	Min. Ø 6x9 (recommended Ø 8x11 L=2.5m max)
Vacuum tube dimensions	Min. Ø 6x9 (recommended Ø 8x11 L=1.8m max)
Storage conditions	-40 °C ÷ +60 °C
Working conditions	-18 °C ÷ +50 °C

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Copyright©