Adherence[®] DC

Dual Cure Permanent Resin Cement with Fluoride

INDICATIONS:

Adherence[®] DC is indicated for the following applications:

- Final cementation of implant prostheses, crowns and bridges, inlays and onlays, Maryland bridges and veneers.
- Final cementation for pins and posts.
- Periodontal splinting.

PROPERTIES:

- *Adherence DC* is radiopaque and provides the high compressive and flexural strength and low solubility required for these types of materials.
- *Adherence DC* is designed to meet requirements per the ISO 4049:2000 standard (Type 2, Class 3 product).

FORMULATION:

- **Catalyst:** BisGMA based resin system with glass filler (68wt%) and peroxide initiator.
- **Base:** BisGMA based resin system with glass filler (61 wt%) with amine catalyst and camphorquinone photosensitizer.

DIRECTIONS FOR USE:

IMPORTANT: Please double bleed cartridge before each use. (Failure to do so may prevent proper setting)

- Using an auto-mixing tip, equal parts of catalyst and base are automatically mixed. Initial self cure set occurs in 3 - 5 minutes, but cement clean up can be started in 1.5 to 2.5 minutes when material becomes rubbery. Final set occurs in 5 - 9 minutes. Light exposure accelerates cure time.

I. Ceramic or composite indirect restoration to enamel/ dentin

- After try-in, place etched ceramic restoration in acetone or ethanol in ultrasonic cleaner for 1-2 minutes.
- Place a silane-coupling agent on internal surface of restoration, following the manufacturer's instructions, and allow to dry completely.
- Clean the tooth with non-fluoridated prophylaxis paste or flour of pumice, wash and dry as indicated with the bonding system used.
- For veneers, use a bonding agent of your choice. Etch and apply the bonding system according to instructions. **Use a dual-curing bonding agent for inlays, onlays and crowns. Avoid the use of self etch-primer adhesive.**
- Select the appropriate shade of *Adherence DC* and place directly onto the silanated surface of the restoration using the auto-mixing tip.
- Seat the restoration in place on the tooth. When you are satisfied that the restoration is properly seated, light cure for 1-2 seconds (3-4 seconds with white opaque shade) and then remove all excess cement. Continue light curing all surfaces for a minimum of 40 seconds per surface (minimum radiance, 400 mW/cm²). **Note:** For thicker veneers and darker shades a longer light curing time is recommended; see inlay placement procedure.
- Finish and polish margins with appropriate materials.

II. Metal to enamel/dentin

Metal restorations to teeth - crowns, inlays, bridges, and posts, including Maryland bridges:

- After try-in, sandblast the tooth side of restoration, wash and dry.
- If desired, use a metal priming system of your choice to treat the metal to be exposed to the cement.
- Clean the tooth with non-fluoridated prophylaxis paste or flour of pumice; wash and dry.

D. Use a dual-curing bonding agent. Avoid the use of self etch-primer adhesive. Etch and apply the bonding system according to instructions.

- Select an appropriate shade of *Adherence DC* and place directly onto the etched/primed surfaces of the restoration using the auto-mixing tip.
- Gently seat the restoration in place on the tooth. When you are satisfied that the restoration is properly seated, light cure exposed margins for 1-2 seconds (3-4 seconds with white opaque shade) and then remove all excess cement. Continue light curing any exposed margins a minimum of 40 seconds per surface (minimum radiance, 400 mW/cm²). The cement without light exposure will give 3 minutes of working time and completely set in 9 minutes.
- Finish and polish margins with appropriate materials.

III. Metal to Metal

Metal to metal cementation (implant prosthesis or crown on post):

- Use an air-abrasive system (sandblasting) on all restoration metal surfaces to be exposed to cement. Intraoral metal should be roughened with a bur or micro-etcher; then wash and dry.
- If desired, apply a metal priming system of your choice to prepare the metal restoration and the intraoral metal for the cement.
- Clean posts or abutments with non-fluoridated pumice prior to cementation to remove any plaque or other debris from the surface, wash and dry.
- Select an appropriate shade of *Adherence DC* and place directly onto the etched/primed surface of the restoration using the auto-mixing tip.
- Seat the restoration in place on the tooth. When you are satisfied that the restoration is properly seated, light cure exposed margins for 1-2 seconds (3-4 seconds with white opaque shade) and then remove all excess cement. Light cure any exposed margins a minimum of 40 seconds per surface (minimum radiance, 400 mW/cm²). The cement without light exposure will give 3 minutes of working time and completely set in 9 minutes.
- Polish margins with appropriate materials.

WARNINGS:

- Contains methacrylate resins. Avoid use of this product on patients with known methacrylate allergies. To reduce the risk of allergic response, avoid exposure to uncured resin.
- Do not use eugenol-containing materials for pulp protection and temporary sealing since they can retard the curing process.
- Do not allow saliva or water contamination of etched tooth surface or unset material.
- Any non-sterilizable item that is handled in the delivery of the dental service should be disinfected by standard dental office hygiene procedures.
- Wear gloves.

STORAGE:

- 59°F - 75°F / 15°C - 24°C
- Storage in a refrigerator ensures longest possible shelf life. Allow cement to come to room temperature before use.
- Do not store in proximity to eugenol-containing products.
- Do not expose the cement to elevated temperatures.

For professional use only in the practice of dentistry

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