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# Anclacem attachment

Editing date: 29.11.2011

## Dual curing attachment luting material.

### 1. Product description:

- **ANCLACEM ATTACHMENT** is a gingival-coloured, dual curing attachment luting material for bonding attachments in acrylic dentures.
- **ANCLACEM ATTACHMENT** is appropriate for use in most of dental attachments into an overdenture or partial denture, above all for ball and Locator® elements.

### 2. Product specifications:

- First setting range: 40 seconds (light-curing).
- Total setting time (auto-curing): 7 - 9 minutes.

### 3. Components:

- 1 syringe ANCLACEM ATTACHMENT 8g (Dual resin)
- 1 bottle 3g Resin primer
- 1 bottle 5ml Metal primer
- 8 cannulas Anclacem attachment

### 4. Composition:

- Monomer matrix: Diurethandimethacrylate; 1,4 Butandioldimethacrylate
- Total filler content: 75 % weight

### 5. Instructions for use:

#### A) Preparation of denture base:

1. Place attachment on implants, check fit of attachments
2. Mark position of matrix on denture base (e.g: using felt-tip pen or occluspray)
3. Use an adequate bur to make the hole in the denture where the attachment must be positioned.  
NOTE: The hole in the denture must be enough deep to avoid the attachment doesn't come into contact with the denture base and also to leave enough space for ANCLACEM ATTACHMENT.
4. Make a control accuracy of fit of matrix in denture (classical silicone check)

#### B) Insulating:

1. Insulate with vaseline all areas that are not to be bonded, e.g: implants or denture base surface. Block out undercuts and gingival margin around the implants with a suitable silicone o wax.
2. Alternatively, cover these areas with Teflon tape to prevent material flowing into undercuts around the implant.

#### C) Preparation of denture base:

1. Apply a thin layer of RESIN PREPARATION PRIMER on the parts of denture base to be bonded. Light-curing (e.g 30 seconds LED lamp) the attachment with ANCLACEM ATTACHMENT.

#### D) Preparation of matrix:

1. After checking a correct positioning of the attachment, clean with alcohol. Apply METAL PREPARATION PRIMER in the attachment and let dry for approx 1 min.  
NOTE: Do not touch the surface after cleaning.

#### E) Application of ANCLACEM ATTACHMENT:

1. Apply ANCLACEM ATTACHMENT in the hole maked to positioning the attachment.
2. Place denture in the mouth and let the patient bite with moderate pressure. Check fit of denture.
3. Light-curing during 40sec with a lamp to accelerate the cured. After wait 2-3 minutes approx before removing the denture from the patient's mouth.  
NOTE: Remove the cannulas each time after use and replace cap of syringe (light protection).

#### F) Finishing and final check:

1. After cementation attachment, remove block-out material and/or teflón tape.
2. Finish denture surface with a carbide bur and polish with an adequate polisher.
3. Subsequently, control fit of denture and occlusion in the mouth of the patient.

### 6. Attention!

- The required physical properties can only be achieved if light sources are working properly. We strongly recommend to regularly check light intensity according to the manufacturers indications.
- **Contra-indications/interactions:**  
If a patient has known hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. The dentist should consider known interactions and crossreactions of the product with other materials already in the patient's mouth before using the product.
- **Side effects:**  
If processed and used properly undesired side effects of this medical device are not to be expected but immune responses (e. g. allergies) or local discomfort cannot be ruled out completely. In case of any undesired side effects – even in cases of doubt – please contact us.
- Use only dentist and prosthetic (professional use).

# Ancladén

Ancladén, S.L. | Rubinstein, 2 Entlo 1ª | 08022 Barcelona  
93 417 29 84 | ancladen@ancladen.com | www.ancladen.com