
Orthocryl® LC

Gebrauchsanweisung | **Lichthärtender Kunststoff**
Instructions for use | **Light-curing acrylic**
Mode d'emploi | **Résine photopolymérisable**
Modo de empleo | **Acrílico fotopolimerizable**
Modalità d'uso | **Resina fotopolimerizzante**

CE 0483

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Dear Customer

Thank you for choosing this high-quality product from Dentaureum.

In order to achieve the best results with this product, it is important to carefully study and follow these instructions for use. The written instructions cannot cover all possible situations that can occur during use. For this reason, our hotline is available to answer any other questions or concerns that may arise.

Due to constant developments, we recommend, even when using the same products regularly, that you study the enclosed current instructions for use or refer to our website at www.dentaureum.de.

1. Manufacturer

Dentaureum GmbH & Co. KG
Turnstraße 31 | 75228 Ispringen | Germany

2. General product description

Orthocryl® LC is a light-curing acrylic for removable and fixed orthodontic appliances. In addition, Orthocryl® LC can be used to fabricate bite plates, anti-snoring devices and drill templates for implantology. Orthocryl® LC does not contain methyl methacrylate or dibenzoyl peroxide!

3. Form of delivery

Orthocryl® LC is delivered in opaque cartridges that fit the injector (REF 160-420-00) (Fig. 1). The Orthocryl® LC bite paste is provided in opaque jars.

4. Preparation

- Remove small blue protective cap on Orthocryl® LC cartridge by turning it counterclockwise (Fig. 2).
- Remove large blue protective cap on Orthocryl® LC cartridge by pushing up protruding flap (Fig. 3).
- Shorten a replaceable cartridge tip (REF 160-421-00) to the correct length and rotate onto the top of the cartridge as far as it will go (Fig. 4).
- The button on the upper side of the injector frees the forward and back feed of the injector plunger. With the button pressed, bring the plunger of the injector into initial position.
- Stick the cartridge onto the injector and lock it with a 90° turn (Fig. 5).
- With the button on the upper side pressed, push forward the plunger of the injector up to the stop.
- Apply emerging Orthocryl® LC on wires, screws and models with care.

5. Application of Orthocryl® LC

Example: expansion plate

Apply Orthocryl® LC on a dry, not soaked, coated plaster model. Mark the position of the expansion screws on the model and drill a hole or a groove with a conical bur (e.g. REF 123-598-00) for the screw holder. Fill the hole on the model with heated wax (e.g. REF 119-150-00) and insert the expansion screw with the lower, basal part of the holder. Remove expansion screw after wax has hardened.

Fig. 1


Fig. 2

Fig. 3

Fig. 4


Fig. 5

Fix bent wire elements to the model with sticky wax (REF 122-300-00).

Application of Orthocryl® separating medium (REF 162-800-00)

- Apply separating medium from the dosing bottle (REF 162-100-00) directly on the dry model and spread evenly using a paintbrush.
- Stand the model onto its back so that excess fluid can drain away.
- After approx. 30 seconds, carefully blow away separating medium residues or droplets with compressed air and, if necessary, remove residues from screws and wire retentions with a dry paintbrush.

Start application of acrylic immediately afterwards!

Apply Orthocryl® LC to the removed expansion screw so that the screw is completely surrounded by the material (Fig. 6). Please observe that the expansion screw is completely embedded in acrylic without any air bubbles. Afterwards, insert the expansion screw back into the drilled hole with the wax coating. Due to the wax, the expansion screw will regain the same position as determined before.

With wire retentions, apply Orthocryl® LC to the space between wire and plaster model first and then embed the complete wire (Fig. 7).

Next, apply Orthocryl® LC until the desired form and thickness of the appliance (maxilla approx. 3 mm / mandible approx. 4 – 5 mm) is achieved (Fig. 8). Should this process require more than a single step, slightly dip the tip of the cartridge into the already applied material, when starting the following step. This prevents air inclusions between the layers.

Eliminating air bubbles

If air bubbles should build up during application, they can be opened and filled up with a pointy instrument, e.g. a probe. You can also apply additional Orthocryl® LC below the air bubble so that it moves towards the surface and disintegrates.

Application of colored Orthocryl® LC

In order to guarantee secure hardening of the color layers, we recommend applying a layer of clear Orthocryl® LC first.

Orthocryl® Disco glitter

Spray a bit of Orthocryl® Disco glitter on the first layer of Orthocryl® LC directly from the spray bottle. Then, cover it with another layer of Orthocryl® LC.

Application of novel designs

Place the cut novel design on the first Orthocryl® LC layer, at about the height of the wire retentions. Then, cover it with another layer of Orthocryl® LC.

5.1. Application of the Orthocryl® LC bite paste

Example: activator

After the upper and lower model have been placed in the fixator by means of an activating bite*, both parts of the fixator are separated and models are prepared as described in section 5.

Orthocryl® LC is then applied to the palatal side of the upper jaw and to the lingual side of the lower jaw. Thereafter the Orthocryl® LC bite paste is taken from the jar with a spatula or a wax knife and applied onto the occlusal or the incisal part of the lower model. To ensure the material is entering all areas of the structure (bite ramp) and to avoid air bubbles, make sure you apply an excess amount of bite paste.

Once the bite ramps have been recovered with the material, both parts of the fixator are put together again, the excess material is removed, if necessary, and the closed fixator is placed in the light-curing unit.

* Models of an average size placed in the fixator 072-004-00 in a space-saving manner (close to telescopic bars) have enough space in the Solidilite V unit of Shofu.

Fig. 6


Fig. 7


Fig. 8


Fig. 9


Fig. 10

6. Polymerization

The specifications on polymerization times refer to the recommended light-curing unit Solidilite V from SHOFU Dental GmbH. Please observe the instructions for use of the unit. The rotary plate in the curing chamber is height-adjustable.

During polymerization of the described expansion plate, the rotary plate should be at the highest position possible according to the instructions for use of the light-curing unit Solidilite V.

The polymerization of Orthocryl® LC in the light-curing unit takes place at a wave length of 400 – 550 nm and a light output of 4 x 150 watt (600 watt). The time of polymerization for a single jaw appliance is 3 min. (6 min. for bimaxillary appliances) (Fig. 9). Afterwards, carefully remove the appliance from the model (Fig. 10) and cure the appliance for another minute (3 min in case of bimaxillary appliances) with the basal side up.

Remove the air-inhibition layer (smear layer) on the polymerized appliance, see point 7.

7. Removing inhibition layer

Completely remove the air-inhibition layer (smear layer). Apply Orthocryl® LC cleaning liquid (REF 160-430-00 / REF 160-431-00) on a piece of fabric and remove the smear layer.

Caution!

Please make sure to completely remove the air-inhibition layer (smear layer) to eliminate patient risks (see safety data sheet).

A direct application in the mouth is not permitted!

8. Finishing

When finishing Orthocryl® LC appliances, proceed as with Orthocryl® cold polymer appliances. However, do not use cross cut hard metal burs. Hard metal burs with transversal tothing or burs for soft acrylics are better suited.

Prepolish appliance with silicone polishers or sand paper, and then with pumice powder and polishing brushes. Achieve a high polish finish with the acrylic polishing material "Edelweiß" (REF 190-100-00 or REF 190-200-00) and a cloth polishing brush.

Note:

Due to stabilizing fillers, Orthocryl® LC appliances exhibit a slight milky haze.

9. Information on the constituents

See safety data sheet.

10. Safety instructions

For curing Orthocryl® LC, use light-curing units that emit light at a wavelength of 400 – 550 nm and have a light output of 4 x 150 watt (600 watt).

In order to avoid direct skin contact with Orthocryl® LC during processing, we recommend wearing protective gloves (e.g. Dermatril® disposable nitrile gloves, REF 162-907-00, REF 162-908-00, REF 162-909-00), see safety data sheet.

A direct application of Orthocryl® LC in the mouth is not permitted!

11. Storage and shelf life

Since Orthocryl® LC is a light-curing acrylic, close cartridges after use and store away from light.

The shelf life in the closed original package is 30 months from the date of fabrication:

- Store in a cool environment (not above 25 °C/ 77°F)
- Protect from contamination and sunlight as they may cause premature polymerization.

12. Cleaning and disinfecting Orthocryl® LC appliances

Appliances should be cleaned regularly and stored in water when not in use to prevent problems with hygiene. This procedure prevents crack corrosion or brown discoloration of the wire components.

For professional cleaning in practices, we recommend using the activeblue® high performance cleaning unit (REF 098-000-00). In combination with activefluid® cleaning liquid (REF 155-020-00) specifically developed for orthodontic appliances, fast and effective cleaning of plaque and tartar is ensured. Specially in areas that cannot be cleaned with a normal tooth brush.

13. Delivery

Orthocryl® LC Starter Set

REF 160-400-00

Content: Orthocryl® LC, 30 g of each color (clear, transparent pink, red, blue, yellow, green), 1 injector, 50 cartridge tips and step-by-step instructions

Orthocryl® LC Single cartridges and accessories

Designation	Quantity		Included in starter set
Orthocryl® LC, clear	30 g	REF 160-401-00	•
Orthocryl® LC, transparent pink	30 g	REF 160-402-00	•
Orthocryl® LC, red	30 g	REF 160-403-00	•
Orthocryl® LC, blue	30 g	REF 160-404-00	•
Orthocryl® LC, yellow	30 g	REF 160-405-00	•
Orthocryl® LC, green	30 g	REF 160-406-00	•
Orthocryl® LC bite paste, transparent	50 g	REF 160-410-00	
Orthocryl® LC injector	1	REF 160-420-00	•
Orthocryl® LC cartridge tip	50	REF 160-421-00	•
Separating medium	1000 ml	REF 162-800-00	
Orthocryl® LC LC cleaning liquid in spray bottle	100 ml	REF 160-430-00	
Orthocryl® LC cleaning liquid in bottle	1000 ml	REF 160-431-00	

14. With reference to our quality standards

Dentaurum ensures a faultless quality of the products that we manufacture. These recommendations are based upon our own experiences. The user himself is responsible for the processing of the products. Responsibility for failures cannot be taken, as we have no influence on the processing on site.

15. Disposal

Please observe the requirements for the disposal of non-polymerized Orthocryl® LC components stated in the safety data sheet.

Please refer to the label. Additional information can be found at www.dentaurum.de (Explanation of the symbols REF 989-313-00).

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