**SUPRATHEL<sup>®</sup>CW**

THE **TEMPORARY** SECOND SKIN

**Gebrauchsinformation**

**DE**

# SUPRATHEL<sup>®</sup>CW

THE TEMPORARY SECOND SKIN

## Produktbeschreibung

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THE TEMPORARY SECOND SKIN

## Product description

SUPRATHEL<sup>®</sup>CW is an absorbable, microporous membrane and approved as alloplastic skin substitute for the treatment of epidermal and dermal wounds.

## Indications

SUPRATHEL<sup>®</sup>CW is used for chronic wounds, such as venous and arterial ulcers, as well as for diabetic wounds.

## Composition

Copolymer of polylactide, trimethylene carbonate and  $\epsilon$ -caprolactone (Lactocapromer).

## Properties and mechanism

SUPRATHEL<sup>®</sup>CW is an alloplastic, absorbable skin substitute that is highly permeable to oxygen and water vapour, providing a particularly favourable environment for wound healing.

Clinical studies confirmed significant pain reduction when using SUPRATHEL<sup>®</sup>CW.

SUPRATHEL<sup>®</sup>CW exhibited good adhesion to the wound base by absorbing wound exudate and prevents mechanical alteration.

SUPRATHEL<sup>®</sup>CW is elastic/plastic and is therefore optimally conformable to all parts of the body, including the hands, feet, and arms. The extremities can be mobilized without loss of the wound covering.

SUPRATHEL<sup>®</sup>CW protects against bacterial and viral infections. Local infections may be treated with antibiotics through the permeable membrane without removal.

## Application

### Preparation of the wound

Prior to application, the wound should be sufficiently cleaned. For necrotic or sloughy wounds, thorough debridement with assured removal of these areas is indicated. Before applying SUPRATHEL<sup>®</sup>CW, bleeding should be stemmed and any excessive disinfectant solution dried off.

### SUPRATHEL<sup>®</sup>CW-application

Please open the package and retrieve under sterile conditions the white membrane from the green paper protection sheet.

For chronic wounds it is particularly important that the intact skin is not covered by the membrane. It is therefore recommended to cut SUPRATHEL<sup>®</sup>CW to a size exactly matching the wound, place the membrane smooth on the wound area and remove any air bubbles under the membrane by stroking laterally from the centre outwards. There is therefore no maceration of the surrounding skin. The anti-inflammatory action of SUPRATHEL<sup>®</sup>CW can therefore exert itself fully on the actual wound.

### Wound dressing

In chronic wounds, the subsequent dressing should be selected in such a way as to support an optimal wound-healing environment, in accordance with guidelines.

A combination with a tulle gras dressing is highly recommended, which should also be changed in the same rhythm as the SUPRATHEL<sup>®</sup>CW dressing.

# SUPRATHEL<sup>®</sup>CW

THE TEMPORARY SECOND SKIN

## Wound healing/Removal of SUPRATHEL<sup>®</sup>CW

Regular, generally once-weekly changing of the SUPRATHEL<sup>®</sup>CW membrane is indicated for chronic wounds. For heavily exuding wounds, the membrane can be changed more frequently as appropriate. In this case, appropriate exudate management is recommended

## Additional information

If local infection occurs, it may initially be treated with antibiotics through the permeable membrane.

Even in the unlikely case of an in-growth of the membrane into the wound base an operative removal of the membrane is not required as SUPRATHEL<sup>®</sup>CW degrades within 4 to 6 weeks and will be resorbed subsequently.

## Adverse reactions

So far no adverse side effects which could be attributed to the product were observed.

In rare cases very mild inflammatory response typical for foreign body reaction may be observed. This usually declines during degradation of the membrane.

## Contraindications

SUPRATHEL<sup>®</sup>CW should not be applied to infected areas.

SUPRATHEL<sup>®</sup>CW should not be applied to chronic dry wounds.

SUPRATHEL<sup>®</sup>CW should not be applied on severe bleeding wounds without additional hemostatic treatment.

## Additional informations

SUPRATHEL<sup>®</sup>CW may only be stored in original packaging dry and cool (22°C max.). In this case storage duration is max. 3 years.

Content is sterile unless packaging is undamaged. In case of any damage of the packaging the sterility of the product is not ensured anymore. The unused contents of opened or damaged packages are to be discarded.

Do not reuse and do not resterilize, this may change the material properties.

SUPRATHEL<sup>®</sup>CW is not toxic, not allergenic.

SUPRATHEL<sup>®</sup>CW is gamma sterilised.

Protect the product from children.

# SUPRATHEL<sup>®</sup>CW

THE TEMPORARY SECOND SKIN

## Explanation of symbols

 CE mark with identification number of the notified body, certified in accordance with ISO 13485; Directive 93/42/ECC with addendum 2007/47/EC

 Do not reuse

 Sterilization by radiation

 Temperature limit

 Use by

 Attention! Reference to accompanying document

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PolyMedics Innovations GmbH