DentAurum Vest

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Dear customer

Thank you for purchasing this product from Dentaaurum. You have chosen a product that will provide you with consistently high quality, reliability and service.

Optimum performance of our dental technology products is based on working precisely to the instructions issued by our Research and Development Department. We have therefore included a clear, detailed summary of each working stage in these instructions for use.

If you have any further questions, our team of dental technicians will be pleased to offer advice.

Simply call: Customer Support Hotline Tel. no. +49 72 31/803-410

Information, instructions for use and safety data sheets for Dentaaurum products are also available on the internet under: www.dentaaurum.de.

Warranty

Dentaaurum guarantees faultless quality of the products manufactured by us. The information in these instructions for use is based upon our own experiences. The user himself is responsible for the correct processing or use of the products. We cannot be held responsible for failures as we have no influence on the processing. However, any claim for compensation only applies to the commercial value of our products.

Instructions for use – DentAurum Vest

DentAurum Vest is an investment material specially developed for precious metal alloys and can also be used universally. The amount of expansion required can be exactly controlled, this guarantees an excellent fit in all areas of crown and bridgework. In order to benefit from these advantages it is essential to follow the instructions for use carefully.

Safety tips

When using speed preheating, the door of the furnace must not be opened within the first 15 minutes after the ring has been inserted. This precaution is taken in order to avoid burns.

Investments contain quartz. Avoid inhalation of dust as this may cause lung diseases (Silicosis/ Lung Cancer). Wearing of face mask type FFP 2 – EN 149:2001 is recommended. Open pouch with scissors and avoid dust development when pouring into mixing bowl. Rinse empty pouch with water before disposal.

Remove dust from working area only with a wet cloth.

To avoid dust when devesting place ring into water until it is completely soaked after ring has reached room temperature.

When sandblasting use suction with fine dust filter.

Form of supply

DentAurum Vest	40 x 150 g	6 kg	REF 106-650-00
DentAurum Vest mixing liquid	1 l		REF 106-651-00

Shelf life

Powder – 36 months with dry storage

Liquid – 18 months (frost-sensitive)

Values of materials according to DIN EN ISO 15912 (Liquid 100 %)

Type **1**, class **1** **2**

Start of solidification:	7 min
Resistance to pressure:	4.5 MPa
Fluidity:	140 mm
Thermal expansion:	1.16 %

Instructions for use – DentAurum Vest

	<p>Storage</p> <p>Room temperature 20 – 22 °C/68 – 71.6 °F.</p>	<p>It is not usually necessary to cool the liquid! Keep the working temperature between 20 °C and 22 °C/71.6 °F! The liquid must be cooled should temperatures increase!</p>	<p>Tips</p> <p>Temperature influences the amount of expansion. The liquid must not be exposed to frost! Careful with winter deliveries! Storage in thermo cupboard at 20 °C/68 °F.</p>												
<p>Wax surface tension reduction</p>	<p>Not necessary when mixed with a good vacuum mixing machine.</p>	<p>If a surface tension reducing agent is required use "Lubrofilm®". Allow to dry! If no surface tension reducing agent is used, make sure the vacuum mixing machine is in good condition!</p>	<p>Avoid alcohol residues, apply a thin coat and blow dry.</p>												
<p>Casting ring system</p>	<p>A metal ring is essential for speed casting!</p>	<p>For an accurately fitting castings in metal casting rings: For sizes 1 and 3 use one layer of dry liner. For sizes 6 and 9 use two layers of dry liner.</p>	<p>As an alternative silicon and rubber casting rings may be used. If so, make sure they are soft and elastic.</p>												
<p>Mixing ratio</p>	<p>150 g : 38 ml.</p>	<p>60 g : 15 ml; 300 g : 76 ml; 450 g : 114 ml</p>	<p>The mixing ratio must be followed strictly. Ensure that the mixing beakers are clean!</p>												
<p>Liquid concentration</p>	<p>Approx. 60 % concentration for standard gold alloys. Approx. 70 % concentration for precious metal bonding alloys. Do not use with non-precious alloys.</p>	<table border="1"> <thead> <tr> <th data-bbox="314 247 385 347">Low melting temp, precious metal alloys</th> <th data-bbox="314 347 385 549">Precious metal bonding alloys</th> <th data-bbox="314 549 385 772">Palladium based alloys</th> </tr> </thead> <tbody> <tr> <td data-bbox="385 247 463 347">Inlay/Onlay</td> <td data-bbox="385 347 463 549"></td> <td data-bbox="385 549 463 772"></td> </tr> <tr> <td data-bbox="463 247 542 347">Crowns and bridges</td> <td data-bbox="463 347 542 549">70 %</td> <td data-bbox="463 549 542 772">80 %</td> </tr> <tr> <td data-bbox="542 247 620 347">Secondary telescope in acrylic</td> <td data-bbox="542 347 620 549">90 %</td> <td data-bbox="542 549 620 772">100 %</td> </tr> </tbody> </table>	Low melting temp, precious metal alloys	Precious metal bonding alloys	Palladium based alloys	Inlay/Onlay			Crowns and bridges	70 %	80 %	Secondary telescope in acrylic	90 %	100 %	<p>For an exact mixing ratio please see table. Speed heating technique requires a little higher liquid concentration!</p>
Low melting temp, precious metal alloys	Precious metal bonding alloys	Palladium based alloys													
Inlay/Onlay															
Crowns and bridges	70 %	80 %													
Secondary telescope in acrylic	90 %	100 %													
<p>Mixing</p>	<p>Mix under vacuum for 60 s.</p>	<p>First liquid then powder into the mixing beaker. Mix thoroughly with a spatula. Ensure the mixing beaker is clean! Avoid plaster residues!</p>	<p>The mixing speed and shape of the mixing machine's blade have an effect on the mixture consistency and subsequently the amount expansion. Ideal mixing speed approx. 360 min⁻¹.</p>												

Working time	4 to 5 min.	Measured at a room temperature of 20 – 22 °C/68 – 71.6 °F.	Liquid must be cooled if temperatures exceed this.
Investing	Gentle vibration.	Set the vibrator to a low frequency and do not vibrate for too long.	Pour the investment into the casting ring up to the top of the crowns. Then carefully fill the crown itself and continue to fill the ring up to the upper edge of the asbestos liner. Only fill the ring up to the liner.
Investing under pressure	Max 5 bar recommended. Only with standard heating techniques.	Use this method only if the vacuum mixing machine is insufficient and if the pressure pot has an accurately fitting lid.	
Conventional heating technique	Setting time 40 min. Holding time 250 °C / 482 °F/1h.	Place the casting ring into the furnace with the opening at the bottom! Ideal heating up time 5 °C/41 °F/min. Holding time at end temperature 40 – 60 min according to the size of the ring. Take steps to prevent the ring from drying out if left over the weekend!	The base plate inside the furnace should be rilled. If the casting rings are to be left over the weekend, seal the investment material with wax!
Speed heating technique	Pull the asbestos liner above the casting ring periphery. After initial setting of the investment, scratch the surface with a knife. Place into the furnace after 20 min! Do not re-open the furnace door until 15 min have elapsed! – Danger of explosion!	Use only metal casting rings! Avoid contact of investment material with the metal casting ring. Place into the furnace after 20 min. Tolerance: + - 3 min (the larger the ring, the more important accurate timing is).	Danger when using the plaster trimmer as it does not run smoothly. Always cover acrylic patterns with wax!
Cooling down, devesting	Allow the ring to cool down to room temperature slowly.	Moisten the investment material and carefully devest the castings. Do not strike with a hammer. Sandblast using glass beads or aluminiumoxid 50 µm.	

Expansion control Recommended liquid concentration

Precious metal alloys	Crowns and bridges		Inlays				Telescopic crowns				Conical crowns (6°)				Implants									
	Wax		Covering one surface		more than one surface		Acrylic soft		Acrylic hard		Acrylic soft		Acrylic hard		Acrylic		Acrylic							
	Metal ring	without ring	Metal ring	without ring	Metal ring	without ring	Metal ring	without ring	Metal ring	without ring	Metal ring	without ring	Metal ring	without ring	Metal ring	without ring	Metal ring	without ring						
High gold casting alloy	65 %	70 %	60 %	70 %	80 %	80 %	65 %	70 %	85 %	95 %	80 %	85 %	75 %	75 %	60 %	70 %	55 %	65 %	50 %	50 %	70 %	80 %	65 %	
DentAurum Classic																								
Casting alloy reduced gold	65 %	70 %	75 %	85 %	70 %	70 %	85 %	70 %	90 %	100 %	85 %	90 %	80 %	95 %	65 %	75 %	60 %	70 %	60 %	70 %	55 %	75 %	85 %	70 %
High gold bonding alloy, high melting temp.	70 %	75 %	65 %																					
DentAurum Bio																								
High gold bonding alloy, low melting temp.	65 %	70 %	75 %	85 %	70 %	70 %	85 %	70 %	90 %	100 %	85 %	90 %	80 %	95 %	65 %	75 %	60 %	70 %	60 %	70 %	55 %	75 %	85 %	70 %
DentAurum LFC																								
Bonding alloy reduced gold high melting temp.	70 %	80 %	65 %																					
Palladium based alloy	80 %	90 %	75 %	100 %	85 %	100 %	100 %	100 %	100 %	100 %	100 %	100 %	100 %	100 %	80 %	90 %	75 %	85 %	70 %	85 %	75 %	90 %	100 %	85 %

This information is based upon in-house examinations. Different working conditions can require alterations to this data!

Various die sizes (anterior crowns, molars) can also result in deviations!

Acrylic, soft = adapted acrylic foil, light cure wax

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