# Chélasolv

Solución de irrigación canalar

## Composición

Edetato de sodio 15 %, cetrimida, excipientes.

## Propiedades

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15 % sodium edetate, cetrimide, excipients.

### Properties

Root canal preparation is a prerequisite for a successful root canal filling and most often requires the use of an irrigation solution combined with mechanical reaming.

It is indeed essential to remove the smear layer produced during root canal instrumentation and all debris which cover the dentinal tubule orifices.

Thanks to the chelating effect of sodium edetate, CHELASOLV dissolves the mineral content of root canal dentin and thus facilitates the removal of debris and smear layer.

The presence of cetrimide (quaternary ammonium compounds) reduces the surface tension of the product, allowing it to easily flow into the narrowest dentinal tubules.

### Indication

Root canal irrigation during root canal preparation.

### Contraindication

Allergy to any constituent.

### Instructions for use

Prepare the canal using traditional techniques. Use the syringe to irrigate the canal and allow the solution to stand for 1 to 3 minutes, at the end of preparation.

Rinse the root canal with distilled water and dry it thoroughly using paper points.

### Warnings and precautions for use

Check apical closure to avoid any risk of overfilling beyond the apex. The presence of hydrophobic agents (eugenol) inhibits the activity of CHELASOLV.

### Presentation

Carton containing 20 6-ml syringes.

*For professional use only in the practice of dentistry.*