# Alert™

Condensable Composite

English

INSTRUCTIONS FOR USE

Français

MODE D'EMPLOI

French

Deutsch

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German

Italiano

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Italian

Español

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Spanish

Čeština

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Czech

English

INSTRUCTIONS FOR USE:

Alert™ Condensable Composite has a "feel" similar to amalgam. However, it is neither amalgam nor conventional composite. Due to its unique fiber matrix, Alert Condensable Composite will seem to fray when cut with an instrument. This is normal. Condensing any frayed edges into the bulk of the material is easy and will enhance clinical performance. For easiest handling when creating anatomical features, use the instruments provided in the kit and minimize any "smearing" or "dragging" actions.

**Indications:** Tooth colored posterior restorative material. For usage in a similar situation that would be applicable to an amalgam restoration.

### Step 1 Treating the Cavity Preparation (Bond-1™ Primer/ Adhesive)

1. Develop cavity preparation and apply matrix band.
2. Apply 37% Phosphoric Acid Etching Gel\* to enamel and dentin surfaces for 15 seconds. Rinse thoroughly for 15 seconds with water. It is important not to splash etch into eyes, onto soft tissues or skin.
3. Dry with clean, oil-free, air without desiccating the dentin. Do not allow etched surface to become contaminated with saliva or blood during the bonding procedure. NOTE: Never use excessive force on syringe plunger to dispense etch which may lead to over dispensing of the material. If resistance is felt replace current tip with a new tip and check flow extraorally.
4. Apply two consecutive coats of Bond-1 Primer/Adhesive (Ensure a shiny surface. If shine disappears, repeat step 3).
5. Air dry gently for 10 seconds.
6. Light cure for 10 seconds.

### Step 2 Wetting the Surface (Flow-It™ ALC™ Flowable Composite)

1. Apply 0.5 to 1.0 mm Flow-It ALC Flowable Composite to walls and floor of preparation.
2. Light cure for 10 seconds.

### Step 3 Restoring the Tooth (Alert Condensable Composite)

1. Choose desired shade of Alert composite.
2. Place Alert Condensable Composite into the cavity preparation and condense.
3. Create occlusal anatomic form by pushing (not pulling) the instrument.

4. Light cure the restoration for 40 seconds. Remove matrix band. **NOTE: If preparation exceeds a depth of 5.0mm, it will be necessary to build the material incrementally. In this case, fill the cavity half way and light cure 20 seconds.**
5. Finish and polish restoration with conventional polishing materials and techniques.

#### **Step 4 Sealing the Surface (Protect-It™ Composite Surface Sealant)**

1. Use 37% Phosphoric Acid Etching Gel\* on the composite surfaces and extend out to cover approximately 2mm of enamel adjacent to restoration margin. Rinse thoroughly for 15 seconds with water. It is important not to splash etch into eyes, onto soft tissues or skin.
2. Dry with clean, oil-free, air without desiccating the dentin. Do not allow etched surface to become contaminated with saliva or blood during the bonding procedure. **NOTE: Never use excessive force on syringe plunger to dispense etch which may lead to over dispensing of the material. If resistance is felt replace current tip with a new tip and check flow extraorally.**
3. Apply Protect-It Composite Surface Sealant. Blow thin.
4. Light cure for 15 seconds.

Alert Condensable Composite is offered in syringe and single dose delivery.

**NOTE:** Single dose delivery, tips and applicator brushes are for single patient use only, to prevent cross-contamination between patients.

**NOTE:** Refer to outer package for expiration date. Do not use after expiration date.

**\*CAUTION:** 37% Phosphoric Acid Etching Gel contains phosphoric acid. Avoid contact with skin, eyes, and soft tissue. In case of contact with skin or eyes, flush immediately with water; get medical attention for eyes. Do not take internally.

**CAUTION:** Bond-1 Primer/ Adhesive contains acetone and is flammable and volatile. Keep away from sources of ignition and excessive heat.

**CAUTION:** Uncured methacrylate resin may cause contact dermatitis and damage the pulp. Avoid contact with skin, eyes and soft tissue. Wash thoroughly with water after contact.

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Store at ambient temperature.

## LIMITED WARRANTY —

Pentron's technical advice, whether verbal or in writing, is designed to assist dentists in using Pentron's product. The dentist assumes all risk and liability for damages arising out of the improper use of Pentron's product.

In the event of a defect in material or workmanship, Pentron will replace the defective product or part thereof, or reimburse the actual cost of the defective product. In order to take advantage of this limited warranty, the defective product must be returned to Pentron.

## PENTRON

Manufactured by Pentron Clinical

1717 West Collins Ave. Orange, CA 92867 USA

800-551-0283 / 714-516-7557 / [www.pentron.com](http://www.pentron.com)

**EC REP** European Union Representative

SpofaDental a.s.

Markova 238

CZ-506 01 Jičín

Phone: +420 493 583 111