1

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WHAT IS THE PELVIC  
FLOOR?

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The pelvic floor is made up of muscular-aponeurotic elements controlled by the nervous system. Its functions include holding the viscera at rest, cushioning abdominal pressure during exertion, contributing to continence and excretion mechanisms (urinary and faecal) and sexual function. The correct execution of all these functions will largely depend on the integrity and strength of the pelvic floor muscles. It is highly important to be aware and prevent pelvic floor dysfunction, which has such a major effect on a woman's quality of life.

The major muscular fibres in the pelvic floor are type-I fibres, which give muscle tone, and, to a lesser extent, type-II fibres, which help the sphincters to automatically close in urgent situations.

When the pelvic floor muscles lose their tone we can find ourselves facing a range of pathologies, such as pelvic organ prolapse (bladder, uterus and rectum), incontinence (urinary and faecal), and sexual

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dysfunction; for this reason, we must exercise the Pelvic Floor muscles and keep them well-toned.

The main risk factors causing damage or weakness to the perineal muscles are:

• **Pregnancy:**

The weight of the uterus, abdominal distension and the relaxing effect of hormones can all weaken the pelvic floor.

• **Childbirth:**

This can cause perineal neurological and muscle-aponeurotic lesions during the second stage of labour, or expulsion (pushing, Kristeller manoeuvre, forceps/ventouse, episiotomy and tearing, etc.).

• **Menopause:**

The lack of oestrogen leads to loss of tone and flaccidity in the perineal muscles.

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**• Sports:**

Impact sports such as aerobics, jogging, traditional sit-ups, tennis, basketball, etc. increase intra-abdominal pressure, and if practiced regularly, the pelvic floor muscles can become relaxed.

**• High-risk professions:**

Sports men and women, singers, wind instrument musicians...

**• Other risk factors:**

Obesity, constipation (straining), chronic cough, straining while urinating, pelvic surgery, hereditary factors, lack of knowledge or awareness of the perineal region... Apart from treatment with Boltex Inertial, it is very important to remove any risk factors that can be avoided.

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# 2

## **BOLTEX INERTIAL, TONIFICA Y FORTALECE TU SUELO PÉLVICO**

**BOLTEX INERTIAL,  
TONES AND STRENGTHENS  
YOUR PELVIC FLOOR**

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
Boltex Inertial is a medical device designed to tone and recondition women's pelvic floor muscles. It has an **exclusive and patented design** based on a set of hexagonal cells on the internal surface of the sphere, which give it extraordinary properties that help to provide **effective stimulation of the pelvic floor muscles**, regardless of its position after having been inserted into the vagina and of the user's anatomical and morphological characteristics.

The efficiency of an inertial device lies in its capacity to tone the PFM. The accelerometry graph shows the comparison between Boltex Inertial, in which the interior sphere knocks against the shell (Top graph) and an inertial device in which the interior sphere rolls instead of knocks (bottom graph, other models on the market). After knocking, a vibration is produced, which leads to a reflex toning action in the perineal muscles.

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
Gráfica superior (Boltex Inertial) / Top graph (Boltex Inertial) / Graphique supérieur (Boltex Inertial)

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
Likewise, the mechanical stimulation with Boltex, combined with voluntary contraction exercises (Kegel) make the toning more effective.

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The first part of the graph shows the progressive toning of the pelvic floor muscles with mechanical vibratory stimulation; in the second part, the mechanical stimulus is stopped; and in the third part, mechanical stimulation is used again, along with voluntary contraction exercises (Kegel).

Boltex Inertial also helps you to know the tone level of your pelvic muscles. Should the device drop during use and need to be moved into the vaginal passage again, the degree of hypotonia is very mild. If it is pushed out shortly after being placed, the loss of tone is severe. The more you use it, the better tone you will achieve, and so it will be easier to hold it in.

During menstruation and ovulation, the device can have a tendency to drop more easily, due to hormonal changes. Be constant and persistent with the treatment.

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
Boltex Inertial is a medical device that meets existing health regulations. The materials used have been tried and tested, and meet all biocompatibility requirements. Thanks to its features, it is ideal for assisting insertion and removal, ensuring it is easy and comfortable for the user and, in particular, guaranteeing contact between the surface of the sphere and the vaginal walls. It is totally hermetic with a smooth outer wall, preventing any build-up of dirt and bacteria and making it easier to clean. It is also comfortable and safe to use.

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### 1. INTERIOR BALL

Coated with silicone to prevent noise.

### 2. THE SPHEROID

It is made from PC/ABS, as this material is a good transmitter of vibration.

### 3. THE STRING

Makes it easier to extract from the vagina. It is made from TPE, which guarantees its medical quality and health and hygiene conditions.

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HOW TO EXERCISE  
WITH BOLTEX INERTIAL

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# How to Use

**Boltex**  
INERTIAL

Follow the instructions given by your health professional.

This device has been designed for use and application in women of any age and build.

As a general rule, follow the recommendations for use listed below:

- 1) **Empty your bladder**, if necessary.
- 2) **Clean Boltex** with warm water and mild soap.
- 3) **Lubrication:** Use a small amount of lubricating gel in case of vaginal dryness.
- 4) **Insertion:** Relax your pelvic muscles both for inserting and for removing the device.

Bending one or both legs (squatting) can help insertion. Place your index or middle finger on the flat area of the sphere (situated between the removal strings) and push the sphere **into the vaginal passage**. Part of the string will remain outside the body to enable the device to be removed.

5) **Exercising:** Exercise your muscles by walking for an hour or the time specified by your doctor. If you do Kegel exercises as well as using Boltex, the benefits will be greater.

As a general recommendation, start with 15 min/day exercises for 2 weeks. After 2 weeks, you can continue with 30-60 min. exercise sessions.

During treatment, combine fast voluntary contractions (close anus tightly) with slow contractions lasting 15-20 seconds.

Boltex Inertial can be used every day. Improvement in muscle tone is usually noticeable after 6-8 weeks of use. Most women recover their muscle tone after 3-6 months, depending on their degree of initial low muscle tone.

*Note: At first you may experience some muscular discomfort (cramps) from using Boltex Inertial. Should you experience any severe discomfort, see your doctor before*

*continuing use.*

6) **Removal:** Relax your muscles and use the string to pull out the device.

Women suffering from stress urinary incontinence are advised to remove the device in the knee-chest position and with their muscles relaxed, as they have urinary bladder descent and this position makes it easier for them to remove the device in the first few weeks.

Knee-chest position

7) **Cleaning:** Clean the device before and after use with warm water and pH-neutral soap. Store it in its original case or bag to keep it in good condition and prevent contact with dust or other materials.

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

# Contraindicaciones

## Counter-indications / Contre-indications

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### You should not use Boltex Inertial in the following cases:

- If you have a vaginal and/or urinary tract infection.
- During menstruation, with absorbent tampons.
- During sex when there is penetration.

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- If you are, or think you may be pregnant, please see your doctor before using the device.
- See your doctor if you feel any pain and/or inflammation after use.
- Do not use for longer than the time recommended by your health professional.
- Carefully follow health & hygiene guidelines before, during and after using the Boltex Inertial.
- Boltex is for personal use only. Do not share.

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# Pelvic System Care S.L.

**Boltex**  
INERTIAL

Pelvic System Care S.L. / Pelvic System Care S.L.

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It is a Research & Development centre for new medical devices concerning Pelvic Floor Pathology Prevention and Treatment. Science and research have enabled it to design the Boltex Inertial device as part of its work aimed at women.

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Boltex Inertial is patented:

"Inertial device for stimulating pelvic floor muscles with energy build-up"  
**(ES201100595)**

For further information, please visit our website:  
[www.pelvicsystemcare.com](http://www.pelvicsystemcare.com)

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